



US009265429B2

(12) **United States Patent**
St. Pierre et al.

(10) **Patent No.:** **US 9,265,429 B2**
(45) **Date of Patent:** **Feb. 23, 2016**

(54) **PHYSIOLOGICAL PARAMETER
MEASURING PLATFORM DEVICE
SUPPORTING MULTIPLE WORKFLOWS**

(75) Inventors: **Shawn C. St. Pierre**, Syracuse, NY (US); **Esther N. Ahn**, San Francisco, CA (US); **Eric J. Bailey**, Oakland, CA (US); **Michael D. Garrant**, Marcellus, NY (US); **Roger W. Graves**, San Francisco, CA (US); **Jennifer M. Kilian**, New York, NY (US); **Thomas A. Myers**, Syracuse, NY (US)

(73) Assignee: **Welch Allyn, Inc.**, Skaneateles Falls, NY (US)

(*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 1150 days.

(21) Appl. No.: **12/751,579**

(22) Filed: **Mar. 31, 2010**

(65) **Prior Publication Data**

US 2011/0071420 A1 Mar. 24, 2011

Related U.S. Application Data

(60) Provisional application No. 61/243,872, filed on Sep. 18, 2009.

(51) **Int. Cl.**
G08B 1/08 (2006.01)
A61B 5/0205 (2006.01)
(Continued)

(52) **U.S. Cl.**
CPC **A61B 5/02055** (2013.01); **A61B 5/7445** (2013.01); **A61B 5/7475** (2013.01);
(Continued)

(58) **Field of Classification Search**
CPC .. A61B 5/0205; A61B 5/024; A61B 5/02438;
A61B 5/14551; A61B 5/6843; A61B 5/6887;
A61B 5/742
See application file for complete search history.

(56) **References Cited**

U.S. PATENT DOCUMENTS

5,140,519 A 8/1992 Friesdorf et al.
D366,460 S 1/1996 Jorgenson
(Continued)

FOREIGN PATENT DOCUMENTS

CN 1356783 A 7/2002
CN 1649538 A 8/2005
(Continued)

OTHER PUBLICATIONS

mCare 300 Vital Signs Monitor, Spacelabs Medical, Inc., Copyright 2006, 4 pages.

(Continued)


Primary Examiner — Travis Hunnings

(74) *Attorney, Agent, or Firm* — Merchant & Gould P.C.

(57) **ABSTRACT**

A device obtains a series of measurements of a physiological parameter of a monitored patient when the device is operating within a monitoring workflow. The device displays a monitoring workflow home screen when the device is operating within the monitoring workflow. The monitoring workflow home screen contains a representation of the physiological parameter of the monitored patient. In addition, the device obtains a measurement of the physiological parameter of each patient in a series of patients when the device is operating within a non-monitoring workflow. The device displays a non-monitoring workflow home screen when the device is operating within the non-monitoring workflow. The non-monitoring workflow home screen contains a representation of the physiological parameter of a given patient in the series of patients. The monitoring workflow home screen is different than the non-monitoring workflow home screen.

17 Claims, 27 Drawing Sheets



(51)	Int. Cl.		8,038,593 B2	10/2011	Friedman
	<i>A61B 5/00</i>	(2006.01)	8,046,705 B2	10/2011	Hunleth et al.
	<i>G06F 19/00</i>	(2011.01)	8,055,514 B2	11/2011	Elsholz
	<i>G06Q 10/06</i>	(2012.01)	8,078,983 B2	12/2011	Davis
	<i>G06Q 50/22</i>	(2012.01)	D652,051 S	1/2012	Judy
(52)	U.S. Cl.		D656,153 S	3/2012	Imamura
	CPC	<i>G06F19/3406</i> (2013.01); <i>G06F 19/3487</i>	D656,157 S	3/2012	Khan
		(2013.01); <i>G06Q 10/06</i> (2013.01); <i>G06Q 50/22</i>	D656,946 S	4/2012	Judy
		(2013.01); <i>A61B 5/746</i> (2013.01)	D657,368 S	4/2012	Magee et al.
			D658,196 S	4/2012	Wood et al.
			D658,667 S	5/2012	Cho
			D664,971 S	8/2012	Lee et al.
			D664,984 S	8/2012	Lee et al.
			D666,625 S	9/2012	Gilmore et al.
			D667,837 S	9/2012	Magee
			D667,838 S	9/2012	Magee et al.
			D674,401 S	1/2013	Trumble et al.
			D675,218 S	1/2013	Arnold et al.
			D676,863 S	2/2013	Ho Kushner et al.
			D676,864 S	2/2013	Velasco et al.
			D688,685 S	8/2013	Rhee et al.
			D689,899 S	9/2013	Lee et al.
			D695,781 S	12/2013	Edwards et al.
			8,732,604 B2	5/2014	Okamoto et al.
			2002/0054141 A1	5/2002	Yen et al.
			2002/0078097 A1	6/2002	Chen et al.
			2002/0126137 A1	9/2002	Kaestner, Jr.
			2002/0177758 A1	11/2002	Schoenberg et al.
			2003/0060727 A1	3/2003	Kline
			2003/0135087 A1 *	7/2003	Hickle et al. 600/26
			2004/0002874 A1	1/2004	Shaffer et al.
			2004/0088199 A1	5/2004	Childress et al.
			2004/0111045 A1	6/2004	Sullivan et al.
			2005/0114374 A1	5/2005	Juszkiewicz et al.
			2005/0229110 A1	10/2005	Gegner et al.
			2005/0288571 A1	12/2005	Perkins et al.
			2006/0149597 A1	7/2006	Powell et al.
			2006/0200009 A1	9/2006	Wekell et al.
			2006/0228096 A1	10/2006	Hoshino et al.
			2006/0229557 A1	10/2006	Fathallah
			2006/0247948 A1	11/2006	Ellis et al.
			2007/0000494 A1	1/2007	Banner et al.
			2007/0067005 A1	3/2007	Schatz et al.
			2007/0150810 A1	6/2007	Katz et al.
			2007/0156456 A1	7/2007	McGillin et al.
			2007/0167173 A1	7/2007	Halcrow et al.
			2007/0167694 A1	7/2007	Causevic
			2007/0174079 A1	7/2007	Kraus
			2007/0185390 A1	8/2007	Perkins
			2007/0197878 A1	8/2007	Shklarski
			2007/0203744 A1	8/2007	Scholl
			2007/0215157 A1	9/2007	Straw
			2007/0276261 A1	11/2007	Banet
			2008/0012833 A1 *	1/2008	Beck et al. 345/173
			2008/0021379 A1	1/2008	Hickle
			2008/0055074 A1	3/2008	Gao
			2008/0058609 A1	3/2008	Garibaldi et al.
			2008/0058614 A1	3/2008	Banet
			2008/0072896 A1	3/2008	Setzer et al.
			2008/0077026 A1	3/2008	Banet
			2008/0086691 A1	4/2008	Hopermann et al.
			2008/0103375 A1	5/2008	Kiani
			2008/0115081 A1	5/2008	Sankaravivelu et al.
			2008/0155406 A1	6/2008	Naka
			2008/0208812 A1	8/2008	Quoc et al.
			2008/0229248 A1	9/2008	Fagans et al.
			2008/0249377 A1	10/2008	Molducci
			2008/0249801 A1	10/2008	Zaleski
			2008/0281168 A1	11/2008	Gibson
			2008/0281637 A1	11/2008	Matz
			2008/0318529 A1	12/2008	Harb
			2009/0005651 A1	1/2009	Ward et al.
			2009/0054743 A1	2/2009	Stewart
			2009/0054798 A1	2/2009	Varney et al.
			2009/0112630 A1	4/2009	Collins, Jr. et al.
			2009/0131761 A1	5/2009	Moroney
			2009/0132588 A1	5/2009	Mahesh et al.
			2009/0143652 A1	6/2009	Warburton
			2009/0149927 A1	6/2009	Kneuer

(56) **References Cited**

U.S. PATENT DOCUMENTS

5,701,894 A	12/1997	Cherry et al.	
D427,574 S	7/2000	Sawada et al.	
6,219,046 B1	4/2001	Thomas et al.	
D454,139 S	3/2002	Feldcamp	
6,364,834 B1	4/2002	Reuss et al.	
6,398,727 B1	6/2002	Bui et al.	
6,416,471 B1	7/2002	Kumar et al.	
6,469,717 B1	10/2002	Wineke et al.	
D465,226 S	11/2002	Friedman	
D468,322 S	1/2003	Walker et al.	
6,535,714 B2	3/2003	Melker et al.	
6,658,276 B2 *	12/2003	Kianl et al. 600/322	
6,675,031 B1	1/2004	Porges et al.	
6,707,476 B1	3/2004	Hochstedler	
6,783,573 B2	8/2004	Richardson	
6,828,910 B2	12/2004	VanRyzin et al.	
D510,582 S	10/2005	Hoang	
D523,440 S	6/2006	Hernandez et al.	
D525,982 S	8/2006	Suzuki	
D527,011 S	8/2006	Bixler	
7,124,366 B2	10/2006	Foreman et al.	
7,134,994 B2	11/2006	Alpert	
D545,829 S	7/2007	Fletcher	
7,336,187 B2	2/2008	Hubbard	
D575,296 S	8/2008	Fairfield	
7,409,399 B2	8/2008	Miyamoto	
D576,634 S	9/2008	Clark et al.	
D579,456 S	10/2008	Chen et al.	
7,433,827 B2	10/2008	Rosenfeld et al.	
D586,818 S	2/2009	Luck	
D590,413 S	4/2009	Bhat et al.	
D590,414 S	4/2009	Bhat et al.	
D592,156 S	5/2009	Drews et al.	
D592,675 S	5/2009	Bhat et al.	
7,565,616 B2	7/2009	Buchmann	
D598,923 S	8/2009	Chen et al.	
D598,929 S	8/2009	Bhat et al.	
D599,358 S	9/2009	Hoefnagels et al.	
D599,398 S	9/2009	Laidlaw et al.	
D603,416 S	11/2009	Poling et al.	
D607,463 S	1/2010	Krieter et al.	
D608,366 S	1/2010	Matas	
7,662,106 B2	2/2010	Daniels et al.	
D612,860 S	3/2010	Tarara et al.	
RE41,236 E	4/2010	Seely	
D614,634 S	4/2010	Nilsen	
7,765,479 B2	7/2010	Goodwin et al.	
7,774,060 B2	8/2010	Westenskow et al.	
7,782,302 B2	8/2010	Lee et al.	
7,787,946 B2	8/2010	Stahmann et al.	
7,810,497 B2	10/2010	Pittman et al.	
D632,699 S	2/2011	Judy et al.	
7,895,527 B2	2/2011	Zaleski et al.	
D635,150 S	3/2011	Sykes et al.	
D637,603 S	5/2011	Godgart	
7,937,134 B2	5/2011	Uber et al.	
D640,264 S	6/2011	Fujii et al.	
7,967,759 B2	6/2011	Couvillon	
D643,043 S	8/2011	Loken	
8,001,967 B2	8/2011	Wallace et al.	
8,020,558 B2	9/2011	Christopher	
D646,689 S	10/2011	Ulliot	
8,028,701 B2	10/2011	Al-Ali	

(56)

References Cited

U.S. PATENT DOCUMENTS

2009/0158415	A1	6/2009	Dillon
2009/0240116	A1	9/2009	Bluth
2009/0240193	A1	9/2009	Mensingher et al.
2009/0241958	A1	10/2009	Baker, Jr.
2009/0275805	A1	11/2009	Lane et al.
2009/0275810	A1	11/2009	Ayers
2009/0282340	A1	11/2009	Akaikie et al.
2009/0306482	A1	12/2009	Davis et al.
2009/0306488	A1	12/2009	Al-Ali
2009/0312648	A1	12/2009	Zhang et al.
2010/0050075	A1	2/2010	Thorson et al.
2010/0069004	A1	3/2010	Bioebaum
2010/0094096	A1	4/2010	Petruzzelli et al.
2010/0097380	A1	4/2010	Daniels et al.
2010/0218766	A1	9/2010	Milne
2010/0274098	A1	10/2010	Belford et al.
2010/0324380	A1	12/2010	Perkins et al.
2011/0010621	A1	1/2011	Wallaert et al.
2011/0015502	A1	1/2011	Peyser
2011/0071420	A1	3/2011	St. Pierre et al.
2011/0169644	A1	7/2011	Muhsin et al.
2011/0190600	A1	8/2011	McKenna et al.
2011/0205577	A1	8/2011	Mori et al.
2011/0209707	A1	9/2011	Terhark
2011/0227739	A1	9/2011	Gilham et al.
2011/0245630	A1	10/2011	St. Pierre et al.
2011/0246565	A1	10/2011	Irwin et al.
2011/0276338	A1	11/2011	Warner et al.
2011/0290250	A1	12/2011	Olson et al.
2011/0313301	A1	12/2011	Lane et al.
2012/0000470	A1	1/2012	Milne et al.
2012/0016251	A1	1/2012	Zhang et al.
2012/0048090	A1	3/2012	Etter et al.
2012/0095778	A1	4/2012	Gross et al.
2012/0096367	A1	4/2012	DelloStritto et al.
2012/0110444	A1	5/2012	Li et al.
2012/0117099	A1	5/2012	Gross
2012/0215075	A1	8/2012	Surace et al.
2012/0296183	A1	11/2012	Kinsley et al.
2013/0151285	A1	6/2013	McLaren et al.
2013/0265327	A1	10/2013	Vann et al.
2013/0267792	A1	10/2013	Petersen et al.
2013/0267793	A1	10/2013	Meador et al.
2013/0267861	A1	10/2013	Vassallo et al.
2013/0268283	A1	10/2013	Vann et al.
2013/0293373	A1	11/2013	Gegner et al.
2013/0311926	A1	11/2013	Keegan et al.
2014/0040429	A1	2/2014	Irwin et al.
2014/0098209	A1	4/2014	Neff
2014/0108041	A1	4/2014	Bechtel et al.
2014/0135588	A1	5/2014	Al-Ali et al.

FOREIGN PATENT DOCUMENTS

CN	101069638	A	11/2007
EP	0 707 824	A2	4/1996
EP	2 093 683	A2	8/2009
GB	2409951	A	7/2005
JP	10-91687	A	4/1998
JP	3682617	B2	8/2005
JP	2010-15193	A	1/2010
WO	0126021		4/2001
WO	0189362		5/2001
WO	2006/076498	A2	7/2006
WO	WO 2006/076498	*	7/2006
WO	2010/102069	A2	9/2010
WO	2011/001302	A1	1/2011

OTHER PUBLICATIONS

Ultraview DM3—Dual Mode Vital Signs Monitor, Spacelabs Healthcare, Copyright 2010, 4 pages.

Colin Prodigy Press-Mate Prodigy II® Portable Vital Signs Monitors, DRE, Copyright 2009, accessed at: http://www.dremed.com/catalog/product_info.php/products_id/1181; 5 pages.

Fingertip Pulse Oximeter SPO₂ Monitor Oxigen Oximeter; Copyright 1995-2010, accessed at: [http://74.125.45.132/search?q=cache:ANfHne9je7gJ:cgi.ebay.com.sg/ws/eBayISAPI.dll%3FViewItem%26item%3D260412542031+\"Fingertip+Pulse+Oximeter+SPO2+Monitor+Oxigen+Oximeter\"&cd=2&hl=en&ct=clnk&gl=us](http://74.125.45.132/search?q=cache:ANfHne9je7gJ:cgi.ebay.com.sg/ws/eBayISAPI.dll%3FViewItem%26item%3D260412542031+\); 6 pages.

Handheld Pulse Oximeter PM-60A, Contec Medical System Co. Ltd., Jul. 31, 2009; accessed at: <http://www.tradeindia.com/selloffer/185893/Handheld-Pulse-Oximeter-PM-60A.html>; 5 pages.

Multiple Vital Signs from One Non-invasive Sensor, Starr™ Life Sciences Corp., Copyright 2009, accessed at: <http://www.starlifesciences.com/mouseox.html>; 1 page.

NTIA Handheld Pulse Oximeter, Frontline Systems; Dec. 31, 2009, accessed at: <http://www.tradeindia.com/fp385774/NTIA-Handheld-Pulse-Oximeter.html>; 3 pages.

Portable Patient Monitors, Welch Allyn, Copyright 2005-2006 Med-Direct, Copyright 2005-2006, accessed at: <http://www.meddirect.co.nz/Product.aspx?ProductID=3476>; 1 page.

International Search Report and Written Opinion in PCT/US2010/048450 dated Apr. 12, 2011, 9 pages.

Nag et al., Wireless E-jacket ofr Multiparameter Biophysical Monitoring and Telemedicine Applications, Medial Devices and Biosensors, 2006, 3rd IEEE/EMBS International Summer School, Sep. 4-6, 2006, pp. 40-44.

O'Donoghue et al., Design and Implementation of a Framework for Monitoring Patients in Hospitals Using Wireless Sensors in Ad Hoc configuration, Engineering in Medicine and Biology Society, EMBS '06, 28th Annual International Conference of the IEEE, vol., No., Aug. 30, 2006-Sep. 3, 2006, pp. 6449-6452.

Adams, A.P., Breathing System Disconnections, British Journal of Anaesthesia 1994; 73: pp. 46-54.

Capnography: An Objective Tool for Assessing Respiratory Status, Physio-Control, Inc., 2008, pp. 1-4.

Al-Qutayri et al.; Framework for Secure Wireless Health monitoring and remote Access System, Inderscience Enterprises Ltd. copyright 2010, 19 pages.

International Search Report and Written Opinion in PCT/US2013/031342 mailed Jun. 24, 2013, 12 pages.

International Search Report and Written Opinion in PCT/US2013/031458 mailed Jun. 28, 2013, 10 pages.

International Search Report and Written Opinion in PCT/US2013/031486 mailed Jun. 28, 2013, 10 pages.

International Search Report and Written Opinion in PCT/US2013/031582 mailed Jun. 28, 2013, 11 pages.

Kozlovsky et al.; Network and Service Management and Diagnostics Solution of a Remote Patient Monitoring System, IEEE copyright 2011, 4 pages.

Lamberti et al.; Ubiquitous Real-Time Monitoring of Critical-Care Patients in Intensive Care Units, IEEE copyright 2003, 4 pages.

Panorama™ Central Station, Surgical Product Guide copyright 2011, 2 pages.

Station—Dictionary.com, [online], retrieved on Sep. 2, 2014, Retrived from, <URL: <http://dictionary.reference.com/browse/station>>, 4 pages.

U.S. Appl. No. 29/417,592, filed Apr. 5, 2012.

U.S. Appl. No. 29/417,611, filed Apr. 5, 2012.

Zhao et al.; a Portable, Low-Cost, Batter-Powered Wireless Monitoring System for Obtaining Varying Physiologic Parameters from Multiple Subjects, IEEE copyright 2006, 4 pages.

Extended European Search Report for PCT/US2010/048450, Nov. 18, 2014, 7 pages.

* cited by examiner




FIG. 1




FIG. 2A




FIG. 2B




FIG. 3A




FIG. 3B




FIG. 3C

314C →







FIG. 3D

380

 Alarm limit exceeded. Pulse rate HIGH.

NIBP ••• 120/80 SYS/DIA mmHg (MAP 93) INTERVALS: Transfusion	 START 0:14:39	SYS 160 90 DIA 90 50	PULSE RATE 122 ♥/MIN SOURCE = NIBP
SpO2 •• 97% (PI 19.3) MODE: Fast	 25	100 95	TEMPERATURE 101.5° °F (38.6° C)
PATIENT Barker, David A. 206A A 05/14/1945	Adult	HEIGHT 72 in WEIGHT 197 lb PAIN 3 RR 15 bpm	Save

Home **Patients** **Alarms** **Review** **Settings**

FIG. 3E

380

(:89) Alarm limit exceeded. Pulse rate high.

NIBP ••• 120/80 SYS/DIA mmHg (MAP 93) INTERVALS: Transfusion	START 0:14:39	SYS 160 90 DIA 90 50	PULSE RATE 122 ✓/MIN SOURCE = NIBP	120 50
SpO2 •• 97% (PI 19.3) MODE: Fast	25 	100 95	TEMPERATURE 101.5° °F (38.6° C)	102.1 88.3
PATIENT Barker, David A. 206A A 05/14/1945	Adult	HEIGHT 72 in WEIGHT 197 lb PAIN 3 RR 15 bpm		Save
Home	Patients	Alarms	Review	Settings

FIG. 3F




FIG. 4




FIG. 5




FIG. 6A




FIG. 6B




FIG. 6C




FIG. 6D




FIG. 6E




FIG. 6F




FIG. 6G




FIG. 7A




FIG. 7B




FIG. 7C




FIG. 7D




FIG. 7E




FIG. 8

800





FIG. 9






FIG. 11

**PHYSIOLOGICAL PARAMETER
MEASURING PLATFORM DEVICE
SUPPORTING MULTIPLE WORKFLOWS**

RELATED APPLICATION

The present application claims priority to U.S. Provisional Patent Application Ser. No. 61/243,872, filed Sep. 18, 2009, the entirety of which is hereby incorporated by reference.

BACKGROUND

Health care practitioners, such as nurses and physicians, use various types of health-care equipment to assist with the task of providing health care to a patient, also referred to herein as a health-care recipient. Some health-care equipment, referred to as single function equipment, is designed to perform a particular function, such as temperature measurement. Some health-care equipment, referred to as multi-function equipment, is designed to implement the performance of more than one function, such as temperature measurement and blood pressure measurement. Such multi-function equipment may impose excess bulk and/or weight upon a user if such multi-function equipment is used for only one function or a subset of the functions implemented by the multi-function equipment.

SUMMARY

In one aspect, a physiological measuring platform (PMP) device comprises a central processing unit (CPU) that is configured to control operation of the PMP device. The PMP device also comprises a display screen. In addition, the PMP device comprises one or more computer readable data storage media storing software instructions that, when executed by the CPU, cause the PMP device to obtain a series of measurements of a physiological parameter of a monitored patient when the PMP device is operating within a monitoring workflow. The software instructions, when executed by the CPU, further cause the PMP device to display, on the display screen, a monitoring workflow home screen when the PMP device is operating within the monitoring workflow. The monitoring workflow home screen contains a first representation of the physiological parameter of the monitored patient. The first representation is based on a measurement in the series of measurements. In addition, the software instructions, when executed by the CPU, further cause the PMP device to obtain a measurement of the physiological parameter of each patient in a series of patients when the PMP device is operating within a non-monitoring workflow. In addition, the software instructions, when executed by the CPU, further cause the PMP device to display, on the display screen, a non-monitoring workflow home screen when the PMP device is operating within the non-monitoring workflow. The non-monitoring workflow home screen contains a second representation of the physiological parameter of a given patient in the series of patients. The second representation is based on the measurement of the physiological parameter of the given patient in the series of patients. The monitoring workflow home screen is different than the non-monitoring workflow home screen.

Another aspect is a method for performing physiological parameter measurements for a plurality of patients. The method comprises, displaying, by a physiological measurement platform (PMP) device, a settings screen on a display screen, the settings screen enabling a user to select a workflow within which the PMP device is to operate. The method

further comprises, when the PMP device is operating within a monitoring workflow: obtaining, by the PMP device, a series of measurements of a physiological parameter of a monitored patient; and displaying, by the PMP device, a monitoring workflow home screen on the display screen. The monitoring workflow home screen contains a first representation of the physiological parameter of the given patient, the first representation based on a measurement in the series of measurements. The method also comprises when the PMP device is operating within a non-monitoring workflow: obtaining, by the PMP device, a measurement of the physiological parameter of each patient in a series of patients; and displaying, by the PMP device, a non-monitoring workflow home screen on the display screen. The non-monitoring workflow home screen contains a second representation of the physiological parameter of a given patient in the series of patients, the second representation based on the measurement of the physiological parameter of the given patient in the series of patients.

Yet another aspect is a method comprising obtaining, by a physiological measuring platform (PMP) device, a measurement of a physiological parameter of a given patient. The method also comprises displaying an alarms screen on a display screen of the PMP device. The alarms screen enables a user of the PMP device to configure global alarm settings and parameter-specific alarm settings. The alarms screen includes a compressed parameter reporting area. The compressed parameter reporting area contains a compressed parameter frame containing a representation of the physiological parameter of the given patient.

Yet another aspect is a method comprising obtaining, by a physiological measuring platform (PMP) device, a measurement of a physiological parameter of a patient. The method also comprises displaying, on a display screen of the PMP device, a parameter reporting area containing a representation of the physiological parameter of the patient. In addition, the method comprises displaying, on the display screen of the PMP device, an alarm message describing a first alarm when the first alarm is active. The first alarm is active when the measurement of the physiological parameter is outside an alarm range for the physiological parameter. The method also comprises modifying, on the display screen, the alarm message to describe a second alarm when the second alarm is active and when a user selects the alarm message.

Yet another aspect is a method comprising displaying, on a display screen of a physiological measuring platform (PMP) device a settings screen. The settings screen containing a set of interval program selection controls. The method also comprises receiving, by the PMP device, a selection by a user of one of the interval program selection controls. Furthermore, the method comprises obtaining, by the PMP device, measurements of a physiological parameter of a patient. In addition, the method comprises recording, by the PMP device, measurements of the physiological parameter at intervals specified by an interval program for a duration specified by the interval program. The interval program corresponds to the selected one of the interval program selection controls.

Yet another aspect is a computer-readable storage medium comprising software instructions that, when executed, cause a physiological measurement platform (PMP) device to display, on a display screen, a settings screen that enables a user to select a workflow within which the PMP device is to operate. When the user selects a monitoring workflow from the settings screen, the software instructions, when executed, cause the PMP device to:

obtain a series of measurements of a physiological parameter of a monitored patient;

display a monitoring workflow home screen on the display screen, the monitoring workflow home screen containing a first representation of the physiological parameter, the first representation of the physiological parameter based on at least one measurement in the series of measurements; and

display, on the display screen, an alarms screen that enables the user to configure global alarm settings and parameter-specific alarm settings. The alarms screen includes a compressed parameter reporting area. The compressed parameter reporting area contains a compressed parameter frame containing a second representation of the physiological parameter of the monitored patient.

When the user selects a spot check workflow from the settings screen, the software instructions, when executed, cause the PMP device to:

obtain a measurement of the physiological parameter of each patient in a series of previously identified patients; and

display a spot check workflow home screen on the display screen. The spot check workflow home screen contains a third representation of the physiological parameter of a given patient in the series of previously identified patients. The third representation of the physiological parameter is based on the measurement of the physiological parameter of the given patient in the series of previously identified patients.

When the user selects a triage workflow from the settings screen, the software instructions, when executed, cause the PMP device to:

obtain a measurement of the physiological parameter of each patient in a series of unidentified patients; and display a triage workflow home screen on the display screen. The triage workflow home screen contains a fourth representation of the physiological parameter of a given patient in the series of unidentified patients, the fourth representation of the physiological parameter based on the measurement of the given patient in the series of unidentified patients.

The monitoring workflow home screen, the spot check workflow home screen, and the triage workflow home screen are each different from one another. The software instructions, when executed, also cause the PMP device to display, on the display screen, an alarm message describing a first alarm when the first alarm is active. The first alarm is active when a given measurement of the physiological parameter is outside an alarm range for the physiological parameter. The software instructions, when executed, also cause the PMP device to modify, on the display screen, the alarm message to describe a second alarm when the second alarm is active and when the user selects the alarm message. The software instructions, when executed, also cause the PMP device to emit an alarm sound when the first alarm is active. The software instructions, when executed, also cause the PMP device to suspend emitting the alarm sound for a given time period in response to user input. The software instructions, when executed, also cause the PMP device to extend the given time period when the PMP device receives additional user input before the given time period expires. The software instructions, when executed, also cause the PMP device to resume emitting the alarm sound when the given time period expires.

BRIEF DESCRIPTION OF THE DRAWINGS

The objects and features of the invention can be better understood with reference to the claims and drawings described below. The drawings are not necessarily to scale,

and the emphasis is instead generally being placed upon illustrating the principles of the invention. Within the drawings, like reference numbers are used to indicate like parts throughout the various views. Differences between like parts may cause those like parts to be each indicated by different reference numbers. Unlike parts are indicated by different reference numbers.

FIG. 1 is a block diagram illustrating an example system for collecting measurements of physiological parameters of patients.

FIG. 2A illustrates a view of an example physiological parameter measuring platform (PMP) device.

FIG. 2B illustrates an example user interface displayed on a user interface display of the PMP device of FIG. 2A.

FIG. 3A illustrates an example monitoring workflow home screen.

FIG. 3B illustrates an example spot check workflow home screen.

FIG. 3C illustrates an example triage workflow home screen.

FIG. 3D illustrates an example alternative view of a patient's SpO2 level.

FIG. 3E illustrates the example monitoring workflow home screen when an alarm is active.

FIG. 3F illustrates the example monitoring workflow home screen when an alarm is active and when an alarm sound for the alarm has been temporarily silenced.

FIG. 4 illustrates an example patient selection screen.

FIG. 5 illustrates an example review screen.

FIG. 6A illustrates an example intervals pane of a settings screen.

FIG. 6B illustrates an example device pane of the settings screen.

FIG. 6C illustrates an example clinician pane of the settings screen.

FIGS. 6D-6F illustrate example profile panes of the settings screen.

FIG. 6G illustrates an example advanced pane of the settings screen.

FIG. 7A illustrates an example global pane of an alarms screen.

FIG. 7B illustrates an example temperature pane of the alarms screen.

FIG. 7C illustrates an example NIBP pane of the alarms screen.

FIG. 7D illustrates an example SpO2 pane of the alarms screen.

FIG. 7E illustrates an example pulse rate pane of the alarms screen.

FIG. 8 is a flowchart illustrating an example operation performed by the PMP device.

FIG. 9 is a flowchart representing a continuation of the operation in which the workflow selection input indicates the triage workflow.

FIG. 10 is a flowchart illustrating an operation performed by the PMP device when the PMP device displays a workflow home screen.

FIG. 11 illustrates example physical components of the PMP device.

DETAILED DESCRIPTION

Embodiments of the present invention are directed to a physiological parameter measuring platform (PMP) device having a user interface configured to operate within and transition between each of a monitoring workflow and a non-monitoring workflow. In the following detailed description,

references are made to the accompanying drawings that form a part hereof, and in which are shown by way of illustrations specific embodiments or examples. These embodiments may be combined, other embodiments may be utilized, and structural changes may be made without departing from the spirit or scope of the present invention. The following detailed description is therefore not to be taken in a limiting sense and the scope of the present invention is defined by the appended claims and their equivalents.

Referring now to the drawings, in which like numerals refer to like elements through the several figures, aspects of the present invention and an exemplary computing operating environment will be described.

FIG. 1 is a block diagram illustrating an example system 100 for collecting measurements of physiological parameters of patients. As illustrated in the example of FIG. 1, the system 100 comprises an Electronic Medical Records (EMR) system 102, an interface system 104, a set of client devices 106A-106N (collectively, "client devices 106"), and a network 108.

The network 108 is an electronic communication network that facilitates communication between the client devices 106 and the between the client devices 106 and the interface system 104. An electronic communication network is a set of computing devices and links between the computing devices. The computing devices in the network use the links to enable communication among the computing devices in the network. The network 108 can include routers, switches, mobile access points, bridges, hubs, intrusion detection devices, storage devices, standalone server devices, blade server devices, sensors, desktop computers, firewall devices, laptop computers, handheld computers, mobile telephones, and other types of computing devices. In various embodiments, the network 108 includes various types of links. For example, the network 108 can include wired and/or wireless links. Furthermore, in various embodiments, the network 108 is implemented at various scales. For example, the network 108 can be implemented as one or more local area networks (LANs), metropolitan area networks, subnets, wide area networks (such as the Internet), or can be implemented at another scale.

The EMR system 102 is a computing system that allows storage, retrieval, and manipulation of electronic medical records. As used herein, a computing system is a system of one or more computing devices. A computing device is a physical, tangible device that processes data. Example types of computing devices include personal computers, standalone server computers, blade server computers, mainframe computers, handheld computers, smart phones, special purpose computing devices, and other types of devices that process data.

Each client device in the set of client devices 106 is a computing device. The client devices 106 can provide various types of functionality. For example, the set of client devices 106 can include one or more PMP devices (such as the PMP device 200). In addition, the set of client devices 106 can include one or more wall-mounted devices. Such wall-mounted devices can have similar functionality to the PMP device 200 but are stationary instead of portable. In addition, the set of client devices 106 can include one or more monitor devices. Such monitor devices can display representations of physiological parameters, but do not directly obtain measurements of the physiological parameters from patients. A monitor device could, for example, be used by a clinician to monitor the physiological parameters of multiple patients at one time.

The client devices 106 can communicate with each other through the network 108. In various embodiments, the client devices 106 can communicate various types of data with each

other through the network 108. For example, in embodiments where the set of client devices 106 includes a set of PMP devices and a monitor device, each of the PMP devices can send data representing measurements of physiological parameters of patients to the monitor device. In this way, the monitor device can display representations of physiological parameters to a clinician.

The interface system 104 is a computing system that acts as an interface between the EMR system 102 and the client devices 106. In some embodiments, the interface system 104 is a Connex system. Different EMR systems have different software interfaces. For example, the EMR system used by two different hospitals can have two different software interfaces. The interface system 104 provides a single software interface to each of the client devices 106. The client devices 106 send requests to software interface provided by the interface system 104. When the interface system 104 receives a request from one of the client devices 106, the interface system 104 translates the request into a request that works with the software interface provided by the EMR system 102. The interface system 104 then provides the translated request to the software interface provided by the EMR system 102. When the interface system 104 receives a response from the EMR system 102, the interface system 104 translates the response from a format of the EMR system 102 to a system understood by the client devices 106. The interface system 104 then forwards the translated response to an appropriate one of the client devices 106.

The client devices 106 can send various types of data to the interface system 104 for storage in the EMR system 102 and can receive various types of data from the EMR system 102 through the interface system 104. For example, in some embodiments, the client devices 106 can send measurements of physiological parameters to the interface system 104 for storage in the EMR system 102. In another example, a monitor device can retrieve past measurements of physiological parameters of patients from the EMR system 102 through the interface system 104.

FIG. 2A illustrates a view of a PMP device 200. The PMP device 200 is classified and referred to as a portable monitor platform device. The PMP device 200 includes multiple health care equipment (HCE) modules. Each of the HCE modules is configured to measure one or more physiological parameters of a health-care recipient, also referred to herein as a patient.

A temperature measurement module 212 is accessible from the front side of the PMP device 200. A SpO2 module 214 and a non-invasive blood pressure (NIBP) module 216 are accessible from a left hand side of the PMP device 200. An upper handle portion 220 enables the PMP device 200 to be carried by hand.

A front side of the PMP device 200 includes a display screen 218 and an outer surface of the temperature measurement module 212. The temperature measurement module 212 is designed to measure the body temperature of a patient. As used in this document, a "module" is a combination of a physical module structure which typically resides within the PMP device 200 and optional peripheral components (not shown) that typically attach to and reside outside of the PMP device 200.

The temperature measurement module 212 includes a front panel 212a. The front panel 212a has an outer surface that is accessible from the front side of the PMP device 200. The front panel 212a provides access to a wall (not shown) storing a removable probe (not shown), also referred to as a temperature probe, that is attached to a probe handle 212b. The probe and its attached probe handle 212b are tethered to the tem-

perature measurement module 212 via an insulated conductor 212c. The probe is designed to make physical contact with a patient in order to sense a body temperature of the patient.

A left hand side of the PMP device 200 includes an outer surface of the SpO2 module 214 and an outer surface of the NIBP module 216. The SpO2 module 214 is a HCE module designed to measure oxygen content within the blood of a patient. The NIBP module 216 is a HCE module designed to measure blood pressure of a patient.

As shown, the SpO2 module 214 includes a front panel 214a. The front panel 214a includes an outer surface that is accessible from the left side of the PMP device 200. The front panel 214a includes a connector 214b that enables a connection between one or more peripheral SpO2 components (not shown) and a portion of the SpO2 module 214 residing inside the PMP device 200. The peripheral SpO2 components reside external to the PMP device 200. The peripheral SpO2 components are configured to interoperate with the SpO2 module 214 when connected to the SpO2 module 214 via the connector 214b. In some embodiments, the peripheral SpO2 components include a clip that attaches to an appendage of a patient, such as a finger. The clip is designed to detect and measure a pulse and an oxygen content of blood flowing within the patient.

As shown, the NIBP module 216 includes a front panel 216a having an outer surface that is accessible from the left side of the PMP device 200. The front panel 216a includes a connector 216b that enables a connection between one or more peripheral NIBP components (not shown) and a portion of the NIBP module 216 residing inside the PMP device 200. The peripheral NIBP components reside external to the PMP device 200. The peripheral NIBP components are configured to interoperate with the NIBP module 216 when connected to the NIBP module 216 via the connector 216b. In some embodiments, the peripheral NIBP components include an inflatable cuff that attaches to an appendage of a patient, such as an upper arm of the patient. The inflatable cuff is designed to measure the systolic and diastolic blood pressure of the patient, the mean arterial pressure (MAP) of the patient, and the pulse rate of blood flowing within the patient.

The PMP device 200 is able to operate within one or more workflows. A workflow is a series of one or more tasks that a user of the PMP device 200 performs. When the PMP device 200 operates within a workflow, the PMP device 200 provides functionality suitable for assisting the user in performing the workflow. When the PMP device 200 operates within different workflows, the PMP device 200 provides different functionality.

When the PMP device 200 is manufactured, the PMP device 200 is configured to be able to operate within one or more workflows. After the PMP device 200 is manufactured, the PMP device 200 can be reconfigured to operate within one or more additional workflows. In this way, a user can adapt the PMP device 200 for use in different workflows as needed.

In various embodiments, the PMP device 200 operates within various workflows. For example, in some embodiments, the PMP device 200 can operate within a monitoring workflow or a non-monitoring workflow. Example types of non-monitoring workflows include, but are not limited to, a spot check workflow and a triage workflow.

In example embodiments, the names for the workflows can be defined by the user. For example, the user can rename a "triage workflow" as "ED 3 North" or any other nomenclature as desired to provide more context to the user.

When the PMP device 200 is operating within the monitoring workflow, the PMP device 200 obtains a series of measurements of one or more physiological parameters of a

single monitored patient over a period of time. In addition, the PMP device 200 displays, on the display screen 218, a monitoring workflow home screen. The monitoring workflow home screen contains a representation of a physiological parameter of the monitored patient. The representation is based on at least one measurement in the series of measurements. A representation of a physiological parameter is a visible image conveying information about the physiological parameter. For example, when the PMP device 200 is operating within the monitoring workflow, the PMP device 200 can obtain a temperature measurement of a single patient once every ten minutes for six hours. In this example, the PMP device 200 displays a monitoring workflow home screen that contains a representation of the patient's body temperature based on a most recent one of the temperature measurements. In this way, a user of the PMP device 200 can monitor the status of the patient.

When the PMP device 200 is operating within a non-monitoring workflow, the PMP device 200 obtains a measurement of one or more physiological parameters from each patient in a series of patients. In addition, the PMP device 200 displays a non-monitoring workflow home screen on the display screen 218. The non-monitoring workflow home screen contains a representation of the physiological parameter of a given patient in the series of patients. The representation is based on the measurement of the physiological parameter of the given patient. In one example, when the PMP device 200 is operating within a spot check workflow, the PMP device 200 obtains blood pressure measurements from a series of previously-identified patients. In this other example, the PMP device 200 displays a spot check workflow home screen containing a blood pressure measurement of a given patient in the series of previously-identified patients. In this way, a user of the PMP device 200 can perform spot checks on the blood pressures of patients who have already been admitted to a hospital. As used in this document, a patient is a previously identified patient when the PMP device 200 stores information regarding the identity of the patient. In another example, when the PMP device 200 is operating within a triage workflow, the PMP device 200 can obtain a single blood pressure measurement from each patient in a series of unidentified patients as the patients arrive at a hospital. In this example, the PMP device 200 displays a triage workflow home screen containing a representation of the patients' blood pressure based on the single blood pressure measurements of the patients. In this way, a user of the PMP device 200 can perform triage on the series of unidentified patients as they arrive. As used in this document, a patient is an unidentified patient when the PMP device 200 does not store information regarding the identity of the patient.

The monitoring workflow home screen is different than the non-monitoring workflow home screen. In various embodiments, the monitoring workflow home screen is different than the non-monitoring workflow home screen in various ways. For example, in some embodiments, the monitoring workflow home screen includes at least one user-selectable control that is not included in the non-monitoring workflow home screen. In other embodiments, a representation of a physiological parameter in the monitoring workflow home screen has a different size than a representation of the same physiological parameter in the non-monitoring workflow home screen.

FIG. 2B illustrates an example user interface displayed on the display screen 218 of FIG. 2A. The PMP device 200 outputs and displays user interfaces discussed in this document on the display screen 218.

FIG. 3A illustrates an example monitoring workflow home screen 300. The PMP device 200 displays the monitoring workflow home screen 300 while the PMP device 200 is operating within a monitoring workflow. The monitoring workflow is designed for obtaining a series of physiological measurements associated with an identified patient over a period of time.

The PMP device 200 is functionally connected to one or more sensors that enable monitoring of at least one physiological parameter that is associated with a patient. Typically, each sensor is physically attached to the patient while the PMP device 200 is operating within the monitoring workflow. These sensors include a temperature probe, a SpO2 clip, and a NIBP blood pressure cuff that are each attachable to the PMP device 200 as described above.

As shown in the example of FIG. 3A, the monitoring workflow home screen 300 includes a device status area 312, a navigation area 318, and a content area 320. The content area 320 is divided into a parameter reporting area 314 and a patient attribute area 316.

The device status area 312 contains data regarding a status of the PMP device 200. In the example of FIG. 3A, the device status area 312 includes text that identifies a clinician ("Patricia Jones") and a health care facility location ("West 4A"). A current time of day value ("03:00") is located towards the center of the device status area 312. A date value ("12/29/2009") is located to the right side of the time of day value. A remaining time of a battery ("1:10") value is located at the right side of the device status area 312.

The navigation area 318 includes a home tab 319a, a patients tab 319b, an alarms tab 319c, a review tab 319d, and a settings tab 319e. Collectively, the home tab 319a, the patients tab 319b, the alarms tab 319c, the review tab 319d, and the settings tab 319e are referred to herein collectively as the screen tabs 319. Selection of screen tabs 319b-319e causes substitution of the monitoring workflow home screen 300 with another screen associated with the screen tabs 319b-319e. For example, the PMP device 200 displays a patient screen when a user selects the patients tab 319b. FIG. 4A, discussed elsewhere in this document, illustrates an example patient selection screen. Furthermore, the PMP device 200 displays an alarms screen when a user selects the alarms tab 319c. FIGS. 7A-7E, discussed elsewhere in this document, illustrate an example alarms screen. Furthermore, the PMP device 200 displays a review screen when a user selects the review tab 319d. FIG. 5, discussed elsewhere in this document, illustrates an example review screen. Furthermore, the PMP device 200 displays a settings screen when a user selects the settings tab 319e. FIGS. 6A-6G, discussed elsewhere in this document, illustrate an example settings screen. When the PMP device 200 displays a screen other than the monitoring workflow home screen 300 and a user selects the home tab 319a, the PMP device 200 displays the monitoring workflow home screen 300.

The parameter reporting area 314 includes one or more parameter reporting frames. Each of the parameter reporting frames contains a representation of a different physiological parameter of a patient. The representations are based on one or more measurements of the physiological parameters of a monitored patient. In addition, each of the parameter reporting frames contains an alarm reporting area. The alarm reporting areas specify upper alarm limits and lower alarm limits for the physiological parameters. The upper alarm limits and the lower alarm limits define the alarm ranges for the physiological parameters. Alarms associated with the physi-

ological parameters are active when measurements of the physiological parameters are outside the alarm range for the physiological parameters.

In the example of FIG. 3A, the parameter reporting area 314 contains a NIBP frame 314a, a pulse rate frame 314b, a SpO2 frame 314c, and a temperature frame 314d. The NIBP frame 314a is located within an upper left portion of the parameter reporting area 314. The pulse rate frame 314b is located within an upper right portion of the parameter reporting area 314. The SpO2 frame 314c is located within a lower left portion of the parameter reporting area 314. The temperature frame 314d is located within a lower right portion of the parameter reporting area 314.

The NIBP frame 314a contains a representation of the blood pressure of the patient. The representation of the blood pressure of the patient is based on one or more measurements of the blood pressure of the patient. In various embodiments, the NIBP frame 314a contains various representations of the blood pressure of the patient. In the example of FIG. 3A, the NIBP frame 314a includes enlarged numerical text that represents a systolic blood pressure value ("120") and a diastolic blood pressure value ("80"), separated from each other via a slash '/' text character. The systolic blood pressure value is located at the left side of the NIBP frame 314a and the diastolic blood pressure is located to the right side of the systolic blood pressure value.

An NIBP alarm status area 322c is located at the right side of the NIBP frame 314a. The NIBP alarm status area 322c specifies an upper alarm limit and a lower alarm limit for the patient's systolic blood pressure and an upper alarm limit and a lower alarm limit for the patient's diastolic blood pressure. The upper alarm limit and the lower alarm limit for the patient's systolic blood pressure define a systolic blood pressure alarm range. The upper alarm limit and the lower alarm limit for the patient's diastolic blood pressure define a diastolic blood pressure alarm range. An alarm associated with the patient's blood pressure is active when the patient's systolic blood pressure is outside the systolic blood pressure alarm range or when the patient's diastolic blood pressure is outside the diastolic blood pressure alarm range.

The NIBP frame 314a also contains a NIBP cuff inflation stop button 322a. The NIBP cuff inflation stop button 322a is labeled with the text "Stop." The NIBP frame 314a also contains a NIBP automatic interval timer 322b. The NIBP automatic interval timer 322b is located between the diastolic blood pressure value and the NIBP alarm status area 322c. Selection of the NIBP cuff inflation stop button 322a ceases inflation of the NIBP cuff and toggles the label of the NIBP cuff inflation stop button 322a to display the ("Start") text. As used herein, a user selects a button or control when the user provides input to the PMP device 200 that specifies the control. For example, a user can select a control by pressing the control, by pressing another button while the control is highlighted, or by another means. Selection of the NIBP cuff inflation stop button 322a (now labeled the "Start" button) restarts inflation of the NIBP cuff and toggles the label of the NIBP cuff inflation stop button 322a to display the ("Stop") text. The NIBP automatic interval timer 322b indicates an amount of time remaining before the next scheduled inflation of the NIBP cuff. Additionally, a user can determine the age of the current NIBP reading on the NIBP frame 314a by subtracting the remaining time on the NIBP automatic interval timer 322b from the original interval duration. In the example of FIG. 3A, the NIBP frame 314a also displays a MAP value in an extended label field 322d.

The pulse rate frame 314b contains a representation of the patient's pulse rate. The representation of the patient's pulse

rate is based on one or more measurements of the patient's pulse rate. In different embodiments, the pulse rate frame **314b** contains different representations of the patient's pulse rate. In the example of FIG. 3A, the pulse rate frame **314b** includes enlarged numerical text that represents a pulse rate value ("122"). The pulse rate value ("122") is located at the left side of the pulse rate frame **314b**. A pulse rate alarm status area **324a** is located at the right side of the pulse rate frame **314b**. The pulse rate frame **314b** also indicates a source of the pulse rate in an extended label field **324b**.

The pulse rate alarm status area **324a** specifies an upper alarm limit and a lower alarm limit. The upper alarm limit and the lower alarm limit define a pulse rate alarm range. An alarm associated with the patient's pulse rate is active when the patient's pulse rate is outside the pulse rate alarm range.

The SpO2 frame **314c** contains a representation of the patient's SpO2 level. The representation of the patient's SpO2 level is based on one or more measurements of the patient's SpO2 level. In different embodiments, the SpO2 frame **314c** contains different representations of the patient's SpO2 level. In the example of FIG. 3A, the SpO2 frame **314c** includes enlarged numerical text that represents an SpO2 value ("97%"). The SpO2 value ("97%") is located at the left side of the SpO2 frame **314c** and is accompanied by a '%' text character on the right side of the SpO2 value. A SpO2 alarm status area **326a** is located at the right side of the SpO2 frame **314c**. An SpO2 alarm parameter **326d**, appearing as a circle adjacent to the text ("25"), indicates a duration of time. The SpO2 alarm status area **326a** specifies an upper alarm limit and a lower alarm limit. The upper alarm limit and the lower alarm limit define a SpO2 alarm range. An alarm associated with the patient's SpO2 level is active when the patient's SpO2 level is outside the SpO2 alarm range for the duration of time indicated by the SpO2 alarm parameter **326d**. The SpO2 frame **314c** also includes a pulse amplitude blip bar **326b** which indicates pulse beat and shows the relative pulse amplitude. As the detected pulse becomes stronger, more bars in the pulse amplitude blip bar **326b** light up with each pulse. The SpO2 frame **314c** also includes an SpO2 response time control button **326c** that is configured for a user to control the SpO2 alarm parameter **326d**.

The temperature frame **314d** contains a representation of the patient's body temperature. The representation of the patient's body temperature is based on one or more measurements of the patient's body temperature. In different embodiments, the temperature frame **314d** contains different representations of the patient's body temperature. In the example of FIG. 3A, the temperature frame **314d** includes enlarged numerical text that represents a temperature value ("101.5"). The temperature value ("101.5") is located at the left side of the temperature frame **314d** and is accompanied by a Fahrenheit degree indicating symbol on the right side of the temperature value. A temperature alarm status area **328a** is located at the right side of the temperature frame **314d**. The temperature alarm status area **328a** specifies an upper alarm limit and a lower alarm limit. The upper alarm limit and the lower alarm limit define a temperature alarm range. An alarm associated with the patient's temperature is active when the patient's temperature level is outside the temperature alarm range.

In some embodiments, the PMP device **200** can measure the patient's temperature in either a predictive mode or in a direct mode. When the PMP device **200** measures the patient's temperature in the predictive mode, the PMP device **200** predicts the patient's current temperature based on periodic readings of the patient's temperature. When the PMP

device **200** measures the patient's temperature in the direct mode, the PMP device **200** continually measures the patient's temperature.

The temperature value in the temperature frame **314d** is based on measurements received from a thermometer attached to a patient. When the PMP device **200** measures the patient's temperature in the predictive mode, the thermometer can be located at various places on the patient's body. Example locations on the patient's body where the thermometer can be located include in the patient's mouth, on the patient's thigh, in the patient's armpit, in the patient's rectum, and other locations. The temperature frame **314d** includes a thermometry location control **328b**. The thermometry location control **328b** indicates a location on the patient's body where the thermometer is located. In the example of FIG. 3A, the thermometry location control **328b** indicates that the thermometer is located in the patient's mouth.

When a user selects the thermometry location control **328b**, the PMP device **200** updates the thermometry location control **328b** such that the thermometry location control **328b** indicates a different location on the patient's body or whether the PMP device **200** is to obtain measurements of the patient's temperature in direct mode. The user can continue selecting the thermometry location control **328b** until the thermometry location control **328b** indicates a location where the thermometer is located on the patient's body or until the thermometry location control **328b** indicates that measurements are to be obtained in direct mode. For example, in some embodiments, the PMP device **200** accepts readings from a thermometer when the thermometer is located in the patient's mouth, in an adult patient's armpit, or in a pediatric patient's armpit. In this example, the thermometry location control **328b** initially indicates that the thermometer is in a patient's mouth. When the user selects the thermometry location control **328b** one time, the thermometry location control **328b** indicates that the thermometer is located in an adult patient's armpit. When the user selects the thermometry location control **328b** again, the thermometry location control **328b** indicates that the thermometer is located in a pediatric patient's armpit. When the user selects the thermometry location control **328b** yet again, the thermometry location control **328b** indicates that the measurements of the patient's temperature are to be obtained in direct mode. When the user selects the thermometry location control **328b** yet again, the thermometry location control **328b** again indicates that the thermometer is located in a patient's mouth.

The display screen **218** enables a user to select the parameter reporting frames **314a-314d** in order to change how the physiological parameters are represented in the parameter reporting frames **314a-314d**. In other words, each of the parameter reporting frames **314a-314d** contains an initial representation of a physiological parameter. The parameter reporting frame displays an alternate representation of the physiological parameter instead of the initial representation of the physiological parameter when a user selects the parameter reporting frame. For example, selecting the temperature frame **314d** toggles the temperature value between being expressed in Fahrenheit or Centigrade. This feature is referred to as "tap to toggle." In another example, when a user selects the pulse rate frame **314b**, the PMP device **200** displays a waveform in the pulse rate frame **314b** instead of a number representing the patient's current pulse rate. The waveform represents a patient's pulse over time. In this example, when the user selects the pulse rate frame **314b** again, the PMP device **200** displays a number in the pulse rate frame **314b** representing the patient's current pulse rate. In yet another example, when the user selects the SpO2 frame **314c**, the

PMP device 200 displays a plethysmographic waveform view in the SpO2 frame 314c. FIG. 3D illustrates an example alternate representation of a patient's SpO2 level. In the example of FIG. 3D, the SpO2 frame 314c contains a plethysmographic waveform view 325.

The patient attribute area 316 contains data that specify attributes of a patient. In various embodiments, the patient attribute area 316 contains data that specify various attributes of a patient. For example, in some embodiments, text that identifies a patient is located at the left side of the patient attribute area 316. In the example of FIG. 3A, the patient attribute area 316 contains the value "83645211" to identify the patient.

Furthermore, in the example of FIG. 3A, the patient attribute area 316 contains a patient type button 316b that is labeled with the text "Adult." The patient type button 316b is located towards the center of the patient attribute area 316. Patient-related attribute values that are labeled with the text ("HEIGHT"), ("WEIGHT"), ("PAIN") and ("RR") are located to the right side of the patient type button 316b. A save button 316g that is labeled with the text ("Save") is located at the right side of the patient attribute area 316.

The patient type button 316b indicates a value of a patient type parameter associated with the current patient. The patient type parameter can store a value of "Adult", "Pediatric", or "Neonatal" patient type. The patient type parameter controls the amount of air pressure applied to the NIBP cuff. Furthermore, in some embodiments, the patient type parameter controls the default alarm limits for the patient's physiological parameters. Thus, by using the patient type button 316b to change the patient type parameter, the user can automatically change the alarm limits for the patient's physiological parameters. As described elsewhere in this document, the user can also manually set the alarm limits for physiological parameters. When a user selects the patient type button 316b, the monitoring workflow home screen 300 is replaced by the patient selection screen as if the patients tab 319b was selected.

When the save button 316g is selected while the PMP device is operating in the monitoring workflow, the PMP device 200 saves a patient reading to local non-volatile storage within the PMP device 200. The patient reading is a set of data that includes measurements of the physiological parameters of the patient. If applicable, the patient reading can also include data indicating attributes of the patient. A user can use the review screen to review saved patient reading. Furthermore, in some embodiments, the PMP device 200 automatically attempts to transmit the patient reading to another computing node. A user can use the settings screen to specify the other computing node. In some embodiments, the other computing node is an electronic medical records system. Optionally, the other computing node is the interface system 104. The PMP device 200 does not clear the monitoring workflow home screen 300 when the save button 316g is selected.

The monitoring workflow is designed for obtaining a series of measurements of one or more physiological parameters of an identified patient over a period of time. A user uses the patient selection screen to select the identified patient. The PMP device 200 is programmable via an intervals pane of the settings screen to periodically record measurements of one or more physiological parameters over time from the identified patient.

A user may select the save button 316g one or more times. Each time the user selects the save button 316g, the PMP device 200 locally saves the patient reading and attempts to send the patient reading to another computing node. The sending is dependent on the workflow in which the PMP

device 200 is operating. In the spot check workflow (see FIG. 3B), the PMP device 200 automatically sends the patient reading because a clinician is present when the measurements are obtained. When the PMP device 200 is operating in the monitoring workflow, a clinician may not be always present. Thus, the patient reading is not sent automatically. Instead, a clinician navigates to the review screen (See FIG. 5) and manually selects the patient reading that the user wants to send.

When the PMP device 200 is operating in the monitoring workflow, the representations of the physiological parameters and patient attributes remain displayed on the monitoring workflow home screen 300, regardless of whether or not the patient reading is saved locally and/or transmitted to another computing node in response to a selection of the save button 316g. Unlike the spot check and triage workflows, the PMP device 200 does not clear the monitoring workflow home screen 300 when a user selects the save button 316g.

In some embodiments, the home screens for each workflow supported by the PMP device 200 contain at least one common user interface element. The common user interface element has the same appearance in each of the home screens, but does not necessarily perform the same function in the home screens of the different workflows. In the example of FIG. 3A, the patient type button 316b and the save button 316g are user interface elements that are included within and common to the monitoring workflow home screen 300 and to home screens for the spot check workflow (FIG. 3B) and the triage workflow (FIG. 3C). Further, the patient type button 316b and the save button 316g have a same appearance but also have different functions depending upon which workflow the PMP device 200 is operating within.

When the PMP device 200 is operating in the monitoring workflow, the PMP device 200 clears the monitoring workflow home screen 300 upon discharge of the current (first) identified patient for which physiological parameters are being obtained. When the PMP device 200 clears the monitoring workflow home screen 300, the PMP device 200 modifies the monitoring workflow home screen 300 such that the monitoring workflow home screen 300 no longer contains data representing physiological parameters of a patient and attributes of the patient. The monitored patient is discharged by the selection of another (second) identified patient for which to obtain physiological data, via the patient selection screen 400 or by power-cycling the PMP device 200. Upon selection of the second patient, the PMP device 200 locally stores the patient reading for the first patient and transmits the patient reading for the first patient to another computing node when a user selects a send patient reading control 448a on the review screen (see FIG. 5). Upon attaching the sensors to the second identified patient, the PMP device 200 obtain measurements of a set of one or more physiological parameter from the second identified patient, periodically over time.

Unlike the home screens for non-monitoring workflows (e.g., the spot check workflow and the triage workflow), the monitoring workflow home screen 300 provides direct access to the alarm screen via the alarm status areas 322c, 324a, 326a, 328a located within the parameter reporting frames 314a-314d. Each alarm status area 322c, 324a, 326a, 328a (including the bell-shaped symbols) indicates high and low alarm limits and provides a visual indication of when an alarm is active.

FIG. 3E illustrates the monitoring workflow home screen 300 when an alarm is active. When an alarm is active with respect to a particular physiological parameter, a perimeter around the parameter reporting frame for the particular physiological parameter transitions from a gray color to a red color.

Changing the color of the perimeter provides a visual indication that the alarm is active. In some embodiments, the perimeter around the parameter reporting frame also flashes, thereby providing another visual indication that the alarm is active. Furthermore, the bell-shaped symbol within the alarm status area transitions from a white color to a red or yellow color, depending on a priority of the alarm. This provides another visual indication that the alarm is active. When the alarm is resolved, the bell-shaped symbol within the alarm status area transitions from the red or yellow color to the white color. In the example of FIG. 3E, an alarm associated with the patient's pulse rate is active. Accordingly, the perimeter around the pulse rate frame 314b is red instead of gray.

Furthermore, as illustrated in the example of FIG. 3E, when an alarm associated with a particular physiological parameter is active, the PMP device 200 causes the device status area 312 to display an alarm message 380. The alarm message 380 visually indicates that an alarm is active and indicates a brief description of the alert. In the example FIG. 3E, the alarm message 380 indicates that the reason for the alert is that the patient's heart rate is too high.

Furthermore, when an alarm associated with a particular physiological parameter is active, the PMP device 200 emits an alarm sound. The PMP device 200 continues to emit the alarm sound until the alarm is deactivated or until a user temporarily silences the alarm sound. When the user temporarily silences the alarm sound, the PMP device 200 suspends emitting the alarm sound for a given time period. In various embodiments, the user is able to temporarily silence the alarm sound in various ways. For example, in some embodiments, the user temporarily silences the alarm sound by touching the parameter reporting frame corresponding to the alarm. In other embodiments, the user temporarily silences the alarm sound by touching the alarm message 380. In yet other embodiments, the user temporarily silences the alarm sound by selecting a physical button on the PMP device 200.

When the user temporarily silences the alarm sound, the PMP device 200 resumes emitting the alarm sound after a given time period expires. For example, in some embodiments, the PMP device 200 resumes emitting the alarm sound after 30 seconds. When the user temporarily silences the alarm sound, the alarm message 380 indicates a time remaining before the PMP device 200 resumes emitting the alarm sound.

FIG. 3F illustrates the monitoring workflow home screen 300 when an alarm is active and when the alarm sound for the alarm has been temporarily silenced. In the example of FIG. 3E, the alarm message 380 indicates that eighty-nine seconds remain before the PMP device 200 resumes emitting the alarm sound.

In some embodiments, the user is able to temporarily silence the alarm sound for various lengths of time by selecting a button on the PMP device 200 or a control displayed by the PMP device 200 multiple times. For example, when the user selects the alarm message 380 one time, the PMP device 200 resumes emitting the alarm sound after sixty seconds. In this example, each time the user selects the alarm message 380, the PMP device 200 adds thirty seconds to the length of time before the PMP device 200 resumes emitting the alarm sound. In some embodiments, the PMP device 200 prevents the user from temporarily silencing the alarm sound for more than a given amount of time. For example, in some embodiments, the PMP device 200 prevents the user from temporarily silencing the alarm sound for more than five minutes.

When the user temporarily silences an alarm sound for an alarm, the monitoring workflow home screen 300 visually indicates that the alarm sound has been temporarily silenced.

In various embodiments, the monitoring workflow home screen 300 visually indicates that the alarm sound has been temporarily silenced in various ways. For example, in some embodiments, the alarm message 380 contains a bell-shaped icon when the alarm sound has been temporarily silenced. In this example, the bell-shaped icon has an X-shaped mark over a bell. The lines of the X-shaped mark are dashed. In other example embodiments, the alarm message 380 or other parts of the monitoring workflow home screen 300 contain differently shaped icons or visual indicators.

The alarm status areas 322c, 324a, 326a, 328a act as navigational short cuts to appropriate panes within the alarms screen for each respective physiological parameter (See FIG. 7A). In this way, selecting an alarm status area while no alarm is active results in the PMP device 200 displaying the same user interface as if a user selects the alarms tab 319c and then selects an appropriate pane of the alarms screen. Navigation to the appropriate pane of the alarms screen happens only when the PMP device 200 is not emitting an alarm sound. In other words, the first selection of this area (or anywhere in the parameter reporting frame) would cause the alarm to silence the audio. A subsequent selection of the alarm status area performs the navigation.

Various alarms have various priority levels. For example, an alarm associated with a patient's pulse rate has a high priority level and an alarm associated with detachment of a SpO2 clip has a medium priority level. The PMP device 200 visually and/or sonically indicates alarms having different priority levels in different ways. For example, in some embodiments, when an alarm having a high priority level is active, the PMP device 200 displays visual indications of the alarm in a given color, such as red. When an alarm having a medium or low priority level is active, the PMP device 200 displays visual indications of the alarm in another color, such as yellow. Furthermore, in some embodiments, when an alarm having a high priority level is active, the PMP device 200 emits a given alarm sound, such as a continuous tone. When an alarm having a medium or low priority level is active, the PMP device 200 emits another alarm sound, such as a periodic beep.

When the PMP device 200 is operating in different workflows, the PMP device 200 can activate different alarms and/or emit different alarm sounds. For example, the PMP device 200 activates an alarm when a SpO2 clip detachment event occurs while the PMP device 200 operating within the monitoring workflow, but does not activate the alarm while in a non-operating in the monitoring workflow. This feature can be convenient considering that the PMP device 200 and the patient are more likely to be unattended by a user of the PMP device 200 for periods of time when the PMP device 200 is operating in the monitoring workflow, but less likely when the PMP device 200 is operating in the spot check or the triage workflows. In one embodiment, monitoring of SpO2 may be prevented by the PMP device 200 beyond ten minutes if the PMP device 200 is operating in the spot check or the triage workflows.

Under certain circumstances, multiple alarms can be active concurrently. For example, one alarm can be active because a patient's pulse rate is too high and another alarm can be active because a SpO2 clip has been detached from the patient. When multiple alarms are active concurrently, the PMP device 200 indicates that multiple alarms are active. In various embodiments, the PMP device 200 indicates that multiple alarms are active in various ways. For example, in some embodiments, the alarm message 380 includes an icon that indicates a number of alarms that are currently active. In this

example, icon can be a triangle containing the number of alarms that are currently active.

In some embodiments, the PMP device 200 displays different visual indicators depending on the priority levels of concurrently active alarms. For example, when two or more high level alarms are active concurrently, the PMP device 200 displays visual indications for each of the high level alarms. In another example, when an alarm having a medium priority level (“the medium level alarm”) and an alarm having a high priority level (“the high level alarm”) are active concurrently, the PMP device 200 displays visual indications for the high level alarm and does not display visual indications for the medium level alarm. If the user of the PMP device 200 temporarily silences the alarm sound for the high level alarm, the PMP device 200 does not emit the alarm sound for the high level alarm or the alarm sound for the medium level alarm. When the high level alarm is resolved before the medium level alarm is resolved, the PMP device 200 displays visual indications for the medium level alarm and the PMP device 200 resumes emitting the alarm sound for the medium level alarm.

In some embodiments, the user of the PMP device 200 is able to toggle between multiple concurrently active alarms by selecting one or more controls in the user interface of the PMP device 200 or by selecting one or more buttons on the PMP device 200. In this way, the user causes the PMP device 200 to display the visual indications for each of the multiple concurrently active alarms. For example, a first alarm and a second alarm can be concurrently active. In this example, the PMP device 200 displays an alarm message describing a first alarm. In this example, the PMP device 200 modifies the alarm message to describe a second alarm when a user selects the alarm message.

In some embodiments, the PMP device 200 enables a user to customize the content within the monitoring workflow home screen 300. For example, the PMP device 200 enables a user to adjust the relative sizes of the parameter reporting frames 314a-314d within the monitoring workflow home screen 300. In another example, the PMP device 200 enables the user to add or remove parameter reporting frames from the monitoring workflow home screen 300. Furthermore, in some such embodiments, the PMP device 200 includes one or more predefined templates for the monitoring workflow home screen 300. Each of the predefined templates specifies a predefined set of content within the monitoring workflow home screen 300. For example, one predefined template specifies that the monitoring workflow home screen 300 includes a large parameter reporting frame for the patient’s blood pressure and three smaller parameter reporting frames for the patient’s pulse rate, SpO2 level, and body temperature. In this example, another predefined template specifies that the monitoring workflow home screen 300 includes only a frame for the patient’s SpO2 level and the patient’s blood pressure. The user can customize the content within the monitoring workflow home screen 300 by selecting one of these predefined templates or by modifying one of these predefined templates.

FIG. 3B illustrates an example spot check workflow home screen 330. The PMP device 200 displays the spot check workflow home screen 330 when the PMP device 200 is operating in a spot check workflow. The spot check workflow home screen 330 is referred to as the “Home” screen for the spot check workflow. Like when the PMP device 200 is operating in the monitoring workflow, each sensor is physically attached to an identified patient when the PMP device 200 is operating in the spot check workflow. These sensors include a temperature probe, a SpO2 clip, and a NIBP blood pressure cuff. As described elsewhere in this document, the tempera-

ture probe, the SpO2 clip and the NIBP blood pressure cuff are peripheral to the PMP device 200.

Also, like the monitoring workflow home screen 300 illustrated in the example of FIG. 3A, the spot check workflow home screen 330 includes a device status area 312, a content area 320, and a navigation area 318. The content area 320 includes a parameter reporting area 314 and a patient attribute area 316. The parameter reporting area 314 of the spot check workflow home screen 330 includes an NIBP frame 334a, a pulse rate frame 314b, a SpO2 frame 334c and a temperature frame 334d.

The NIBP frame 334a contains a representation of a patient’s blood pressure. The representation is based on one or more measurements of the blood pressure of a previously identified patient. In the example of FIG. 3B, the NIBP frame 334a includes text representing the patient’s systolic and diastolic blood pressure. In the spot check workflow home screen 330, the text representing the patient’s systolic and diastolic blood pressure is larger than the text in the monitoring workflow home screen 300 representing the patient’s systolic and diastolic blood pressure. The NIBP frame 334a of the spot check workflow home screen 330 does not include the NIBP alarm status area 322c and does not include the NIBP automatic interval timer 322b of the monitoring workflow home screen 300.

The pulse rate frame 334b contains a representation of the patient’s pulse rate. The representation is based on one or more measurements of the pulse rate of the patient. In the example of FIG. 3B, the pulse rate frame 334b includes text representing the patient’s pulse rate. In the spot check workflow home screen 330, the text representing the patient’s pulse rate is larger than the text in the monitoring workflow home screen 300 representing the patient’s pulse rate. The pulse rate frame 334b of the spot check workflow home screen 330 does not include the pulse rate alarm status area 324a included within the pulse rate frame 314b of the monitoring workflow home screen 300.

The SpO2 frame 334c contains a representation of the patient’s SpO2 level. The representation is based on one or more measurements of the SpO2 level of the patient. In the example of FIG. 3B, the SpO2 frame 334c includes text representing the SpO2 value. In the spot check workflow home screen 330, the text representing the patient’s SpO2 level is larger than the text in the monitoring workflow home screen 300 representing the patient’s SpO2 level. The SpO2 frame 334c does not include the alarm status area 326a, the SpO2 response time control button 326c, or the SpO2 alarm parameter 326d of the SpO2 frame 314c of the monitoring workflow home screen 300.

The temperature frame 334d contains a representation of the patient’s body temperature. The representation is based on one or more measurements of the body temperature of the patient. In the example of FIG. 3B, the temperature frame 334d includes text representing the patient’s body temperature. In the spot check workflow home screen 330, the text representing the patient’s body temperature is larger than the text in the monitoring workflow home screen 300 representing the patient’s body temperature. The temperature frame 334d of the spot check workflow home screen 330 does not include the temperature alarm status area 328a included within the temperature frame 314d of the monitoring workflow home screen 300.

The patient attribute area 316 of the spot check workflow home screen 330 includes text that identifies the patient by name, initials, numerical identifier, or location. In the example of FIG. 3B, the patient attribute area 316 includes text that identifies the patient as “Bar, D.” A user is able to

the settings screen to configure the PMP device 200 to identify the patient by name or by number. The user is also able to use the settings screen to configure the PMP device 200 not to save or send patient readings when the PMP device 200 does not store information regarding the identity of the patient.

The navigation area 318 of the spot check workflow home screen 330 includes the home tab 319a, the patients tab 319b, the review tab 319d, and the settings tab 319e. The navigation area 318 excludes the alarms tab 319c included in the navigation area 318 of the monitoring workflow home screen 300. As a result, the spot check workflow home screen 330 does not provide direct navigation to the alarms screen as provided by the monitoring workflow home screen 300.

When the PMP device 200 is operating in the spot check workflow, the PMP device 200 locally saves a patient reading and attempts to send the patient reading to another computing node when a user selects the save button 316g. The patient reading includes measurements of the physiological parameters of the patient and data identifying the patient. Unlike when the PMP device 200 is operating in the monitoring workflow, the PMP device 200 clears the spot check workflow home screen 330 when the user selects the save button 316g. When the PMP device 200 clears the spot check workflow home screen 330, the PMP device 200 modifies the spot check workflow home screen 330 such that the spot check workflow home screen 330 no longer contains representations of the physiological parameters of the patient and attributes of the patient.

Because the PMP device 200 clears the spot check workflow home screen 330 when a user selects the save button 316g of the spot check workflow home screen 330, a first identified patient is essentially discharged when the user selects the save button 316g. After selecting the save button 316g, the user selects a second identified patient via the patient selection screen. Upon attaching the sensors to the second identified patient, the PMP device 200 obtains a set of one or more physiological parameter values from the second identified patient.

The spot check workflow is designed for obtaining measurements of physiological parameters from each patient in a series of identified patients. A user selects each patient in the series of identified patients from the patient selection screen. Alternatively, a user can select each patient in the series of identified patient by scanning barcodes of the patients, thereby bypassing the patient selection screen. Further, a user can identify a patient at any step in the workflow prior to saving the patient reading of the spot check workflow home screen 330. Upon attaching the sensors a patient in the series of identified patients, a user uses the PMP device 200 to obtain measurements of one or more physiological parameters of that patient. Typically, the PMP device 200 is attached to a first patient for no more time than is required to obtain one measurement for each of the NIBP, the pulse rate, the SpO₂, and the body temperature of the first patient. The user then detaches the sensors from the first patient and attaches the sensors to a second patient that is next in the series of patients.

In one use scenario, a user can use the spot check workflow “making rounds” within a health care facility. For example, a clinician can use the PMP device 200 to obtain one set of measurements of physiological parameters for each patient in a group of twelve patients within a health care facility. Optionally, the clinician can transport the PMP device 200 sequentially to each patient in the group of patients within one “round” of obtaining measurements of physiological parameters. Each “round” of obtaining measurements of physiological parameters can be obtained for each patient in the patients of the group every hour during a working shift.

In some embodiments, the PMP device 200 enables a user to customize the content within the spot check workflow home screen 330. For example, the PMP device 200 enables a user to adjust the relative sizes of the parameter reporting frames 334a-334d within the spot check workflow home screen 330. In another example, the PMP device 200 enables the user to add or remove parameter reporting frames from the spot check workflow home screen 330. Furthermore, in some such embodiments, the PMP device 200 includes one or more predefined templates for the spot check workflow home screen 330. Each of the predefined templates specifies a predefined set of content within the spot check workflow home screen 330. The user can customize the content within the spot check workflow home screen 330 by selecting one of these predefined templates or by modifying one of these predefined templates.

FIG. 3C illustrates an example triage workflow home screen 360. The PMP device 200 displays the triage workflow home screen 360 when the PMP device 200 is in a triage workflow.

Like when the PMP device 200 is operating in the monitoring workflow or the spot check workflow, a user physically attaches each sensor to an unidentified patient when the PMP device 200 is operating in the triage workflow. These sensors include the temperature probe, the SpO₂ clip, and the blood pressure cuff that are each attachable to the PMP device 200 as described elsewhere in this document.

The triage workflow home screen 360 includes a device status area 312, a content area 320, and a navigation area 318. The content area 320 of the triage workflow home screen 360 includes a parameter reporting area 314 and a patient attribute area 316. The parameter reporting area 314 of the triage workflow home screen 360 includes an NIBP frame 364a, a pulse rate frame 364b, a SpO₂ frame 364c and a temperature frame 364d.

The NIBP frame 364a contains a representation of the systolic and diastolic blood pressure of the patient. The representation is based on one or more measurements of the blood pressure of the patient. In the example of FIG. 3C, the NIBP frame 364a includes text representing the systolic and diastolic blood pressure of the patient. The size of the text in the NIBP frame 364a is similar to the size of the text in the NIBP frame 334a of the spot check workflow home screen 330 and larger than the size of the text in the NIBP frame 314a of the monitoring workflow home screen 300. The NIBP frame 364a does not include the NIBP alarm status area 322c and does not include the NIBP automatic interval timer 322b included in the NIBP frame 314a of the monitoring workflow home screen 300.

The pulse rate frame 364b contains a representation of the pulse rate of the patient. The representation is based on one or more measurements of the pulse rate of the patient. In the example of FIG. 3C, the pulse rate frame 364b contains text representing the pulse rate value of the patient. The size of the text in the pulse rate frame 364b is similar to the size of the text in the pulse rate frame 334b of the spot check workflow home screen 330 and larger than the size of the text in the pulse rate frame 314b of the monitoring workflow home screen 300. The pulse rate frame 364b does not include the pulse rate alarm status area 324a that is included within the pulse rate frame 314b of the monitoring workflow home screen 300.

The SpO₂ frame 364c contains a representation of the SpO₂ level of the patient. The representation is based on one or more measurements of the SpO₂ level of the patient. In the example of FIG. 3C, the SpO₂ frame 364c includes text representing the SpO₂ value of the patient. The size of the text

in the SpO2 frame 364c is similar to the size of the text in the SpO2 frame 334c of the spot check workflow home screen 330 and larger than the size of the text in the SpO2 frame 314c of the monitoring workflow home screen 300. The SpO2 frame 334c does not include the alarm status area 326a, the SpO2 response time control button 326b, or the SpO2 alarm parameter 326b of the SpO2 frame 314c of the monitoring workflow home screen 300.

The temperature frame 364d contains a representation of the body temperature of the patient. The representation is based on one or more measurements of the body temperature of the patient. In the example of FIG. 3C, the temperature frame 364d includes text representing the body temperature of the patient. The size of the text in the temperature frame 364d is similar to the size of the text in the temperature frame 334d of the spot check workflow home screen 330 and larger than the size of the text in the temperature frame 314d of the monitoring workflow home screen 300. The temperature frame 364d does not include the temperature alarm status area 328a that is included within the temperature frame 314d of the monitoring workflow home screen 300.

The navigation area 318 of the triage workflow home screen 360 includes the home tab 319a, the review tab 319d, and the settings tab 319e. The navigation area 318 of the triage workflow home screen 360 excludes the alarms tab 319c of the monitoring workflow home screen 300. As a result, the triage workflow home screen 360 does not provide direct navigation to the alarms screen as provided by the monitoring workflow home screen 300. Furthermore, the navigation area 318 of the triage workflow home screen 360 excludes the patients tab 319b of the monitoring workflow home screen 300 and the spot check workflow home screen 330. As a result, the triage workflow home screen 360 does not provide direct navigation to the patient selection screen as provided by the monitoring workflow home screen 300 and the spot check workflow home screen 330.

The patient attribute area 316 of the triage workflow home screen 360 includes a patient type button 316b and a save button 316g. The patient attribute area 316 does not include text that identifies the patient.

When the PMP device 200 is operating within the triage workflow, the PMP device 200 locally saves a patient reading of the triage workflow home screen 360. The home screen data of the triage workflow home screen 360 includes measurements of the physiological parameters of the unidentified patient. Furthermore, a user can configure the PMP device 200 to transmit the patient reading of the triage workflow home screen 360 to another computing node in response to a selection of the save button 316g. Unlike the monitoring workflow home screen 300 and like the spot check workflow home screen 330, the PMP device 200 clears the triage workflow home screen 360 when a user selects the save button 316g. When the PMP device 200 clears the triage workflow home screen 360, the PMP device 200 modifies the triage workflow home screen 360 such that the triage workflow home screen 360 no longer contains representations of the physiological parameters of the patient.

When the PMP device 200 is operating in the triage workflow, selection of the patient type button 316b toggles the text label of the patient type button 316b and toggles a value of a patient type parameter within the triage workflow home screen 360. The text label and data value associated with the patient type button 316b, toggles between the label/values of (“Adult”), (“Pediatric”) and (“Neonatal”).

When the user selects the save button 316g, the patient reading transmitted to another computing node includes a

patient type parameter value equal to (“Adult”), (“Pediatric”) and (“Neonatal”). The patient type data substitutes for patient identification data that is absent from the patient reading while operating within the triage workflow.

As a result, when the PMP device 200 is operating in the triage workflow, a first unidentified patient is essentially discharged when a user selects the save button 316g. Upon attaching the sensors to a second unidentified patient, the user uses the PMP device 200 to obtain a set of one or more measurements of physiological parameter from the second unidentified patient.

The triage workflow is designed for obtaining measurements of physiological parameters from each of a series of unidentified patients. A user does not select each patient in the series of unidentified patients from the patient selection screen. Instead, the user attaches the sensors to each unidentified patient in sequence. Upon attaching the sensors to a patient in the series of the unidentified patients, the user uses the PMP device 200 to obtain measurements of one or more physiological parameters of that patient at that time. Typically, the user attaches the sensors to the patient for no more time that is required to obtain one measurement for each of the NIBP, the pulse rate, the SpO2 level, and the body temperature of the patient. The user then detaches the sensors are from the patient and attaches the sensors to another patient that is next in the series of patients.

In one use scenario, the triage workflow can be used for obtaining measurements of physiological parameters from unidentified health care recipients. The recipients may or may not be patients of a health care facility. The other computing node is a personal computer which receives the patient reading from the PMP device 200 with no patient names attached.

In some embodiments, the PMP device 200 enables a user to customize the content within the triage workflow home screen 360. For example, the PMP device 200 enables a user to adjust the relative sizes of the parameter reporting frames 364 within the triage workflow home screen 360. In another example, the PMP device 200 enables the user to add or remove parameter reporting frames from the triage workflow home screen 360. Furthermore, in some such embodiments, the PMP device 200 includes one or more predefined templates for the triage workflow home screen 360. Each of the predefined templates specifies a predefined set of content within the triage workflow home screen 360. The user can customize the content within the triage workflow home screen 360 by selecting one of these predefined templates or by modifying one of these predefined templates.

FIG. 4 illustrates an example patient selection screen 400. A user can navigate to the patient selection screen 400 by selecting the patients tab 319b. The patient selection screen 400 includes a device status area 312, a pane selection tab area 310 and a navigation area 318. The pane selection tab area 310 includes four pane selection tabs that are labeled “List”, “Summary”, “Modifiers” and “Manual.” The list pane 418 includes a column label area 412, a patient listing area 414, and a screen control area 416. The device status area 312 appears like the device status area 312 of the monitoring workflow home screen 300, the device status area 312 of the spot check workflow home screen 330, and the device status area 312 of the triage workflow home screen 360. The navigation area 318 appears like the navigation area 318 of the monitoring workflow home screen 300.

The patient listing area 414 includes lines 415a-415e (collectively, “lines 415”). Each of the lines 415 is associated with a different patient. For example, a patient having a name “Barker, D.” is listed on the line 415a, a patient having a name

“Connor, W.” is listed on the line **415c**, and a patient having a name “Davison, A” is listed on the line **415e**.

The screen control area **416** includes three buttons. A first (left) button that is labeled “Add” is utilized to add a patient to the patient listing area **414**. A second (middle) button that is labeled “Delete” is utilized to delete a patient from the patient listing area **414**. A third (right) button that is labeled “Select” is utilized to highlight a patient that listed within the patient listing area **414**.

FIG. 5 illustrates an example review screen **500**. A user can navigate to the review screen **500** using the review tab **319d**. The review screen **500** includes a device status area **312**, a review table **541**, a screen control area **548**, and a navigation area **318**. The device status area **312** appears like the device status area **312** of the monitoring workflow home screen **300**, the device status area **312** of the spot check workflow home screen **330**, and the device status area **312** of the triage workflow home screen **360**. The navigation area **318** of the review screen **500** appears like the navigation area **318** of the monitoring workflow home screen **300**.

The review table **541** includes a column label line **542**. In addition, the review table **541** includes patient readings **546a-546f** and columns **552a-552i**. Each of the patient readings **546a-546f** comprises a different set of measurements of physiological parameters of a patient. Each of the columns **552b-552h** is associated with a different physiological parameter. Check marks in the column **552a** indicate whether particular ones of the patient readings **546a-546f** are selected. In the example of FIG. 5, the patient readings **546a** and **546e** are selected.

For example, the patient reading **546a** includes a patient name (“Barker, D.”) in the column **552b**, a date and time value in the column **552c**, a systolic and diastolic measurement value in the column **552d**, a mean arterial pressure (MAP) measurement value (not shown), a temperature measurement value in the column **552f**, a pulse measurement value in the column **552g**, an SpO2 measurement value in the column **552h**, and a patient attributes set of values in the column **552i**.

The screen control area **548** includes a send patient reading button **548a**, a print patient reading button **548b**, a delete patient reading button **548c** and a view list box **548d**. The print patient reading button **548b** is utilized to send a selected patient reading. The print patient reading button **548b** is utilized to print selected patient readings. The delete patient reading button **548c** is utilized to delete selected patient readings. The view list box **548d** enables filtering of the review table **541**. A user can use the view list box **548d** to select to view all available patient readings or only particular patient readings based on sent, unsent or alarm status.

FIG. 6A illustrates an example intervals pane **610** of a settings screen **600**. Direct navigation to the settings screen **600** is provided from the monitoring workflow home screen **300**, from the spot check workflow home screen **330**, and from the triage workflow home screen **360**. For example, in some embodiments, a user of the PMP device **200** navigates to the settings screen **600** by selecting the settings tab **319e**. As illustrated in the example of FIG. 6A, the settings screen **600** includes a device status area **312**, a pane selection tab area **604**, the intervals pane **610** and a navigation area **318**.

The device status area **312** of the settings screen **600** appears like the device status area **312** of the monitoring workflow home screen **300**, the device status area **312** of the spot check workflow home screen **330**, and the device status area **312** of the triage workflow home screen **360**. The pane selection tab area **604** includes five pane selection tabs that are labeled “Intervals,” “Device,” “Clinician,” “Profiles,” and

“Advanced.” The navigation area **318** of the settings screen **600** appears like the navigation area **318** of the monitoring workflow home screen **300**.

A user can use the intervals pane **610** to program the PMP device **200** to save NIBP measurements, pulse rate measurements, SpO2 measurements, and temperature measurements for a patient at given intervals over an given time period. Other combinations of physiological parameters can be obtained at other fixed or non-fixed intervals.

The intervals pane **610** includes an automatic control **611**, a program control **612**, a stat control **613**, an off control (not shown), a print control **615**, and a start control **616**. When the automatic control **611** is selected, the settings screen **600** contains one or more controls for entering a fixed time period to right of control. When the program control **612** is selected, the settings screen **600** contains controls for defining and selecting interval programs. When the stat control **613** is selected, the settings screen **600** contains one or more controls for capturing and saving as many measurements of physiological parameters as possible within a time period, for example a five minute period. When the off control is selected, intervals are turned off. When the print control **615** is selected, a captured measurement at each interval may be printed. The start control **616** is a command button to start the selected interval workflow: Automatic, Program, or Stat. The PMP device **200** automatically initiates an NIBP reading and initiates the interval program once the start control **616** is selected.

In the example of FIG. 6A, the program control **612** is selected. Consequently, the intervals pane **610** includes controls for defining and selecting interval programs. An interval program is a set of parameters that governs the intervals at which the PMP device **200** records measurements of one or more physiological parameters of a patient and for lengths of time for which the PMP device **200** records the measurements of the one or more physiological parameters of the patient.

In the example of FIG. 6A, the intervals pane **610** contains interval program selection controls **660a**, **660b** (collectively, interval program selection controls **660**). More or fewer controls **660** can be provided. The interval program selection controls **660** specify the names of interval programs. In some examples, the interval program selection controls **660** specify names associated with the controls, such as “Transplant,” “Transfusion,” “Cardiac,” which represent names of interval programs.

A user selects an interval program by selecting one of the interval program selection controls **660** that corresponds to the interval program. After selecting the interval program, the user starts the interval program by selecting the start control **616**. When the user starts the interval program, the PMP device **200** records measurements of physiological parameters of a patient at the intervals specified by the interval program for the duration specified by the interval program.

In some embodiments, the monitoring workflow home screen **300** contains remaining time controls. The remaining time controls indicate amounts of time remaining before the PMP device **200** records measurements of the physiological parameters. For example, the pulse rate frame **314b** of the monitoring workflow home screen **300** can include a remaining time control that indicates that twenty seconds remain before the PMP device **200** records a measurement of the patient’s pulse rate. Furthermore, in some embodiments, the PMP device **200** displays the intervals pane **610** when the user selects the remaining time controls. In this way, the user of the PMP device **200** can easily access the intervals pane **610**.

Furthermore, in the example of FIG. 6A, the intervals pane **610** contains a program name control **662**, duration controls

664a-664e (collectively, “duration controls 664”), and interval controls 666a-666e (collectively, “interval controls 666”). The duration controls 664 and the interval controls 666 are organized into rows. Each row corresponds to a different physiological parameter of the patient. For example, a first row corresponds to the patient’s blood pressure, a second row corresponds to the patient’s pulse rate, and so on.

When the user selects an interval program using one of the interval program selection controls 660, the program name control 662 contains the name of the interval program. In the example of FIG. 6A, the program name control 662 indicates that the name of the interval program is “Transfusion.” The duration controls 664 specify lengths of time during which the PMP device 200 is to record measurements of the physiological parameters. The interval controls 666 specify lengths of time that the PMP device 200 waits between recording measurements of the physiological parameters. In the example of FIG. 6A, the duration control 664a indicates 1 hour and the interval control 666a indicates 15 minutes. Consequently, the PMP device 200 records measurements of a physiological parameter at 15 minute intervals for one hour after the “Transfusion” interval program starts.

When the user selects an interval program using one of the interval program selection controls 660, the user is able to change the values in the program name control 662, the duration controls 664, and the interval controls 666. For example, the user can change the value of the duration control 664b from one hour to two hours.

FIG. 6B illustrates an example device pane 620 of the settings screen 600. The settings screen 600 includes a device status area 312, a pane selection tab area 604, a device pane 520 and a navigation area 318. The device status area 312, the pane selection tab area 604, and the navigation area 318 appear like that of FIG. 6A.

The device pane 620 provides controls and information for a user to select or view device characteristics. The device pane 620 includes a set date and time control 622, a touch-screen lock-out control 628, and a command button 629. The set date and time control 622 is configured to set up a date and time for the PMP device 200. A user may select a date display format through the date format selection (not shown), such as MM/DD/YYYY. A user may also select a time zone offset from UTC based on where the PMP device 200 is located. The touch-screen lock-out control 628 is configured to allow a user to lock out the touch screen. When the command button 629 is selected, the touch screen is locked out immediately.

FIG. 6C illustrates an example clinician pane of the settings screen 600. The settings screen 600 includes a device status area 312, a pane selection tab area 604, a clinician pane 630 and a navigation area 318. The device status area 312, the pane selection tab area 604 and the navigation area 318 appear like that of FIG. 6A. The clinician pane 630 provides for entry of clinician identification. The clinician identified can include a first name, a last name, a middle initial, and an identification number.

FIGS. 6D-6F illustrate an example profiles pane 640 of the settings screen 600. The settings screen 600 includes the device status area 312, a pane selection tab area 604, the profiles pane 640 and a navigation area 318. The device status area 312, the pane selection tab area 604 and the navigation area 318 appear like that of FIG. 6A. The profiles pane 640 allows a user to select a desired workflow within which the PMP device 200 is to operate. In the examples of FIG. 6D-6F, the profiles pane 640 includes a “Med/Surge” control 642, a “Spot Check” control 644, and an “ED Triage” control 646. A user selects the “Med/Surge” control 642 to select a “Med/Surge” workflow as the workflow in which the PMP device

200 is to operate. The user selects the “Spot Check” control 644 to select the spot check workflow as the workflow in which the PMP device 200 is to operate. The user selects the “ED triage” control 646 to select the triage workflow as the workflow in which the PMP device 200 is to operate. As shown in the examples of FIGS. 6D-6F, the profiles pane 640 shows a preview of the home screen of the selected workflow.

FIG. 6G illustrates an example advanced pane 650 of the settings screen 600. The settings screen 600 includes a device status area 312, a pane selection tab area 604, the advanced pane 650 and a navigation area 318. The device status area 312, the pane selection tab area 504 and the navigation area 318 appear like that of FIG. 6A. The advanced pane 650 enables an administrator to access configuration settings that are likely to be set once or less frequently than those required for a normal use. When a proper access code is entered, the PMP device 200 switches to an administration workflow and the navigation tabs switch to advanced settings tabs. The advanced pane 650 also provides read-only information about the PMP device 200.

FIG. 7A illustrates an example global pane 710 of an alarms screen 700. A user of the PMP device 200 navigates to the alarms screen 700 by selecting the alarms tab 319c. The alarms screen 700 includes controls that enable a user to configure global alarm settings and parameter-specific alarm settings.

As illustrated in the example of FIG. 7A, the alarms screen 700 includes the device status area 312, a compressed parameter reporting area 702, a pane selection tab area 704, the global pane 710, and the navigation area 318. The device status area 312 appears like the device status areas of the monitoring workflow home screen 300, the spot check workflow home screen 330 and the triage workflow home screen 360. The navigation area 318 of the alarms screen 700 appears like the navigation area 318 of the monitoring workflow home screen 300.

The compressed parameter reporting area 702 contains one or more compressed parameter frames containing data representing live physiological parameters of a patient. In the example of FIG. 7A, the compressed parameter reporting area 702 contains a compressed NIBP frame 703a, a compressed SpO2 frame 703b, a compressed pulse rate frame 703c, and a compressed temperature frame 703d. The compressed NIBP frame 703a, the compressed SpO2 frame 703b, the compressed pulse rate frame 703c, and the compressed temperature frame 703d are referred to herein collectively as the compressed frames 703. The compressed parameter reporting area 702 has frame dimensions that are smaller than the frame dimensions of the parameter reporting area 314 of the monitoring workflow home screen 300. Because the compressed parameter reporting area 702 is smaller than the parameter reporting area 314, the compressed parameter reporting area 702 appears “squished” relative to the parameter reporting area 314.

The compressed NIBP frame 703a contains a representation of the systolic and the diastolic blood pressure of the patient. The compressed SpO2 frame 703b contains a representation of a SpO2 percentage of the patient, a pulse amplitude blip bar, and alarm limits for the SpO2 parameter. The compressed pulse rate frame 703c contains a representation of a pulse rate of the patient. The compressed temperature frame 703d contains a representation of the body temperature of the patient. In the example of FIG. 7A, the compressed frames 703 do not contain user-selectable controls. Consequently, the user is unable to toggle views of the physiological parameters from the compressed frames 703.

When an alarm associated with one of the physiological parameters is active, the corresponding one of the compressed frames 703 provides a visual indication of the alarm. For example, when the patient's blood pressure is too high, a perimeter of the compressed NIBP frame 703a changes from one color (e.g., gray) to a different color (e.g., red). Furthermore, when an alarm is active, the device status area 312 of the alarms screen 700 is replaced by an alarm message describing the alarm. In some embodiments, selecting the alarm message causes the PMP device 200 to temporarily stop emitting an alarm sound associated with the alarm.

The pane selection tab area 704 includes five pane selection tabs that are labeled "Global", "NIBP", "SpO2", "Pulse" and "Temperature." When the user selects one of the pane selection tabs in the pane selection tab area 704, the PMP device 200 updates the alarms screen 700 to contain a pane associated with the selected pane selection tab.

The global pane 710 is located below the compressed parameter reporting area 702. The global pane 710 provides controls for configuring global alarm settings for the PMP device 200. The global alarm settings are settings that apply to all alarms provided by the PMP device 200. In the example of FIG. 7A, the global pane 710 includes a reset button 712a, a display alarm limit parameter button 712b, a set of alarm audio buttons 712c, a set of volume buttons 712d, and a silence/reset button 712e. A user can use the reset button 712a to reset alarm limits. The user can use the display alarm limit parameter button 712b to configure the PMP device 200 to display or not to display alarm limits. The user can use the set of alarm audio buttons 712c to configure the PMP device 200 to emit alarm sounds. The user can use the set of volume buttons 712d to set a volume of the alarm sounds to high, medium or low. Also, the user can use the silence/reset button 712e to silence or reset alarms.

FIG. 7B illustrates an example temperature pane 720 of the alarms screen 700. The alarms screen 700 includes the device status area 312, the compressed parameter reporting area 702, the pane selection tab area 704, the temperature pane 720, and a navigation area 318. The device status area 312, the compressed parameter reporting area 702, the pane selection tab area 704 and the navigation area 318 appear like that of FIG. 7A.

The temperature pane 720 is located below the compressed parameter reporting area 702. The temperature pane 720 enables a user of the PMP device 200 to configure parameter-specific alarm settings for the temperature alarm. In the example of FIG. 7B, the temperature pane 720 provides controls 722a, 722b and 722c. The control 722a sets the temperature alarm to ON or OFF. The control 722b sets an upper temperature limit for the temperature alarm. The control 722c sets a lower temperature limit for the temperature alarm.

FIG. 7C illustrates an example NIBP pane 730 of the alarms screen 700. The alarms screen 700 includes a device status area 312, a compressed parameter reporting area 702, a pane selection tab area 704, the NIBP pane 730, and a navigation area 318. The device status area 312, the compressed parameter reporting area 702, the pane selection tab area 704 and the navigation area 318 appear like that of FIG. 7A.

The NIBP pane 730 enables a user of the PMP device 200 to configure parameter-specific alarm settings for the NIBP alarm. In the example of FIG. 7C, the NIBP pane 730 contains a control 732a to set the systolic and diastolic alarm ON or OFF. The NIBP pane 730 also contains controls 732b-732e to set the systolic and diastolic upper and lower alarm limits for the systolic and diastolic alarm. Further, the NIBP pane 730 provides controls 732f, 732g and 732h. The control 732f sets a MAP alarm to ON or OFF. The control 732g sets set an

upper MAP limit for the MAP alarm. The control 732h sets a lower MAP limit for the MAP alarm.

FIG. 7D illustrates an example SpO2 pane 740 of the alarms screen 700. The alarms screen 700 includes the device status area 312, the compressed parameter reporting area 702, the pane selection tab area 704, the SpO2 pane 740, and the navigation area 318. The device status area 312, the compressed parameter reporting area 702, the pane selection tab area 704, and the navigation area 318 appear like that of FIG. 7A.

The SpO2 pane 740 enables a user of the PMP device 200 to configure parameter-specific alarm settings for the SpO2 alarm. In the example of FIG. 7D, the SpO2 pane 740 provides controls 742a, 742b and 742c. The control 742a sets the SpO2 alarm to ON or OFF. The controls 742b and 742c set an upper SpO2 limit and a lower SpO2 limit for the SpO2 alarm. The SpO2 pane 740 also includes a SatSeconds™ control 743d. When a numeric value is selected, an SpO2 alarm is not activated until a length of time specified by the SatSeconds control 743d has expired. The length of time specified by the SatSeconds control 743d is based on a duration of a low SpO2 event multiplied by a number of percentage points that the patient's SpO2 falls outside the lower SpO2 alarm threshold. Also, the SatSeconds control 743d appears next to the controls 742b and 742c when the control 742a is set to ON. However, when the control 742a is set to OFF, the SatSeconds control 743d is not displayed and the alarm condition is not delayed per a SatSeconds algorithm.

FIG. 7E illustrates an example pulse rate pane 750 of the alarms screen 700. The alarms screen 700 includes the device status area 312, the compressed parameter reporting area 702, the pane selection tab area 704, the pulse rate pane 750 and the navigation area 318. The device status area 312, the compressed parameter reporting area 702, the pane selection tab area 704 and the navigation area 318 appear like that of FIG. 7A.

The pulse rate pane 750 enables a user of the PMP device 200 to configure parameter-specific alarm settings for the pulse rate alarm. In the example of FIG. 7E, the pulse rate pane 750 provides an on/off control 752a, an upper limit control 752b and a lower limit control 752c. The on/off control 752a allows a user to turn the pulse rate alarm on or off. The upper limit control 752b allows a user to set an upper pulse rate limit for the pulse rate alarm. The lower limit control 752c allows a user to set a lower pulse rate limit for the pulse rate alarm.

The upper limit control 752b includes a numerical portion, an up button and a down button. The numerical portion of the upper limit control 752b specifies an upper pulse rate limit for the pulse rate alarm. Selecting the up button of the upper limit control 752b incrementally increases the upper pulse rate limit. Selecting the down button of the upper limit control 752b incrementally decreases the upper pulse rate limit. Selecting the numerical portion of the upper limit control 752b causes the PMP device 200 to display a numerical keypad. The user is able to use the numerical keypad to manually enter a value for the upper pulse rate limit.

The lower limit control 752c includes a numerical portion, an up button, and a down button. The numerical portion of the lower limit control 752c specifies the lower pulse rate limit. Selecting the up button of the lower limit control 752c incrementally increases the lower pulse rate limit. Selecting the down button of the lower limit control 752c incrementally decreases the lower pulse rate limit. Selecting the numerical portion of the lower limit control 752c causes the PMP device

200 to display a numerical keypad. The user is able to use the numerical keypad to manually enter a value for the lower pulse rate limit.

In the examples shown, the values that are entered by the user are validated for accuracy. For example, upper and lower thresholds for each parameter can be set. If the user selects a value outside of the upper and lower thresholds, the system can prevent such a selection and/or alert the user. For example, if the user types in 300 minutes on the keypad for a parameter, and the upper threshold set for that parameter is 240 minutes, the system will not allow the user to increase the parameter to 300 minutes.

In some embodiments, various panes of the alarms screen 700 can include other upper limit controls and other lower limit controls like the upper limit control 752b and the lower limit control 752c. The other upper limit controls and the other lower limit controls allow a user to set upper limits and lower limits for alarms for various physiological parameters of patients.

FIG. 8 is a flowchart illustrating an example operation 800 performed by the PMP device 200. As illustrated in the example of FIG. 8, the operation 800 begins when the PMP device 200 displays the settings screen 600 (802). The PMP device 200 then receives workflow selection input from a clinician via the settings screen 600 (804). In response to receiving the workflow selection input, the PMP device 200 determines whether the workflow selection input indicates the monitoring workflow, the triage workflow, or the spot check workflow (806). If the PMP device 200 determines that the workflow selection input indicates the triage workflow, the PMP device 200 performs the steps illustrated in FIG. 9.

If the PMP device 200 determines that the workflow selection input indicates the monitoring workflow, the PMP device 200 displays the patient selection screen 400 (808). The PMP device 200 then receives patient selection input from the clinician via the patient selection screen 400 (810). The patient selection input indicates a patient that the clinician intends to monitor.

Next, the PMP device 200 displays the monitoring workflow home screen 300 (812). In some embodiments, the PMP device 200 displays the monitoring workflow home screen 300 automatically after the clinician selects the monitored patient via the patient selection screen 400. In other embodiments, the PMP device 200 displays the monitoring workflow home screen 300 in response to the clinician selecting the home tab 319a.

Subsequently, the PMP device 200 obtains a measurement of a physiological parameter of the monitored patient (814). The PMP device 200 then displays a representation of the physiological parameter on the monitoring workflow home screen 300 (816). The representation of the physiological parameter is based, at least in part, on the measurement of the physiological parameter. The PMP device 200 then obtains another measurement of the physiological parameter from the same monitored patient (814). In this way, the PMP device 200 enables the clinician to monitor the physiological parameters of the same monitored patient over a period of time.

If, however, the PMP device 200 determines that the workflow selection input indicates the spot check workflow, the PMP device 200 displays the patient selection screen 400 (818). In various embodiments, the PMP device 200 displays the patient selection screen 400 in response to various events. For example, in some embodiments, the PMP device 200 automatically displays the patient selection screen 400 when the PMP device 200 determines that the workflow selection input indicates the spot check workflow. In another example,

the PMP device 200 displays the patient selection screen 400 in response to receiving a selection of the patients tab 319b.

The PMP device 200 then receives patient selection input from the clinician via the patient selection screen 400 (820). The patient selection input indicates a selected patient. The selected patient is a patient on whom the clinician intends to perform a spot check.

Next, the PMP device 200 displays the spot check workflow home screen 330 (822). In various embodiments, the PMP device 200 displays the spot check workflow home screen 330 in response to various events. For example, in some embodiments, the PMP device 200 automatically displays the spot check workflow home screen 330 when the PMP device 200 receives the patient selection input from the clinician. In other embodiments, the PMP device 200 displays the spot check workflow home screen 330 in response to receiving a selection of the home tab 319a.

The PMP device 200 then obtains a measurement of a physiological parameter of the selected patient (824). In response to obtaining the measurement of the physiological parameter of the selected patient, the PMP device 200 displays a representation of the physiological parameter of the selected patient (826). The representation of the physiological parameter is based, at least in part, on the measurement of the physiological parameter.

Subsequently, the PMP device 200 clears the home screen data from the spot check workflow home screen 330 (828). The PMP device 200 then displays the patient selection screen 400 again (818) and receives patient selection input indicating another patient (818), and so on. In this way, the PMP device 200 displays representations of physiological parameters of each patient in a series of patients.

FIG. 9 is a flowchart representing a continuation of the operation 800 in which the workflow selection input indicates the triage workflow. As illustrated in the example of FIG. 9, if the PMP device 200 determines that the workflow selection input indicates the triage workflow, the PMP device 200 displays the triage workflow home screen 360 (900). In various embodiments, the PMP device 200 displays the triage workflow home screen 360 in response to various events. For example, in some embodiments, the PMP device 200 automatically displays the triage workflow home screen 360 when the PMP device 200 determines that the workflow selection input indicates the triage workflow. In another example, the PMP device 200 displays the triage workflow home screen 360 when the PMP device 200 receives a selection of the home tab 319a.

Next, the PMP device 200 obtains a measurement of a physiological parameter of an unidentified patient in a series of patients (902). After obtaining the measurement, the PMP device 200 displays a representation of the physiological parameter on the triage workflow home screen 360 (904). The representation of the physiological parameter is based at least in part on the measurement of the physiological parameter.

Subsequently, the PMP device 200 clears the home screen data from the triage workflow home screen 360 (906). The PMP device 200 then obtains a measurement of the physiological parameter of a different patient in the series of unidentified patients (902), displays a representation of the physiological parameter (904), and so on. In this way, the PMP device 200 displays representations of physiological parameters of each patient in the series of unidentified patients.

FIG. 10 is a flowchart illustrating an operation 1000 performed by the PMP device 200 when the PMP device 200 displays a workflow home screen. In some embodiments, the PMP device 200 can perform the operation 1000 when the

PMP device **200** displays either the monitoring workflow home screen **300**, the spot check workflow home screen **330**, or the triage workflow home screen **360**.

As illustrated in the example of FIG. **10**, the PMP device **200** initially waits to receive an event (**1002**). An event is an occurrence or happening of significance to a task or program.

The PMP device **200** receives a measurement event when the PMP device **200** receives a measurement of a physiological parameter of a patient. If the PMP device **200** receives a measurement event (“YES” of **1004**), the PMP device **200** determines whether the measurement of the physiological parameter is within an alarm range for the physiological parameter (**1006**). If the measurement of the physiological parameter is not within (i.e., is outside) the alarm range for the physiological parameter (“NO” of **1006**), the PMP device **200** displays an alarm message on the workflow home screen (**1008**). In addition, the PMP device **200** emits an alarm sound (**1010**).

The PMP device **200** also generates a representation of the physiological parameter based on the measurement of the physiological parameter (**1012**). The PMP device **200** generates the representation of the physiological parameter regardless of whether the measurement is within the alarm range of the physiological parameter. The PMP device **200** then displays the representation of the physiological parameter on the workflow home screen (**1014**). The PMP device **200** then waits for another event (**1002**).

The PMP device **200** can receive an alarm message selection event when the PMP device **200** displays an alarm message, there are two or more active alarms, and a clinician selects the alarm message. If the PMP device **200** receives an alarm message selection event (“YES” of **1016**), the PMP device **200** modifies the alarm message such that the alarm message describes a next one of the active alarms (**1018**). The PMP device **200** then waits for another event (**1002**).

The PMP device **200** can receive an alarm silence event when an alarm is currently active, the PMP device **200** is emitting an alarm sound, and a clinician selects a control to suspend the alarm sound. If the PMP device **200** receives an alarm silence event (“YES” of **1020**), the PMP device **200** suspends the alarm sound (**1022**). In addition, the PMP device **200** modifies the alarm message to indicate that the alarm has been silenced (**1024**). The PMP device **200** then waits for another event (**1002**).

The PMP device **200** can receive a representation selection event when a clinician selects a representation of a physiological parameter. If the PMP device **200** receives a representation selection event (“YES” of **1026**), the PMP device **200** modifies a parameter reporting area associated with the physiological parameter such that the parameter reporting area contains a different representation of the physiological parameter (**1028**). The PMP device **200** then waits for another event (**1002**).

The PMP device **200** can receive a save selection event when a clinician selects a save button (e.g., save button **316g**) on the workflow home screen. If the PMP device **200** receives a save selection event (“YES” of **1030**), the PMP device **200** saves the home screen data (**1032**). The PMP device **200** then waits for another event (**1002**).

FIG. **11** illustrates example physical components of the PMP device **200**. As illustrated in the example of FIG. **11**, the PMP device **200** include at least one central processing unit (“CPU”) **1108**, a system memory **1112**, and a system bus **1110** that couples the system memory **1112** to the CPU **1108**. The system memory **1112** includes a random access memory (“RAM”) **1118** and a read-only memory (“ROM”) **1120**. A basic input/output system containing the basic routines that

help to transfer information between elements within the PMP device **200**, such as during startup, is stored in the ROM **1120**. The PMP device **200** further includes a mass storage device **1114**. The mass storage device **1114** is able to store software instructions and data.

The mass storage device **1114** is connected to the CPU **1108** through a mass storage controller (not shown) connected to the bus **1110**. The mass storage device **1114** and its associated computer-readable data storage media provide non-volatile, non-transitory storage for the PMP device **200**. Although the description of computer-readable data storage media contained herein refers to a mass storage device, such as a hard disk or CD-ROM drive, it should be appreciated by those skilled in the art that computer-readable data storage media can be any available non-transitory, physical device or article of manufacture from which the PMP device **200** can read data and/or instructions.

Computer-readable data storage media include volatile and non-volatile, removable and non-removable media implemented in any method or technology for storage of information such as computer-readable software instructions, data structures, program modules or other data. Example types of computer-readable data storage media include, but are not limited to, RAM, ROM, EPROM, EEPROM, flash memory or other solid state memory technology, CD-ROMs, digital versatile discs (“DVDs”), other optical storage media, magnetic cassettes, magnetic tape, magnetic disk storage or other magnetic storage devices, or any other medium which can be used to store the desired information and which can be accessed by the PMP device **200**.

According to various embodiments of the invention, the PMP device **200** may operate in a networked environment using logical connections to remote network devices through the network **108**, such as a local network, the Internet, or another type of network. The PMP device **200** connects to the network **108** through a network interface unit **1116** connected to the bus **1110**. It should be appreciated that the network interface unit **1116** may also be utilized to connect to other types of networks and remote computing systems. The PMP device **200** also includes an input/output controller **1122** for receiving and processing input from a number of other devices, including a keyboard, a mouse, a touch user interface display screen, or another type of input device. Similarly, the input/output controller **1122** may provide output to a touch user interface display screen, a printer, or other type of output device.

As mentioned briefly above, the mass storage device **1114** and the RAM **1118** of the PMP device **200** can store software instructions and data. The software instructions include an operating system **1132** suitable for controlling the operation of the PMP device **200**. The mass storage device **1114** and/or the RAM **1118** also store software instructions, that when executed by the CPU **1108**, cause the PMP device **200** to provide the functionality of the PMP device **200** discussed in this document. For example, the mass storage device **1114** and/or the RAM **1118** can store software instructions that, when executed by the CPU **1108**, cause the PMP device to display the monitoring workflow home screen **300** and other screens.

It should be appreciated that various embodiments can be implemented (1) as a sequence of computer implemented acts or program modules running on a computing system and/or (2) as interconnected machine logic circuits or circuit modules within the computing system. The implementation is a matter of choice dependent on the performance requirements of the computing system implementing the invention. Accordingly, logical operations including related algorithms

can be referred to variously as operations, structural devices, acts or modules. It will be recognized by one skilled in the art that these operations, structural devices, acts and modules may be implemented in software, firmware, special purpose digital logic, and any combination thereof without deviating from the spirit and scope of the present invention as recited within the claims set forth herein.

Although the invention has been described in connection with various embodiments, those of ordinary skill in the art will understand that many modifications may be made thereto within the scope of the claims that follow. For example, it should be appreciated that the screens illustrated in this document are merely examples and that in other embodiments equivalent screens can have different contents and appearances. Accordingly, it is not intended that the scope of the invention in any way be limited by the above description, but instead be determined entirely by reference to the claims that follow.

What is claimed is:

1. A physiological measuring platform (PMP) device comprising:

a central processing unit (CPU) that is configured to control operation of the PMP device;

a display screen; and

a set of one or more computer readable data storage media storing software instructions that, when executed by the CPU, cause the PMP device to:

obtain a series of measurements of a physiological parameter of a monitored patient when the PMP device is operating within a monitoring workflow;

display, on the display screen, a monitoring workflow home screen when the PMP device is operating within the monitoring workflow, the monitoring workflow home screen containing at least one non-compressed monitoring parameter frame including:

a first representation of the physiological parameter of the monitored patient and a first representation label, the first representation based on a measurement in the series of measurements; and

at least one user-selectable control;

obtain a measurement of the physiological parameter of each patient in a series of patients when the PMP device is operating within a non-monitoring workflow;

display, on the display screen, a non-monitoring workflow home screen when the PMP device is operating within the non-monitoring workflow, the non-monitoring workflow home screen containing at least one non-compressed non-monitoring parameter frame including a second representation of the physiological parameter of a given patient in the series of patients and a second representation label, the second representation based on the measurement of the physiological parameter of the given patient in the series of patients, the monitoring workflow home screen being different than the non-monitoring workflow home screen;

display, on the display screen, an alarms screen when the PMP device is operating within the monitoring workflow, the alarms screen enabling a user to configure global alarm settings and parameter-specific alarm settings, wherein the alarms screen includes a compressed parameter reporting area, the compressed parameter reporting area containing a compressed parameter frame containing a third representation of live physiological parameters of the given patient;

display, on the display screen, an alarm message describing a first alarm when the first alarm is active, the first alarm being active when a given measurement of the physiological parameter is outside an alarm range for the physiological parameter; and

when a second alarm is active at the same time as the first alarm, modify, on the display screen, the alarm message to describe the second alarm when a user selects the alarm message, the second alarm being active when a second given measurement of a second physiological parameter is outside an alarm range for the second physiological parameter,

wherein the user-selectable control is not included in the compressed parameter frame in the alarms screen.

2. The PMP device of claim 1, wherein the software instructions, when executed by the CPU, cause the PMP device to display a settings screen on the display screen, the settings screen enabling the user to select a workflow within which the PMP device is to operate.

3. The PMP device of claim 1,

wherein the alarms screen includes an upper limit control and a lower limit control;

wherein the upper limit control includes a numerical portion, an up button, and a down button, the numerical portion of the upper limit control specifying an upper limit for an alarm associated with the physiological parameter; and

wherein the lower limit control includes a numerical portion, an up button, and a down button, the numerical portion of the lower limit control specifying a lower limit for the alarm;

wherein the software instructions, when executed by the CPU, cause the PMP device to:

incrementally increase the upper limit for the alarm when the user selects the up button of the upper limit control;

incrementally decrease the upper limit for the alarm when the user selects the down button of the upper limit control;

enable the user to manually enter a value for the upper limit for the alarm when the user selects the numerical portion of the upper limit control;

incrementally increase the lower limit for the alarm when the user selects the up button of the lower limit control;

incrementally decrease the lower limit for the alarm when the user selects the down button of the lower limit control; and

enable the user to manually enter a value for the lower limit for the alarm when the user selects the numerical portion of the lower limit control;

wherein values for the upper and lower limits selected by the user are validated so that values falling outside of predefined thresholds are disallowed.

4. The PMP device of claim 1, wherein the software instructions, when executed by the CPU, further cause the PMP device to:

emit an alarm sound when an alarm is active, the alarm being active when a given measurement of the physiological parameter is outside the alarm range for the physiological parameter;

suspend emitting the alarm sound for a given time period in response to user input;

extend the given time period when the PMP device receives additional user input before the given time period expires; and

resume emitting the alarm sound when the given time period expires.

5. The PMP device of claim 1,
wherein the software instructions, when executed by the CPU, cause the PMP device to display an alternate representation of the physiological parameter in a parameter reporting frame instead of the initial representation of the physiological parameter when the user selects the parameter reporting frame.

6. The PMP device of claim 1,
wherein the physiological parameter is body temperature; wherein the monitoring workflow home screen includes a temperature frame that contains data representing the physiological parameter and also contains a thermometry location control that indicates a location of a thermometer on the body of the monitored patient, the series of measurements of the physiological parameter based on data received from the thermometer; and
wherein the software instructions, when executed by the CPU, cause the PMP device to update the thermometry location control such that the thermometry location control indicates a different location on the body of the monitored patient when the user selects the thermometry location control.

7. The PMP device of claim 1, wherein said non-monitoring workflow is a spot check workflow.

8. The PMP device of claim 1,
wherein said non-monitoring workflow is a triage workflow,
wherein when the PMP device is operating in the triage workflow, each patient in the series of patients is unidentified,
wherein the software instructions, when executed by the CPU, further cause the PMP device to display, on the display screen, a spot check workflow home screen when the PMP device is operating within a spot check workflow, and
wherein when the PMP device is operating within the spot check workflow, the PMP device obtains a measurement of the physiological parameter from a series of one or more previously-identified patients.

9. The PMP device of claim 1,
wherein the monitoring workflow home screen and the non-monitoring workflow home screen contain at least one common user interface element, the common user interface element having a same appearance but not performing a same function in the monitoring workflow home screen and the non-monitoring workflow home screen.

10. The PMP device of claim 9,
wherein the common user interface element is a save button, and
wherein when the PMP device is operating in the monitoring workflow, the software instructions, when executed by the CPU, cause the PMP device to save a first patient reading in response to selection of the save button; and
wherein when the PMP device is operating in the non-monitoring workflow, the software instructions, when executed by the CPU, cause the PMP device to perform the following in response to selection of the save button: save a second patient reading; and clear the non-monitoring workflow home screen.

11. The PMP device of claim 10, wherein selection of the save button communicates the first patient reading to another computing node when the PMP device is operating within the monitoring workflow.

12. The PMP device of claim 1, wherein the software instructions, when executed by the CPU, cause the PMP device to:
display a settings screen, the settings screen containing a set of interval program selection controls; and
record measurements of the physiological parameter at intervals specified by an interval program for a duration specified by the interval program, the interval program corresponding to a selected one of the interval program selection controls;
wherein the device is configured to allow the user to manually enter interval values for the interval program.

13. The PMP device of claim 1, wherein the device is programmed to allow the user to rename at least one of the monitoring workflow and the non-monitoring workflow for display on the device.

14. A method comprising:
obtaining, by a physiological measuring platform (PMP) device, a measurement of a physiological parameter of a given patient;
displaying, on a display screen, a monitoring workflow home screen when the PMP device is operating within a monitoring workflow, the monitoring workflow home screen containing at least one parameter reporting area containing a representation of the physiological parameter of the given patient;
displaying an alarms screen on a display screen of the PMP device, the alarms screen enabling a user of the PMP device to configure global alarm settings and parameter-specific alarm settings
displaying, on the display screen, an alarm message describing a first alarm when the first alarm is active, the first alarm being active when a given measurement of the physiological parameter is outside an alarm range for the physiological parameter;
when a second alarm is active at the same time as the first alarm, modifying, on the display screen, the alarm message to describe the second alarm when the second alarm is active and when the user selects the alarm message, the second alarm being active when a second given measurement of a second physiological parameter is outside an alarm range for the second physiological parameter, and
wherein the alarms screen includes at least one compressed parameter reporting area, the at least one compressed parameter reporting area containing a compressed parameter frame containing another representation of the physiological parameter of the given patient, the compressed parameter reporting area appearing compressed relative to the parameter reporting area,
wherein the compressed parameter frame includes alarm limits associated with the physiological parameter, and
wherein at least one of the at least one parameter reporting areas includes a user-selectable control configured to initiate or stop a measurement by a peripheral device operatively connected to the PMP device, and
wherein the user-selectable control is not included in the compressed parameter reporting area in the alarms screen.

15. The method of claim 14, wherein the parameter-specific alarm settings include alarm limits.

16. The method of claim 14,
wherein the alarms screen includes an upper limit control and a lower limit control;
wherein the upper limit control includes a numerical portion, an up button, and a down button, the numerical

portion of the upper limit control specifying an upper limit for an alarm associated with the physiological parameter; and
 wherein the lower limit control includes a numerical portion, an up button, and a down button, the numerical portion of the lower limit control specifying a lower limit for the alarm;
 wherein the method further comprises:
 incrementally increasing the upper limit for the alarm when the user selects the up button of the upper limit control;
 incrementally decreasing the upper limit for the alarm when the user selects the down button of the upper limit control;
 enabling the user to manually enter a value for the upper limit for the alarm when the user selects the numerical portion of the upper limit control;
 incrementally increasing the lower limit for the alarm when the user selects the up button of the lower limit control;
 incrementally decreasing the lower limit for the alarm when the user selects the down button of the lower limit control;
 enabling the user to manually enter a value for the lower limit for the alarm when the user selects the numerical portion of the lower limit control.

17. A non-transitory computer-readable storage medium comprising software instructions that, when executed, cause a physiological measurement platform (PMP) device to:
 display, on a display screen, a settings screen that enables a user to select a workflow within which the PMP device is to operate;
 when the user selects a monitoring workflow from the settings screen:
 obtain a series of measurements of a physiological parameter of a monitored patient;
 display a monitoring workflow home screen on the display screen, the monitoring workflow home screen containing a first representation of the physiological parameter of the monitored patient and a user-selectable control configured to control a peripheral device operatively connected to the PMP device, the first representation of the physiological parameter based on at least one measurement in the series of measurements;
 display, on the display screen, an alarms screen that enables the user to configure global alarm settings and parameter-specific alarm settings,
 wherein the alarms screen includes a compressed parameter reporting area, the compressed parameter reporting area containing a compressed parameter frame containing a second representation of the physiological parameter of the monitored patient, the second representation of the physiological parameter being smaller than the first representation of the

physiological parameter, and the second representation of the physiological parameter including alarm limits associated with the physiological parameter, and
 when the user selects a spot check workflow from the settings screen:
 obtain a measurement of the physiological parameter of each patient in a series of previously identified patients; and
 display a spot check workflow home screen on the display screen, the spot check workflow home screen containing a third representation of the physiological parameter of a given patient in the series of previously identified patients, the third representation of the physiological parameter based on the measurement of the physiological parameter of the given patient in the series of previously identified patients;
 when the user selects a triage workflow from the settings screen:
 obtain a measurement of the physiological parameter of each patient in a series of unidentified patients; and
 display a triage workflow home screen on the display screen, the triage workflow home screen containing a fourth representation of the physiological parameter of the given patient in the series of unidentified patients, the fourth representation of the physiological parameter based on the measurement of the given patient in the series of unidentified patients,
 wherein the monitoring workflow home screen, the spot check workflow home screen, and the triage workflow home screen are each different from one another; and
 wherein the user-selectable control is not included in the compressed parameter frame in the alarms screen,
 display, on the display screen, an alarm message describing a first alarm when the first alarm is active, the first alarm being active when a given measurement of the physiological parameter is outside an alarm range for the physiological parameter;
 when a second alarm is active at the same time as the first alarm, modify, on the display screen, the alarm message to describe a second alarm when the second alarm is active and when the user selects the alarm message, the second alarm being active when a second given measurement of a second physiological parameter is outside an alarm range for the second physiological parameter;
 emit an alarm sound when the first alarm is active;
 suspend emitting the alarm sound for a given time period in response to user input;
 extend the given time period when the PMP device receives additional user input before the given time period expires; and
 resume emitting the alarm sound when the given time period expires.

* * * * *

专利名称(译)	支持多种工作流程的生理参数测量平台装置		
公开(公告)号	US9265429	公开(公告)日	2016-02-23
申请号	US12/751579	申请日	2010-03-31
[标]申请(专利权)人(译)	ST PIERRE SHAWNÇ 安埃丝特 贝利埃里克J GARRANT MICHAELD 罗杰W的坟墓 Kilian珍妮佛M MYERS托马斯		
申请(专利权)人(译)	ST. PIERRE , SHAWN C. AHN , ESTHER N. BAILEY , ERIC J. GARRANT , MICHAEL D. 格雷夫斯ROGER W. KILIAN , JENNIFER M. MYERS , 托马斯A.		
当前申请(专利权)人(译)	伟伦 , INC.		
[标]发明人	ST PIERRE SHAWN C AHN ESTHER N BAILEY ERIC J GARRANT MICHAEL D GRAVES ROGER W KILIAN JENNIFER M MYERS THOMAS A		
发明人	ST. PIERRE, SHAWN C. AHN, ESTHER N. BAILEY, ERIC J. GARRANT, MICHAEL D. GRAVES, ROGER W. KILIAN, JENNIFER M. MYERS, THOMAS A.		
IPC分类号	G08B1/08 A61B5/205 A61B5/00 G06Q10/06 G06Q50/22 G06F19/00 A61B5/0205		
CPC分类号	A61B5/02055 A61B5/7445 A61B5/7475 G06F19/3406 G06F19/3487 G06Q10/06 G06Q50/22 A61B5/746 G06F19/00 G16H15/00 G16H40/63		
代理机构(译)	商户德律P.C.		
优先权	61/243872 2009-09-18 US		
其他公开文献	US20110071420A1		
外部链接	Espacenet USPTO		

摘要(译)

当设备在监测工作流程内操作时，设备获得被监测患者的生理参数的一系列测量。当设备在监控工作流程内运行时，设备会显示监控工作流程主屏幕。监视工作流程主屏幕包含被监视患者的生理参数的表示。此外，当设备在非监测工作流程内操作时，该设备获

得一系列患者中每个患者的生理参数的测量值。当设备在非监控工作流程内运行时，设备会显示非监控工作流程主屏幕。非监测工作流程主屏幕包含一系列患者中给定患者的生理参数的表示。监视工作流程主屏幕与非监视工作流程主屏幕不同。

200

