



(12) **EUROPEAN PATENT SPECIFICATION**

(45) Date of publication and mention of the grant of the patent:
09.06.2004 Bulletin 2004/24

(51) Int Cl.⁷: **G01N 33/68**

(21) Application number: **02721912.0**

(86) International application number:
PCT/CA2002/000563

(22) Date of filing: **24.04.2002**

(87) International publication number:
WO 2002/088706 (07.11.2002 Gazette 2002/45)

(54) **PROCESS FOR DIFFERENTIAL DIAGNOSIS OF ALZHEIMER'S DEMENTIA AND DEVICE THEREFOR**

PROZESS UND MITTEL FÜR DIE DIFFERENTIALDIAGNOSE DER DEMENZ VOM ALZHEIMER TYP

PROCEDE DE DIAGNOSTIC DIFFERENTIEL D'UNE DEMENCE DE TYPE ALZHEIMER ET DISPOSITIF ASSOCIE

(84) Designated Contracting States:
AT BE CH CY DE DK ES FI FR GB GR IE IT LI LU MC NL PT SE TR
Designated Extension States:
AL LT LV MK RO SI

(30) Priority: **25.04.2001 US 842079**
04.10.2001 US 971740

(43) Date of publication of application:
21.01.2004 Bulletin 2004/04

(73) Proprietor: **Syn.X Pharma, Inc.**
Toronto, Ontario M9W 1E7 (CA)

(72) Inventors:
• **JACKOWSKI, George**
Kettleby, Ontario L0G 1J0 (CA)
• **TAKAHASHI, Miyoko**
North York, Ontario M2N 1B8 (CA)

(74) Representative: **Jones, Helen Marjorie Meredith Gill Jennings & Every,**
Broadgate House,
7 Eldon Street
London EC2M 7LH (GB)

(56) References cited:
WO-A-00/52476 **US-A- 5 445 937**

- **TAKAHASHI MIYOKO ET AL: "Immunoassay for serum glutamine synthetase in serum: Development, reference values, and preliminary study in dementias." CLINICAL CHEMISTRY, vol. 48, no. 2, February 2002 (2002-02), pages 375-378, XP001121876 ISSN: 0009-9147**
- **TUMANI HAYRETTIN ET AL: "Glutamine synthetase in cerebrospinal fluid, serum, and brain: A diagnostic marker for Alzheimer disease?" ARCHIVES OF NEUROLOGY, vol. 56, no. 10, October 1999 (1999-10), pages 1241-1246, XP002227900 ISSN: 0003-9942 cited in the application**
- **BOKSHA IRINA S ET AL: "Glutamine synthetase and glutamine synthetase-like protein from human brain: Purification and comparative characterization." JOURNAL OF NEUROCHEMISTRY, vol. 75, no. 6, December 2000 (2000-12), pages 2574-2582, XP002227901 ISSN: 0022-3042**
- **GUNNERSEN D ET AL: "DETECTION OF GLUTAMINE SYNTHETASE IN THE CEREBROSPINAL FLUID OF ALZHEIMER DISEASED PATIENTS: A POTENTIAL DIAGNOSTIC BIOCHEMICAL MARKER" PROCEEDINGS OF THE NATIONAL ACADEMY OF SCIENCES OF USA, NATIONAL ACADEMY OF SCIENCE. WASHINGTON, US, vol. 89, no. 24, 15 December 1992 (1992-12-15), pages 11949-11953, XP000384420 ISSN: 0027-8424 cited in the application**

Note: Within nine months from the publication of the mention of the grant of the European patent, any person may give notice to the European Patent Office of opposition to the European patent granted. Notice of opposition shall be filed in a written reasoned statement. It shall not be deemed to have been filed until the opposition fee has been paid. (Art. 99(1) European Patent Convention).

- **TUMANI H ET AL: "Purification and immunocharacterization of human brain glutamine synthetase and its detection in cerebrospinal fluid and serum by a sandwich enzyme immunoassay" JOURNAL OF IMMUNOLOGICAL METHODS, ELSEVIER SCIENCE PUBLISHERS B.V.,AMSTERDAM, NL, vol. 188, no. 1, 15 December 1995 (1995-12-15), pages 155-163, XP004020947 ISSN: 0022-1759**
- **LISTRON CHAD D ET AL: "Expression, purification, and characterization of recombinant human glutamine synthetase." BIOCHEMICAL JOURNAL, vol. 328, no. 1, 15 November 1997 (1997-11-15), pages 159-163, XP002227902 ISSN: 0264-6021**
- **KATO K ET AL: "Concentrations of several proteins characteristic of nervous tissue in cerebral cortex of patients with Alzheimer's disease." JOURNAL OF MOLECULAR NEUROSCIENCE: MN. UNITED STATES 1991, vol. 3, no. 2, 1991, pages 95-99, XP009004174 ISSN: 0895-8696 cited in the application**

Description

FIELD OF THE INVENTION

[0001] The present invention relates to a method for the diagnosis of Alzheimer' dementia (AD). The invention particularly relates to a process for quantifying the presence of at least one biochemical marker associated with Alzheimer's dementia. More particularly, the invention relates to a point-of-care immunoassay which utilizes unique antibodies to enable the differential diagnosis of Alzheimer's versus non-Alzheimer's forms of dementia.

BACKGROUND OF THE INVENTION

[0002] Alzheimer's disease, also referred to as Alzheimer's dementia or AD is a progressive neurodegenerative disorder that causes memory loss and serious mental deterioration. Diagnosticians have long sought a means to definitively identify AD during the lifetime of demented patients, as opposed to histopathological examination of brain tissue, which is the only present means available for rendering an ultimate diagnosis of AD. AD is the most common form of dementia, accounting for more than half of all dementias and affecting as many as 4 million Americans and nearly 15 million people worldwide. Dementia may start with slight memory loss and confusion, but advances with time reaching severe impairment of intellectual and social abilities. At age 65, the community prevalence of AD is between 1-2 %. By age 75, the figure rises to 7 %, and by age 85 it is 18 %. The prevalence of dementia in all individuals over age 65 is 8 %. Of those residing in institutions, the prevalence is about 50 %, at any age.

[0003] The social impact of this disease is enormous, caused by the burden placed on caregivers, particularly in the latter stages of the disease. The substantial economic costs are largely related to supportive care and institutional admission. The rapidly increasing proportion of elderly people in society means that the number of individuals affected with AD will grow dramatically, therefore finding an early accurate diagnosis and a cure for AD is becoming an issue of major importance world wide.

[0004] When an individual is suspected of AD, several recommended tests are performed : (1) Mini Mental State Examination (MMSE) - an office-based psychometric test in the form of a Functional Assessment Questionnaire (FAQ) to examine the scale for functional autonomy, (2) Laboratory tests - complete blood count, measurement of thyroid stimulating hormone, serum electrolytes, serum calcium and serum glucose levels, (3) Neuroimaging - most commonly used is computed tomography (CT) which has a role in detecting certain causes of dementia such as vascular dementia (VaD), tumor, normal pressure hydrocephalus or subdural hematoma. However, neuroimaging is less effective in dis-

tinguishing AD or other cortical dementias from normal aging. In primary care settings, some suggest that CT could be limited to atypical cases, but others recommend routine scanning. Magnetic resonance imaging (MRI) currently offers no advantage over CT in most cases of dementia.

[0005] While Alzheimer's is the most common form of dementia, accounting for at least 60 % of cases, diagnostic procedures for determining the exact cause of dementia, among more than 80 different species, is difficult at best. Furthermore, the currently performed tests are inadequate in differentiating AD from other types of dementia.

[0006] In comparison to other disease areas, the field of dementia raises questions concerning the value of diagnosis, since there is currently no cure or effective therapy available. In dementia, as in all other branches of medicine, the certainty of a diagnosis has an important impact on the management of the patient. While AD cannot be cured at present time, there is symptomatic treatment available and the first drugs (acetylcholinesterase inhibitors) for the temporary improvement of cognition and behavior are now licensed by the U.S. Food and Drug Administration. Other drugs are at different stages of clinical trials: (1) Drugs to prevent decline in AD - DESFERRIOXAMINE, ALCAR, anti-inflammatory drugs, antioxidants, estrogen, (2) Neurotrophic Factors : NGF, (3) Vaccine : the recent most exciting report by Schenk et al. (Nature 1999;400:173-7) raises the hope of a vaccine for AD.

[0007] The specificity of the various therapies thus require sophisticated diagnostic methodologies, having a high degree of sensitivity for AD, in order to insure their success.

[0008] Currently there are a multitude of tests available which aid in the diagnosis of AD. However, the only true existing diagnosis is made by pathologic examination of postmortem brain tissue in conjunction with a clinical history of dementia. This diagnosis is based on the presence in brain tissue of neurofibrillary tangles and of neuritic (senile) plaques, which have been correlated with clinical dementia. Neuritic plaques are made up of a normally harmless protein called amyloid-beta. Before neurons begin to die and symptoms develop, plaque deposits form between neurons early on in the disease process. The neurofibrillary tangles are interneuronal aggregates composed of normal and paired helical filaments and presumably consist of several different proteins. The internal support structure for brain neurons depends on the normal functioning of a protein called tau. In Alzheimer's disease, threads of tau protein undergo alterations that cause them to become twisted. The neurohistopathologic identification and counting of neuritic plaques and neurofibrillary tangles requires staining and microscopic examination of several brain sections. However, the results of this methodology can widely vary and is time-consuming and labor-intensive.

[0009] Given the ability of both current and prospec-

tive pharmacological therapies to forestall and/or reverse the onset and/or progress of Alzheimer's dementia, an early diagnosis of AD will assist to better manage the care of patients. There are many cases where non-AD dementia could be confused with AD dementia. Such examples include small, undetected strokes which temporarily interrupt blood flow to the brain. Clinically depressed patients or those with Parkinson's disease can also experience lapses in memory. Many older people are on a variety of medications which as a side effect may, alone or in conjunction, impair their ability to perform cognitive tasks.

[0010] Thus, if diagnostic techniques for the early differentiation of AD could be provided, physician's would achieve an enhanced ability to prescribe appropriate therapeutic intervention at an early stage in the pathogenesis of this disease.

[0011] Various biochemical markers for AD are known and analytical techniques for the determination of such markers have been described in the art. As used herein the term "marker" "biochemical marker" or "marker protein" refers to any enzyme, protein, polypeptide, peptide, isomeric form thereof, immunologically detectable fragments thereof, or other molecule that is released from the brain during the course of AD pathogenesis. Such markers include, but are not limited to, any unique proteins or isoforms thereof that are particularly associated with the brain.

[0012] Glutamine synthetase (GS) is recognized as an astrocyte-specific enzyme involved in the regulation of ammonia and glutamate metabolism that is over-expressed following brain injury (Norenberg and Martinex-Hernandez, *Brain Res* 1979;161:303). A few studies on the clinical role of glutamine synthetase have been reported: Gunnensen and Haley (*Proc Natl Acad Sci USA* 1992;89:11949) found monomeric GS protein in 38 of 39 AD cerebrospinal fluid (CSF) samples, Tumani et al. (*Arch Neurol* 1999;56(10):1241) describe that the concentration of GS in lumbar CSF of patients with AD is increased significantly but nonspecifically (i.e. also increased in VaD, schizophrenia and ALS). On p.1244, the left-hand column, Tumani states that GS was not found in serum.

[0013] Neuron-specific gamma-enolase (NSE_{γγ}) and S100B proteins, abundant in the brain, are also useful markers for assessing the extent of brain damage: NSE_{γγ} for neuronal damage and S100B for astrocyte damage. Concentrations of NSE and S100B proteins from cerebrocortical regions have been examined by means of enzyme linked immunosorbent assay (ELISA). The levels of these proteins in frontal cortex of AD patients were found to be significantly elevated (Kato et al. *J Mol Neurosci*, 1991;3(2):95). Activated astrocytes over-expressing S100B have been intimately associated with the neuritic β-amyloid plaques of AD (Sheng et al. *J Neurosci Res*, 1994;39:398, Mrak et al. *JNeuropathol Exp Neurol* 1996;55:273).

[0014] There are a number of different potential uses

for biomarkers in AD evaluation, and each use could involve a different marker or set of markers. Such uses may include, but are not limited to, the use of a marker to distinguish AD from other causes of dementia; distinguishing dementia from the non-pathological effects of aging; monitoring the progress of the disease after clinical symptoms become apparent; utilization of a surrogate to monitor the efficacy of the forthcoming therapies for AD; and isolating markers which have utility as risk assessment factors for AD; and identifying both the earliest biological changes occurring in the brain and other changes that occur as the disease progresses. Ideally, it would be preferable to isolate a single marker to fulfill all requirements with a high degree of sensitivity and specificity, however this may be an unreasonable goal. Any individual marker needs to be assessed by sensitivity, specificity, reliability and validity for the type of clinical situation to which it is meant to apply. A marker which is poor at distinguishing AD from other causes of dementia, could nevertheless be an excellent marker for monitoring the progression of the disease process or the response to therapy.

[0015] With regard to diagnostic devices, the clinical evaluation and use of point-of-care tests utilizing biological markers are valuable tools for evaluating risk, monitoring disease progression and guiding therapeutic interventions. The advantages which flow from the use of biological markers as diagnostic tools include strengthening the certainty of the clinical diagnosis, distinguishing AD from other causes of dementia, and quantifying the severity of the disease and rate of progression. In addition, tests using biological markers should be rapid, non-invasive, simple to perform and inexpensive.

[0016] What is lacking in the art is a relatively non-invasive method and device therefore effective for definitively diagnosing Alzheimer's dementia in living patients. Additionally, a definitive method of assessing the risk of developing AD is greatly needed.

DESCRIPTION OF THE PRIOR ART

[0017] U.S. Patent No. 5,445,937 to Haley teaches a method for the diagnosis of Alzheimer's disease, as well as a means for the diagnosis and differentiation of other diseases. This is done by use of a disease-specific biochemical marker, glutamine synthetase (GS) and its respective photoaffinity label, or labeled antibody specific for GS, at GS's binding site. The '937 patent focuses on examining cerebrospinal fluid (CSF) to detect the presence of a photoaffinity labeled, or labeled antibody, nucleotide binding protein and correlates the subsequent level with the presence of AD. Haley teaches a variety of immunoassay techniques to accomplish such a method. While Haley hypothesizes regarding the prophetic utility of diagnostic methods which use blood as a sample and further suggests that a monoclonal and/or polyclonal antibody immunoassay could be developed, he nevertheless fails to reduce either to practice. Thus, the

'937 patent is only useful in teaching a diagnostic test which utilizes cerebrospinal fluid. Obtaining a sample of cerebrospinal fluid entails invasive techniques quite uncomfortable to a patient and requires a lengthy period of time to accomplish. In addition, the only polyclonal and/or monoclonal antibodies suggested by Haley are those having a specificity to sheep brain GS as opposed to a human recombinant form of GS as instantly disclosed herein.

[0018] In U.S. Pat. No. 5,508,167, Roses et al. describe methods for diagnosing AD involving the detection of an apolipoprotein E type 4 (ApoE4) isoform or DNA encoding ApoE4. The methods can use blood samples and are analyzed by an immunochemical assay. The blood sample is optionally combined with a reducing agent to reduce the disulfide bond in cysteine residues to the corresponding reactive sulfhydryl groups. Roses et al. further describes a kit for detection of the ApoE4 isoform. The test is based on the differences in the amino acid sequences of the three major ApoE isoforms. The test is not specific for human GS nor does it have sensitivity in differentially diagnosing AD versus non-AD dementia.

[0019] Tumani et al. (Arch. Neurol., (1999) 56, pp1241-1246) examine the levels of GS in CSF and the examination of serum in order to determine whether GS is a useful biochemical marker in the diagnosis of AD. Analysis is by an ELISA utilizing a biotin-labeled monoclonal antibody directed against sheep brain GS. Normal ranges of GS concentration are reported as 4 pg/mL in human CSF and 36 pg/mL in human serum. The samples of CSF in AD patients are elevated with a mean level of GS concentration of 20 ± 12 pg/mL, with ALS patients at 13 ± 13 pg/mL, and vascular dementia (VaD) patients at a mean elevated level 13 ± 7 pg/mL. Vascular dementia and ALS patients show a slightly lower increase. Patients with AD are measured with mean levels of 11 ± 53 pg/mL in serum. However, patients with amyotrophic lateral sclerosis (ALS) and vascular dementia also present with mean elevated levels at 116 ± 62 pg/mL and 72 ± 59 pg/mL in serum, respectively. Thus, no definitive diagnosis regarding AD dementia or differential diagnosis between AD versus non-AD dementia could be elucidated from these assays.

[0020] Gunnarsen and Haley (Proc. Natl. Acad. Sci. (1992) 89, pp 11949-11953) provides evidence of GS detected in CSF of patients with AD but not in that of healthy control subjects or controls with other diseases. The other diseases under consideration are epilepsy, ALS, and Parkinson's. Patients with ALS or Pick's disease in addition to AD do show positive results, meanwhile ALS patients do not show positive results, indicating GS is specific to AD. As with other publications, antibodies raised against non-human GS are utilized for detection of GS.

[0021] Generally, most scientific papers tend to focus on the peptide, β -amyloid, since it is postulated to be a major determinant of AD. This is supported by the ob-

ervation that certain forms of familial AD mutations result in the over production of β -amyloid, particularly the longer form (1-42) which aggregates more readily than the shorter form. Hensley et al. (Proc. Natl. Acad. Sci., (1994), 91, pp3270-3274) examine the neurotoxicity based on free radical generation by the peptide β -amyloid in its aggregation state. Several synthetic fragments of the peptide are tested for resulting neurotoxicity. Based on the fact that oxygen seems to be a requirement for radical generation and glutamate synthetase and creatine kinase enzymes are oxidation-sensitive biomarkers, the inactivation of these enzymes are utilized as indicators of active attack on biological molecules by these fragmented β -amyloid aggregates.

SUMMARY OF THE INVENTION

[0022] The present invention relates to a method for the diagnosis of Alzheimer's dementia (AD), particularly to a method for diagnosing and differentiating Alzheimer's dementia from other forms of dementia by testing for the presence of specific biochemical markers for Alzheimer's disease in bodily fluids, particularly in blood, blood products, urine, saliva and the like. The invention further relates to a process for quantifying the presence of at least one biochemical marker associated with Alzheimer's dementia. More particularly, the invention relates to a point-of-care immunoassay which utilizes unique antibodies to enable the differential diagnosis of Alzheimer's versus non-Alzheimer's forms of dementia.

[0023] According to the invention there is provided a new diagnostic kit for diagnosing or monitoring the progression of Alzheimer's dementia comprising:

at least one monoclonal antibody which is specific for monomeric brain associated human glutamine synthetase having a weight of about 44 kDa, wherein said antibody is selected from the group consisting of the antibody produced by the clone deposited with the ATCC as Accession Number PTA-3339 and the antibody produced by the clone deposited with the ATCC as Accession Number PTA-3340.

[0024] There is also provided a method of diagnosing or monitoring the progression of Alzheimer's dementia comprising:

contacting a sample of blood, a blood product or cerebrospinal fluid from a patient suspected of suffering from Alzheimer's dementia with monoclonal antibody, wherein said antibody is specific for monomeric brain associated human glutamine synthetase having a weight of about 44 kDa and is selected from the group consisting of the antibody produced by the clone deposited with the ATCC as Accession Number PTA-3339 and the antibody produced by the clone deposited with the ATCC as Accession Number PTA-3340.;

detecting specific binding of the antibody to the sample;

wherein detection of an elevated level of said glutamine synthetase in the sample from the patient relative to control samples indicates a diagnosis of Alzheimer's dementia, or wherein the extent of elevation of said glutamine synthetase in the sample from the patient relative to control samples correlates to the progression of Alzheimer's dementia in the patient.

[0025] There is also provided an isolated monoclonal antibody produced by the clone deposited with the ATCC as Accession Number PTA-3339 and an isolated monoclonal antibody produced by the clone deposited with the ATCC as Accession Number PTA-3340 and the clones themselves.

[0026] The present invention relates to methods and ELISA system for diagnosing, subtyping and monitoring Alzheimer's disease. The invention is based on the discovery that GS protein is released from the brain and can be detected in body fluids outside the brain.

[0027] The generation and purification of recombinant human GS are described. These GS proteins may be used to generate monoclonal or polyclonal antibodies that, in turn, can be used in immunoassays, wherein they enter into an immunoreaction which can be monitored and/or quantified to detect circulating GS protein in suspected individuals. Alternatively, the GS protein themselves may be used in immunoassays to detect circulating autoantibodies in such individuals. The occurrence of Alzheimer's dementia is characterised by the recognition of levels of a particular biochemical marker in bodily fluid, said levels correlating to the manifestation of Alzheimer's dementia symptoms as quantified by MMSE testing. As a risk assessment test, the recognition of levels of such markers which are indicative of the development of Alzheimer's dementia further augments the diagnostic capability afforded to the skilled practitioner.

[0028] The instant invention provides a relatively non-invasive and highly sensitive method for the definitive diagnosis of Alzheimer's disease.

[0029] The invention also provides a method which includes analysis of at least one body fluid of a patient to determine the presence of at least one marker indicative of AD vs non-AD dementia.

[0030] The instant invention also provides antibodies specific to neuronal related proteins, as identified by the method of the current invention.

[0031] The immunoassay is effective for the recognition of neuronal specific proteins in one or more human bodily fluids. The instant invention utilises a purified monoclonal antibody specific for human glutamine synthetase. The test kit for the diagnosis of AD is suitable as a non-invasive point-of-care test which can be performed utilizing a sample comprising blood or any blood product.

[0032] Other objects and advantages of this invention

will become apparent from the following description taken in conjunction with the accompanying figures wherein are set forth, by way of illustration and example, certain embodiments of this invention. The figures constitute a part of this specification and include exemplary embodiments of the present invention and illustrate various objects and features thereof.

BRIEF DESCRIPTION OF THE FIGURES

[0033]

Figure 1 is a comparison of statistically significant values for GS, NSE and S100 in the blood of a cohort of patients clinically evaluated as suffering from a form of Alzheimer's dementia;

Figure 2 is a comparison of statistically significant values for GS, NSE and S100 in the blood of a cohort of patients clinically evaluated as suffering from a form of non-Alzheimer's dementia;

Figure 3 depict box whisker plots of GS relative to age in healthy individuals.

DETAILS DESCRIPTION OF THE INVENTION

[0034] The marker which is analyzed according to the method of the invention is released into the circulation and may be present in the blood or in any blood product, for example plasma, serum, cytolized blood, e.g. by treatment with hypotonic buffer or detergents and dilutions and preparations thereof, and other body fluids, e.g. CSF, saliva, urine, lymph, and the like. In another preferred embodiment the concentration of the markers in CSF may be measured.

[0035] The terms "above normal" and "above cut-off" are used to refer to a level of a marker that is greater than the level of the marker observed in normal individuals who are not undergoing a cerebral event, i.e. brain degeneration or atrophy. For some markers, very low levels of the marker may be present normally in an individual's blood. For other markers analyzed according to the invention, detectable levels may be present normally in blood. Therefore, these terms indicate a level that is significantly above the normal level found in healthy individuals. The term "significantly" or "statistically significant" refers to statistical significance and generally means a two standard deviation (SD) above normal, or higher, concentration of the marker. The assay method by which the analysis for the marker protein is carried out must be sensitive to be able to detect the level of the marker which is present over the concentration range of interest and also must be highly specific.

[0036] The marker and two other markers have been compared to evaluate their value as a diagnostic tool are proteins which are released by specific cells in the brain as these cells become damaged during a cerebral event. These markers, e.g. proteins, can be either in their native form or they may be immunologically detectable

fragments of the proteins resulting, for example, by proteolytic breakdown. By "immunologically detectable" is meant that the marker fragments contain an epitope which is specifically recognized by antibody reagents used in the assay.

[0037] The markers are cell type specific. The enzyme enolase catalyzes the interconversion of 2-phosphoglycerate and phosphoenolpyruvate in the glycolytic pathway. The enzyme exists in three isoproteins each the product of a separate gene. The gene loci have been designated ENO1, ENO2 and ENO3. The gene product of ENO1 is the non-neuronal enolase (NNE or α), which is widely distributed in various mammalian tissues. The gene product of ENO2 is the muscle specific enolase (MSE or β) which is localized mainly in the cardiac and striated muscle, while the product of the ENO3 gene is the neuronal specific enolase (NSE or γ) which is largely found in the neurons and neuroendocrine cells. The native enzymes are found as homo- or heterodimeric isoforms composed of three immunologically distinct subunits, α , β and γ . Each subunit (α , β and γ) has a molecular weight of 16 kDa, 44 kDa and 46 kDa respectively.

[0038] The $\alpha\gamma$ and $\gamma\gamma$ enolase isoforms, which have been designated as neuron specific enolase each have a molecular weight of $\sim 80,000$ Dalton. It has been shown that NSE concentration in CSF and blood increases after brain injury (e.g. stroke, head trauma), and the release of NSE from damaged cerebral tissue into the CSF and blood circulation reflects the extent of damage to cerebral tissue. NSE has a biological half-life of ~ 48 hours.

[0039] The S100 protein is a cytoplasmic acidic calcium-binding protein found predominantly in the gray matter of the brain, primarily in astroglial and Schwann cells. The protein exists in several homodimeric or heterodimeric isoforms consisting of two immunologically distinct subunits, A1 (MW 10,400 Dalton) and B (MW 10,500 Dalton). In the central nervous system (CNS), homodimer S 100 B-B (21,000 Dalton) and heterodimer S100 A1-B (20,900 Dalton) constitute over 95 % of the total S100 (Isobe et al. *Biochem Int* 1983;6:419, Zimmer et al. *Brain Res Bull* 1995;37:417). Since a high percentage of S100B is found in the brain, a number of studies have examined this protein as a marker of cerebral injury. The biological half life of S100B is 113 minutes (Usui et al. *Clin Chem* 1989;35:1942). Repeated measurements of S100B serum levels are useful to monitor the course of neurologic damage.

[0040] Glutamine synthetase (GS) is a ubiquitous enzyme that catalyzes the ATP-dependent conversion of glutamate to glutamine using ammonia as the nitrogen source. GS is present at high concentrations in liver, muscle, kidney and brain (De Groot et al. *Biochim Biophys Acta* 1987;908:231). GS in the human brain is an astrocyte-specific enzyme involved in protecting neurons via converting the potentially neurotoxic glutamate and ammonia into glutamine.

[0041] The divalent cation site of GS renders it ex-

tremely sensitive to oxidation.

[0042] Senile plaque-dense regions of the brain of patients with AD represent environments of elevated oxidative stress and that protein in the brain of patients with AD is more oxidized than that of controls. Reactive microglia extensively present with senile plaque regions have been proposed as a source of oxyradicals in the brain.

[0043] Three distinct types of GS have been known: GSI, GSII and GSIII. Genes for GSI and GSIII have been found in bacteria. The human GS gene belongs to the GSII type (Brown et al. *J Mol Evol* 1994;38:566). GS in the brain is thought to exist as an octameric structure with a molecular weight of 360,000 - 400,000 Daltons (Tumani et al. *J Immunol Meth* 1995;188:155). However, in blood circulation the protein is thought to be in the monomeric form with MW 44 ± 1 kDa (Boksha et al. *J Neurochem* 2000;75:2574). High concentrations of GS were reported in lumbar CSF from AD patients (Gunnerson and Haley, *Proc Natl Acad Sci USA* 1992;89:11949, Tumani et al. *Arch Neurol* 1999;56:1241).

[0044] The exact mechanisms leading to increased GS concentration in CSF remain unknown. Over expression of GS in reactive astrocytes and subsequent release into the extracellular space was proposed (Tumani et al, *Arch Neurol* 1999;56:1241).

[0045] In a further contemplated embodiment of the invention, body fluid samples may be taken from a patient at one point in time or at different points in time for ongoing analysis. Typically, a first sample is taken from a patient upon presentation with possible symptoms of AD and analyzed according to the invention. Subsequently, some period of time after presentation, for example, about 3 - 6 months after the first presentation, a second sample is taken and analyzed according to the invention. The data can be used to diagnose AD, rule out AD, distinguish between AD and non-AD dementia.

[0046] The level of any one or all of the three specific markers of interest found in the patient's body fluid were measured from one single sample, although one or more individual markers can be measured in one sample. By "sample" is meant a body fluid such as blood or CSF. All the markers can be measured with one assay device or by using a separate assay device for each marker in which case aliquots of the same sample can be used. It is preferred to measure each of up to the three markers in the same single sample, irrespective of whether the analyses are carried out in a single analytical device or in separate devices so that the level of each marker simultaneously present in a single sample can be used to provide meaningful data.

[0047] The presence of each marker is determined using antibodies specific for each of the markers and detecting specific binding of each antibody to its respective marker. Any suitable direct or indirect assay method may be used, including those which are commercially available to determine the level of each of the specific markers measured according to the invention. The as-

says may be competitive assays, sandwich assays, and the label may be selected from the group of well-known labels such as radioimmunoassay, fluorescent or chemiluminescence immunoassay, or immunoPCR technology. Extensive discussion of the known immunoassay techniques is not required here since these are known to those of skilled in the art. See Takahashi et al. (Clin Chem 1999;45(8):1307) for S100B assay.

[0048] Although not limited thereto, the immunoassay method used in the instant examples comprised a double antibody or sandwich ELISA for measuring the level of the marker proteins in the sample. According to this method, one of the antibodies is a "capture" antibody which is immobilized onto a solid-phase, and the other is a "detector" antibody which is labeled with, for example, an enzyme. The detector antibody binds to marker protein bound to the capture antibody to form a sandwich structure. For each assay of the three markers, marker protein standards are used to prepare a standard or calibration curve of absorbance vs marker protein concentration. This method is of importance since, in addition to the advantages of the preferred embodiment, the added markers, NSE and S100B, assist in indicating an on-going destruction of neurons and monitoring acute events in the brain, respectively.

[0049] The assay methods used to measure the marker proteins should exhibit sufficient sensitivity to be able to measure each protein over a concentration range from normal values found in healthy persons to elevated levels in sick people, i.e. 2SD above normal (= cut-off) and higher. For the GS protein, the cut-off= 0.022 ng/ml, NSE = 8.34 ng/ml and S100B = 0.02 ng/ml.

[0050] The assay may be carried out in various formats, including a microtiter plate format which is preferred for carrying out the assays in a batch mode. The assays may also be carried out in automated analyzers which are well known in the art. Another assay format which can be used according to the invention is a rapid manual test which can be administered at the point-of-care at any location. Typically, such devices will provide a result which is above or below a cut-off, i.e. a semi-quantitative result.

Expression of rhGS and isolation of monoclonal antibody:

[0051] In order to isolate an antibody with specific binding properties for human Glutamine Synthetase (human recombinant GS or rhGS) cDNA of the rhGS was purchased from ATCC. The full length of the rhGS open reading frame (ORF) was obtained by PCR and sub-cloning in pET28a (*Ndel/Shol*). The construct included a poly-histidine tag at the N-termini of the rhGS ORF and no extra sequence at the C-termini. The protein was expressed into *Escherichia coli* BL21 (DE3) following techniques described by Listrom et al. (Biochem J 1997;328:159). The preparation of crude cellular extract as well as solubilization of insoluble expression

products consisting in urea/alkali treatment was achieved following the method of Moreno et al. (J Comp Neurol 1994;350:260). Affinity purification was performed by Ni-NTA chromatography following the supplier's recommendations.

Preparation of Monoclonal Antibody:

[0052] The monoclonal antibody of the present invention was produced by the polyethylene glycol (PEG) mediated cell fusion method.

Preparation of immunocytes: Balb/c mice, a strain with H-2^d haplotype from Charles River Canada, St. Constant, Quebec, Canada, female, 7-9 week old, were immunized with the rhGS emulsified in an equal volume of either Freund's complete adjuvant (FCA) for the first injection and then in Freund's incomplete adjuvant (FIA) for subsequent injections at 2-4 weeks intervals with 25-100 µg of GS. Immunized mice were sacrificed 3-4 days after the final immunization, given either intravenously and/or intraperitoneally, in phosphate buffered saline buffer (PBS), pH 7.4.

Preparation of hybridoma: Spleen cells from the mice immunized with the GS protein and the Sp2/0 cells were fused in the presence of 42 % PEG according to the method described by Fuller, SA, Takahashi, M and Hurrell, JGR (Preparation of Monoclonal Antibodies: In: Subbel F, Brent B, Kingston R., et al., eds. Current Protocols in Molecular Biology. New York: Greene Publishing Associates, 1987: Unit 11). The resulting fused cells were suspended in the HAT selection medium and plated onto four 96-well plates which were pre-seeded with feeder cells, peritoneal exudate cells (PEC), as described by Fuller et al. (see above reference).

Screening of GS-specific antibody-secreting hybridomas :

[0053] The screening of hybridoma cultures were carried out in two methods. (1) Solid-phase ELISA - Screening ELISA -1: Confluent hybridoma culture supernatants were added to 96-well microtiter plates (NUNC MaxiSorp (trade mark), GIBCO BRL) coated with rhGS at 2 µg/ml in 100 mM carbonate buffer, pH 9.6. The excess binding sites were blocked by bovine serum albumin (BSA) in PBS, pH 7.4. After washing the plate with PBS containing 0.05 % Tween 20 (trade mark) (WB), 100 µL culture supernatants containing the monoclonal antibodies were incubated with the immobilized antigen for 1 hour at 37°C. After washing, peroxidase conjugated goat anti-mouse IgG (Jackson ImmunoResearch Lab, Inc., West Grove, Pa.) was added and incubated for 30 minutes at room temperature (RT) on a orbital shaker. After washing, TMB substrate solution (Sigma) was added. After 5 minutes incubation at RT in dark, the reaction was stopped with 1 M H₂SO₄ and optical density was read at 450_{nm}. Selected positive cultures were subjected to cloning by the limiting dilution

method as described by Fuller et al. (see above reference). ELISA screening and cloning procedures were repeated until culture stability and clonality were obtained. (2) Screening ELISA 2 : In the second method, monoclonal antibodies in the hybridoma culture supernatant were captured via goat anti-mouse IgG_{FC} (Jackson immunoResearch) immobilized on the solid-phase of ELISA plate. After 1 hour incubation at 37°C, the plate was washed as for method-1. Then biotinylated GS (prepared using Biotin Labeling Kit from Boehringer Mannheim following the manufacturer's recommendation) diluted 1/2000 dilutions in PBS with 0.5 % BSA was added to each well. After 30 minutes incubation at RT on a shaker, plate was washed and HRP-conjugated streptavidin (Boehringer Mannheim) at 1/10,000 was added and incubated for 30 minutes at RT. After washing, TMB substrate solution was added and reaction was read as for method-1.

[0054] To develop an ELISA assay, two clones designated as 1G3 and 5G4 were selected. These clones were deposited, in accordance with the Budapest Treaty, with the American Type Culture Collection, 10801 University Blvd., Manassas, VA 20110-2209 on April 25, 2001 under Accession Number PTA-3339 and PTA-3340, respectively. In accordance with 37 CFR 1.808, the depositors assure that all restrictions imposed on the availability to the public of the deposited materials will be irrevocable removed upon the granting of a patent. Production of monoclonal antibody:

[0055] The GS-specific monoclonal antibody was produced using ascites. Ascites were produced in Balb/c mice previously treated with 0.25 ml pristane by injecting intraperitoneally with 1 - 5 x 10⁶ hybridoma cells in 0.25 - 0.5 ml PBS, pH 7.4. Ten to 14 days later, ascites were collected. The monoclonal antibody from ascites was purified on affinity column (Protein G, AVIDAL, trade mark), using known procedures. The purified monoclonal antibody was used for immunochemical studies.

Polyclonal antibody production and purification:

[0056] Goats were immunized with biweekly 250 - 500 µg intramuscular and/or subcutaneous injections of purified rhGS emulsified in Freund's adjuvant Titer was monitored routinely by screening serum by a half-sandwich ELISA. Polyclonal antibodies (PAb) specific for GS were subsequently purified from goat serum by affinity purification using cyanogen bromide activated sepharose-4B (trade mark, Pharmacia) coupled to recombinant GS. The purified polyclonal antibodies were dialyzed against 10 mM PBS, pH 7.4, concentrated by ultrafiltration and stored at - 20°C.

Diagnostic Assays and Detection of GS in Biological Samples:

[0057] Another purpose of the present invention is to provide reagents for use in diagnostic assays for the de-

tection of GS from individuals suffering from Alzheimer's disease.

[0058] In one mode of this embodiment, GS of the present invention may be used as antigens in immunoassays for the detection of those individuals suffering from AD. The protein, GS of the present invention may be used in any immunoassay system known in the art including, but not limited to : radioimmunoassay, enzyme-linked immunosorbent assay (ELISA), "sandwich" assays, precipitin reactions, gel diffusion immunodiffusion assay, agglutination assay, fluorescent immunoassays, protein A or G immunoassays and immunoelectrophoresis assays. According to the present invention, monoclonal or polyclonal antibodies produced against human GS are useful in an immunoassay on samples of blood or blood products such as serum, plasma or the like, spinal fluid or other body fluid, e.g. saliva, urine, lymph, and the like, to diagnose patients with AD. Alzheimer's dementia may be determined by utilizing a single monoclonal antibody or a plurality of monoclonal antibodies, singly or combined, which are specific against human glutamine synthetase.

[0059] The antibodies can be used in any type of immunoassay. This includes both the two-site sandwich assay and the single site immunoassay of the non-competitive type, as well as in traditional competitive binding assays.

[0060] Particularly preferred, for ease and simplicity of detection, and its quantitative nature, is the sandwich or double antibody assay of which a number of variations exist, all of which are contemplated by the present invention. For example, in a typical sandwich assay, unlabeled antibody is immobilized on a solid phase, e.g. microtiter plate, and the sample to be tested is added. After a certain period of incubation to allow formation of an antibody-antigen complex, a second antibody, labeled with a reporter molecule capable of inducing a detectable signal, is added and incubation is continued to allow sufficient time for binding with the antigen at a different site, resulting with a formation of a complex of antibody-antigen-labeled antibody. The presence of the antigen is determined by observation of a signal which may be quantitated by comparison with control samples containing known amounts of antigen.

Clinical Studies:

[0061] A prospective observational pilot study was carried out at three geriatric clinics. The study evaluated 38 patients who came to the clinic, wherein Mini Mental State Examination (MMSE) examination and other routine tests were given. Of those, 24 were diagnosed as AD and 14 with other types of dementia other than AD. The mean age of the patients presenting with Alzheimer's dementia was approximately 79 years with an age range of from 54 to 87 years. Mini Mental State Examination score (MMSE) was recorded. Blood sample was obtained, and after clotting, the sample was centrifuged

and aliquots of serum were frozen and stored at 70 °C until analysis for S100, NSE and GS was performed.

[0062] Control subjects included 153 healthy blood donors (age range from 18-87 years: median age 71.03 ± 9.95 years) whose blood samples were used to determine reference values for concentrations of S100, NSE and GS.

[0063] Alzheimer's disease is recognized as a progressive disease process which begins in the basal neocortex, spreads to the hippocampus, and eventually invades all cortical areas. There is no remission in the disease course (Braak and Braak, *Neurobiol Aging* 1997; 18(4):351). This indicates that AD is a disease process and not simply a product of aging. AD pathogenesis involves steady and progressive destruction of neurons. Although NSE does not have specificity for AD, the marker contributes as an indicator of on-going destruction of neurons. Utilizing the monoclonal antibody instantly disclosed for human GS, an elegant, sensitive and specific assay for AD is made possible when the protein is analyzed in non-CSF bodily fluids. This is unlike prior art assays, wherein the markers were specific to non-human GS, and no definitive data could be derived from sera or any other non-CSF bodily fluid. S 100 does not have the sensitivity, but it is a useful marker to monitor acute events in the brain (e.g. TIA, stroke, hypoxia leading to ischemic events, etc) which are common events among older individuals who are the target population of AD.

[0064] All reference values are reported as mean + 2SD. The reference values for the three markers are: GS = 0.022 ng/ml, NSE = 8.34 ng/ml and S 100 = 0.02 ng/ml.

[0065] The levels of S100 and NSE were analyzed using SMART S100 and SMART NSE ELISA assay kit, respectively, available from Syn-X Phanna Inc., Mississauga, Ont, CANADA. For the GS assay, antibody reagents and calibrator (i.e. recombinant human GS) were produced at Syn-X Pharma, and ELISA assay was developed as described above.

[0066] The box whisker plots in Fig. 3 depict the analysis of serum concentrations of GS protein in healthy individuals categorized on the basis of age. Interestingly, the individuals in their 60's show significantly higher levels of GS in blood in comparison to the other younger age groups. However, the mean GS level in blood of individuals 70+ go down similarly to those observed in their 20's and 40's. A similar pattern is observed with S100 protein. On the other hand, age distribution of NSE serum concentration is different from the other two markers in which significantly higher concentrations were observed in the age group of 14 to 40 years. There was no relationship of serum levels of these three proteins to sex. While not wishing to be bound by any particular theories, it was hypothesized that those in their 60's with elevated OS (and/or S100) are indeed already developing brain deterioration, and such individuals in their 70's become classified as "patients" and remaining

so-called "healthy" individuals are well fit health-wise and show low levels of these protein markers. Thus, it was concluded that focusing upon GS levels in the 60's age group is especially valuable as a diagnostic for risk-assessment studies.

[0067] Of the 24 MMSE-based AD serum samples, 20 were mild AD cases, and 4 were moderate (See Figure 1). The sensitivity of the three markers, i.e. GS, NSE and S100B were 100%, 33% and 8% respectively. The GS level correlates well with the severity of AD, i.e. MMSE score, while similar correlation is not observed with NSE. The sensitivity of S 100 is very low, however, when this marker is elevated it may be an indication of ongoing destruction of astrocytes in the brain due to acute events.

[0068] Of the 14 non-AD dementia samples, only one sample shows an elevated level in both GS and S100, while 7 samples exhibit a concentration of NSE above the cut-off level (Figure 2). This indicate GS and S100 to be highly specific markers for AD.

[0069] All patents and publications mentioned in this specification are indicative of the levels of those skilled in the art to which the invention pertains. All patents and publications are herein incorporated by reference to the same extent as if each individual publication was specifically and individually indicated to be incorporated by reference.

[0070] One skilled in the art will readily appreciate that the present invention is well adapted to carry out the objects and obtain the ends and advantages mentioned, as well as those inherent therein. The oligonucleotides, peptides, polypeptides, biologically related compounds, methods, procedures and techniques described herein are presently representative of the preferred embodiments, are intended to be exemplary and are not intended as limitations on the scope. Changes therein and other uses will occur to those skilled in the art which are encompassed within the spirit of the invention and are defined by the scope of the appended claims. Although the invention has been described in connection with specific preferred embodiments, it should be understood that the invention as claimed should not be unduly limited to such specific embodiments.

Claims

1. A diagnostic kit for diagnosing or monitoring the progression of Alzheimer's dementia comprising:

at least one monoclonal antibody which is specific for monomeric brain associated human glutamine synthetase having a weight of about 44 kDa, wherein said antibody is selected from the group consisting of the antibody produced by the clone deposited with the ATCC as Accession Number PTA-3339 and the antibody produced by the clone deposited with the ATCC

as Accession Number PTA-3340.

2. The diagnostic kit of claim 1 in which the antibody is immobilised on a solid support and, further comprising at least one second antibody which is labelled, wherein said second antibody specifically binds to said monomeric brain associated human glutamine synthetase.
3. The diagnostic kit of claim 1, wherein said first antibody is the antibody encoded by the clone deposited with the ATCC as Accession Number PTA-3339.
4. The diagnostic kit of claim 1, wherein said first antibody is the antibody encoded by the clone deposited with the ATCC as Accession Number PTA-3340.
5. An isolated monoclonal antibody produced by the clone deposited with the ATCC as Accession Number PTA-3339.
6. An isolated monoclonal antibody produced by the clone deposited with the ATCC as Accession Number PTA-3340.
7. The isolated clone deposited with the ATCC as Accession Number PTA-3339.
8. The isolated clone deposited with the ATCC as Accession Number PTA-3340.
9. A method of diagnosing or monitoring the progression of Alzheimer's dementia comprising:

contacting a sample of blood, a blood product or cerebrospinal fluid from a patient suspected of suffering from Alzheimer's dementia with monoclonal antibody, wherein said antibody is specific for monomeric brain associated human glutamine synthetase having a weight of about 44 kDa and is selected from the group consisting of the antibody produced by the clone deposited with the ATCC as Accession Number PTA-3339 and the antibody produced by the clone deposited with the ATCC as Accession Number PTA-3340;
detecting specific binding of the antibody to the sample;

wherein detection of an elevated level of said glutamine synthetase in the sample from the patient relative to control samples indicates a diagnosis of Alzheimer's dementia, or wherein the extent of elevation of said glutamine synthetase in the sample from the patient relative to control samples correlates to the progression of Alzheimer's dementia in

the patient.

10. The method of claim 9, wherein said antibody is the antibody produced by the clone deposited with the ATCC as Accession Number PTA-3339.
11. The method of claim 9, wherein said antibody is the antibody produced by the clone deposited with the ATCC as Accession Number PTA-3340.
12. The method of claim 9 wherein said step of detecting includes at least one direct or indirect immunoassay selected from the group consisting of a competitive binding assay, a non-competitive binding assay, a radioimmunoassay, an enzyme-linked immunosorbent assay (ELISA), a sandwich assay, a precipitin reaction, a gel diffusion immunodiffusion assay, an agglutination assay, a fluorescent immunoassay, chemiluminescence immunoassay, immunoPCR immunoassay, a protein A or protein G immunoassay, and an immunoelectrophoresis assay.

Patentansprüche

1. Diagnostischer Testsatz für die Diagnose und Kontrolle des Fortschreitens von Alzheimer-Demenz, umfassend
mindestens einen monoklonalen Antikörper, der spezifisch ist für monomere Hirn-assoziierte humane Glutamin-Synthetase mit einem Gewicht von etwa 44 kDa, wobei der Antikörper ausgewählt ist aus der Gruppe bestehend aus dem durch den Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3339 hinterlegt ist, produzierten Antikörper und dem durch den Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3340 hinterlegt ist, produzierten Antikörper.
2. Diagnostischer Testsatz nach Anspruch 1, bei dem der Antikörper auf einem festen Träger immobilisiert ist, und der einen zweiten Antikörper umfasst, der markiert ist, wobei der zweite Antikörper spezifisch an die monomere Hirn-assoziierte humane Glutamin-Synthetase bindet.
3. Diagnostischer Testsatz nach Anspruch 1, wobei der erste Antikörper der von dem Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3339 hinterlegt ist, kodierte Antikörper ist.
4. Diagnostischer Testsatz nach Anspruch 1, wobei der erste Antikörper der von dem Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3340 hinterlegt ist, kodierte Antikörper ist.
5. Isolierter monoklonaler Antikörper, produziert von

dem Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3339 hinterlegt ist.

6. Isolierter monoklonaler Antikörper, produziert von dem Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3340 hinterlegt ist. 5
7. Isolierter Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3339 hinterlegt ist. 10
8. Isolierter Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3340 hinterlegt ist.
9. Methode zur Diagnose oder Kontrolle des Fortschreitens von Alzheimer-Demenz, umfassend: 15

Kontaktieren einer Blutprobe, eines Blutprodukts oder einer Cerebrospinalflüssigkeit eines Patienten, von dem angenommen wird, dass er an Alzheimer-Demenz leidet, mit einem monoklonalen Antikörper, wobei der Antikörper spezifisch ist für monomere Hirn-assoziierte humane Glutamin-Synthetase mit einem Gewicht von etwa 44 kDa und ausgewählt ist aus der Gruppe bestehend aus dem durch den Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3339 hinterlegt ist, produzierten Antikörper und dem durch den Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3340 hinterlegt ist, produzierten Antikörper; 20
Nachweis von spezifischer Bindung des Antikörpers an die Probe; 25

wobei der Nachweis einer erhöhten Konzentration der Glutamin-Synthetase in der Probe des Patienten im Vergleich zu Kontrollproben eine Diagnose von Alzheimer-Demenz anzeigt oder wobei das Ausmaß des Anstiegs von Glutamin-Synthetase in der Probe des Patienten im Vergleich zu Kontrollproben mit dem Fortschreiten der Alzheimer-Demenz in dem Patienten korreliert. 30

10. Methode nach Anspruch 9, wobei der Antikörper der von dem Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3339 hinterlegt ist, produzierte Antikörper ist. 35
11. Methode nach Anspruch 9, wobei der Antikörper der von dem Klon, der bei ATCC unter der Hinterlegungsnummer PTA-3340 hinterlegt ist, produzierte Antikörper ist. 40
12. Methode nach Anspruch 9, wobei der Nachweis-schritt mindestens einen direkten oder indirekten Immunassay umfasst, ausgewählt aus der Gruppe bestehend aus einem kompetitiven Bindungsassay, einem nicht-kompetitiven Bindungsassay, einem Radioimmunassay, einem Enzym-verknüpften Im- 45

munosorbent-Assay (ELISA), einem Sandwich-Assay, einer Präzipitin-Reaktion, einem Geldiffusions-Immunodiffusionsassay, einem Agglutinationsassay, einem Fluoreszenzimmunassay, einem Chemilumineszenzimmunassay, einem Immuno-PCR-Immunassay, einem Protein A- oder Protein G-Immunassay und einem Immunelektrophoreseassay.

Revendications

1. Trousse de diagnostic pour le diagnostic ou la surveillance de la progression de la démence d'Alzheimer, qui comprend :

au moins un anticorps monoclonal qui est spécifique à la glutamine synthétase humaine monomère associée au cerveau, d'un poids d'environ 44 kDa, ledit anticorps ayant été sélectionné dans le groupe constitué de l'anticorps produit par le clone déposé auprès de l'ATCC sous le n° d'accès PTA-3339 et de l'anticorps produit par le clone déposé auprès de l'ATCC sous le n° d'accès PTA-3340.

2. Trousse de diagnostic selon la revendication 1, dans laquelle l'anticorps est immobilisé sur un support solide, et qui comprend en outre au moins un deuxième anticorps qui est marqué, ledit deuxième anticorps se liant spécifiquement à ladite glutamine synthétase humaine monomère associée au cerveau.
3. Trousse de diagnostic selon la revendication 1, dans laquelle ledit premier anticorps est l'anticorps codé par le clone déposé auprès de l'ATCC sous le n° d'accès PTA-3339.
4. Trousse de diagnostic selon la revendication 1, dans laquelle ledit premier anticorps est l'anticorps codé par le clone déposé auprès de l'ATCC sous le n° d'accès PTA-3340.
5. Anticorps monoclonal isolé produit par le clone déposé auprès de l'ATCC sous le n° d'accès PTA-3339.
6. Anticorps monoclonal isolé produit par le clone déposé auprès de l'ATCC sous le n° d'accès PTA-3340.
7. Clone isolé déposé auprès de l'ATCC sous le n° d'accès PTA-3339.
8. Clone isolé déposé auprès de l'ATCC sous le n° d'accès PTA-3340.

9. Procédé de diagnostic ou de surveillance de la progression de la démence d'Alzheimer, qui comprend les étapes qui consistent à :

mettre en contact avec un anticorps monoclonal un échantillon de sang, un produit sanguin ou du fluide cérébrospinal d'un patient soupçonné d'être atteint par la démence d'Alzheimer, ledit anticorps étant spécifique à la glutamine synthétase humaine monomère associée au cerveau, d'un poids d'environ 44 kDa, et ayant été sélectionné dans le groupe constitué de l'anticorps produit par le clone déposé auprès de l'ATCC sous le n° d'accès PTA-3339 et de l'anticorps produit par le clone déposé auprès de l'ATCC sous le n° d'accès PTA-3340, détecter la liaison spécifique de l'anticorps à l'échantillon, la détection d'un niveau de ladite glutamine synthétase dans l'échantillon provenant du patient, plus élevé que dans des échantillons-témoins, indiquant un diagnostic de la démence d'Alzheimer, ou la valeur du niveau de ladite glutamine synthétase dans l'échantillon provenant du patient, plus élevé que dans des échantillons-témoins, étant corrélée à la progression de la démence d'Alzheimer chez le patient.

10. Procédé selon la revendication 9, dans lequel ledit anticorps est l'anticorps codé par le clone déposé auprès de l'ATCC sous le n° d'accès PTA-3339.
11. Procédé selon la revendication 9, dans lequel ledit anticorps est l'anticorps codé par le clone déposé auprès de l'ATCC sous le n° d'accès PTA-3340.
12. Procédé selon la revendication 9, dans lequel ladite étape de détection comprend au moins une immunodétermination directe ou indirecte sélectionnée dans le groupe constitué d'une détermination par liaison compétitive, d'une détermination par liaison non compétitive, d'une radioimmunodétermination, d'une détermination par immunosorption liée à enzyme (ELISA), d'une détermination sur sandwich, d'une réaction à la précipitine, d'une détermination immunologique par diffusion dans un gel, d'une détermination par agglutination, d'une immunodétermination par fluorescence, d'une immunodétermination par chimioluminescence, d'une immunodétermination par immunoPCR (réaction en chaîne de polymérase), d'une immunodétermination par protéine A ou par protéine G et d'une immunodétermination par électrophorèse.

55

FIGURE 1

AD samples (Serum)

Syn-X CODE	Age	Sex	S100B (ng/ml) Cut-off=0.02	NSE (ng/ml) Cut-off=8.34	GS (ng/ml) Cut-off=0.022	History, Diagnosis : MMSE	20-30-mild 10-19 moderate 0-9 severe
ADH-011	66	M	0.010	9.065	0.058	22 mild AD	
ADH-012	79	M	0.008	4.253	0.023	14 moderate AD	
ADH-013	76	M	0.005	6.166	0.068	17 moderate AD	
ADH-014	85	F	0.009	4.784	0.162	17 moderate AD	
ADH-015	81	M	0.011	8.608	0.058	25 mild AD	
ADH-016	83	F	0.009	14.679	0.077	22 mild AD	
ADH-017	84	F	0.015	10.724	0.092	22 mild AD	
ADH-018	87	F	0.011	4.651	0.042	21 mild AD	
ADH-019	74	F	0.010	6.211	0.131	22 mild AD	
ADH-020	74	M	0.004	3.504	0.035	25 mild AD	
ADH-021	77	M	0.006	1.888	0.045	20 mild AD	
ADH-022	84	M	0.010	4.784	0.049	23 mild AD	
ADH-023	79	M	0.009	7.698	0.052	27 mild AD	
ADH-024	87	M	0.014	8.699	0.071	20 moderate AD	
ADH-025	75	F	0.009	5.853	0.036	21 mild AD	
ADH-026	79	F	0.006	6.750	0.043	24 mild AD	
ADH-027	81	M	0.005	3.284	0.048	23 mild AD	
ADH-028	54	M	0.006	4.121	0.022	21 mild AD / poss VaD	
ADH-029	82	F	0.004	12.969	0.080	23 mild AD	
ADH-030	81	F	0.009	4.563	0.048	24 mild AD	
ADH-031	83	F	0.009	16.316	0.067	20 mild AD	
ADH-032	84	F	0.017	3.460	0.029	24 mild AD	
ADH-033	87	F	0.020	2.628	0.067	21 mild AD	
ADH-034	74	F	0.022	9.111	0.070	22 mild AD	
Sensitivity			8 %	33 %	100 %		


FIGURE 2

Non-AD Dementia Controls

DC-1		Age	Sex	S100B (ng/ml) Cut-off=0.02	NSE (ng/ml) Cut-off=8.34	GS (ng/ml) Cut-off=0.022	History, Diagnosis
Syn-X CODE							
DC-001	86	M	0	3.816	0.019	mild dementia	
DC-002	69	M	0	3.816	0.008	mild dementia (AS), Angina	
DC-003	74	F	0	28.097	0.010	severe dementia (AS)	
DC-004	52	M	0	5.831	0.003	VaD	
DC-005	61	F	0.007	6.220	0.016	FTD, progressive speech problem	
DC-006	75	M	0.015	5.261	0.011	VaD, Hypothyroidism	
ADH-009	86	F	0.009	4.386	0.02	Dementia - type unknown	
ADH-037	78	M	0.008	10.585	0.042	LBD	
AD-042	78	M	0.001	8.146	0.015	16 moderate VaD	
AD-069	86	F	0.002	15.731	0.015	25 mild VaD	
AD-078	56	M	0.004	17.589	0.010	21 mild FTD	
AD-082	72	F	0.017	24.562	0.019	26 mild MCI	
AD-084	51	F	0.011	10.515	0.006	4 severe FTD	
AD-086	67	F	0.025	32.101	0.013	20 mild CBD	
Specificity			93 %	50 %	93 %		

MCI : mild cognitive impairment, CBD : corticobasal degeneration, FTD : fronto-temporal dementia, VaD : vascular dementia
 LBD or DLBD : dementia with Lewy bodies

FIGURE 3
GLUTAMINE SYNTHETASE LEVELS WITH AGE



专利名称(译)	阿尔茨海默氏痴呆的鉴别诊断方法及其装置		
公开(公告)号	EP1381866B1	公开(公告)日	2004-06-09
申请号	EP2002721912	申请日	2002-04-24
申请(专利权)人(译)	SYN.X PHARMA , INC.		
当前申请(专利权)人(译)	SYN.X PHARMA , INC.		
[标]发明人	JACKOWSKI GEORGE TAKAHASHI MIYOKO		
发明人	JACKOWSKI, GEORGE TAKAHASHI, MIYOKO		
IPC分类号	G01N33/53 C07K16/40 C12N5/10 C12P21/08 G01N33/48 G01N33/566 G01N33/573 G01N33/577 G01N33/68		
CPC分类号	G01N33/6896 C07K16/40 C07K2317/20 G01N33/573 G01N2333/4709 G01N2333/4727 G01N2333 /9015 G01N2333/988 G01N2800/2821		
优先权	09/842079 2001-04-25 US 09/971740 2001-10-04 US		
其他公开文献	EP1381866A2		
外部链接	Espacenet		

摘要(译)

公开了一种诊断阿尔茨海默病 (AD) 的方法。该方法涉及直接检测体液，优选血液或血液制品中生化标志物，特别是人谷氨酰胺合成酶的存在。通过包含对人谷氨酰胺合成酶特异的抗体的免疫测定进行检测。另外，公开了一种用于区分AD和非AD痴呆的方法。

FIGURE 1

Syn-X CODE	Age	Sex	S100B		CS	NSE		CS	History, Diagnosis		20-30 mild 16-19 moderate 0-15 severe
			(ng/ml) Ct=ad=0.02	(ng/ml) Ct=ad=0.02		(ng/ml) Ct=ad=0.02	(ng/ml) Ct=ad=0.02		(ng/ml) Ct=ad=0.02	(ng/ml) Ct=ad=0.02	
ADH-011	66	M	0.010	9.065	0.058						22 mild AD
ADH-012	79	M	0.008	4.255	0.023						14 moderate AD
ADH-013	76	M	0.005	6.166	0.068						17 moderate AD
ADH-014	85	F	0.009	4.794	0.162						17 moderate AD
ADH-015	81	M	0.011	8.608	0.058						25 mild AD
ADH-016	83	F	0.009	14.679	0.077						22 mild AD
ADH-017	84	F	0.015	10.774	0.092						22 mild AD
ADH-018	87	F	0.011	4.651	0.048						21 mild AD
ADH-019	74	F	0.010	6.211	0.131						22 mild AD
ADH-020	74	M	0.004	3.594	0.055						25 mild AD
ADH-021	77	M	0.006	1.888	0.045						20 mild AD
ADH-022	84	M	0.010	4.794	0.049						23 mild AD
ADH-023	79	M	0.009	7.688	0.052						27 mild AD
ADH-024	87	M	0.014	8.669	0.071						20 moderate AD
ADH-025	75	F	0.009	5.835	0.056						21 mild AD
ADH-026	79	F	0.006	6.550	0.043						24 mild AD
ADH-027	81	M	0.005	3.284	0.048						22 mild AD
ADH-028	54	M	0.006	4.121	0.022						21 mild AD, poss VAD
ADH-029	82	F	0.004	12.969	0.080						23 mild AD
ADH-030	81	F	0.009	4.553	0.048						24 mild AD
ADH-031	83	F	0.009	16.316	0.067						20 mild AD
ADH-032	84	F	0.017	3.460	0.029						24 mild AD
ADH-033	87	F	0.020	2.628	0.067						21 mild AD
ADH-034	74	F	0.022	9.111	0.070						22 mild AD
Sensitivity											
8 % 33 % 100 %											