



US 20030152948A1

(19) **United States**

(12) **Patent Application Publication**
Shah et al.

(10) **Pub. No.: US 2003/0152948 A1**

(43) **Pub. Date: Aug. 14, 2003**

(54) **METHODS FOR THE SIMULTANEOUS DETECTION OF HCV ANTIGENS AND HCV ANTIBODIES**

(76) **Inventors: Dinesh O. Shah, Libertyville, IL (US); George J. Dawson, Libertyville, IL (US); A. Scott Muerhoff, Kenosha, WI (US); Lily Jiang, Mundelein, IL (US); Robin A. Gutierrez, Gurnee, IL (US); Thomas P. Leary, Kenosha, WI (US); Suresh Desai, Libertyville, IL (US); James L. Stewart, Libertyville, IL (US)**

Correspondence Address:
STEVEN F. WEINSTOCK
ABBOTT LABORATORIES
100 ABBOTT PARK ROAD
DEPT. 377/AP6A
ABBOTT PARK, IL 60064-6008 (US)

(21) **Appl. No.: 10/173,480**

(22) **Filed: Jun. 17, 2002**

Related U.S. Application Data

(63) Continuation-in-part of application No. 09/891,983, filed on Jun. 26, 2001.


Publication Classification

(51) **Int. Cl.⁷ C12Q 1/70; C12Q 1/68; G01N 33/53; C12N 5/06; C12N 5/16**
(52) **U.S. Cl. 435/6; 435/5; 435/7.1; 435/345**

(57) **ABSTRACT**

The subject invention relates to methods for the simultaneous detection of Hepatitis C Virus (HCV) antigens as well as antibodies produced in response to HCV antigens. Furthermore, the subject invention allows one to detect antigens in the early, acute stage of infection, even prior to the development of antibodies, thereby allowing for early detection of infected blood and blood products, thus improving the safety of the blood supply.

A PROTOTYPE ABBOTT PRISM[™] HCV Ag/Ab COMBO ASSAY FORMAT



A PROTOTYPE ABBOTT PRISM™ HCV Ag/Ab COMBO
ASSAY FORMAT





FIG.1


PRISM HCV CORE ANTIGEN ASSAY FORM

STEP 1



ANTI-HCV CORE
MONOCLONAL ANTIBODIES

STEP 2



ACRIDINYLATED
ANTI-HCV CORE
MONOCLONAL ANTIBODIES

FIG.2

ABBOTT PRISM™ HCV
ASSAY FORMAT




FIG.3

SEQ ID NO:1

5'HO- ata gaa ttc cat gca gaa aaa aaa caa acg taa cac caa c
-3'OH

SEQ ID NO:2

5'HO- cgg ctg aga acg ttc aga ggt ttt aac gat ctg acc acc acc cgg
g -3'OH

SEQ ID NO:3

5'HO- aaa acc tct gaa cgt tct cag ccg -3'OH

SEQ ID NO:4

5'HO- tat gga tcc tta tta cgg aga cag cag cca acc agc
-3'OH

SEQ ID NO:5

1 GAATTCATG CAGAAAAAAAA ACAAACGTAA CACCAACCGT CGTCCGCAGG
51 ACGTTAAATT CCCGGGTGGT GGTCAGATCG TTAAACCTC TGAACGTTCT
101 CAGCCGCGTG GGCCTCGTCA GCCGATCCCG AAAGCTCGTC GTCCGGAAGG
151 TCGTACCTGG GCTCAGCCGG GTTACCCGTG GCCGCTGTAC GGTAACGAAG
201 GTTGCGGTTG GGCAGGTTGG CTGCTGTCTC CGTAATAAGG ATCC

SEQ ID NO:6

1 MQKKNKRNTN RRPQDVKFPG GGQIVKTSER SQPRGRRQPI PKARRPEGRT
51 WAQPGYPWPL YGNEGCGWAG WLLSP**

SEQ ID NO:7

1 GAATTCATG CAGAAAAAAAA ACAAACGTAA CACCAACCGT CGTCCGCAGG

FIG.4A

51 ACGTTAAATT CCCGGGTGGT GGTCAGATCG TTGGTCTGCT GCCGCGTCGT
 101 GGTCCGCGTC TGGGTCGTAA AACCTCTGAA CGTTCTCAGC CGCGTGGGCG
 151 TCGTCAGCCG ATCCCGAAAG CTCGTCGTCC GGAAGGTCGT ACCTGGGCTC
 201 AGCCGGGTTA CCCGTGGCCG CTGTACGGTA ACGAAGGTTG CGGTTGGGCT
 251 GGTGGCTGC TGTCTCCGTA ATAAGGATCC

SEQ ID NO:8

1 MQKKNKRNTN RRPQDVKFPG GGQIVGLLPR RGPRLGRKTS ERSQPRGRRQ
 51 PIPKARRPEG RTWAQPGYPW PLYGNEGCGW AGWLLSP**

SEQ ID NO:9

5'HO- acc cag acg cgg acc acg acg cgg cag cag acc aac gat ctg acc
 acc acc c -3'OH

SEQ ID NO:10

5'HO- ccg cgt cgt ggt ccg cgt ctg ggt cgt aaa acc tct gaa cgt tct
 cag -3'OH

SEQ ID NO:11

1 GAATTCCATG CAGAAAAAAA ACAAACGTAA CACCAACCGT CGTCCGCAGG
 51 ACGTTAAATT CCCGGGTGGT GGTCAGATCG TTGGTGGTGT TTACGTTCTG
 101 CCGCGTCGTG GTCCGCGTCT GGGTGTCTG GCTACGCGTA AAACCTCTGA
 151 ACGTTCTCAG CCGCGTGGGC GTCGTCAGCC GATCCCGAAA GTCGTCGTG
 201 CGGAAGGTCG TACCTGGGCT CAGCCGGGTT ACCCGTGGCC GCTGTACGGT
 251 AACGAAGGTT GCGGTTGGGC TGGTGGCTG CTGTCTCCGT AATAAGGATC
 301 C

FIG.4B

SEQ ID NO:12

1 MQKKNKRNTN RRPQDVKFPG GGQIVGGVYV LPRRGPR LGV LATRKTSERS
51 QPRGRRQPIP KARRPEGRTW AQP GYPWPLY GNEGCGWAGW LLSP**

SEQ ID NO:13

5'H0- acc cag acg cgg acc acg acg cgg cag aac gta aac acc acc aac
-3'OH

SEQ ID NO:14

5'H0- ccg cgt cgt ggt ccg cgt ctg ggt gtt ctg gct acg cgt aaa acc
-3'OH

SEQ ID NO:15

1 GAATTCCATG CAGAAAAAAA ACAAACGTAA CACCAACCGT CGTCCGCAGG
51 ACGTTAAATT CCCGGGTGGT GGTCAGATCG TTGGTGGTGT TTACCTGCTN
101 CCGCGTCGTG GTCCGCGTCT GGGTGTTCGT GCTACGCGTA AAACCTCTGA
151 ACGTTCTCAG CCGCGTGGGC GTCGTCAGCC GATCCGAAAG CTCGTCGTCC
201 GGAAGGTCGT ACCTGGGCTC AGCCGGGTTA CCCGTGGCCG CTGTACGGTA
251 ACGAAGGTTG CGGTTGGGCT GGTTGGCTGC TGTCTCCGTA ATAAGGATCC

SEQ ID NO:16

1 MQKKNKRNTN RRPQDVKFPG GGQIVGGVYL LPRRGPR LGV RATRKTSERS
51 QPRGRRQPIP KARRPEGRTW AQP GYPWPLY GNEGCGWAGW LLSP**

SEQ ID NO:17

1 GAATTCCATG GCTGTTGACT TTATCCCGGT TGAAAATCTC GAGACTACTA
51 TGC GTTCTCC GGTTTTC ACT GACA ACTCTT CTCCGCCGGT TGTTCCGCAG
101 TCTTTCCAGG TTGCTCACCT GCATGCTCCG ACTGGTTCTG GTAAATCTAC
151 TAAAGTTCCA GCTGCTTACG CTGCTCAGGG TTACAAAGTT CTGGTTCTGA
201 ACCCGTCTGT TGCTGCTACT CTGGGTTTCG GCGCCTACAT GTCTAAAGCT

FIG.4C

251 CACGGTATCG ACCCGAACAT TCGTACTGGT GTACGTACTA TCACTACTGG
301 TTCTCCGATC ACTTACTCTA CTTACGGTAA ATTCCTGGCT GACGGTGGTT
351 GCTCTGGTGG TGCTTACGAT ATCATCATCT GCGACGAATG CCACTCTACT
401 GACGCTACTT CTATCCTGGG TATCGGTACC GTTCTGGACC AGGCTGAAAC
451 TGCAGGTGCT CGTCTGGTTG TTCTGGCTAC TGCTACTCCG CCGGGTTCTG
501 TTA CTGTTCC GCACCCGAAC ATCGAAGAAG TTGCTCTGTC GACTACTGGT
551 GAAATCCCGT TCTACGGTAA AGCTATCCCG CTCGAGGTTA TCAAAGGTGG
601 TCGTCACCTG ATTTTCTGCC ACTCTAAAAA AAAATGCGAC GAACTGGCTG
651 CTAAGCTTGT TGCTCTGGGT ATCAACGCTG TTGCTTACTA CCGTGGTCTG
701 GACGTTTCTG TTATCCCGAC TTCTGGTGAC GTTGTGTTG TGGCCACTGA
751 CGCTCTGATG ACTGGTTACA CTGGTGACTT CGACTCTGTT ATCGATTGCA
801 ACACTTGCAA TTCCATGTCT ACCAACCCGA AACCGCAGAA AAAAAACAAA
851 CGTAACACCA ACCGTCGTCC GCAGGACGTT AAATTCCCGG GTGGTGGTCA
901 GATCGTTAAA ACCTCTGAAC GTTCTCAGCC GCGTGGGCGT CGTCAGCCGA
951 TCCCGAAAGC TCGTCGTCCG GAAGGTCGTA CCTGGGCTCA GCCGGGTTAC
1001 CCGTGGCCGC TGTACGGTAA CGAAGGTTGC GGTGGGCTG GTTGGCTGCT
1051 GTCTCCGTAA TAAGGATCC

SEQ ID NO:18

1 MAVDFIPVEN LETTMRSPVF TDNSSPPVVP QSFQVAHLHA PTGSGKSTKV
51 PAAYAAQGYK VLVLNPSVAA TLGFGAYMSK AHGIDPNIRT GVRTITTTGSP
101 ITYSTYKFL ADGGCSGGAY DIIICDECHS TDATSILGIG TVLDQAETAG
151 ARLVVLATAT PPGSVTVPHP NIEEVALSTT GEIPFYKAI PLEVIKGRH
201 LIFCHSKKCC DELAAKLVAL GINAVAYYRG LDVSVIPTSG DVVVVATDAL
251 MTGYTGDFDS VIDCNTNSM STNPKPQKKN KRNTNRRPQD VKFPGGGQIV
301 KTSERSQPRG RRQPIPKARR PEGRTWAQPG YPWPLYGNEG CGWAGWLLSP

FIG.4D

SEQ ID NO:19

1 GAATTCCATG GCTGTTGACT TTATCCCGGT TGAAAATCTC GAGACTACTA
51 TGC GTTCTCC GGTTTTC ACT GACA ACTCTT CTCCGCCGGT TGTCCGCAG
101 TCTTCCAGG TTGCTCACCT GCATGCTCCG ACTGGTCTG GTAAATCTAC
151 TAAAGTTCCA GCTGCTTACG CTGCTCAGGG TTACAAAGTT CTGGTTCTGA
201 ACCCGTCTGT TGCTGCTACT CTGGGTTTCG GCGCCTACAT GTCTAAAGCT
251 CACGGTATCG ACCCGAACAT TCGTACTGGT GTACGACTA TCACTACTGG
301 TTCTCCGATC ACTTACTCTA CTTACGGTAA ATTCCTGGCT GACGGTGGTT
351 GCTCTGGTGG TGCTTACGAT ATCATCATCT GCGACGAATG CCACTCTACT
401 GACGCTACTT CTATCCTGGG TATCGGTACC GTTCTGGACC AGGCTGAAAC
451 TGCAGGTGCT CGTCTGGTTG TTCTGGCTAC TGCTACTCCG CCGGGTCTG
501 TTA CTGTTCC GCACCCGAAC ATCGAAGAAG TTGCTCTGTC GACTACTGGT
551 GAAATCCCGT TCTACGGTAA AGCTATCCCG CTCGAGGTTA TCAAAGGTGG
601 TCGTCACCTG ATTTTCTGCC ACTCTAAAA AAAATGCGAC GAACTGGCTG
651 CTAAGCTTGT TGCTCTGGGT ATCAACGCTG TTGCTTACTA CCGTGGTCTG
701 GACGTTTCTG TTATCCCGAC TTCTGGTGAC GTTGTTGTTG TGGCCACTGA
751 CGCTCTGATG ACTGGTTACA CTGGTGACTT CGACTCTGTT ATCGATTGCA
801 ACACTTGCAA TTCCATGCAG AAAAAAACA AACGTAACAC CAACCGTCGT
851 CCGCAGGACG TTAAATTCCC GGGTGGTGGT CAGATCGTTA AAACCTCTGA
901 ACGTTCTCAG CCGCGTGGGC GTCGTCAGCC GATCCCGAAA GCTCGTCGTC
951 CGGAAGGTCG TACCTGGGCT CAGCCGGGT ACCCGTGGCC GCTGTACGGT
1001 AACGAAGGTT GCGGTTGGGC TGGTTGGCTG CTGTCTCCGT AATAAGGATC
1051 C

SEQ ID NO:20

1 MAVDFIPVEN LETTMRSPVF TDNSSPPVVP QSFQVAHLHA PTGSGKSTKV
51 PAAYAAQGYK VLVLNPSVAA TLGFGAYMSK AHGIDPNIRT GVRTITTGSP
101 ITYSTYKFL ADGGCSGGAY DIIICDECHS TDATSILGIG TVLDQAETAG
151 ARLVVLATAT PPGSVTPHP NIEEVALSTT GEIPFYGKAI PLEVIKGRH
201 LIFCHSKKCC DELAAKLVAL GINAVAYYRG LDVSVIPTSG DVVVVATDAL
251 MTGYTGDFDS VIDCNTNSM QKKNKRNTNR RPQDVKFP GG GQIVKTSERS
301 QPRGRROPIP KARRPEGRW AQPYPWPLY GNEGCGWAGW LLSP

FIG.4E

SEQ ID NO:21

5' - TATAGAATTC CATGGCTGTT GACTTTATCC -3'

SEQ ID NO:22

5' - GGAATTGCAA GTGTTGCAAT CGATAAC -3'

SEQ ID NO:23

5' - GTTATCGATT GCAACACTTG CAATTCCATG CAGAAAAAAAA ACAAACGTAA C -3'

SEQ ID NO:24

1 GAATTCCATG GCTGTTGACT TTATCCCGGT TGAAAATCTC GAGACTACTA
 51 TCGTCTCTCC GGTTTTCACT GACAACTCTT CTCCGCCGGT TGTTCGCAG
 101 TCTTTCCAGG TTGCTCACCT GCATGCTCCG ACTGGTCTG GTAATCTAC
 151 TAAAGTTCCA GCTGCTTACG CTGCTCAGGG TTACAAAGTT CTGGTCTGA
 201 ACCCGTCTGT TGCTGCTACT CTGGGTTTCG GCGCCTACAT GTCTAAAGCT
 251 CACGGTATCG ACCCGAACAT TCGTACTGGT GTACGACTA TCACTACTGG
 301 TTCTCCGATC ACTTACTCTA CTTACGGTAA ATTCCTGGCT GACGGTGGT
 351 GCTCTGGTGG TGCTTACGAT ATCATCATCT GCGACGAATG CCACTCTACT
 401 GACGCTACTT CTATCCTGGG TATCGGTACC GTTCTGGACC AGGCTGAAAC
 451 TGCAGGTGCT CGTCTGGTTG TTCTGGCTAC TGCTACTCCG CCGGGTCTG
 501 TTAAGTTTCC GCACCCGAAC ATCGAAGAAG TTGCTCTGTC GACTACTGGT
 551 GAAATCCCGT TCTACGGTAA AGCTATCCCG CTCGAGGTTA TCAAAGGTGG
 601 TCGTCACCTG ATTTTCTGCC ACTCTAAAAA AAAATGCGAC GAACTGGCTG
 651 CTAAGCTTGT TGCTCTGGGT ATCAACGCTG TTGCTTACTA CCGTGGTCTG
 701 GACGTTTCTG TTATCCCGAC TTCTGGTGAC GTTGTGTTG TGGCCACTGA
 751 CGCTCTGATG ACTGGTTACA CTGGTACTT CGACTCTGTT ATCGATTGCA
 801 ACACTTGCAA TTCCGGTGGT GGTGGTTCTA TGCAGAAAAA AAACAAACGT
 851 AACACCAACC GTCGTCCGCA GGACGTTAAA TTCCCGGGTG GTGGTCAGAT
 901 CGTTAAAACC TCTGAACGTT CTCAGCCGCG TGGCGTCGT CAGCCGATCC
 951 CGAAAGCTCG TCGTCCGAA GGTCGTACCT GGGCTCAGCC GGGTTACCCG
 1001 TGGCCGCTGT ACGGTAACGA AGGTTGCGGT TGGGCTGGTT GGCTGCTGTC
 1051 TCCGTAATAA GGATCC

FIG.4F

SEQ ID NO:25

1 MAVDFIPVEN LETTMRSPVF TDNSSPPVVP QSFQVAHLHA PTGSGKSTKV
51 PAAYAAQGYK VLVLNPSVAA TLGFGAYMSK AHGIDPNIRT GVRTITTGSP
101 ITYSTYGFKL ADGGCSGGAY DIIICDECHS TDATSILGIG TVLDQAETAG
151 ARLVVLATAT PPGSVTVPHP NIEEVALSTT GEIPFYGKAI PLEVIKGRH
201 LIFCHSKKCC DELAAKLVAL GINAVAYYRG LDVSVIPTSG DVVVVATDAL
251 MTGYTGDFDS VIDCNTCNSG GGGSMQKKNK RNTNRRPQDV KFPGGGQIVK
301 TSERSQPRGR RQPIPKARRP EGRTWAQPGY PWPLYGNEGC GWAGWLLSP

SEQ ID NO:26

5' - GTTATCGATT GCAACACTTG CAATTCCGGT GGTGGTGGTT CTATGCAGAA
AAAAAACAAA CGTAAC -3'

SEQ ID NO:27

1 GAATTCCATG GCTGTTGACT TTATCCCGGT TGAAAATCTC GAGACTACTA
51 TGC GTTCTCC GGTTTTCACT GACAACTCTT CTCCGCCGGT TGTCCGCAG
101 TCTTTCCAGG TTGCTCACCT GCATGCTCCG ACTGGTCTG GTAAATCTAC
151 TAAAGTTCCA GCTGCTTACG CTGCTCAGGG TTACAAAGTT CTGGTTCTGA
201 ACCCGTCTGT TGCTGCTACT CTGGGTTTTG GCGCCTACAT GTCTAAAGCT
251 CACGGTATCG ACCCGAACAT TCGTACTGGT GTACGTAATA TCACTACTGG
301 TTCTCCGATC ACTTACTCTA CTTACGGTAA ATTCCTGGCT GACGGTGGTT
351 GCTCTGGTGG TGCTTACGAT ATCATCATCT GCGACGAATG CCACTCTACT
401 GACGCTACTT CTATCCTGGG TATCGGTACC GTTCTGGACC AGGCTGAAAC
451 TGCAGGTGCT CGTCTGGTTG TTCTGGCTAC TGCTACTCCG CCGGGTCTG
501 TTA CTGTTCC GCACCCGAAC ATCGAAGAAG TTGCTCTGTC GACTACTGGT
551 GAAATCCCGT TCTACGGTAA AGCTATCCCG CTCGAGGTTA TCAAAGGTGG
601 TCGTCACCTG ATTTTCTGCC ACTCTAAAAA AAAATGCGAC GAACTGGCTG

FIG.4G

651 CTAAGCTTGT TGCTCTGGGT ATCAACGCTG TTGCTTACTA CCGTGGTCTG
701 GACGTTTCTG TTATCCCGAC TTCTGGTGAC GTTGTGTTG TGGCCACTGA
751 CGCTCTGATG ACTGGTTACA CTGGTGACTT C GACTCTGTT ATCGATTGCA
801 ACACTTGCAA TTCCGGTGGT GGTGGTTCTA TGTCTACCAA CCCGAAACCG
851 CAGAAAAAAA ACAAACGTAA CACCAACCGT CGTCCGCAGG ACGTTAAATT
901 CCCGGGTGGT GGTCAGATCG TTGGTGGTGT TTACCTGCTG CCGCGTCGTG
951 GTCCGCGTCT GGGTGTTCGT GCTACGCGTA AACCTCTGA ACGTTCTCAG
1001 CCGCGTGGGC GTCGTCAGCC GATCCCGAAA GCTCGTCGTC CGGAAGGTCG
1051 TACCTGGGCT CAGCCGGGTT ACCCGTGGCC GCTGTACGGT AACGAAGGTT
1101 GCGGTTGGGC TGGTTGGCTG CTGTCTCCGC GTGGATCTCG TCCGTCTTGG
1151 GGTCCGACCG ACCCGCGTCG TCGTTCTCGT AACCTTGGTA AAGTTATCGA
1201 TACCCTGACC TCGGTTTCG CTGACCTGAT GGGTTACATA CCGCTGGTTG
1251 GAGCTCCGCT GGGTGGTGCT GCTCGTGCTT AACCCATGGA TCC

SEQ ID NO:28

1 MAVDFIPVEN LETTMRSPVF TDNSSPPVVP QSFQVAHLHA PTGSGKSTKV
51 PAAYAAQGYK VLVLNPSVAA TLGFGAYMSK AHGIDPNIRT GVRTITTGSP
101 ITYSTYGKFL ADGGCSGGAY DIIICDECHS TDATSILGIG TVLDQAETAG
151 ARLVVLATAT PPGSVTVPHP NIEEVALSTT GEIPFYGKAI PLEVIKGRH
201 LIFCHSKKCC DELAAKLVAL GINAVAYYRG LDVSVIPTSG DVVVVATDAL
251 MTGYTGDFDS VIDCNTCNSG GGGSMSTNPK PQKKNKRNTN RRPQDVKFPG
301 GGQIVGGVYL LPRRGPRLGV RATRKTSEERS QPRGRRQPIP KARRPEGRTW
351 AQPGYPWPLY GNEGCGWAGW LLSPRGSRPS WGPTDPRRRS RNLGKVIDTL
401 TCGFADLMGY IPLVGAPLGG AARA

FIG.4H

SEQ ID NO:29

5' - GTTATCGATT GCAACACTTG CAATTCGGT GGTGGTGGTT CTATGTCTAC
CAACCCGAAA CCGCAG -3'

SEQ ID NO:30

5' - TATAGGATCC ATGGGTAAAG CACGAGC-3'

SEQ ID NO:31

1 GAATTCATG GCTGTTGACT TTATCCCGT TGAAAATCTC GAGACTACTA
51 TCGTTCCTCC GGTTTCACT GACAACTCTT CTCCGCCGGT TGTTCGCAG
101 TCTTCCAGG TTGCTCACCT GCATGCTCCG ACTGGTTCTG GTAAATCTAC
151 TAAAGTTCCA GCTGCTTACG CTGCTCAGGG TTACAAAGTT CTGGTTCTGA
201 ACCCGTCTGT TGCTGCTACT CTGGGTTTCG GCGCCTACAT GTCTAAAGCT
251 CACGGTATCG ACCCGAACAT TCGTACTGGT GTACGTACTA TCACTACTGG
301 TTCTCCGATC ACTTACTCTA CTTACGGTAA ATTCCTGGCT GACGGTGGTT
351 GCTCTGGTGG TGCTTACGAT ATCATCATCT GCGACGAATG CCACTCTACT
401 GACGCTACTT CTATCCTGGG TATCGGTACC GTTCTGGACC AGGCTGAAAC
451 TGCAGGTGCT CGTCTGGTTG TTCTGGCTAC TGCTACTCCG CCGGGTTCTG
501 TTAAGTTCC GCACCCGAAC ATCGAAGAAG TTGCTCTGTC GACTACTGGT
551 GAAATCCCGT TCTACGGTAA AGCTATCCCG CTCGAGGTTA TCAAAGGTGG
601 TCGTCACCTG ATTTTCTGCC ACTCTAAAAA AAAATGCGAC GAACTGGCTG
651 CTAAGCTTGT TGCTCTGGGT ATCAACGCTG TTGCTTACTA CCGTGGTCTG
701 GACGTTTCTG TTATCCCGAC TTCTGGTGAC GTTGTGTTG TGGCCACTGA
751 CGCTCTGATG ACTGGTTACA CTGGTACTT CGACTCTGTT ATCGATTGCA
801 ACACTTGCAA TTCCATGTCT ACCAACCCGA AACCGCAGAA AAAAAACAAA

FIG.4I

851 CGTAACACCA ACCGTCGTCC GCAGGACGTT AAATTCCTGG GTGGTGGTCA
901 GATCGTTTAC CTGCTGCCGC GTCGTGGTCC GCGTCTGGGT GTTACGCGTA
951 AAACCTCTGA ACGTTCTCAG CCGCGTGGGC GTCGTCAGCC GATCCCGAAA
1001 GCTCGTCGTC CGGAAGGTCG TACCTGGGCT CAGCCGGGTT ACCCGTGGCC
1051 GCTGTACGGT AACGAAGGTT GCGGTTGGGC TGGTTGGCTA CTGTCTCCGT
1101 AATAAGGATC C

SEQ ID NO:32

1 MAVDFIPVEN LETTMRSPVF TDNSSPPVVP QSFQVAHLHA PTGSGKSTKV
51 PAAYAAQGYK VLVLNPSVAA TLGFGAYMSK AHGIDPNIRT GVRTITTGSP
101 ITYSTYGFLL ADGGCSGGAY DIIICDECHS TDATSILGIG TVLDQAETAG
151 ARLVVLATAT PPGSVTVPHP NIEEVALSTT GEIPFYGKAI PLEVIKGRH
201 LIFCHSKKCC DELAAKLVAL GINAVAYYRG LDVSVIPTSG DVVVVATDAL
251 MTGYTGDFDS VIDCNTNSM STNPKPQKKN KRNTNRRPQD VKFPGGGQIV
301 YLLPRRGPRL GVTRKTSERS QPRGRRQPIP KARRPEGRTW AQPYPWPLY
351 GNEGCGWAGW LLSP

SEQ ID NO:33

5' - ACCCAGACGC GGACCACGAC GCGGCAGCAG GTAAACGATC TGACCACCAC
CC -3'

SEQ ID NO:34

5' - CCGCGTCGTG GTCCGCGTCT GGGTGTACG CGTAAACCT CTGAACGTTCC
TCAG -3'

SEQ ID NO:35

KTKRNTNRRPQDVKFPGGGQIVYLLPRRGPRLGVTRKTS

FIG.4J

SEQ ID NO:36

KTKRNTNRRPQDVKFPGGGQIVGGVYLLPRRGPRLGVRATRKTS

SEQ ID NO:37

KTKRNTNRRPQDVKFPGGGQIVKTS

FIG.4K

METHODS FOR THE SIMULTANEOUS DETECTION OF HCV ANTIGENS AND HCV ANTIBODIES

[0001] The subject application is a Continuation-In-Part of pending U.S. patent application Ser. No. 09/891,983, filed on Jun. 26, 2001, hereby incorporated in its entirety by reference.

BACKGROUND OF THE INVENTION

[0002] 1. Technical Field

[0003] The subject invention relates to methods for the simultaneous detection of Hepatitis C Virus (HCV) antigens as well as antibodies produced in response to HCV antigens. Furthermore, the subject invention allows one to detect antigens in the early, acute stage of infection, even prior to the development of antibodies, thereby allowing for early detection of infected blood and blood products, and thus improving the safety of the blood supply.

[0004] 2. Background Information

[0005] Recent epidemiological studies indicate that HCV infects more than 170 million people worldwide and that, in more than 50% of the cases, the infection is chronic. In the United States, there are approximately 4 million people infected, and 30,000 new infections are estimated to occur annually (NIH Conference, *Hepatology Suppl* 1:2S (1997)). In addition, HCV is responsible for 8,000-10,000 deaths annually in the United States and is the leading indicator for liver transplantation.

[0006] The HCV genome is a single-stranded RNA molecule of positive polarity that is approximately 9400-9500 nucleotides in length. The organization of the coding regions resembles that of other flaviviruses [Major et al., *Hepatology* 25:1527 (1997)] as well as the more recently discovered GB viruses [Muerhoff A S, et al., *J Virol* 69:5621 (1995)]. The HCV genome possesses a large open reading frame (ORF) encoding a polyprotein precursor of 3010 to 3033 amino acids depending on the particular isolate [Choo et al., *Proc Natl Acad Sci USA* 88:2451 (1991); Grakoui et al., *J Virol* 67:1385 (1993)]. HCV structural genes (core and envelope) are encoded near the 5'-end of the genome, followed by the proteases and helicase, the helicase cofactor and the replicase. Noncoding regions (NCR), thought to be important in replication, are found at each end of the genome.

[0007] HCV infection occurs primarily through parenteral exposure, i.e., through shared needles, by tattooing, or through transfusion of contaminated blood or blood products. Following exposure, the virus enters a susceptible hepatocyte and viral replication occurs. There is an eclipse phase period of approximately 10 days during which time there is no evidence of viral presence (i.e., viral RNA cannot be detected), serum transaminase levels are within normal limits, and there is no evidence of an immune response to HCV [Busch et al., *Transfusion* 40:143 (2000)]. Typically, about 10 days following exposure, HCV RNA can be detected, often with viral loads between 100,000-120,000,000 HCV RNA copies per ml of serum. Several weeks later, there is typically an increase in ALT levels indicating inflammation of the liver; antibodies are detected an average of about 70 days after exposure.

[0008] One of the preventive measures employed to limit the spread of HCV infections is to screen blood for exposure

to HCV, either by the detection of antibodies to HCV or by the detection of viral-specific molecules (e.g., HCV RNA or HCV core proteins) in serum/plasma. Blood or blood products derived from individuals identified as having been exposed to HCV, by these tests, are removed from the blood supply and are not utilized for distribution to recipients of blood products (see, e.g., U.S. Pat. No. 6,172,189). These tests may also be utilized in the clinical setting to diagnose liver disease attributable to HCV infection.

[0009] Due to the unavailability of native, intact HCV virions, serologic antibody tests have relied on recombinant antigens or synthetic peptides, representing selected fragments of the viral polyprotein. The first generation anti-HCV screening tests were based on detection of antibodies directed against a recombinant protein (HCV genotype 1a) originating from sequences located in the nonstructural NS-4 protein (C100-3) [Choo et al., *Science* 244:359 (1989); Kuo et al., *Science* 244:362 (1989)]. The first generation assays failed to detect antibodies in approximately 10% of individuals having chronic HCV infection and up to 10-30% of individuals presenting with acute HCV infection. The second generation anti-HCV assays have incorporated recombinant proteins from three different regions of the HCV genome (HCV genotype 1a), including amino acid sequences from the core, NS3, and NS4 protein [Mimms et al., *Lancet* 336:1590 (1990); Bresters et al., *Vox Sang* 62:213 (1992)], allowing a marked improvement over the first generation tests in identifying HCV infected blood donors [Aach et al., *N Engl J Med* 325:1325 (1991); Kleinman et al., *Transfusion* 32:805 (1992)]. The second generation assays detect antibodies in close to 100% of chronic HCV cases [Hino K., *Intervirology* 37:77 (1994)] and in nearly 100% of the acute cases by 12 weeks post infection [Alter et al., *N Engl J Med* 327:1899 (1992); Bresters et al., *Vox Sang* 62:213 (1992)]. The third generation test includes a recombinant protein expressing amino acid sequences from the NS5 region, as well as antigens from the core, NS3 and NS4. Some studies have indicated a slight improvement in sensitivity in comparing the third generation tests to second generation tests [Lee et al., *Transfusion* 35:845 (1995); Courouce et al. *Transfusion* 34:790-795 (1994)], but this improvement is largely attributed to changes in the NS3 protein rather than the inclusion of NS5 [Courouce et al., *Lancet* 343:853 (1994)].

[0010] In general, the second and third generation HCV antibody tests detect exposure to HCV about 70 days after exposure. Since HCV establishes persistent, and in many cases lifelong infection, the detection of antibodies to HCV represents a very efficient method for determining exposure to HCV. However, antibody testing alone will frequently fail to detect HCV infected individuals during the first 70 days after exposure.

[0011] The existing HCV antigen tests rely on detecting the presence of the HCV core antigen in serum or plasma. The core (or nucleocapsid) protein comprises the first 191 amino acids of the polyprotein. Two different types of serologic assays have been developed which permit detection of HCV core antigens in serum. One assay format detects HCV core antigens in subjects prior to seroconversion and is utilized in screening blood donors, while the other assay format detects core antigens only in hepatitis C patients, regardless of their HCV antibody status and is

utilized in clinical laboratories to diagnose exposure to HCV or to monitor antiviral therapy.

[0012] Recent data on samples obtained during the pre-seroconversion period indicate that the HCV antigen test detects exposure to HCV significantly earlier than antibody testing [Aoyagi et al., *J Clin Microbiol* 37:1802 (1999); Peterson et al., *Vox Sang* 78:80(2000); Dawson et al., *Transfusion*, SD161, 40(2000); Muerhoff et al., 7th *International Meeting on Hepatitis C virus and related viruses*, Dec. 3-7, 2000], and represents an alternative to nucleic acid testing for detecting exposure to HCV during the pre-seroconversion period. The advantages of HCV antigen detection are that the test is rapid, simple, may not require sample extraction or other pretreatment, and is not as prone to handling errors (e.g., contamination) as may occur in the HCV RNA tests.

[0013] In clinical laboratories, the HCV antigen test has comparable sensitivity to the HCV DNA tests in detecting exposure to HCV in patients infected with different HCV genotypes [Dickson et al., *Transplantation* 68:1512 (1999)] and in monitoring antiviral therapy [Tanaka et al., *Hepatology* 32:388 (2000); Tanaka et al., *J Hepatol* 23:742 (1995)]. Thus, HCV core antigen tests present a practical alternative to HCV RNA for screening blood donors or for monitoring antiviral therapy.

[0014] The uniqueness of the current invention lies in its ability to detect HCV antibodies and HCV antigens simultaneously (see also International Application No. PCT/JP99/04129). This combination test or "combo" assay utilizes antigen detection to identify exposure to HCV during the pre-seroconversion "window period" and antibody detection to identify exposure to HCV after seroconversion.

[0015] All U.S. patents and publications referred to herein are hereby incorporated in their entirety by reference.

SUMMARY OF THE INVENTION

[0016] The subject invention encompasses a method of simultaneously detecting at least one Hepatitis C Virus (HCV) antigen and at least one HCV antibody in a test sample comprising the steps of: a) contacting the test sample with: 1) at least one HCV viral antigen or portion thereof coated on a solid phase (e.g., a microparticle), for a time and under conditions sufficient for the formation of antibody/antigen complexes and 2) at least one antibody to HCV or portion thereof coated on the solid phase, for a time and under conditions sufficient for the formation of antigen/antibody complexes; b) detecting the antibody/antigen complexes, presence of the complexes indicating presence of at least one HCV antigen in the test sample; and c) detecting the antigen/antibody complexes, presence of the complexes indicating presence of at least one HCV antibody in the test sample. The at least one HCV antigen coated on the solid phase may be, for example, core antigen, NS3, NS4, NS5, and portions (or fragments) thereof. The at least one antibody coated on the solid phase may be, for example, a monoclonal antibody selected from the group consisting of 13-959-270, 14-1269-281, 14-1287-252, 14-153-234, 14-153-462, 14-1705-225, 14-1708-269, 14-1708-403, 14-178-125, 14-188-104, 14-283-112, 14-635-225, 14-726-217, 14-886-216, 14-947-104, 14-945-218, 13-975-157, 14-1350-210, 107-35-54, 110-81-17, 13-975-157, 14-1350-210, C11-3, C11-7, C11-10, C11-14 and C11-15. Further, the at least one monoclonal

antibody coated on the solid phase preferably is not reactive with the at least one antigen coated on the solid phase. In particular, the at least one monoclonal antibody may be a HCV anti-core monoclonal antibody and the at least one antigen may be a recombinant HCV core protein. The recombinant core protein does not contain epitopes to which the anti-core monoclonal antibody binds.

[0017] Additionally, the present invention includes a method for simultaneously detecting the presence of at least one HCV antigen and at least one HCV antibody in a test sample comprising the steps of: a) contacting the test sample with: 1) at least one HCV viral antigen or portion thereof coated on a solid phase, wherein the solid phase is, for example, a microparticle, for a time and under conditions sufficient for the formation of antibody/antigen complexes and 2) at least one HCV antibody or portion thereof coated on the solid phase, for a time and under conditions sufficient for the formation of antigen/antibody complexes; b) adding a first conjugate to the resulting antibody/antigen complexes for a time and under conditions sufficient to allow the conjugate to bind to the bound antibody in (a)(1), wherein the conjugate comprises a second antibody (e.g., mouse anti-human IgG) attached to a label (for example, a chemiluminescent compound) capable of generating a detectable signal and simultaneously adding a second conjugate to the resulting antigen/antibody complexes for a time and under conditions sufficient to allow said second conjugate to bind to the bound antigen in (a)(2), wherein said second conjugate comprises a third antibody (e.g., a monoclonal antibody to anti-HCV core antigen such as C11-10) attached to the label, for example, chemiluminescent compound, capable of generating a detectable signal; and b) detecting the presence of the generated signal, presence of the signal indicating the presence of at least one HCV antigen or at least one HCV antigen in the test sample. Again, the at least one HCV antigen coated on the solid phase may be selected from the group consisting of core antigen, NS3, NS4, NS5, and portions thereof. Further, the at least one antibody coated on the solid phase may be a monoclonal antibody selected from the group consisting of, for example, 13-959-270, 14-1269-281, 14-1287-252, 14-153-234, 14-153-462, 14-1705-225, 14-1708-269, 14-1708-403, 14-178-125, 14-188-104, 14-283-112, 14-635-225, 14-726-217, 14-886-216, 14-947-104, 14-945-218, 13-975-157, 14-1350-210, 107-35-54, 110-81-17, C11-3, C11-7, C11-10, C11-14 and C11-15. The at least one monoclonal antibody coated on the solid phase is preferably not reactive with the at least one antigen coated on the solid phase.

[0018] Also, the present invention encompasses a kit comprising: a) a container containing at least one HCV antigen coated on a solid phase, wherein the solid phase is, for example, a microparticle; and b) a container containing at least one HCV antibody coated on a solid phase, wherein the solid phase is preferably a microparticle.

[0019] The present invention also includes a kit comprising: a container containing: 1) at least one HCV antigen coated on a solid phase, wherein the solid phase is preferably a microparticle, and 2) at least one HCV antibody, coated on the solid phase. The kit may further comprise at least one conjugate comprising a signal-generating compound attached to a HCV antigen or HCV antibody. The signal-generating compound may be, for example, acridinium or an acridinium-containing compound.

[0020] Additionally, the present invention includes a method of detecting HCV antigen in a test sample comprising the steps of: a) contacting the test sample with at least one HCV antibody (e.g., monoclonal) coated on a solid phase, wherein the solid phase is a microparticle, for a time and under conditions sufficient for the formation of antibody/antigen complexes; and b) detecting the presence of antibody/antigen complexes, presence of the complexes indicating presence of antigen in the test sample.

[0021] The invention also encompasses a method of detecting HCV antigen in a test sample comprising the steps of: a) contacting the test sample with at least one HCV antibody (e.g., monoclonal) coated on a solid phase, wherein the solid phase is, preferably, a microparticle, for a time and under conditions sufficient for the formation of antibody/antigen complexes; b) adding a conjugate to the resulting antibody/antigen complexes for a time and under conditions sufficient to allow the conjugate to bind to the bound at least one antibody, wherein the conjugate comprises a second antibody attached to a label, for example, a chemiluminescent compound capable of generating a detectable signal; and c) detecting the signal generated by the label, for example, chemiluminescent compound, a signal generated by the label indicating the presence of antigen in the test sample.

[0022] Also, the present invention includes a recombinant protein comprising an amino acid sequence selected from the group consisting of, for example, SEQ ID NO: 6, SEQ ID NO: 8, SEQ ID NO: 12 and SEQ ID NO: 16 as well as an amino acid sequence comprising conservative amino acid substitutions of these sequences. (A conservative substitution is defined as one or more amino acid substitutions in a sequence which do not change the function of the sequence.) The present invention also includes a recombinant protein comprising an amino acid sequence encoded by a nucleotide sequence selected from the group consisting of, for example, SEQ ID NO: 5, SEQ ID NO: 7, SEQ ID NO: 11 and SEQ ID NO: 15. (Substitutions, deletions and additions within the sequences which do not affect functionally affect the protein encoded by the sequence are also considered to be within the scope of the present invention.)

[0023] Additionally, the present invention includes a vector or construct comprising a nucleotide sequence selected from the group consisting of, for example, SEQ ID NO: 5, SEQ ID NO: 7, SEQ ID NO: 11 and SEQ ID NO: 15. The invention also includes a host cell comprising the vector or construct.

[0024] Furthermore, the present invention includes an immunoassay which may simultaneously detect at least one HCV antigen or at least one HCV antibody in a test sample.

BRIEF DESCRIPTION OF THE DRAWINGS

[0025] FIG. 1 illustrates the Abbott PRISM® HCV Ab assay format. The assay uses a 2-step format that consists of microparticles coated with recombinant HCV antigens from the core, NS3, NS4 and NS5 regions of the HCV genome. These microparticles, when combined with the donor specimen, a diluent, and a complex of goat polyclonal anti-human F(ab')₂ fragment/murine monoclonal anti-biotin:Acridinium conjugate yield an amount of photons representing a qualitative measurement of anti-HCV antibodies in the specimen when triggered with the PRISM® Activator solution.

[0026] FIG. 2 illustrates the HCV Ag assay format. This assay also uses a 2-step format. When microparticles coated with HCV Mab (e.g., c11-14) are combined with the donor specimen, a diluent and acridinium-labeled Mab (e.g., acridinium labeled c11-10), an amount of photons representing a qualitative measurement of anti-HCV antigens in the specimen will result. The measured amount of photons indicates the amount of HCV antigens in the specimen when triggered with the PRISM® Activator solution.

[0027] FIG. 3 illustrates the HCV Ag/Ab combo assay format. The assay uses a 2-step format. When HCV recombinant antigen and monoclonal antibody blended microparticles (e.g., HCV peptide from the core, and recombinant antigens from the NS3, NS4 and NS5 regions of the HCV genome blended with microparticles coated with c11-14) are combined with the donor specimen, a diluent and blended Acridinium-labeled Mabs (e.g., acridinium labeled c11-10 and acridinylated mouse-anti-human IgG), an amount of photons representing a qualitative measurement of anti-HCV antigens or anti-HCV antibodies or both in the specimen will result when triggered with PRISM® Activator solution.

[0028] FIG. 4 lists all of the nucleotide and amino acid sequences referred to herein as well as the corresponding sequence identifier numbers.

DETAILED DESCRIPTION OF THE INVENTION

[0029] The subject invention relates to various methods which may be utilized in order to simultaneously detect antigens of HCV and antibodies to HCV in a biological sample. Thus, if an individual has either developed specific antibodies to HCV and/or has HCV specific antigens in the biological sample tested, the methods of the present invention will yield a positive result. Such results may be used, for example, to diagnose the patient in terms of presence and status of infection (i.e., acute or chronic) as well as to determine the suitability of a donor blood or blood product sample for transfusion.

[0030] Also, the present invention overcomes the problems associated with the "window period" (i.e., 50-60 days post infection) wherein an individual may be infected with HCV but may not have developed antibodies yet. Such individuals may transmit HCV to others during this period. Thus, by detecting HCV during this "window period", the present invention allows for a quick diagnosis of HCV, as opposed to waiting for the development of antibodies, and prevents contamination of the blood supply.

[0031] In one embodiment of the present invention, HCV viral antigens (e.g., core, N3, N4 and N5), or portions thereof, are coated on a solid phase (or are in a liquid phase). The test or biological sample (e.g., serum, plasma, urine, etc.) is then contacted with the solid phase. If antibodies are present in the sample, such antibodies bind to the antigens on the solid phase and are then detected by either a direct or indirect method. The direct method comprises simply detecting presence of the complex itself and thus presence of the antibodies. In the indirect method, a conjugate is added to the bound antibody. The conjugate comprises a second antibody, which binds to the first bound antibody, attached to a signal-generating compound or label. Should the second antibody bind to a bound first antibody, the signal-generating

compound generates a measurable signal. Such signal then indicates presence of the first antibody in the test sample.

[0032] Examples of solid phases used in diagnostic immunoassays are porous and non-porous materials, latex particles, magnetic particles, microparticles (see U.S. Pat. No. 5,705,330), beads, membranes, microtiter wells and plastic tubes. The choice of solid phase material and method of labeling the antigen or antibody present in the conjugate, if desired, are determined based upon desired assay format performance characteristics.

[0033] As noted above, the conjugate (or indicator reagent) will comprise an antibody (or perhaps anti-antibody, depending upon the assay), attached to a signal-generating compound or label. This signal-generating compound or "label" is itself detectable or may be reacted with one or more additional compounds to generate a detectable product. Examples of signal-generating compounds include chromogens, radioisotopes (e.g., 125I, 131I, 32P, 3H, 35S and 14C), chemiluminescent compounds (e.g., acridinium), particles (visible or fluorescent), nucleic acids, complexing agents, or catalysts such as enzymes (e.g., alkaline phosphatase, acid phosphatase, horseradish peroxidase, beta-galactosidase and ribonuclease). In the case of enzyme use (e.g., alkaline phosphatase or horseradish peroxidase), addition of a chromo-, fluoro-, or lumo-genic substrate results in generation of a detectable signal. Other detection systems such as time-resolved fluorescence, internal-reflection fluorescence, amplification (e.g., polymerase chain reaction) and Raman spectroscopy are also useful.

[0034] Examples of biological fluids which may be tested by the above immunoassays include plasma, urine, whole blood, dried whole blood, serum, cerebrospinal fluid, saliva, tears, nasal washes or aqueous extracts of tissues and cells.

[0035] At the same time as the antibodies are being detected, HCV antigens are also being detected; thus, the present invention obviates the need for the running of two different tests. This is accomplished by exposing the test sample to a solid phase (or liquid phase) coated with specific antibodies to HCV (e.g., human or animal monoclonal antibodies to core, polyclonal antibodies, chimeric antibodies, etc.). Antigens, if present in the sample, bind to the solid phase and may then be detected by a direct or indirect method as described above. More specifically, the indirect method involves the addition of a conjugate comprising a second antibody (which binds to the bound antigen) attached to a label or signal-generating compound. When the second antibody binds to the bound antigen, a detectable signal is then generated indicating presence of HCV antigen in the test sample.

[0036] The antibodies which are coated on the solid phase as well as the "second antibody" may be, as noted above, monoclonal antibodies or polyclonal antibodies. For example, if one chooses to utilize monoclonal antibodies, they may be selected from Abbott monoclonal antibodies 13-959-270, 14-1269-281, 14-1287-252, 14-153-234, 14-153-462, 14-1705-225, 14-1708-269, 14-1708-403, 14-178-125, 14-188-104, 14-283-112, 14-635-225, 14-726-217, 14-886-216, 14-947-104 and 14-945-218. The following anti-core monoclonal antibodies may also be utilized for purposes of the present invention: 107-35-54, 110-81-17, 13-975-157, 14-1350-210 (see U.S. Pat. No. 5,753,430) and Tonen HCV core monoclonals C11-3, 7, 10, 14 and 15 (see

PCT Application WO 099/06836), all of which are available from the American Type Culture Collection, 10801 University Boulevard, Manassas, Va. 20110-2209. (For a discussion of the manner in which monoclonal antibodies may be created, see Kohler and Milstein, *Nature* (1975) 256:494, and reviewed in *Monoclonal Hybridoma Antibodies: Techniques and Applications*, ed. Hurrell (CRC Press, Inc., 1982); see also J. W. Goding in *Monoclonal Antibodies: Principles and Practice* (Academic Press, N.Y., 1983; see also U.S. Pat. No. 5,753,430).

[0037] It should be noted that HCV core protein may be one possible target of the HCV antigen portion of the assay. More specifically, the detection of the core protein is accomplished by using monoclonal antibodies directed towards epitopes within the core protein. These anti-core monoclonals are placed on the solid phase and facilitate the capture of core antigen proteins from the test sample. For detection of HCV antibodies in the test sample, recombinant HCV core protein is also placed on the solid phase. It should be noted however that there are significant problems associated with the use of a single protein as the target for an antigen test and as the capture reagent for antibody detection, namely there is significant "cross-reactivity" between the core antigen and the anti-core monoclonal antibodies coated onto the solid phase(s). This results in a false positive signal, even in the absence of the test sample, since the monoclonal antibodies will bind to epitopes present on the recombinant protein.

[0038] In order to avoid such cross-reactivity, the core protein used in the antibody detection portion of the assay may be modified such that the ability of the anti-core monoclonals to bind HCV core is eliminated. Such modification may be achieved by use of recombinant DNA technology in which the epitope region (i.e., the short sequence of amino acids needed for monoclonal antibody binding) is eliminated or modified. Thus, use of the modified recombinant core protein would consequently maintain several human epitopes to which antibodies present in the serum of infected individuals would bind; however, the anti-core monoclonal antibodies used for antigen capture would not bind the modified protein. Alternatively, one could replace the HCV core recombinant protein with polypeptides that include sequences known to bind to antibodies present in the serum of most infected individuals, but do not include sequences containing the epitopes recognized by the anti-core monoclonals used to detect HCV core antigens.

[0039] More specifically, as noted above, in order to avoid cross-reactivity, one may use core antigens for antibody detection in the assay. In particular, in the present invention, the solid phase may be coated with nonstructural proteins (NS) 3, 4 and/or 5 (i.e., NS3, NS4 and/or NS5) and/or the core protein. Alternatively, in the present invention, the solid phase may be coated with any of the above-mentioned HCV proteins, or segments or portions thereof, either individually or in combination (for antibody detection). The antigens used for coating the solid phase may be generated as a contiguous recombinant protein, expressed as recombinant proteins, either as a single entity or as discrete entities, or as synthetic peptides designed either as a single entity or discrete entities.

[0040] It should also be noted that one may also detect antibodies to HCV E2 in the combo assay. Thus, using the

present assay described herein, one may replace an assay which detects anti-core antibody. Alternatively, one may supplement such an anti-core antibody assay with the antigen assay portion of the combo assay described herein. (See, e.g., U.S. Pat. No. 6,156,495 relating to detection of HGBV E2 antibody or antigen.)

[0041] With respect to detection of antigens in the present invention, as noted above, the monoclonal or polyclonal antibodies coated on the solid phase must not recognize the core antigens used on the solid phase (for antibody detection). Thus, for example, in the present invention, one may use the full antibody or a fragment thereof. (For purposes of the present invention, a "fragment" or "portion" of an antibody is defined as a subunit of the antibody which reacts in the same manner, functionally, as the full antibody with respect to binding properties.)

[0042] Additionally, it should also be noted that the initial capture antibody (for detecting HCV antigens) used in the immunoassay may be covalently or non-covalently (e.g., ionic, hydrophobic, etc.) attached to the solid phase. Linking agents for covalent attachment are known in the art and may be part of the solid phase or derivatized to it prior to coating.

[0043] The second manner in which to use the solid phase for detecting HCV antibodies involves elimination of the core antigens entirely. For example, the solid phase is coated with NS3, NS4 and/or NS5 and a substitute for the core protein or regions thereof (e.g., E2). In contrast, the antibodies coated on the solid phase for detection of antigen are directed against the core protein of HCV.

[0044] Other assay formats which may be used for purposes of the present invention, in order to simultaneously detect antigens and antibodies include, for example, Dual assay strip blots, a rapid test, a Western blot, as well as the use of paramagnetic particles in, for example, an Architect® assay (Frank Quinn, *The Immunoassay Handbook*, Second edition, edited by David Wild, pages 363-367, 2001). Such formats are known to those of ordinary skill in the art.

[0045] It should also be noted that the assays of the present invention may also be used to solely detect HCV antigens or HCV antibodies, rather than both, if desired. Certainly, if one desires to establish that an infection initially exists, one may simply want to determine the presence of antigen in a test sample such as during the "window period". On the other hand, if one wants to establish the stage of infection (e.g., acute versus chronic), one may wish to look for the presence of antibodies and titer thereof.

[0046] It should also be noted that the elements of the assay described above are particularly suitable for use in the form of a kit. The kit may also comprise one container such as vial, bottles or strip, with each container with a pre-set solid phase, and other containers containing the respective conjugates. These kits may also contain vials or containers of other reagents needed for performing the assay, such as washing, processing and indicator reagents.

[0047] The present invention may be illustrated by the use of the following non-limiting examples:

EXAMPLE I

Mapping of HCV Core Epitopes Recognized by Monoclonal Antibodies

[0048] To determine the region within the HCV core protein to which each of the monoclonal antibodies binds, a

series of overlapping, biotinylated peptides were synthesized (Table I). These peptides were used to develop EIAs, as described below. It should be noted that all monoclonals were able to detect a recombinant HCV core fusion protein using EIA methodologies similar to that described below (data not shown).

TABLE I

HCV-Core Derived Peptides		
Peptide	Sequence	Core Region Represented
A	MSTNPKPQKKNKRNTNRR	1-18
B	NKRNTNRRPQDVKFPGGG	11-28
C	DVKFPGGGQIVGGVYLLP	21-38
D	VGGVYLLPRRGPRLGVRA	31-48
E	GPRLGVRATRKTSESRQP	41-58
F	KTSESRQPRRRQPIPKA	51-68
G	RRQPIPKARRPEGRTWAQ	61-78
H	PEGRTWAQPGYPWPPLYGN	71-88
I	QYPWPPLYGNEGCGWAGWLL	81-98
J	CGWAGWLLSPRGRSPSW	91-107
1	WLLSPRGRSPSWGPTDPRRRSRNLG	96-120
2	SWGPTDPRRRSRNLGKVIDTLTCGF	106-130
3	SRNLGKVIDTLTCGFADLMGYIPLV	116-140
4	LTCGFADLMGYIPLVGAPLGGGAARA	126-150
5	YIPLVGAPLGGGAARALAHGVRVLED	136-160
6	GAARALAHGVRVLEDGVNYATGNLPL	146-170
7	LEDGVNYATGNLPGCSFSIFLLA	158-180
8	LPGCSFSIFLLALLSCLTVPASA	169-191

[0049] Coating of polystyrene beads: One quarter-inch polystyrene beads were used as the solid phase for the peptide EIAs. Prior to coating, beads were washed with 15% isopropanol (in water) at room temperature for 30 minutes without agitation. Isopropanol was removed and the beads were rinsed once with deionized water. The washed beads were then added to a vial containing the peptide diluted to 5 µg/ml in 0.1 M sodium phosphate, pH 7.5 buffer (0.233 ml per bead). Beads were incubated at 56° C. for 2 hours with end-over-end mixing. Beads were then washed three times with PBS and then incubated in PBS containing 0.1% Triton X-100 at 40° C. for 1 hour with end-over-end mixing. They were again washed three times in PBS and then incubated at 40° C. in 5% BSA/PBS for 1 hour end-over-end mixing. Beads were washed four times with PBS and then incubated at room temperature in PBS containing 5% sucrose without mixing for 20 minutes. Sucrose buffer was removed and beads air-dried. Coated beads were stored desiccated at 4° C.

[0050] Bead coating validation: To determine whether the biotinylated peptides were actually coated onto the beads, an assay was performed in which beads were incubated in buffer containing horseradish peroxidase-labeled streptavidin (200-400 ng/ml). The beads were then washed with deionized water and substrate added. Product was detected by absorbance at 492 nm. All peptides in Table I were shown to be coated onto the polystyrene beads by this assay (data not shown).

[0051] HCV peptide EIAs: Monoclonal antibodies generated against a recombinant HCV core protein (see Example I) were tested for their ability to bind to each of peptide-coated beads as follows: monoclonal antibodies were diluted to 50 ng/ml in sample diluent buffer (Tris buffer containing 20% goat serum, 10% calf serum, 0.2% Triton X-100 and sodium azide) of which 0.2 ml was added into a reaction well containing the peptide-coated bead and incubated at room temperature for 2 hours with mixing. Beads

EXAMPLE II

Epitope Mapping of Monoclonal Antibodies

[0053] A. Preparation of HCV Gene Fragment Library.

[0054] A plasmid containing nucleotides 14-5294 of the H strain of HCV (Ogata et al., *Proc. Natl. Acad. Sci. USA* 88:3392-3396 (1991)) in pGEM-9Zf(-) (Promega Corp., Madison, Wis.) was partially digested using DNase I by the following method in order to obtain random epitope-encoding fragments:

[0055] Five μg aliquots of plasmid DNA were incubated at 15° C. for 10 minutes in 0.5 M Tris-HCl, pH 7.6, and 10 mM MnCl₂ in the presence of anywhere from 0.1 to 0.7 units of DNase I. Aliquots from each digestion were analyzed by agarose gel electrophoresis. The two digestion mixtures containing 0.6 and 0.7 units of DNase I were found to give the largest amount of fragments in the 50-200 bp range. These two mixtures were pooled and extracted one time with an equal volume of phenol-chloroform (1:1, v/v) then precipitated by the addition of one tenth volume 3 M sodium acetate and 2.5 volumes 100% ethanol followed by centrifugation at 14,000 \times g for 10 minutes. The ends of the DNA molecules were then made blunt using the PCR Polishing kit (Stratagene, Inc., La Jolla, Calif.) as per manufacturer's directions. The DNA was again extracted and precipitated as described above, followed by ligation to a double-stranded adaptor in a 10 μl reaction volume using a T4 DNA ligase kit (Stratagene, Inc., La Jolla, Calif.) as directed by the manufacturer. The sequence of this double stranded adaptor was:

5'-GATCGCTCGAATTCCTCG-3' (SEQUENCE ID NO:1)

3'-TTCTAGCGAGCTTAAGGAGC-5' (SEQUENCE ID NO:2)

[0056] The sense-strand oligonucleotide of the adaptor (SEQ ID NO: 1) was then used as a primer in a PCR reaction such that all DNAs were amplified independent of their sequence. This method is a modification of that described by Akowitz et al., *Gene* 81:295-306 (1989) and Reyes et al., *Mol. Cell. Probes* 5:473-481 (1991). PCR was performed in the presence of the sense-strand oligonucleotide primer at a final concentration of 1 μM in a reaction volume of 100 μl using the GeneAmp Gold PCR kit (PE Applied Biosystems, Foster City, Calif.) as directed by the manufacturer in a PE-9600 thermocycler. A pre-incubation at 94° C. for 8 min was followed by twenty-five cycles of PCR as follows: denaturation at 94° C. for 20 seconds, annealing at 55° C. for 30 seconds, and extension at 72° C. for 1.0 min. This was followed by a final extension step at 72° C. for 10 min. The PCR product was extracted and precipitated as described above. The entire PCR was run on a 1.2% agarose gel and a gel slice containing DNA fragments between approximately 70 and 250 bp was removed. The DNA was extracted from the gel slice using the QIAEX II kit (QIAGEN, Inc., Valencia, Calif.) as per manufacturer's directions. The DNA was digested with the restriction enzyme EcoRI (New England Biolabs, Beverly, Mass.) as directed by the manufacturer. The digested DNA was then extracted and precipitated as described above.

[0057] T7Select10-3b (Novagen, Inc., Madison, Wis.) was digested with EcoRI and dephosphorylated with calf intes-

tinal alkaline phosphatase (New England Biolabs, Beverly, Mass.) as directed by the manufacturer. Size-selected digested DNA fragments (30 ng) (supra) were ligated with 0.5 μg digested T7Select10-3b in a 5 μl reaction volume at 16° C. overnight. The entire ligate was packaged into phage heads using T7Select packaging extract (Novagen, Inc., Madison, Wis.) and titered as directed by the manufacturer. The resulting unamplified library contained a total of 3.9 \times 10⁶ members (PFU). The packaged phage were amplified by liquid lysate amplification in *E. coli* BLT5403 (20 ml culture) as directed by the *T7 Select System Manual* (Novagen, Inc., Madison, Wis.). The amplified library had a titer of 1.3 \times 10¹¹ PFU per ml.

[0058] B. Biopanning of HCV Gene Fragment Library.

[0059] Each monoclonal antibody (20 μg) that was to be used for biopanning was incubated at 4° C. for 4 hours on an end-over-end rocker in 300 μl blocking buffer (2% BSA, 3% nonfat dry milk, 0.2% Tween 20, 0.02% sodium azide in phosphate-buffered saline). During the incubation of the monoclonal antibody, an aliquot of the amplified HCV gene fragment library (supra) containing approximately 10¹¹ phage was precipitated as follows: 1/10 volume of 5 M NaCl was added to the phage, mixed thoroughly, followed by addition of 1/2 volume polyethylene glycol (MW 8000), mixed thoroughly again, and incubated on ice for 1-2 hours. The phage were centrifuged at 6000 \times g for 10 min at room temp, all supernatant was removed and the phage pellet was vigorously resuspended in 120 μl buffer containing 1 M NaCl, 10 mM Tris-HCl pH 8.0, 1 mM EDTA. The phage were added to the pre-incubated monoclonal antibody and incubated at 4° C. overnight on an end-over-end rocker.

[0060] The next morning, the antibody-phage complexes were captured on paramagnetic particles coupled to goat anti-mouse IgG (Fc specific) as follows. A 0.2 ml aliquot of Goat Anti-Mouse IgG Fc BioMag particles (Polysciences, Inc., Warrington, Pa.) was washed three times with 0.4 ml 0.1% Tween 20, 0.1% BSA, 0.02% sodium azide in phosphate-buffered saline (PBS) by gentle vortexing followed by capture on a magnetic stand for 0.5-1 minute. The supernatant was removed carefully without disturbing the particles. Particles were then resuspended in the IgG-phage from overnight incubation above and incubated at room temperature on an end-over-end rocker for 3 hours. Particles were washed six times as above using 6.0 ml 0.5% Tween 20, 0.1% BSA in PBS per wash. Bound phage were eluted using 0.2 ml 0.1% Tween 20, 0.1% BSA, 1.0% SDS in PBS at room temperature on an end-over-end rocker for 90 minutes. The tube was placed on a magnetic stand for 1-2 minutes, after which the supernatant containing the eluted phage was removed to a clean tube. The sample containing eluted phage was titered as directed in the *T7 Select System Manual*.

[0061] The eluted phage was amplified as follows. Ten ml LB Broth (Gibco BRL, Gaithersburg, Md.) plus 100 $\mu\text{g}/\text{ml}$ ampicillin was inoculated with *E. coli* BLT5403 and incubated at 37° C. overnight with shaking. The following morning, 35 ml LB Broth plus 100 $\mu\text{g}/\text{ml}$ ampicillin, 1X M9 salts, 0.4% glucose, 1 mM MgSO₄ was inoculated with 0.2 ml of the overnight culture and incubated at 37° C. with shaking until the A600 absorbance was 0.5-0.6. Eluted phage (185 μl) from first round biopanning (supra) was added and incubation at 37° C. was continued for 1.5-2

hours, until the A600 absorbance of the culture had dropped to approximately 0.5, indicating lysis. The culture was centrifuged at 8000×g for 10 minutes and the supernatant was removed to a clean tube and stored at 4° C. The culture supernatant was titered as directed in the *T7 Select System Manual*.

[0062] One to two subsequent rounds of biopanning and amplification were performed as above with the following modifications. After pre-blocking the monoclonal antibody for 4 hours at 4° C., 150 μ l amplified phage from the previous round of biopanning was added instead of 10¹¹ PEG-precipitated phage from the starting library. In addition, after biopanning, a 20 ml culture rather than a 35 ml culture was used to amplify the eluted phage, and 100 μ l rather than 185 μ l of eluted phage was added to the culture.

[0063] C. Selection and Sequencing of HCV Core-containing Clones.

[0064] A DNA fragment containing a region of the HCV genome that encodes amino acids 1-173 of the HCV nucleocapsid protein was utilized as a hybridization probe. This region was chosen because all of the monoclonal antibodies analyzed in the biopanning experiments recognize epitopes in the HCV core protein. Phage resulting from 2-3 rounds of biopanning and amplification were plated on *E. coli* BLT5403 and incubated at 37° C. until plaques formed. DNA was transferred onto Hybond-N+ membranes (Amersham Life Sciences, Inc., Arlington Heights, Ill.), denatured, neutralized, and UV cross-linked, as described by the manufacturer. The membranes were pre-hybridized, hybridized with the HCV nucleocapsid gene ³²P-labeled probe, washed and exposed as described and known in the art. Individual hybridizing plaques were isolated and the inserts were amplified by PCR using T7SelectUP and T7SelectDOWN primers (Novagen, Inc., Madison, Wis.) as directed in the *T7Select System Manual*. For each monoclonal antibody, 30-50 independent hybridizing plaques were amplified and then purified using the QIAquick PCR purification kit (Qiagen, Inc., Chatsworth, Calif.). Purified PCR products were sequenced directly on an ABI Model 377 DNA Sequencer using the ABI Big Dye Terminator Cycle Sequencing Ready Reaction kit (Perkin-Elmer) and the T7SelectUP primer. All of the sequences resulting from biopanning with a particular monoclonal antibody were aligned with the HCV nucleocapsid gene sequence and the minimum region of overlap among all clones was identified. This overlap region defined the epitope recognized by the monoclonal antibody. The epitopes recognized by several monoclonal antibodies that were identified using this method are shown in TABLE III.

TABLE III

Monoclonal Antibody	Region of HCV Core Recognized
C11-15	Amino acids 19-27
C11-10	Amino acids 32-36
C11-14	Amino acids 45-50
C11-3	Amino acids 104-110
C11-7	Amino acids 112-124
14-635-225	Amino acids 49-53
14-153-462	Amino acids 50-63
14-726-217	Amino acids 57-63
14-178-125	Amino acids 59-64

TABLE III-continued

Monoclonal Antibody	Region of HCV Core Recognized
14-1269-281	Amino acids 59-64
14-947-104	Amino acids 59-64
14-188-104	Amino acids 59-64
14-1708-269	Amino acids 59-64
107-35-54	Amino acids 102-109

EXAMPLE III

Construction of Recombinant Antigens for Use in an HCV Core Antibody/Antigen Combination Assay

[0065] A. Background.

[0066] The human immune response to Hepatitis C Virus (HCV) core is, for the most part, exclusive to the N-terminal half of the native protein. Multiple epitopes (regions comprising a defined number of amino acids, usually <10) have been identified within the first 115 amino acids of the native protein (Sallberg et al). Therefore, recombinant antigens utilized in assays for the detection of human anti-core antibodies present in the serum of infected individuals need only contain this portion of the native protein. Conversely, in vitro assays for the detection of HCV core protein utilize murine monoclonal antibodies to capture and detect native core protein also present in the serum of infected individuals. Combination assays for the simultaneous detection of both core antigen and human anti-core antibody in a single assay combine the two assay formats. In this case, a recombinant core antigen is necessary that will be recognized by human anti-core antibodies present in the serum, while escaping recognition by the murine monoclonal antibodies used to capture and detect native core antigen also present in serum. Such recombinants can be constructed by eliminating small regions (1-30 or more amino acids) within the core antigen, thus disrupting or eliminating the epitope(s) recognized by the murine monoclonal antibodies while at the same time leaving undisturbed numerous other epitopes that will allow human anti-core antibody detection.

[0067] B. Antigen Construction.

[0068] The first of these antigens constructed contained HCV amino acids 8-100 in which amino acids 32-50 were deleted from the recombinant. In this manner, antibodies C11-10 and C11-14, which bind to epitopes at amino acids 32-36 and 45-50, respectively (Example II), will not bind to the resulting recombinant antigen. A plasmid containing a bacterial codon-optimized version of nucleotides 342-791 from the H strain of HCV (Ogata, supra (1991)) was used as template for the PCR. Oligonucleotide primers designed to this sequence were utilized to amplify two distinct fragments (1=SEQ ID NO: 1 and SEQ ID NO: 2; 2=SEQ ID NO: 3 and SEQ ID NO: 4) of the HCV core antigen, one upstream of the deletion and the other downstream of the deletion. Additionally, the oligonucleotide primers flanking the deletion on each fragment (SEQ ID NO: 2 and SEQ ID NO: 3) contained regions of overlap with one another. PCR was performed in a PE-9600 thermocycler in the presence of 0.5 μ M of each oligonucleotide primer and 3-5 μ g of plasmid template using the TaKaRa LA Taq PCR Kit (Pan Vera,

Corp., Madison, Wis.) as per manufacturer's instructions (50 μ l volume). A pre-incubation at 94° C. for 1 minute was followed by 35 cycles of PCR (denaturation at 94° C. for 20 seconds; annealing at 50° C. for 30 seconds; extension at 72° C. for 30 seconds), which was then followed by a final extension at 72° C. for 10 minutes. Following the independent amplification of the two fragments, each was purified using the QIAquick PCR Purification Kit (QIAGEN, Inc., Valencia, Calif.) and eluted into 25 μ l of water. A second round of PCR was then performed as supra for 10 cycles to tether the two purified fragments (1 μ l each in a 20 μ l volume reaction) to one another at the regions of overlap within the oligonucleotide primers (SEQ ID NO: 2 and SEQ ID NO: 3). Finally, the product of the second PCR was amplified in a third PCR for 35 cycles utilizing 0.5 μ M of the flanking oligonucleotide primers (SEQ ID NO: 1 and SEQ ID NO: 4) and the conditions described supra (100 μ l volume).

[0069] C. Recombinant Expression.

[0070] The product of the third PCR was purified using the QIAquick PCR Purification Kit, then digested with the restriction enzymes EcoRI and BamHI, digestion sites that were incorporated into the flanking oligonucleotide primers (SEQ ID NO: 1 and SEQ ID NO: 4). The fragment encoding the recombinant was ligated into the bacterial expression vector pJO200 that had been similarly digested with EcoRI and BamHI using the pGEM-T Easy Ligation Kit (Promega, Madison, Wis.), then transformed into XL1-Blue competent cells (Stratagene, La Jolla, Calif.). After selection of clones containing the appropriately sized inserts, the nucleotide sequence of the recombinant was confirmed (SEQ ID NO: 5), exhibiting the deduced amino acid sequence (SEQ ID NO: 6). (See U.S. Pat. Nos. 5,322,769 and 6,172,189 for a description of the expression of recombinant proteins.)

[0071] D. Other Recombinants.

[0072] Several other HCV core recombinants have been constructed, cloned and expressed in a manner identical to that described supra. The first of these (see nucleotide sequence SEQ ID NO: 7 and corresponding amino acid sequence SEQ ID NO: 8) contains HCV amino acids 8-100 in which amino acids 33-35 and 46-49 have been deleted. The oligonucleotide primers utilized to amplify the two fragments in the first PCR were SEQ ID NO: 1 and SEQ ID NO: 9, and SEQ ID NO: 10 and SEQ ID NO: 4, respectively, followed by final amplification in the third PCR with SEQ ID NO: 1 and SEQ ID NO: 4. The second of these recombinants (see nucleotide sequence SEQ ID NO: 11 and corresponding amino acid sequence SEQ ID NO: 12) encodes HCV amino acids 8-100 in which the leucine residue at amino acid 36 has been substituted with valine, and the arginine residue at amino acid 47 has been substituted with leucine. The oligonucleotide primers utilized to amplify the two fragments in the first PCR were SEQ ID NO: 1 and SEQ ID NO: 13, and SEQ ID NO: 14 and SEQ ID NO: 4, respectively, followed by final amplification in the third PCR with SEQ ID NO: 1 and SEQ ID NO: 4. Finally, a recombinant (see nucleotide sequence SEQ ID NO: 15 and corresponding amino acid SEQ ID NO: 16) encoding HCV amino acids 8-100 was constructed using oligonucleotide primers SEQ ID NO: 1 and SEQ ID NO: 4 in a single PCR reaction.

EXAMPLE IV

Construction of Additional Recombinant Antigens for Use in an HCV Core Antibody/Antigen Combination Assay

[0073] A. Background.

[0074] Additional recombinant antigens constructed for use in an HCV antigen/antibody combination assay included antigens that contained the 33c region of HCV (amino acids 1192-1457) tethered to a core region of the virus. The template used for such constructions was a plasmid containing a bacterial codon-optimized sequence of amino acids 1192-1457, followed by amino acids 1-150 from the H strain of HCV (Ogata, 1991), with two non-HCV coding amino acids separating the two sequences. This recombinant (HC-43) is routinely used in multiple commercial assays for the detection of HCV. The HC43 recombinant is expressed as a non-fusion protein from the pL promoter of bacterial phage lambda. (See U.S. Pat. Nos. 5,705,330, 5,616,460 and 5,773,212 for a discussion of HC43 and U.S. Pat. Nos. 6,153,377 and 5,859,193 for a discussion of the lambda pL vector system.) The additional recombinants were constructed by one of two methods. First, existing clones encoding distinct, related recombinants were joined by DNA ligation to form a third unique recombinant, or, new unique clones were constructed by tethering PCR described in Example III, Part B.

[0075] B. Initial Antigen Construction.

[0076] The first of the newer antigens constructed (p9MB-18) contained HCV amino acids 1192-1457 (representing a segment of NS3) tethered to amino acids 1-100 (representing a segment of core protein) in which amino acids 32-50 had been deleted. This recombinant was constructed by restricting plasmid pHC43 with the endonucleases, Xma I and Bam HI. Xma I cuts pHC43 near amino acid 24 of the core-encoding region while Bam HI cuts downstream of the translation termination codon. This region of pHC43 was replaced by DNA ligation using the pGEM-T Easy Ligation Kit (Promega Corp., Madison, Wis.), with the Xma I-Bam HI fragment obtained from the pJO200 vector encoding HCV core amino acids 8-100 with 32-50 deleted described in Example III (SEQ ID NO: 5 and SEQ ID NO: 6). The new plasmid was then used to transform XL1-Blue competent cells (Stratagene, La Jolla, Calif.). After selection of clones containing the appropriately sized insert, the nucleotide sequence of the recombinant was confirmed (SEQ ID NO: 17), exhibiting the deduced amino acid sequence in SEQ ID NO: 18.

[0077] C. Other Recombinants Constructed by Tethering PCR.

[0078] Recombinant antigens made by tethering PCR were constructed as detailed in Example III, Part B. The first of these recombinants, p9MB-19 (SEQ ID NO: 19 (nucleotide sequence) and SEQ ID NO: 20 (amino acid sequence)), contains HCV amino acids 1192-1457 followed by amino acids 8-100 in which amino acids 32-50 had been deleted. The oligonucleotide primers used to amplify the initial two fragments in the first PCR were SEQ ID NO: 21 and SEQ ID NO: 22, and SEQ ID NO: 23 and SEQ ID NO: 4, respectively. Final amplification in the third PCR utilized oligonucleotide primers SEQ ID NO: 21 and SEQ ID NO: 4.

[0079] Recombinant p9MB-20 (SEQ ID NO: 24 and SEQ ID NO: 25) contains HCV amino acids 1192-1457, four glycine residues and a serine residue, followed by HCV amino acids 8-100 in which amino acids 32-50 had been deleted. The oligonucleotides primers used to amplify the initial two fragments in the first PCR were SEQ ID NO: 21 and SEQ ID NO: 22, and SEQ ID NO: 26 and SEQ ID NO: 4, respectively. Final amplification in the third PCR utilized oligonucleotides primers SEQ ID NO: 21 and SEQ ID NO: 4.

[0080] Recombinant p9MB-22 (SEQ ID NO: 27 and SEQ ID NO: 28) contains HCV amino acids 1192-1457, four glycine residues and a serine residue, followed by HCV amino acids 1-150. The oligonucleotides primers used to amplify the initial two fragments in the first PCR were SEQ ID NO: 21 and SEQ ID NO: 22, and SEQ ID NO: 29 and SEQ ID NO: 30, respectively. Final amplification in the third PCR utilized oligonucleotides primers SEQ ID NO: 21 and SEQ ID NO: 30.

[0081] Recombinant p9MB-31 (SEQ ID NO: 31 and SEQ ID NO: 32) contains HCV amino acids 1192-1457 followed by amino acids 1-100 in which amino acids 31, 32, 33, 47 and 48 had been deleted. The oligonucleotides primers used to amplify the initial two fragments in the first PCR were SEQ ID NO: 21 and SEQ ID NO: 33, and SEQ ID NO: 34 and SEQ ID NO: 4, respectively. Final amplification in the third PCR utilized oligonucleotides primers SEQ ID NO: 21 and SEQ ID NO: 4.

[0082] D. Other Recombinants Constructed by DNA Ligation.

[0083] Recombinant 9MB-24 (SEQ ID NO: 35 and SEQ ID NO: 36) contains HCV amino acids 1192-1457 followed by amino acids 1-100 in which amino acids 33-35 and 46-49 were deleted from the recombinant. This recombinant was constructed by restricting plasmid pHc43 with the endonucleases, Xma I and Bam HI, and replacing the fragment by DNA ligation, with the Xma I-Bam HI fragment obtained from SEQ ID NO: 7.

[0084] Recombinant 9MB-25 (SEQ ID NO: 37 and SEQ ID NO: 38) contains HCV amino acids 1192-1457, four glycine residues and a serine residue, followed by HCV amino acids 1-100 in which amino acids 33-35 and 46-49 were deleted from the recombinant. This recombinant was constructed by restricting plasmid p9MB-22 with the endonucleases, Xma I and Bam HI, and replacing the fragment by DNA ligation, with the Xma I-Bam HI fragment obtained from SEQ ID NO: 7.

[0085] Recombinant 9MB-25 (SEQ ID NO: 39 and SEQ ID NO: 40) contains HCV amino acids 1192-1457, four glycine residues and a serine residue, followed by HCV amino acids 1-100 in which amino acids 32-50 were deleted from the recombinant. This recombinant was constructed by restricting plasmid p9MB-22 with the endonucleases, Xma I and Bam HI, and replacing the fragment by DNA ligation, with the Xma I-Bam HI fragment obtained from p9MB-18 (SEQ ID NO: 27).

EXAMPLE V

Preparation of Microparticles

[0086] Microparticles, coated with several monoclonal antibodies, were prepared by coating several separate popu-

lations of microparticles with HCV monoclonal antibodies which recognize different regions within the HCV core protein. Similarly, microparticles were coated with recombinant antigens cloned from the NS3 and NS4 regions of HCV. The peptide used for microparticle coating was from the core region of HCV.

[0087] Microparticles for Antibody Assay:

[0088] The following recombinant proteins and peptides were used to coat the microparticles for antibody assays.

[0089] A. Preparation of Recombinant Proteins:

[0090] i. HCV HC43 antigen HCV. HC43 recombinant antigen was obtained from Chiron Corporation, Emeryville, Calif. It contained amino acid sequence 1-150 (corresponding to the core protein) and 1192-1457 (corresponding to amino acid residues within NS3) of HCV-1 (amino acid sequence available from GenBank®, as described hereinabove).

[0091] ii. HCV C-100 antigen. HCV C-100 recombinant antigen was obtained from Chiron Corporation, Emeryville, Calif. It contained amino acid sequence 1569-1961 (corresponding to amino acid residues within NS4) of HCV-1 (available from GenBank®, as described hereinabove).

[0092] iii. HCV NS5 antigen. HCV NS5 recombinant antigen was obtained from Chiron Corporation, Emeryville, Calif. It contained amino acid sequence 2054-2995 of HCV (available from GenBank®, described hereinabove).

[0093] iv. HCV NS3 NS4 *E. Coli* construct CKS-33c-BCD antigen. HCV HC31 recombinant antigen was obtained from Chiron corporation, Emerville, Calif. It contained amino acid Sequence 1192-1457 of HCV, and amino acid sequence 1676-1931 of the NS4 region. In addition, it consists of 239 amino acids of CKS (available from GenBank®, described hereinabove).

[0094] A1. Preparation of R-Antigen Coated Microparticles.

[0095] i. Preparation of HCV HC43/C100 Microparticles. Microparticles coated with both HC43 and c-100 were prepared in the following manner. Briefly, a 500 μ l aliquot of microparticles (10% weight/volume, 0.7-0.9 micron, available from Seradyn, Indianapolis, Ind.) was mixed with 962 μ l of a coating buffer (Phosphate buffer, pH 5.0 with Tween-20) for approximately 1 minute at room temperature. Then, 154 μ l of an HCV C100-3 antigen solution (0.65 mg/ml) and 308 μ l of an HC43 antigen solution (650 μ g/ml) were added to the microparticle solution, mixed and tumbled for 16 hours at room temperature. The microparticles were pelleted at 12,000 \times g for 10 minutes in an Eppendorf microfuge. The suspension was removed, and the microparticles were washed with wash buffer (Phosphate, NaCl, dithiothreitol-DTT, EDTA, sodium dodecyl sulfate-SDS, pH 6.5) and heat stressed at 56° C. for 20 hours. The microparticles were then resuspended in 2.5 ml of microparticle diluent (Phosphate Buffer, pH 6.5, EDTA, DTT, NaCl and SDS, Sucrose, azide) at a final concentration of 2.0%.

[0096] ii. Preparation of HCV NS5 Microparticles. Five hundred and thirty microliters of an HCV NS5 coating buffer (Carbonate, pH 10, SDS) and 200 μ l of a 10% weight/volume 0.7-0.9 micron microparticles (available from Seradyn, Indianapolis, Ind.) were added to 270 μ l of the HCV

NS5 antigen solution (concentration of 650 $\mu\text{g/ml}$). The microparticles were mixed and tumbled for 16 hours at room temperature. The microparticles were pelleted at 12,000 \times g for 10 minutes in an Eppendorf microfuge. The suspension was removed and the microparticles were washed with wash buffer (Phosphate, NaCl, DTT, EDTA, SDS, pH 6.5) and heat stressed at 56° C. for 20 hours. The washed microparticles were then resuspended in 2.5 ml of microparticle diluent (Phosphate Buffer, pH 6.5, EDTA, DTT, NaCl and SDS, Sucrose, azide) at a final concentration of 0.4%.

[0097] iii. Preparation of HCV NS3 NS4E. *E. Coli* Construct CKS-33c-BCD Microparticles. A 100 μl aliquot of microparticles (10% weight/volume, 0.7-0.9 micron, available from Seradyn, Indianapolis, Ind.) was mixed with 452 μl of coating buffer (Phosphate buffer, pH 5.0 with Tween-20) for approximately 10 minutes at room temperature. Then, 200 μg of CKS-33C-BCD Ag was added and mixed for 16 hours at room temperature.

[0098] The microparticles were pelleted at 12,000 \times g for 10 minutes in an Eppendorf microfuge. The prepared microparticles were washed with wash buffer (DTT, EDTA, SDS in PBS, pH 6.5). The supernatant was removed, and the microparticles were resuspended in 1 ml of microparticle diluent (Phosphate Buffer, pH 6.5, EDTA, DTT, NaCl, Sucrose and SDS, Sucrose).

[0099] iv. Blending of HCV HC43/C100 and HCV NS5 Microparticles. Two hundred twenty microliters of HCV HC43/C100 microparticles prepared as described in Example (IV)(A1)(i) and 330 μl of HCV NS5 microparticles prepared as described in Example (IV)(A1) (ii) were mixed together. This mixture was incubated at room temperature for 15 minutes and diluted to 50 ml in microparticle diluent. (Phosphate Buffer, pH 6.5, EDTA, DTT, NaCl, Sucrose and SDS, Sucrose).

[0100] v. Preparation of Biotinylated Core Peptide. HCV core peptide aa 11-28 was biotinylated at N-terminus during synthesis using an automated peptide synthesizer with $\geq 90\%$ purity.

[0101] vi. Preparation of Streptavidin-Coated Microparticles. A four ml aliquot of carboxylated microparticles (10% weight/volume, 0.227 micron, Seradyn, Indianapolis, Ind.) was mixed with 2486 μl of coupling buffer (MES (2-(N-morpholino) ethanesulfonic acid) pH 6.7) for 10 minutes at room temperature. Then, 114.4 μl of EDAC solution (10 mg/ml in coupling buffer) was added to the microparticle solution and mixed for 15 minutes at room temperature. Subsequently, 1 ml of Streptavidin solution (1 mg/ml in PBS) was added to the activated microparticles and tumbled for 16 hours at room temperature. The prepared microparticles were then pelleted at 12,000 \times g for 3 minutes in an Eppendorf microfuge. The supernatant was removed, and the microparticles were resuspended in 4 ml of PBS. The centrifugation process was repeated one more time, and microparticles were stored in 4 ml of PBS to yield a final concentration of approximately 1%.

[0102] vii. Preparation of Core Peptide Coated Microparticles.

[0103] To 1 ml of coated microparticles from Example (IV) (A1)(vi) was added 375 μl of HCV Core peptide from Example (V)(A1)(v) and 11-28 aa at 1 mg/ml in PBS buffer. The mixture was then incubated for 2 hours at room tem-

perature. The prepared microparticles were washed with wash buffer (DTT, EDTA, SDS in PBS, pH 6.5), and the microparticles were resuspended in 1 ml of microparticle diluent (Calf Bovine Serum, Horse IgG, TWEEN 20, BSA, Casein, EDTA, Sucrose and Proclin, pH 6.5) yielding 1% solids final concentration.

[0104] Microparticles for the Antigen Assay:

[0105] B. Preparation of Monoclonal Antibodies:

[0106] The methods for generating monoclonal antibodies are presented in U.S. Pat. No. 5,753,430. Briefly, *E. coli* derived recombinant antigens encoded by HCV sequences, designated as pHCV34 (HCV-core, a.a. 1-150), were employed as immunogens for antibodies to core. Detailed information on the cloning of pHCV34 is disclosed in U.S. patent application Ser. No. 07/572,822, incorporated herein by reference). The protein was prepared for immunization with appropriate adjuvants after purification, as would be performed by those skilled in the art.

[0107] BALB/c mice were injected intraperitoneally with 15 μg of purified pHCV34 with 15 μl each of Trehalose dimycolate (TDM) and *M. phlei* in a buffer emulsion prepared according to the manufacturer's instructions. Subsequent immunizations were performed on day 14, 28 and 42. Mice were bled on days 21 and 49, and the immune response was studied by enzyme-linked immunosorbent assay utilizing pHCV34 coated on polystyrene beads, as detailed in U.S. Pat. No. 5,753,430.

[0108] Upon demonstration of specific anti-HCV antibody present at reasonable titers in the sera of immunized mice, mice were boosted with 40 μl of pHCV34 antigen. The mice were sacrificed and their spleens were removed; the white cells were mixed and fused with SP2/0 cells. The cell mixture was cultured in Biscoe's Modified Dubach's Medium (IMDM) supplemented with 20% fetal calf serum, and the hybridized cells were selected by using a hypoxanthine and thymidine medium. Hybridoma cell lines were established, and all monoclonal antibodies specific for antibodies to core were prepared from ascite fluids of the mice and were purified by chromatography on a protein-A column (Pharmacia, Uppsala, Sweden). The epitopes of the monoclonal antibodies were analyzed by an ELISA test described in Example I.

[0109] B1. Preparation of Monoclonal Antibody-coated Microparticles for Antigen Assay.

[0110] i. Preparation of HCV C11-14 Microparticles. Briefly, a 1 ml aliquot of carboxylated microparticles (10% weight/volume, 0.227 micron, available from Seradyn, Indianapolis, Ind.) was mixed with 9 ml of coupling buffer (MES (2-(N-morpholino) ethanesulfonic acid), pH 6.7) for approximately 10 minutes at room temperature. Then, 150 μl of 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide (EDAC, 10 mg/ml in coupling buffer, Sigma Chemical Company) was added to the microparticle solution and mixed for 15 minutes at room temperature. Eighteen hundred and twenty-two microliters of C11-14 monoclonal antibody solution (2.13 mg/ml) was added to the activated microparticles, mixed and tumbled for 16 hours at room temperature. The microparticles were then pelleted at 12,000 \times g for 3 minutes in an Eppendorf microfuge. The supernatant was removed, and the microparticles were washed with microparticle wash buffer (Phosphate Buffer

Saline-PBS, Tween 20, pH 7.2), followed by washing with microparticle coating buffer (Tris Buffer Saline-TBS, Bovine Serum Albumin-BSA, pH 7.2) and final washing with microparticle final exchange buffer (PBS, Tween 20, pH 7.2). The microparticles were resuspended in 5 ml of final exchange buffer and heat stressed at 45° C. for 72 hours. After heat stress, 5 ml of microparticle diluent (Calf Bovine Serum, Horse IgG, Tween 20, BSA, Casein, Ethylene diamine tetraacetic acid (EDTA), Sucrose and Proclin, pH 6.5) was added to give a final concentration of approximately 1.0%.

[0111] ii. Preparation of HCV A5 (14-635-225) Mab Microparticles. A similar procedure as mentioned in Example (V)(B1)(i) was used, except that instead of C11-14 Mab, A5 14-635-225) Mab was used for coating on microparticles.

[0112] iii. Preparation of HCV C11-3 Microparticles. Six point six microliters of 1 N HCl was added to 300 μ l (1.45 mg/ml) of C11-3 monoclonal antibody to bring the pH to 2.5. The monoclonal was then incubated at this pH for 5 minutes. The pH was then brought to 6.5 by adding 50 mM MES buffer. A 100 μ l aliquot of carboxylated microparticles (10% weight/volume, 0.227 micron, Seradyn, Indianapolis, Ind.) was then mixed with 333 μ l of coupling buffer (MES, pH 6.7) for 10 minutes at room temperature. Then, 15 μ l of EDAC solution (10 mg/ml in coupling buffer) was added to the microparticle solution and mixed for 5 minutes at room temperature. Five hundred and fifty-two microliters of pH shocked C11-3 monoclonal antibody solution (0.725 mg/ml) were added to the activated microparticles, mixed and tumbled for 16 hours at room temperature. The microparticles were then pelleted at 12,000 \times g for 3 minutes in an Eppendorf microfuge. The supernatant was removed, and the microparticles were washed with microparticle wash buffer (Phosphated Buffered Saline (PBS), Tween 20, pH 7.2), followed by a wash with microparticle coating buffer (Tris Buffered Saline (TBS), Bovine Serum Albumin (BSA), pH 7.2) and a final wash with microparticle final exchange buffer (PBS, Tween 20, pH 7.2). The microparticles were resuspended in 0.5 ml of final exchange buffer and heat stressed at 45° C. for 72 hours. After heat stress, 0.5 ml of microparticle diluent (Calf Bovine Serum, Horse IgG, Tween 20, BSA, Casein, ethylene diamine tetraacetic acid (EDTA), Sucrose, and ProClin, pH 6.5) was added to give a final concentration of approximately 1.0%.

[0113] iv. Blending of HCV C11-14 and C11-3 Microparticles.

[0114] Thirty-six microliters of HCV C11-3 microparticles (1% solid) prepared as described in Example (V)(B1)(iii) and 84 μ l HCV C11-14 microparticles (1% solids) prepared as described in Example (V)(B1)(i) were mixed with 880 μ l microparticle diluent (Calf Bovine Serum, Horse IgG, Tween 20, BSA, Casein, EDTA, Sucrose, and Proclin, pH 6.5).

[0115] C. Preparation of Microparticles for Combo Assay:

[0116] For the dual assay, two separate PRISM® channels were used, one for the HCV antibody assay and one for the HCV antigen assay. For the combo assay, both the antibody and antigen assays were performed on a single channel where the reagents for both antigen and antibody assay were combined in one kit.

[0117] i. Blending of C11-14 mAb Coated Microparticles With Core Antigen (Peptide)-Coated Microparticles and HCV HC33 Antigen Coated Microparticles. Three hundred and fifty microliters of core peptide coated microparticles (1% solids stock) prepared as in Example (V)(A1)(vii) and 700 μ l of HCV NS3 NS4 *E. coli* Construct CKS-33C-BCD Ag coated microparticles (1% solids stock) prepared as described in Example (V)(A1)(iii) and 319 μ l HCV C11-14 microparticles (1.0099% solids stock) prepared as described in Example (V)(B1)(i) were mixed with 5631 μ l microparticle diluent (calf Bovine Serum, Horse IgG, Tween 20, BSA, Casein, EDTA, Sucrose and Proclin, pH 6.5).

[0118] D. Preparation of Microparticles for Combo Assay Using p9MB18-Coated Microparticles, c200-Coated Microparticles, and C11-14-Coated Microparticles:

[0119] i. HCV C-200 antigen. HCV C-200 recombinant antigen was obtained from Chiron Corporation, Emeryville, Calif. In particular, the antigen comprises amino acid sequence 1192-1932 of HCV (available from GenBank, as described hereinabove) and is from the NS3 and NS4 regions. The c200 antigen is a chimeric fusion protein, with 154 amino acids of human superoxide dismutase (hSOD).

[0120] ii. Preparation of c200 Microparticles. A 100 μ l aliquot of microparticles (10% weight/volume, 0.7-0.9 micron, available from Seradyn, Indianapolis, Ind.) was mixed with 830.6 μ l of coating buffer (MES buffer, pH 6.5 with SDS) for approximately 10 minutes at room temperature. Then, 69.4 μ l of c200 antigen solution (0.72 mg/ml) was added to the microparticle solution, mixed and tumbled for 16 hours at room temperature. The microparticles were pelleted at 12,000 \times g for 10 minutes in an Eppendorf microfuge. The suspension was removed, and the microparticles were washed with wash buffer (Phosphate, NaCl, EDTA, SDS, pH 6.5) and heat stressed at 56° C. for 16-20 hours. The microparticles were then resuspended in 500 μ l of microparticle diluent (Phosphate Buffer, pH 6.5, EDTA, NaCl, Sucrose, SDS, azide) at a final concentration of 2.0%.

[0121] iii. Preparation of p9MB-18 Microparticles. A 100 μ l aliquot of microparticles (10% weight/volume, 0.7-0.9 micron (available from Seradyn, Indianapolis, Ind.)) was mixed with 788 μ l of coating buffer (MES buffer, pH 6.5 with SDS) for approximately 10 minutes at room temperature. Then, 112 μ l of an HCV p9MB-18 antigen solution (0.89 mg/ml) was added to the microparticle solution, mixed and tumbled for 16 hours at room temperature. The microparticles were pelleted at 12,000 \times g for 10 minutes in an Eppendorf microfuge. The suspension was removed, and the microparticles were washed with wash buffer (phosphate, NaCl, EDTA, SDS, pH 6.5) and heat stressed at 56° C. for 16-20 hours. The microparticles were then resuspended in 500 μ l of microparticle diluent (Phosphate Buffer, pH 6.5, EDTA, NaCl, Sucrose, SDS, azide) at a final concentration of 2.0%.

[0122] iv. Blending of C11-14 mAb Coated Microparticles with HCV p9MB-18-Coated Microparticles and HCV c200 Antigen-Coated Microparticles. Fifty microliters of p9MB-18 coated microparticles (2% solids stock) prepared as in Example (V)(D)(iii), 62.5 μ l of HCV c200 Ag coated microparticles (2% solids stock) prepared as described in Example (V)(D) ii) and 125 μ l HCV C11-14 microparticles (2% solids stock) prepared as described in Example (V)(B1)(i) were mixed with 4762.5 μ l microparticle diluent (EDTA, SDS, Sucrose and Proclin, PBS pH 6.5).

[0123] Sensitivity: The seroconversion sensitivity was 95.8% as compared to nucleic acid testing data. The

PRISM® HCV Ag/Ab Real Combo assay detected 23/24 positive bleeds as reactive. Data is summarized in Table V, below. Overall sensitivity for the seroconversion panels shown in Table IV are comparable between the blended microparticles prepared as described in Example (V)(C)(i) and microparticles prepared as described in Example (V)(C)(iv) using assay format provided in FIG. 3. However, for anti-HCV Core, specific sample P9MB18-coated microparticles showed significant improvement over core specific peptide (11-28 aa, refer to Table VI), when used in combination in combo assay format (Table VI).

[0124] D. Preparation of p9MB31 coated Microparticles:

[0125] i. A 100 µl aliquot of microparticles (10% weight/volume, 0.7-0.9 micron, available from Seradyn, Indianapolis, Ind.) was mixed with 788 µl of coating buffer (MES

TABLE IV

Panel or Sample	p9MB18 coated uPS (V) (D) (iv) S/CO	p9MB31 coated uP (V) (E) (i) S/CO
Panel A	13.03	9.43
NABI #2	21.9	27.76
NABI #15	28.21	33.04

Panel A is mainly anti-Core panel
NABI #2 is anti-NS3 sample.
NABI #15 is an anti-Core and anti-NS3 sample.

[0127]

TABLE V

Geno-type	Ag earlier than Ab by	Day Number	HCV RNA Test Results	HCV Ab S/CO	HCV Ag S/CO	HCV Ab/Ag Real Combo S/CO
1a	38 Days	0	-	0.09	0.49	0.39
		24	+	0.11	10.63	2.47
		27	+	0.10	43.77	4.66
		31	+	0.11	72.92	8.09
		62	+	5.19	44.41	5.61
		64	+	5.22	69.55	5.86
		69	+	5.91	12.92	4.22
		71	+	6.29	7.09	3.77
1b	18 Days	0	+	0.09	76.61	6.21
		4	+	0.09	56.09	4.29
		7	+	0.08	39.63	2.97
		13	+	0.34	32.14	2.44
		18	+	1.53	14.93	1.68
		21	+	3.20	19.97	2.64
		164	Not Tested	5.86	0.61	3.61
1a	23 Days	0	-	0.09	0.45	0.39
		2	-	0.08	0.45	0.42
		17	+	0.07	20.06	1.58
		19	+	0.09	45.84	4.11
		24	+	0.09	81.03	6.26
		26	+	0.07	63.30	5.95
		36	+	0.31	74.78	8.44
		40	+	4.03	49.53	5.74
1a	32 Days	0	-	0.46	0.54	0.29
		22	-	0.42	0.53	0.43
		24	-	0.43	0.47	0.40
		42	+	0.46	8.79	0.99
		46	+	0.44	22.26	2.02
		74	+	4.22	19.82	2.40
		76	+	4.50	23.78	2.99

buffer, pH 6.5 with SDS) for approximately 10 minutes at room temperature. Then, 112 µl of an HCV p9MB-31 antigen solution (0.89 mg/ml) was added to the microparticle solution, mixed, and tumbled for 16 hours at room temperature. The microparticles were pelleted at 12,000xg for 10 minutes in an Eppendorf microfuge. The suspension was removed, and the microparticles were washed with wash buffer (Phosphate, NaCl, EDTA, SDS, pH 6.5) and heat stressed at 56° C. for 16-20 hours. The microparticles were then resuspended in 500 µl of microparticle diluent (Phosphate Buffer, pH 6.5, EDTA, NaCl, Sucrose, SDS, azide) at a final concentration of 2.0%.

[0126] Sensitivity: The sensitivity of the HCV panels was compared between microparticles coated with p9MB-18 and p9MB-31 r-antigen using the assay format provided in FIG. 1. Overall, sensitivity for both of these r-antigens is comparable as shown in Table IV.

[0128]

TABLE VI

Panel or sample	Peptide blended Combo ups (V) (C) (i) S/CO	P9MB-18 blended Combo ups (V) (D) (iv) S/CO
Panel A	0.95	1.25
NABI #7	0.44	2.27
NABI 0141044662	1.15	2.04

*S/CO values ≥ 1.00 are considered reactive
NABI #7 and NABI 0141044662 are anti-HCV Core samples from NABI.
Panel A is mainly anti-Core panel.

EXAMPLE VI

Preparation of Acridinium-labeled Conjugates

[0129] A. Conjugate for Antibody Assay:

[0130] For the antibody assay, either mouse anti-human IgG directly labeled with acridinium or a pre-complex of biotinylated anti-human F(ab')₂ and acridinium anti biotin conjugate was used.

[0131] i. Pre-complex of biotinylated anti-human F(ab')₂ and acridinium anti-biotin conjugate. The labeled anti-biotin antibody was prepared as disclosed in U.S. Pat. No. 5,705,330. The pre-complex of biotinylated Anti-human F(ab')₂ and acridinium anti-biotin conjugate were also prepared as disclosed in U.S. Pat. No. 5,705,330.

[0132] ii. Acridinium labeled Mouse anti-Human IgG. Fifty-three point six microliters of conjugation buffer (CB) containing sodium phosphate, NaCl, 3-(3-chloramidopropyl)-dimethylammonio-1-propane-sulfonate (CHAPS, Sigma Chemical Company, Saint Louis, Mo.), pH 8.0 and 7.2 μ l of N-hydroxysuccinimide ester of 10-(3-sulfopropyl)-N-tosyl-N-(2-carboxyethyl)-9-acridinium carboxamide (4 mg/ml in dimethyl formamide) was added to 131 μ l of Mouse anti-Human IgG (4.59 mg/ml) and 601 μ l of PBS at room temperature. The reaction mixture was mixed with a rotator for 20 minutes at room temperature. The reaction was quenched by loading the reaction mixture onto the HPLC. This was applied to a 300 \times 7.8 mm Bio-Sil SEC-250 gel filtration column (Bio-Rad, Richmond, Calif.) which had been equilibrated with buffer containing CHAPS, NaCl and sodium phosphate, pH 6.3. The column was eluted at 1.0 ml/min with the same buffer using a Beckman 421A controller equipped with a model 114M pump. Fractions of 1 ml were collected, and the absorbance determined at 280 nm and 370 nm with a Beckman DU-7 spectrophotometer. The extent of acridinium incorporation was calculated as described in U.S. Pat. No. 5,705,330. The acridinium to IgG ratio (mole/mole) obtained was approximately 2.5. The conjugate was stored at 4° C.

[0133] B. Conjugate for Antigen Assay:

[0134] i. Acridinylation of c11-10 conjugate. A similar procedure as mentioned in Example (V)(A)(ii) was used except for the following changes. Seven hundred microliters of conjugate buffer, 300 μ l (1 mg/ml) of C11-10 Mab and 2.9 μ l (4mg/ml) of acridinium derivative were mixed for 10 minutes at room temperature. The acridinium to IgG ratio (mole/mole) obtained was approximately 2.0. The conjugate was stored at 4° C.

[0135] C. Conjugate for Combo Assay:

[0136] Blending of Acridinylated Mouse anti-Human IgG and Acridinylated C11-10 Mab Conjugate:

[0137] Fourteen microliters of Acridinylated Mouse Anti-Human IgG (1 μ g/ml) was mixed with 390 μ l of Acridiny-

lated C11-10 Mab (1.79 μ g/ml) conjugate yielding 2 ng/ml Mouse anti-Human IgG with 100 ng/ml C11-10, incubated for 2 hours and filtered before use. Preparation of Mouse anti-Human IgG and Acridinylated C11-10 Mab conjugate are described in Example (V)(A)(ii) and (VI)(B), respectively.

EXAMPLE VII

Detection of HCV Core Protein by Monoclonal Antibodies

[0138] Since a large number of anti-HCV core monoclonal antibodies were available for use in developing an antigen detection assay, it was necessary to determine which combination of monoclonal antibodies would provide the greatest sensitivity. Because the number of combinations possible when using more than one monoclonal on the solid phase (i.e. for capture) and in the liquid phase (i.e. detection) is extremely large, a simplified "screening" method was used to identify the best performing pair of monoclonals. It was assumed that once the most sensitive pair was identified, other monoclonals could be added to improve assay sensitivity, if necessary.

[0139] In order to identify the best pairs, therefore, monoclonal antibodies were coated onto microparticles or conjugated with acridinium as described in Example IV and V. Screening assays used monoclonal antibody-coated microparticles (0.40 μ m diameter) at a working concentration of 0.09-0.15% solids and conjugated monoclonals at a working concentration of 100-125 ng/ml. For all experiments, the same positive and negative control plasma were used (0.1 ml for each assay). The positive control serum was from an HCV-infected individual who tested negative for HCV antibodies but whose plasma had an HCV RNA titer of 19,000,000 copies per ml. The negative control plasma was from a normal blood donor who was negative for HCV antibodies and RNA. Assays were performed using the instrumentation and operation methods as described in Example IX.

[0140] Table VII shows the mean signal-to-negative (S/N) values obtained upon testing the various pairs of monoclonal antibodies for their ability to detect HCV core antigen in the positive control human plasma (nd: not determined). From this data, it is apparent that some pairs of monoclonal antibodies exhibit greater sensitivity than others and that the sensitivity was dependent upon the proper configuration of the assay. For example, when monoclonal antibody A05 was used as the capture reagent and C11-10 was used as the detection reagent, the resulting S/N value was 150.0; however, when the opposite configuration was used, the resulting S/N value was only 6.8

TABLE VII

		Detection of HCV Core Antigen in Human Plasma by Various Pairs of Anti-core Monoclonal Antibodies													
		CONJUGATE													
Mab		107-35-54	14-635-225	14-726-217	13-975-157	13-959-270	14-178-125	14-886-216	110-81-17	14-945-218	C11-3	C11-7	C11-10	C11-14	C11-15
		Mab													
MICRO-PARTICLE		A04	A05	A06	A08	A09	A13	A14	A15	A218	C11-3	C11-7	C11-10	C11-14	C11-15
14-1708-403	A01	0.6	nd	1.3	2.5	1.0	1.3	5.2	1.0	1.4	2.4	0.9	53.0	1.9	3.5
14-153-462	A02	nd	nd	1.4	nd	0.9	1.5	1.5	nd	2.4	0.6	1.0	22.0	1.8	0.2
14-1269-281	A03	0.8	nd	1.3	2.4	1.1	1.3	7.4	0.9	1.4	2.2	2.9	101.6	3.0	6.0
107-35-54	A04	nd	nd	0.9	nd	nd	nd	nd	nd	nd	1.7	1.2	7.6	12.9	1.3
14-635-225	A05	1.2	nd	1.5	2.9	0.9	1.0	8.4	1.0	1.5	4.2	2.1	150.0	1.9	4.0
14-726-217	A06	nd	nd	0.6	nd	1.6	1.8	2.5	nd	2.7	1.7	2.4	9.5	3.8	0.7
13-975-157	A08	nd	nd	0.7	nd	nd	nd	nd	nd	nd	1.4	0.8	4.6	1.6	3.2
13-959-270	A09	nd	nd	1.6	nd	nd	nd	nd	nd	nd	1.5	1.2	2.4	1.5	1.0
14-1705-255	A10	nd	nd	0.7	nd	1.3	1.7	1.2	nd	2.5	1.6	1.3	13.6	4.5	1.2
14-188-104	A11	nd	nd	0.8	nd	1.0	1.2	2.8	nd	1.2	2.8	0.7	16.5	1.8	2.4
14-1287-252	A12	nd	nd	0.9	nd	nd	nd	nd	nd	nd	2.8	1.9	6.3	2.4	6.1
14-886-216	A14	nd	nd	0.9	nd	nd	nd	nd	nd	nd	3.2	3.1	3.0	4.2	2.1
110-81-17	A15	nd	nd	1.5	nd	nd	nd	nd	nd	nd	1.0	1.9	6.3	2.4	2.3
14-947-104	A16	0.5	nd	1.1	6.1	1.5	1.3	6.1	1.1	1.1	1.6	1.4	69.0	4.7	1.2
C11-3	C11-3	0.6	4.7	nd	nd	1.1	2.1	1.0	nd	0.9	2.0	1.1	11.1	4.6	3.6
C11-7	C11-7	1.2	3.9	nd	nd	1.1	0.8	3.6	nd	1.3	2.5	1.8	6.2	3.4	2.0
C11-10	C11-10	0.6	6.8	nd	nd	0.9	1.4	1.1	nd	1.0	8.2	1.9	4.3	4.6	2.8
C11-14	C11-14	0.8	2.0	1.5	nd	0.8	1.1	13.4	nd	1.7	7.9	0.9	208.0	2.2	14.4
C11-15	C11-15	1.6	4.6	nd	nd	1.2	1.4	1.3	nd	1.5	5.1	1.5	4.9	3.8	5.7

EXAMPLE VIII

[0141] HCV Core Antigen Assay Sample Diluent Buffer The HCV core antigen assay for PRISM®, as described in Example XI, utilizes a sample diluent buffer (SDB) for dilution of the human serum or plasma sample to be tested. The monoclonal antibody-coated microparticles are then added to form a reaction mixture. It is possible that the sensitivity and specificity of the antigen detection assay is affected by the composition of the SDB, in terms of the ingredients and their concentration.

[0142] It was hypothesized that, since HCV is believed to be an enveloped virus, it would be necessary to include detergent (surfactant) in the SDB to remove the lipid envelope, thereby exposing the core protein to solution. In addition, it was surmised that addition of chaotropic salts to the SDB might aid in dissolution of the nucleocapsid complex which could enhance detectability of core antigen.

[0143] To investigate the possible effects of SDB composition on the HCV core antigen assay sensitivity, a series of buffers was prepared and tested in an HCV core antigen assay comprised of monoclonal antibody C11-7 or C11-14 coated microparticles (as stated in the table legends) and acridinium labeled monoclonal antibody C11-10 conjugate. The simplest SDB used (also referred to as basal buffer), in terms of number of components, consisted of 0.1 M potassium phosphate, pH 7.2, 10 mM EDTA. This is the buffer to which detergents and salts were added. The performance of the SDBs was determined by examining their effect on the signal-to-negative (S/N) ratio obtained upon testing of a positive control human plasma from an individual who tested negative for HCV antibodies but whose plasma had an HCV RNA titer of 19,000,000 copies per ml. The negative control plasma was from a normal blood donor who was

negative for HCV antibodies and RNA. Screening assays used coated microparticles at a working concentration of 0.09-0.15% solids and C11-10 conjugate at a working concentration of 100-125 ng/ml. For all experiments, the same positive and negative control plasma was used (0.1 ml for each assay). Assays were performed using the instrumentation and operation methods as described in Example VIII.

[0144] As shown in Table VIII, the S/N value obtained varies greatly depending upon the detergent added to the sample diluent buffer and its concentration. Addition of the zwitterionic surfactant SB-12 (lauryl sulfobetaine) gave the highest S/N values. In addition, as shown in Table IX, the highest S/N values were again seen with SB-12 when compared to other zwitterionic detergents of the same class but with different alkyl chain lengths. Titration of the amount of SB-12 added to the basal buffer in the presence of 0.5% or 2% Triton X-100 is shown in Table X. Increasing the SB-12 concentration over 6% diminished S/N values obtained in the core antigen assay significantly.

[0145] Further experiments examined the effect of the addition of salts or different combinations of zwitterionic or nonionic detergents to the sensitivity of the core antigen assay. Results presented in Tables XI and XII suggest a marginal effect on S/N is observed when KCl is substituted for NaCl, the same is true for the addition of urea. The sample diluent buffers containing SB-16 (palmityl sulfobetaine) appear to exhibit enhanced S/N values. The effect of urea was examined by including increasing concentrations in one of the SDBs that gave reasonably high S/N values compared to the others in a previous experiment (Table XIII). In this particular buffer, addition of urea to a final concentration of 2.0-2.5 M appears to have increased S/N values most significantly.

[0146] The effect on S/N values by the addition of various proteins or serum from nonhuman sources to a sample diluent buffer was also examined (Table XIV). The inclusion of bovine serum albumin, with or without other proteinaceous components, only marginally increased the S/N values obtained upon testing of the HCV positive control serum. In contrast, some combinations of protein or sera actually decreased the S/N value relative to that observed for the protein-free buffer.

[0148]

TABLE X

Detergent/Additive (in basal buffer)	S/N with 0.5% Triton X100	S/N with 2% Triton X100
2% SB-12	13.3	14.8
4% SB-12	14.3	16.6

TABLE VIII

Effect of Detergent on HCV Core Antigen Detection			
Detergent/Additive (in basal buffer)	Acronym	S/N @ 0.5%	S/N @ 2%
Dodecyltrimethylammonium bromide	SB-12	2.4	7.9
1-dodecylpyridinium chloride	DPC	0.1	6.9
Sodium dodecylsulfate	SDS	2.7	5.0
Cholamidopropyltrimethylammonium propanesulfonate	CHAPS	1.0	4.7
3 α ,7 α ,12 α -Trihydroxy-5 β -cholanic acid	Cholate	1.2	2.2
t-Octylphenoxyethoxyethanol	Triton X-100	2.0	1.8
Carboxymethyltrimethylammonium	Betaine	1.9	1.5
Taurocholic acid	TCA	1.1	1.3
Dodecyltrimethylammonium bromide	DTAB	0.2	1.1
Mixture of steroids, polysacc., detergents	Saponin	3.1	0.9
N-Dodecanoyl-N-methylglycine (N-lauroyl sarcosine)	NLS	1.6	0.2
Cetyltrimethylammonium bromide	CTAB	1.6	nd
Tetradecyltrimethylammonium bromide	TDTAB	1.6	nd

Assay format: Anti-HCV core monoclonal C11-7 coated microparticles (0.40 μ m) and C11-10 conjugate.

nd: not determined.

[0147]

TABLE IX

Detergent/Additive (in basal buffer)	Acronym	S/N @ 0.5%	S/N @ 2%
(3-[(3-cholamidopropyl) dimethylammonio]-2-hydroxy-1-propanesulfonate	CHAPSO	2.9	11.5
N-dodecyl-N,N-(dimethylammonio)butyrate	DDMAB	1.9	nd
N-dodecyl-N,N-(dimethylammonio)undecanoate	DDMAU	4.1	9.4
N,N-dimethyldodecylamine-N-Oxide	LDAO	3.8	1.8
N-octyl-N,N-dimethyl-3-ammonia-1-propanesulfonate	SB-8	2.5	3.8
N-decyl-N,N-dimethyl-3-ammonia-1-propanesulfonate	SB-10	3.0	5.2
N-dodecyl-N,N-dimethyl-3-ammonia-1-propanesulfonate	SB-12	16.0	38.3
N-tetradecyl-N,N-dimethyl-3-ammonia-1-propanesulfonate	SB-14	4.5	0.8
N-hexadecyl-N,N-dimethyl-3-ammonia-1-propanesulfonate	SB-16, 0.125%	5.0	nd

Assay format: Anti-HCV core monoclonal C11-14 coated microparticles (0.40 μ m) and C11-10 conjugate.

nd: not determined.

TABLE X-continued

Detergent/Additive (in basal buffer)	S/N with 0.5% Triton X100	S/N with 2% Triton X100
6% SB-12	10.7	15.3
8% SB-12	0.8	0.91

Assay format: Anti-HCV core monoclonal C11-14 coated microparticles (0.40 μ m) and C11-10 conjugate.

[0149]

TABLE XI

SDB Components and Final Concentration							
NaCl	SB-12	SB-16	CTAB	Triton X-100	Urea	Exp 1 S/N	Exp 2 S/N
(no buffer added)						3.85	nd
Water						3.41	4.7
0.5 M						2.84	6.7
0.5 M				1.80%		32.17	34.2
0.5 M	2%					93.4	78.2
0.5 M	2%				2.5 M	106.3	82.9
0.5 M	2%		0.10%			97.1	67.4
0.5 M	2%		0.10%		2.5 M	93.9	86.3
0.5 M	2%		0.10%	1.80%		98.8	79.1
0.5 M	2%		0.10%	1.80%	2.5 M	84.8	90.5
0.5 M		2%			2.5 M	105.2	92.9
0.5 M		2%	0.10%			106.8	108.1
0.5 M		2%	0.10%		2.5 M	142.2	102.3
0.5 M		2%	0.10%	1.80%		115.1	101.6
0.5 M		2%	0.10%	1.80%	2.5 M	nd	103.0

All SDBs were prepared in basal buffer containing 0.1 M potassium phosphate, pH 7.2, 10 mM EDTA. Assay format: Anti-HCV core monoclonal C11-14 coated microparticles (0.227 μ m) and C11-10 conjugate.

[0150]

TABLE XII

SDB Components and Final Concentration						
KCl	SB-12	SB-16	CTAB	Triton X-100	Urea	S/N
0.5 M	2%					48.6
0.5 M	2%				2.5 M	86.8
0.5 M	2%		0.10%			52.8
0.5 M	2%		0.10%		2.5 M	77.6
0.5 M	2%		0.10%	1.80%		83.6
0.5 M	2%		0.10%	1.80%	2.5 M	106.8
0.5 M		2%			2.5 M	122.6
0.5 M		2%	0.10%		2.5 M	136
0.5 M		2%	0.10%	1.80%		127.9
0.5 M		2%	0.10%	1.80%	2.5 M	113.7

All SDBs were prepared in basal buffer containing 0.1 M potassium phosphate, pH 7.2, 10 mM EDTA. Assay format: Anti-HCV core monoclonal C11-14 coated microparticles (0.227 μ m) and C11-10 conjugate.

[0151]

TABLE XIII

Effect of Urea on Antigen Assay Sensitivity								
HCV positive control plasma dilution factor	S/N values at various urea concentrations							
	0.0 M	0.5 M	1.0 M	1.5 M	2.0 M	2.5 M	3.0 M	5.0 M
1:2	35.2	31.3	32.9	33.4	43.0	34.6	33.2	1.1
1:4	19.9	15.7	17.7	13.6	21.5	18.9	15.6	1.1
1:8	9.1	8.2	9.3	6.1	10.5	8.1	8.0	0.9
1:16	5.9	4.3	5.7	5.3	6.6	5.6	5.5	1.0
1:32	3.8	3.2	3.6	3.3	3.6	3.5	3.0	1.0

Urea was added to increasing final concentrations in the following buffer: 0.1 M potassium phosphate, pH 7.2, 10 mM EDTA, 0.5 M NaCl, 2.0% SB-16, 0.10 CTAB, 1.8% Triton X-100. HCV positive control plasma was diluted in negative control plasma. Assay format: Anti-HCV core monoclonal C11-14 coated microparticles (0.227 μ m) and C11-10 conjugate.

[0152]

Table XIV

Effect of Protein or Serum on Antigen Assay Sensitivity		
Component(s) Added (final conc., w/v)	S/N PC 1:2	S/N PC 1:16
No additions	69.3	9.6
1% BSA, 2% mouse serum	73.7	10.3
1% BSA, 0.1% casein	70.5	10.7
1% BSA	70.2	10.3
3% horse serum	68.2	9.6
2% mouse serum	65.6	10.9
0.1% casein	51.7	8.1
2% mouse serum, 0.1% casein	50.9	7.9
1% BSA, 3% horse serum	50.6	8.8
0.1% casein, 3% horse serum	40.7	6.9
2% mouse serum, 3% horse serum	34.3	4.9

Protein or sera were added to buffer containing 100 mM potassium phosphate, pH 7.2, 10 mM EDTA, 0.5 M NaCl, 2% SB-16, 1.10% CTAB, 1.8% Triton X-100, 2.5 M urea. Positive control plasma (PC) was diluted 1:2 or 1:16 in negative control plasma.

EXAMPLE IX

PRISM® HCV Ab, PRISM® HCV Ag, and PRISM® HCV Ab/Ag Combo Assays

[0153] The PRISM® antibody assay is described in U.S. Pat. No. 5,705,330, incorporated herein by reference. and the PRISM® antigen and antibody assays are described in Shah and Stewart, *The Immunoassay Handbook*, second edition, edited by David Wild, p 297-303 (2001), also incorporated herein by reference.

[0154] With respect to the present invention, the following procedures were utilized:

[0155] HCV Ab Assay:

[0156] Assay Format is Provided in FIG. 1.

[0157] Generally, at station 1, 50 μ l of control or sample, 50 μ l of specimen diluent buffer (SDB, Phosphate buffer, pH 7.0 containing Tween 20, newborn calf serum, NaCl, superoxide dismutase (SOD), *E. coli* lysate and azide), and 50 μ l of HCV antigen coated microparticles (prepared as described in Example (V)(A1)(iv) above) were dispensed into each incubation well and assay timing was started.

These were mixed by mutual diffusion of each into the other without external agitation or shaking to form a reaction mixture. At station 4, the reaction mixture was transferred to a detection well which contained a fibrous matrix and washed twice with 300 μ l of transfer wash (TW, containing borate buffer, pH 7.0, with NaCl, Tween-20, Glycerol and Proclin 300). After 18 minutes of incubation at room temperature, 50 μ l of a pre-complexed biotinylated F(ab')₂/acridinium labeled anti-biotin, (biotinylated F(ab')₂ fragment of goat anti-human IgG and acridinium labeled anti-biotin antibody), was dispensed into the matrix of the detection well at station 5. The well was incubated for 23 minutes at 37° C., and the fibrous matrix containing the reaction mixture was washed three times with 100 μ l of FW, containing MES (2-[N-morpholino] ethanesulfonic acid), pH 5.7, with NaCl and Proclin 300 at station 8. At station 9, as in all of the assays described below, a chemiluminescence (CL) signal was generated by addition of an alkaline hydrogen peroxide solution, and the photons were measured by a photo multiplier tube. The amount of light emitted is proportional to the amount of antibody in the sample. The presence or absence of antibody in the sample is determined by comparing the number of photons collected from the sample to a cutoff (S/CO) value determined from a calibration performed in the batch. The results are expressed as S/CO (signal to cutoff) in Table XV below. The cutoff value is calculated by the sum of product of the average chemiluminescence counts of the positive control (n=4) times 0.55 plus the average chemiluminescence counts of the negative control (n=4).

[0158] Sensitivity: The seroconversion sensitivity was 100% as compared to the HCV RNA data provided in vendor certificate of analysis data for the selected seroconversion panels tested. Data is summarized in Table XV.

[0159] Specificity: Based on repeat reactive rates, the specificity of the HCV Ab assay was >99% with the population tested (Table XVII).

[0160] HCV Ag Assay:

[0161] Assay Format is Provided in FIG. 2.

[0162] Mab C11-14/Mab C11-10 pair:

[0163] Generally, at station 1, 100 μ l of control or sample, 50 μ l of specimen diluent buffer (SDB, Sodium phosphate, EDTA, Triton X-100, Urea and sodium azide), and 50 μ l of HCV Mab coated microparticles (prepared as described in Example (V)(B1)(i)) were dispensed into each incubation well and the assay timing was started. These were mixed by mutual diffusion of each into the other without external agitation or shaking to form a reaction mixture. At station 4, the reaction mixture was transferred to a detection well which contained a fibrous matrix and washed twice with 300 μ l of transfer wash (TW) (MES, NaCl, Triton X-100, PEG, Antifoam, Proclin 300, pH 5.6) after 18 minutes of incubation at room temperature. At station 5, 50 μ l of acridinylated C11-10 Mab conjugate (as mentioned in Example (VI)(B)) was dispensed into the matrix of the detection well. The contents of the well were incubated for 23 minutes, and the fibrous matrix containing the reaction mixture was washed one time with 200 μ l of final wash (FW) (Tris buffer with LiCl and LDS) followed by three times with 100 μ l of FW. The CL signal was triggered and measured at station 9. The results are expressed as S/CO (signal to cutoff) in Table XV. The cutoff value is 2.2 times the average chemiluminescence count of the negative control (n=5).

[0164] Sensitivity: Two groups of commercially available seroconversion panels containing serially collected samples from individuals who developed antibodies to HCV were tested in the prototype PRISM® HCV antigen test and in the PRISM® HCV antibody test. For the first group of seroconversion samples, the first available bleed date was negative for HCV RNA. In subsequent bleed dates, HCV RNA was detected for one or more bleed dates, followed in all cases by detection of antibodies to HCV. For the second group of seroconversion panels, the first bleed date was already positive for HCV RNA; antibodies to HCV were detected in subsequent bleed dates. For the two groups seroconversion sensitivity was 98.5% as compared to data obtained by HCV RNA testing. The PRISM® HCV Ag detected 67/68 HCV RNA positive bleeds as reactive. Data is summarized in Table XVI. These data indicate that HCV Ag testing detects HCV infection in individuals who have not yet mounted an antibody response.

[0165] Specificity: Based on repeat reactive rates, the specificity of the HCV Ag assay was >99% with the population tested (Table XVIII).

[0166] HCV Ag Assay: Mab A5 (14-635-225)/Mab C11-10 pair:

[0167] An assay procedure, similar to that mentioned for C11-14/C11-10, was used. The only difference was that the test used Abbott A5 (14-635-225) Mab coated microparticles instead of Mab C11-14 coated microparticles.

[0168] Sensitivity: A total of 4 seroconversion panels were evaluated and sensitivity was compared with the data generated using C11-14/C11-10 pair. Both these pairs detected the same number of positive bleeds. Sensitivity data for A5 (14-635-225)/C11-10 pair is summarized in Table XVI.

[0169] Specificity: Based on repeat reactive rates, the specificity of the HCV Ag assay was 100% with the mini population (n=100) tested (Table XIX).

[0170] HCV Ag Assay: Mab C11-14 and C11-3/Mab C11-10 pair:

[0171] An assay procedure, similar to that mentioned for C11-14/C11-10, was used. The only difference was the use of C11-14 and C11-3 blended microparticles (Example (IV)(B1)(iv)) instead of Mab C11-14 coated microparticles.

[0172] Sensitivity: The performance of this pair was assessed by comparing the S/N ratio against panels consisting of recalcified human plasma positive for HCV Core antigen (termed "PC") and a human plasma negative for HCV antigens and antibodies (termed "NC") (Table XX). The S/N was determined by the formula:

$$S/N = \text{Average of PC} / \text{Average of NC}$$

[0173] The average chemiluminescence counts of four specimens were used to determine each average.

[0174] PRISM® HCV Ag/Ab Combo Assay:

[0175] Two different formats (i.e., Dual Combo assay and Real combo assay) were evaluated on the PRISM® HCV Ag/Ab Assay as follows:

[0176] Dual Combo Assay: The HCV Ag/Ab dual combo assay is run simultaneously on PRISM® using two different channels. A total of six channels in PRISM® are used simultaneously to run several assays (HIV, HBcore, HBsAg, HTLV, and HCVAb) five of which are currently in use, while

one channel remains open for new markers (e.g. HCV Ag assay) or can be reserved in case one of the channels become problematic. Thus, by using one channel for an HCV Ag assay and five other channels for five other assays, a reserve channel is not available for use.

[0177] The PRISM® HCV Ab and PRISM® HCV Ag assays were performed individually. The results from both assays were combined to produce a single, final report.

[0178] Real Combo Assay: The PRISM® HCV Ag and HCV Ab assays were combined and performed as a single assay in one of the PRISM® channels.

[0179] PRISM® HCV Ag/Ab Dual Combo Assay:

[0180] Sensitivity: The seroconversion sensitivity of the HCV Ab/Ag Dual Combo assay was 98.5%. Data is summarized in Table XV.

[0181] Specificity: Based on repeat reactive rates, the specificity of the HCV Ab/Ag Dual Combo assay was >99% with the population tested (Table XXI).

[0182] PRISM® HCV Ag/Ab Real Combo Assay:

[0183] Assay format is provided in FIG. 3. The 2 step PRISM® HCV Combo assay was performed as mentioned above for the HCV Ab or Ag assay with the following changes: At station 1, 100 μ l of control or sample, 50 μ l of specimen diluent buffer (Phosphate buffer, pH 7.0 containing Tween 20, newborn calf serum, NaCl, Tween-20, super-

oxide dismutase (SOD), *E. coli* lysate and azide), and 50 μ l of HCV antigen and Mab blended microparticles (prepared as described in Example (V)(C)(i) above) were dispensed into each incubation well and the assay timing was started. At station 4, the reaction mixture was transferred to a detection well which contained a fibrous matrix and washed twice with 300 μ l of transfer wash (MES, NaCl, Triton X-100, PEG, Antifoam, Proclin 300, pH 5.6). After 18 minutes of incubation at 37 degree C., 50 μ l of blended conjugate acridinylated C11-10 and Acridinylated Mouse anti-Human IgG (prepared as described in Example (VI)(C)) was dispensed into the matrix of the detection well at station # 5. The well was incubated for 23 minutes, and the fibrous matrix containing the reaction mixture was washed three times with 100 μ l of final wash (Tris buffer with LiCl and LDS). The CL signal was triggered and measured at station 9. The results are expressed as S/CO (signal to cutoff) in Table XV below. The cutoff value is 2.2 times the average chemiluminescence count of the negative control (n=3).

[0184] Sensitivity : The seroconversion sensitivity was 95.8% as compared to PCR data. The PRISM® HCV Ag/Ab Real Combo assay detected 23/24 positive bleeds as reactive. Data is summarized in Table XVI.

[0185] Specificity : Based on repeat reactive rates the specificity of HCV Ag/Ab Real Combo assay was 100% with the population tested (Table XX):

TABLE XV

Genotype	Ag earlier than Ab by	Day Number	HCV RNA		HCV Ab S/CO	HCV Ag S/CO	HCV Ab/Ag Dual Combo S/CO
			Test Results				
1a	38 Days	0	-		0.09	0.49	NEG
		24	+		0.11	10.63	POS
		27	+		0.10	43.77	POS
		31	+		0.11	72.92	POS
		62	+		5.19	44.41	POS
		64	+		5.22	69.55	POS
		69	+		5.91	12.92	POS
		71	+		6.29	7.09	POS
1a	20 Days	0	+		0.10	107.78	POS
		3	+		0.11	97.81	POS
		10	+		0.15	63.52	POS
		20	+		1.57	Not Tested	POS
1b	18 Days	0	+		0.09	76.61	POS
		4	+		0.09	56.09	POS
		7	+		0.08	39.63	POS
		13	+		0.34	32.14	POS
		18	+		1.53	14.93	POS
		21	+		3.20	19.97	POS
2b	14 Days	0	+	Not Tested	5.86	0.61	POS
		2	+		0.10	8.41	POS
		7	+		0.40	23.62	POS
		9	+		2.62	14.09	POS
		9	+		3.09	16.67	POS
		14	+		3.99	4.74	POS
1a	42 Days	0	-		0.09	0.76	NEG
		2	-		0.08	0.47	NEG
		7	-		0.09	0.41	NEG
		9	-		0.09	0.40	NEG
		15	-		0.08	0.42	NEG
		17	-		0.08	0.49	NEG
		22	-		0.07	0.49	NEG
		24	-		0.08	0.45	NEG
29	-		0.09	0.52	NEG		

TABLE XV-continued

Genotype	Ag earlier than Ab by	Day Number	HCV RNA	HCV Ab S/CO	HCV Ag S/CO	HCV Ab/Ag
			Test Results			Dual Combo S/CO
		31	-	0.08	0.51	NEG
		36	-	0.08	0.57	NEG
		38	-	0.09	0.51	NEG
		43	-	0.08	0.52	NEG
		45	-	0.09	0.43	NEG
		50	-	0.09	0.69	NEG
		52	-	0.08	0.52	NEG
		57	-	0.09	0.49	NEG
		64	-	0.09	0.90	NEG
		67	-	0.08	0.52	NEG
		74	-	0.09	0.51	NEG
		79	-	0.08	0.52	NEG
		84	-	0.09	0.40	NEG
		105	-	0.09	0.45	NEG
		108	-	0.08	0.66	NEG
		112	-	0.09	0.57	NEG
		119	-	0.09	0.45	NEG
		121	-	0.09	0.42	NEG
		140	+	0.12	10.65	POS
		143	+	0.09	3.81	POS
		147	+	0.10	9.30	POS
		150	+	0.09	34.08	POS
		154	+	0.09	58.01	POS
		157	+	0.09	80.90	POS
		161	+	0.08	107.11	POS
		164	+	0.09	114.41	POS
		168	+	0.09	93.29	POS
		171	+	0.10	89.06	POS
		182	+	1.83	63.62	POS
		186	+	4.39	68.72	POS
		189	+	5.20	119.62	POS
1a	21 Days	0	+	0.08	63.86	POS
		4	+	0.07	50.76	POS
		17	+	0.09	73.66	POS
		21	+	1.02	45.58	POS
		25	+	4.03	60.94	POS
		29	+	5.08	47.38	POS
1a	37 Days	0	+	0.08	33.66	POS
		2	+	0.06	30.83	POS
		7	+	0.08	30.10	POS
		9	+	0.07	39.66	POS
		26	+	0.07	25.51	POS
		32	+	0.12	15.29	POS
		37	+	2.43	15.51	POS
		41	+	3.36	3.10	POS
1a	28 Days	0	+	0.09	67.75	POS
		2	+	0.09	87.93	POS
		10	+	0.10	36.53	POS
		12	+	0.10	60.67	POS
		19	+	0.10	39.62	POS
		21	+	0.11	26.25	POS
		28	+	2.78	9.94	POS
		30	+	4.00	17.02	POS
		35	+	4.71	15.26	POS
		37	Not tested	4.84	13.02	POS
1a	25 Days	0	+	0.15	4.73	POS
		2	+	0.40	6.63	POS
		8	+	0.16	7.48	POS
		10	+	0.11	5.20	POS
		16	+	0.17	7.60	POS
		18	+	0.11	7.58	POS
		23	+	0.64	8.66	POS
		25	+	2.11	9.58	POS
		30	+	2.76	6.21	POS
		32	+	3.39	7.84	POS
		49	+	6.12	1.83	POS
		53	+	6.13	1.93	POS
		56	+	6.34	1.63	POS
1a	28 Days	0	-	0.13	1.27	POS
		2	-	0.23	0.52	NEG
		8	-	0.09	0.50	NEG

TABLE XV-continued

Genotype	Ag earlier than Ab by	Day Number	HCV RNA	HCV Ab S/CO	HCV Ag S/CO	HCV Ab/Ag
			Test Results			Dual Combo S/CO
		11	+	0.10	0.54	NEG
		15	+	0.10	1.92	POS
		18	+	0.11	1.90	POS
		28	+	0.12	2.42	POS
		30	+	0.10	7.04	POS
		35	+	0.14	6.01	POS
		37	+	0.98	13.68	POS
		43	+	4.74	10.07	POS
		46	+	5.27	4.91	POS
1a	13 Days	0	+	0.10	1.82	POS
		3	+	0.10	1.72	POS
		5	+	0.15	1.35	POS
		11	+	0.97	1.70	POS
		13	+	1.26	3.63	POS
		19	-	3.70	2.94	POS
		25	-	4.89	2.43	POS
		27	-	5.20	1.61	POS
		32	-	5.61	1.35	POS
		35	-	5.86	1.30	POS
		41	-	6.11	0.88	POS
		45	-	5.69	0.67	POS
		48	-	5.95	1.94	POS
1a	23 Days	0	-	0.09	0.45	NEG
		2	-	0.08	0.45	NEG
		17	+	0.07	20.06	POS
		19	+	0.09	45.84	POS
		24	+	0.09	81.03	POS
		26	+	0.07	63.30	POS
		36	+	0.31	74.78	POS
		40	+	4.03	49.53	POS
1a	33 Days	0	-	0.09	0.54	NEG
		3	-	0.08	0.50	NEG
		7	-	0.08	0.58	NEG
		12	-	0.09	0.54	NEG
		14	-	0.09	0.56	NEG
		19	-	0.10	0.53	NEG
		25	-	0.10	0.51	NEG
		28	-	0.09	0.50	NEG
		32	-	0.09	0.50	NEG
		35	-	0.10	0.35	NEG
		39	-	0.09	0.53	NEG
		45	+	0.10	9.43	POS
		47	+	0.11	42.00	POS
		52	+	0.11	28.05	POS
		56	+	0.11	25.63	POS
		60	+	0.09	78.15	POS
		73	+	0.18	9.54	POS
		78	+	1.83	5.34	POS
		80	+	2.13	3.40	POS
1a	32 Days	0	-	0.46	0.54	NEG
		22	-	0.42	0.53	NEG
		24	-	0.43	0.47	NEG
		42	+	0.46	8.79	POS
		46	+	0.44	22.26	POS
		74	+	4.22	19.82	POS
		76	+	4.50	23.78	POS
3a	141 Days	0	+	0.12	2.81	POS
		4	+	0.49	1.95	POS
		11	+	2.48	1.41	POS
		13	+	2.54	1.41	POS
		44	+	4.29	0.43	POS
		46	+	4.68	0.43	POS
1a	38 Days	0	-	0.09	0.49	0.38
		24	+	0.11	10.63	1.60
		27	+	0.10	43.77	5.32
		31	+	0.11	72.92	11.04
		62	+	5.19	44.41	10.95
		64	+	5.22	69.55	10.16
		69	+	5.91	12.92	4.51
		71	+	6.29	7.09	3.68

TABLE XV-continued

Genotype	Ag earlier than Ab by	Day Number	HCV RNA		HCV Ab/Ag	
			Test Results	HCV Ab S/CO	HCV Ag S/CO	Dual Combo S/CO
1b	18 Days	0	+	0.09	76.61	14.09
		4	+	0.09	56.09	5.59
		7	+	0.08	39.63	4.75
		13	+	0.34	32.14	3.81
		18	+	1.53	14.93	4.16
		21	+	3.20	19.97	3.81
1a	23 Days	164	Not Tested	5.86	0.61	2.72
		0	-	0.09	0.45	0.31
		2	-	0.08	0.45	0.48
		17	+	0.07	20.06	1.36
		19	+	0.09	45.84	2.72
		24	+	0.09	81.03	4.57
		26	+	0.07	63.30	4.15
		36	+	0.31	74.78	6.90
1a	32 Days	40	+	4.03	49.53	5.11
		0	-	0.46	0.54	0.48
		22	-	0.42	0.53	0.31
		24	-	0.43	0.47	0.32
		42	+	0.46	8.79	0.99
		46	+	0.44	22.26	2.11
		74	+	4.22	19.82	3.99
		76	+	4.50	23.78	3.27

* S/CO values ≥ 1.00 are considered reactive.

[0186]

TABLE XVI

Genotype	Ag earlier than Ab by	Day Number	HCV RNA		HCV Ag		HCV Ab/Ag Dual Combo S/CO
			Test Results	HCV Ab S/CO	c11-14/c11-10 S/CO	HCV Ag A5/c11-10 S/CO	
1a	38 Days	0	-	0.09	0.49	0.53	NEG
		24	+	0.11	10.63	6.47	POS
		27	+	0.10	43.77	29.30	POS
		31	+	0.11	72.92	64.97	POS
		62	+	5.19	44.41	37.22	POS
		64	+	5.22	69.55	39.55	POS
		69	+	5.91	12.92	7.54	POS
		71	+	6.29	7.09	4.42	POS
1b	18 Days	0	+	0.09	76.61	68.95	POS
		4	+	0.09	56.09	48.69	POS
		7	+	0.08	39.63	57.02	POS
		13	+	0.34	32.14	34.85	POS
		18	+	1.53	14.93	16.61	POS
		21	+	3.20	19.97	13.78	POS
		164	Not Tested	5.86	0.61	0.62	NEG
1a	23 Days	0	-	0.09	0.45	0.80	NEG
		2	-	0.08	0.45	0.43	NEG
		17	+	0.07	20.06	7.77	POS
		19	+	0.09	45.84	17.28	POS
		24	+	0.09	81.03	30.76	POS
		26	+	0.07	63.30	27.95	POS
		36	+	0.31	74.78	36.66	POS
		40	+	4.03	49.53	21.78	POS
1a	32 Days	0	-	0.46	0.54	0.45	NEG
		22	-	0.42	0.53	0.47	NEG
		24	-	0.43	0.47	0.53	NEG
		42	+	0.46	8.79	5.04	POS
		46	+	0.44	22.26	14.18	POS
		74	+	4.22	19.82	14.15	POS
		76	+	4.50	23.78	8.54	POS

* S/CO values ≥ 1.00 are considered reactive.

[0187]

TABLE XVII

N Tested	989*
RR	9
RRR	0.91
N negative	980
Mean S/CO Neg	0.12
Pop	
SD	0.09
SD to CO	9.76

[0188]

TABLE XVIII

N Tested	989*
RR	1
RRR	0.10
N negative	988
Mean S/CO Neg	0.44
Pop	
SD	0.08
SD to CO	7.15

*Volunteer blood donors

[0189]

TABLE XIX

N Tested	100*
RR	0
RRR	0.00
N negative	100
Mean S/CO Neg	0.56
Pop	
SD	0.15
SD to CO	2.98

*Volunteer blood donor

[0190]

TABLE XX

Panel	C11-10 conjugate:		
	C11-14 microparticle only	C11-3 microparticle only	C11-14 + C11-3 Blended microparticle
PC S/N	213.7	87.8	179.4
NC	78.8	90	93.25
Counts			

[0191]

TABLE XXI

N Tested	989*
RR	9
RRR	0.91
N negative	980
Mean S/CO Neg	0.12
Pop	
SD	0.09
SD to CO	9.76

[0192]

TABLE XXII

N Tested	92*
RR	0
RRR	0
N negative	92
Mean S/CO Neg	0.49
Pop	
SD	0.15
SD to CO	3.4

*Volunteer blood donors

EXAMPLE X

Binding of Monoclonal Antibodies to Three Peptides

[0193] Three new peptides were synthesized, two of which are compatible with an HCV Ab/Ag combo format and one suitable for use as a control. Each peptide was synthesized with a N-terminal biotin for ease of tracking during preparation of the solid phase. Peptide aa 10-53 (ALAM-17) (KTKRNTNRRPQDVKFPGGGQIVG-GVYLLPRRGPRLGVRATRKT) (SEQ ID NO: 37) contains HCV core amino acids 10-53 with no intervening deletions. Peptide aa 10-53³²⁻⁵⁰ (ALAM-18) (KTKRNTNRRPQDVKFPGGGQIVKTS) (SEQ ID NO: 38) contains HCV core amino acids 10-53 with a 19 amino acid deletion encompassing aa 32-50. Peptide aa 10-53³¹⁻³³47-48 contains HCV core amino acids 10 through 53 where amino acids 31-33 and 47-48 were deleted (ALAM16)(KTKRNTNRRPQDVKFPGGGQIVYLLPRRGPRLGVTRKT) (SEQ ID NO: 35). Peptides ALAM-16 and ALAM-18 are both compatible with an HCV Ab/Ag combo format. In the case of ALAM-16, the deletion of amino acids 31-33 prevents monoclonal c11-10 (epitope 32-36) from binding to the antigen, and the deletion of amino acids 47 and 48 prevents binding of the c11-14 monoclonal (epitope 45-50). ALAM-18 contains a deletion that encompasses both the c11-10 and c11-14 binding regions. Data showing the lack of binding to ALAM-16 and ALAM-18 by monoclonals c11-10 and c11-14 are shown in the following table:

TABLE XXIII

Mono-clonal	Epitope (aa)	ALAM-17 Peptide aa 10-53 S/N	ALAM-18 Peptide aa 10-5332-50 S/N	ALAM-16 Peptide aa 10-5331-3347-48 S/N
C11-15	19-27	>125.0	>200.0	>153.8
C11-10	32-36	>125.0	0.9	0.4
C11-14	45-50	>125.0	0.4	0.4
C11-3	104-110	1.0	0.4	0.4
C11-7	112-124	0.9	0.5	0.4

*All monoclonals were run at a concentration of 1 ug/ml.

[0194] In addition, 254 HCV genotyped seropositive specimens, representing HCV genotypes 1, 2, 3, 4 and 6, were tested by peptide aa 10-53³¹⁻³³47-48 (ALAM-16) to determine feasibility of this peptide as an antigenic target in an HCV Ab/Ag assay format. All 254 (100%) specimens were reactive toward this peptide. Thus, the deletions present in this peptide, which serve to eliminate binding by the monoclonals necessary for Ag detection, do not negatively impact reactivity of antibodies toward the remaining core epitopes.

 SEQUENCE LISTING

<160> NUMBER OF SEQ ID NOS: 63
 <210> SEQ ID NO 1
 <211> LENGTH: 40
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

 <400> SEQUENCE: 1

 atagaattcc atgcagaaaa aaaacaaacg taacaccaac 40

 <210> SEQ ID NO 2
 <211> LENGTH: 46
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

 <400> SEQUENCE: 2

 cggctgagaa cgttcagagg ttttaacgat ctgaccacca cccggg 46

 <210> SEQ ID NO 3
 <211> LENGTH: 24
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

 <400> SEQUENCE: 3

 aaaacctctg aacgttctca gccg 24

 <210> SEQ ID NO 4
 <211> LENGTH: 36
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

 <400> SEQUENCE: 4

 tatggatcct tattacggag acagcagcca accagc 36

 <210> SEQ ID NO 5
 <211> LENGTH: 244
 <212> TYPE: DNA
 <213> ORGANISM: Hepatitis C Virus

 <400> SEQUENCE: 5

 gaattccatg cagaaaaaaaa acaaacgtaa caccaaccgt cgtccgcag acgttaaatt 60
 cccgggtggt ggtcagatcg ttaaacctc tgaacgttct cagccgcgtg ggcgtcgtca 120
 gccgatcccg aaagctcgtc gtccggaagg tcgtacctgg gctcagccgg gttaccctgtg 180
 gccgctgtac ggtaacgaag gttgcggttg ggcaggttg ctgctgtctc cgtaataagg 240
 atcc 244

 <210> SEQ ID NO 6
 <211> LENGTH: 75
 <212> TYPE: PRT
 <213> ORGANISM: Hepatitis C Virus

 <400> SEQUENCE: 6

-continued

```

Met  Gln  Lys  Lys  Asn  Lys  Arg  Asn  Thr  Asn  Arg  Arg  Pro  Gln  Asp  Val
 1          5          10          15

Lys  Phe  Pro  Gly  Gly  Gln  Ile  Val  Lys  Thr  Ser  Glu  Arg  Ser  Gln
 20          25          30

Pro  Arg  Gly  Arg  Arg  Gln  Pro  Ile  Pro  Lys  Ala  Arg  Arg  Pro  Glu  Gly
 35          40          45

Arg  Thr  Trp  Ala  Gln  Pro  Gly  Tyr  Pro  Trp  Pro  Leu  Tyr  Gly  Asn  Glu
 50          55          60

Gly  Cys  Gly  Trp  Ala  Gly  Trp  Leu  Leu  Ser  Pro
65          70          75
    
```

```

<210> SEQ ID NO 7
<211> LENGTH: 280
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Variant - HVC-Core Recombinant
    
```

```

<400> SEQUENCE: 7
gaattccatg cagaaaaaaa acaaacgtaa caccaaccgt cgtccgcagg acgttaaatt      60
cccgggtggt ggtcagatcg ttggtctgct gccgcgtcgt ggtccgcgtc tgggtcgtaa      120
aacctctgaa cgttctcagc cgcgtgggcg tcgtcagccg atcccgaag ctgctcgtcc      180
ggaaggtcgt acctgggctc agccgggtta cccgtggccg ctgtacggta acgaaggttg      240
cggttgggct ggttggtcgc tgtctccgta ataaggatcc      280
    
```

```

<210> SEQ ID NO 8
<211> LENGTH: 87
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Variant - HVC-Core Recombinant
    
```

```

<400> SEQUENCE: 8

Met  Gln  Lys  Lys  Asn  Lys  Arg  Asn  Thr  Asn  Arg  Arg  Pro  Gln  Asp  Val
 1          5          10          15

Lys  Phe  Pro  Gly  Gly  Gly  Gln  Ile  Val  Gly  Leu  Leu  Pro  Arg  Arg  Gly
 20          25          30

Pro  Arg  Leu  Gly  Arg  Lys  Thr  Ser  Glu  Arg  Ser  Gln  Pro  Arg  Gly  Arg
 35          40          45

Arg  Gln  Pro  Ile  Pro  Lys  Ala  Arg  Arg  Pro  Glu  Gly  Arg  Thr  Trp  Ala
 50          55          60

Gln  Pro  Gly  Tyr  Pro  Trp  Pro  Leu  Tyr  Gly  Asn  Glu  Gly  Cys  Gly  Trp
65          70          75          80

Ala  Gly  Trp  Leu  Leu  Ser  Pro
 85
    
```

```

<210> SEQ ID NO 9
<211> LENGTH: 52
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer
    
```

```

<400> SEQUENCE: 9
accagagcgc ggaccacgac gcggcagcag accaacgatc tgaccaccac cc      52
    
```

-continued

<210> SEQ ID NO 10
 <211> LENGTH: 48
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer
 <400> SEQUENCE: 10
 ccgcgctcgtg gtccgctct gggtcgtaaa acctctgaac gttctcag 48

<210> SEQ ID NO 11
 <211> LENGTH: 301
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Variant - HVC-Core Recombinant
 <400> SEQUENCE: 11
 gaattccatg cagaaaaaaaa acaaacgtaa caccaaccgt cgtccgcagg acgttaaatt 60
 cccgggtggt ggtcagatcg ttggtggtgt ttacgttctg ccgctcgtg gtcgctcgtc 120
 ggggtgttctg gctacgcgta aaacctctga acgttctcag ccgctgggc gtcgctcagcc 180
 gatccgaaa gctcgtcgtc cggaaggctg tacctgggct cagccgggtt acccgtggcc 240
 gctgtacggt aacgaaggtt gcggttgggc tggttggctg ctgtctccgt aataaggatc 300
 c 301

<210> SEQ ID NO 12
 <211> LENGTH: 94
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Variant - HVC-Core Recombinant
 <400> SEQUENCE: 12
 Met Gln Lys Lys Asn Lys Arg Asn Thr Asn Arg Arg Pro Gln Asp Val
 1 5 10 15
 Lys Phe Pro Gly Gly Gly Gln Ile Val Gly Gly Val Tyr Val Leu Pro
 20 25 30
 Arg Arg Gly Pro Arg Leu Gly Val Leu Ala Thr Arg Lys Thr Ser Glu
 35 40 45
 Arg Ser Gln Pro Arg Gly Arg Arg Gln Pro Ile Pro Lys Ala Arg Arg
 50 55 60
 Pro Glu Gly Arg Thr Trp Ala Gln Pro Gly Tyr Pro Trp Pro Leu Tyr
 65 70 75 80
 Gly Asn Glu Gly Cys Gly Trp Ala Gly Trp Leu Leu Ser Pro
 85 90

<210> SEQ ID NO 13
 <211> LENGTH: 45
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer
 <400> SEQUENCE: 13
 acccagacgc ggaccacgac gcggcagaac gtaaacacca ccaac 45

<210> SEQ ID NO 14

-continued

<211> LENGTH: 45
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

<400> SEQUENCE: 14

ccgcgtcgtg gtccgcgtct ggggtgttctg gctacgcgta aaacc 45

<210> SEQ ID NO 15
 <211> LENGTH: 300
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Variant - HCV-Core Recombinant
 <220> FEATURE:
 <221> NAME/KEY: misc_feature
 <222> LOCATION: (100)...(100)
 <223> OTHER INFORMATION: n = a or g or c or t/u, unknown, or other at position 100

<400> SEQUENCE: 15

gaattccatg cagaaaaaaaa acaaacgtaa caccaaccgt cgtccgcagg acgttaaatt 60
 cccgggtggt ggtcagatcg ttggtggtgt ttacctgctn ccgcgtcgtg gtccgcgtct 120
 ggggtgttctg gctacgcgta aaacctctga acgttctcag ccgcgtgggc gtcgtcagcc 180
 gatccgaaag ctcgctgctc ggaaggtcgt acctgggctc agccgggtta cccgtggccg 240
 ctgtacggta acgaaggttg cggttgggct ggttggtctgc tgtctccgta ataaggatcc 300

<210> SEQ ID NO 16
 <211> LENGTH: 94
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Variant - HCV-Core Recombinant

<400> SEQUENCE: 16

Met Gln Lys Lys Asn Lys Arg Asn Thr Asn Arg Arg Pro Gln Asp Val
 1 5 10 15
 Lys Phe Pro Gly Gly Gly Gln Ile Val Gly Gly Val Tyr Leu Leu Pro
 20 25 30
 Arg Arg Gly Pro Arg Leu Gly Val Arg Ala Thr Arg Lys Thr Ser Glu
 35 40 45
 Arg Ser Gln Pro Arg Gly Arg Arg Gln Pro Ile Pro Lys Ala Arg Arg
 50 55 60
 Pro Glu Gly Arg Thr Trp Ala Gln Pro Gly Tyr Pro Trp Pro Leu Tyr
 65 70 75 80
 Gly Asn Glu Gly Cys Gly Trp Ala Gly Trp Leu Leu Ser Pro
 85 90

<210> SEQ ID NO 17
 <211> LENGTH: 1069
 <212> TYPE: DNA
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: p9MB-18

<400> SEQUENCE: 17

gaattccatg gctggtgact ttatcccggg tgaaaaatctc gagactacta tgcgtttctcc 60
 ggttttcact gacaactctt ctccgccggg tgttccgcag tctttccagg ttgctcacct 120

-continued

```

gcatgctccg actggttctg gtaaattctac taaagttcca gctgcttacg ctgctcaggg 180
ttacaaagtt ctgggtctga acccgtctgt tgctgctact ctgggtttcg ggcctacat 240
gtctaaagct cacggtatcg acccgaacat tcgtactggg gtacgtacta tcactactgg 300
ttctccgac acttactcta cttacggtaa attcctggct gacggtgggt gctctggtg 360
tgcttacgat atcatcatct ggcacgaatg ccaactctact gacgctactt ctatcctggg 420
tatcgggtacc gttctggacc aggctgaaac tgcagggtgct cgtctgggtg ttctggctac 480
tgctactccg ccgggttctg ttactgttcc gcacccgaac atcgaagaag ttgctctgtc 540
gactactggg gaaatcccg tctacggtaa agctatcccg ctcgagggta tcaaagggtg 600
tcgtcacctg attttctgcc actctaaaaa aaaatgacgac gaactggctg ctaagcttgt 660
tgctctgggt atcaacgctg ttgcttacta ccgtggctg gacgtttctg ttatcccgac 720
ttctgggtgac gttgttggtg tggccactga cgctctgatg actggttaca ctgggtgactt 780
cgactctgtt atcgattgca acaacttgcaa ttccatgtct accaaccgga aaccgcagaa 840
aaaaaaciaa cgtaacacca accgtcgtcc gcaggacgtt aaattcccgg gtggtggtca 900
gatcgttaaa acctctgaac gttctcagcc gcgtgggctg cgtcagccga tccgaaagc 960
tcgtcgtccg gaaggtcgta cctgggctca gccgggttac ccgtggccgc tgtacggtaa 1020
cgaaggttgc gttgggctg gttggctgct gtctccgtaa taaggatcc 1069
    
```

```

<210> SEQ ID NO 18
<211> LENGTH: 350
<212> TYPE: PRT
<213> ORGANISM: Hepatitis C Virus
<220> FEATURE:
<223> OTHER INFORMATION: p9MB-18
    
```

<400> SEQUENCE: 18

```

Met Ala Val Asp Phe Ile Pro Val Glu Asn Leu Glu Thr Thr Met Arg
 1             5             10             15
Ser Pro Val Phe Thr Asp Asn Ser Ser Pro Pro Val Val Pro Gln Ser
          20             25             30
Phe Gln Val Ala His Leu His Ala Pro Thr Gly Ser Gly Lys Ser Thr
          35             40             45
Lys Val Pro Ala Ala Tyr Ala Ala Gln Gly Tyr Lys Val Leu Val Leu
 50             55             60
Asn Pro Ser Val Ala Ala Thr Leu Gly Phe Gly Ala Tyr Met Ser Lys
65             70             75             80
Ala His Gly Ile Asp Pro Asn Ile Arg Thr Gly Val Arg Thr Ile Thr
          85             90             95
Thr Gly Ser Pro Ile Thr Tyr Ser Thr Tyr Gly Lys Phe Leu Ala Asp
          100            105            110
Gly Gly Cys Ser Gly Gly Ala Tyr Asp Ile Ile Ile Cys Asp Glu Cys
          115            120            125
His Ser Thr Asp Ala Thr Ser Ile Leu Gly Ile Gly Thr Val Leu Asp
130            135            140
Gln Ala Glu Thr Ala Gly Ala Arg Leu Val Val Leu Ala Thr Ala Thr
145            150            155            160
Pro Pro Gly Ser Val Thr Val Pro His Pro Asn Ile Glu Glu Val Ala
          165            170            175
    
```

-continued

Leu Ser Thr Thr Gly Glu Ile Pro Phe Tyr Gly Lys Ala Ile Pro Leu
 180 185 190
 Glu Val Ile Lys Gly Gly Arg His Leu Ile Phe Cys His Ser Lys Lys
 195 200 205
 Lys Cys Asp Glu Leu Ala Ala Lys Leu Val Ala Leu Gly Ile Asn Ala
 210 215 220
 Val Ala Tyr Tyr Arg Gly Leu Asp Val Ser Val Ile Pro Thr Ser Gly
 225 230 235 240
 Asp Val Val Val Val Ala Thr Asp Ala Leu Met Thr Gly Tyr Thr Gly
 245 250 255
 Asp Phe Asp Ser Val Ile Asp Cys Asn Thr Cys Asn Ser Met Ser Thr
 260 265 270
 Asn Pro Lys Pro Gln Lys Lys Asn Lys Arg Asn Thr Asn Arg Arg Pro
 275 280 285
 Gln Asp Val Lys Phe Pro Gly Gly Gly Gln Ile Val Lys Thr Ser Glu
 290 295 300
 Arg Ser Gln Pro Arg Gly Arg Arg Gln Pro Ile Pro Lys Ala Arg Arg
 305 310 315 320
 Pro Glu Gly Arg Thr Trp Ala Gln Pro Gly Tyr Pro Trp Pro Leu Tyr
 325 330 335
 Gly Asn Glu Gly Cys Gly Trp Ala Gly Trp Leu Leu Ser Pro
 340 345 350

<210> SEQ ID NO 19
 <211> LENGTH: 1051
 <212> TYPE: DNA
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: p9MB-19

<400> SEQUENCE: 19

gaattccatg gctggtgact ttatcccggg tgaaaatctc gagactacta tgcgttctcc 60
 ggttttcact gacaactcct ctccgcccgg tgttccgcag tctttccagg ttgctcacct 120
 gcatgctccg actggttctg gtaaactctac taaagttcca gctgcttacg ctgctcaggg 180
 ttacaaagt ctggttctga acccgtctgt tgctgctact ctgggtttcg gcgctacat 240
 gtctaaagct cacggatcg acccgaacat tcgtactggg gtacgtacta tcaactactg 300
 ttctccgatc acttactcta cttacggtaa attcctggct gacggtggtt gctctgggg 360
 tgcttacgat atcatcatct ggcacgaatg ccaactctact gacgctactt ctatcctggg 420
 tatcgggtacc gttctggacc aggctgaaac tgcagggtgct cgtctgggtt ttctggctac 480
 tgctactccg ccgggttctg ttactgttcc gcacccgaac atcgaagaag ttgctctgtc 540
 gactactggt gaaatcccgt tctacggtaa agctatcccg ctccgaggta tcaaagggtg 600
 tcgtcacctg attttctgcc actctaaaaa aaaatgacac gaactggctg ctaagcttgt 660
 tgctctgggt atcaacgctg ttgcttacta ccgtggctct gacgtttctg ttatcccgc 720
 ttctgggtgac gttggtggtg tggccactga cgctctgatg actggttaca ctgggtgactt 780
 cgactctggt atcgattgca aacttgcaa ttccatgcag aaaaaaaaaa aacgtaaacac 840
 caaccgtcgt ccgcaggacg ttaaattccc ggggtgggtg cagatcggtt aaacctctga 900
 acgttctcag ccgcgtgggc gtcgtcagcc gatcccgaag gctcgtcgtc cggaaggctc 960
 tacctgggct cagccgggtt acccgtggcc gctgtacggg aacgaagggt gcggttgggc 1020

-continued

tggttgctg ctgtctccgt aataaggatc c

1051

<210> SEQ ID NO 20
 <211> LENGTH: 344
 <212> TYPE: PRT
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: p9MB-19

<400> SEQUENCE: 20

Met	Ala	Val	Asp	Phe	Ile	Pro	Val	Glu	Asn	Leu	Glu	Thr	Thr	Met	Arg
1				5					10					15	
Ser	Pro	Val	Phe	Thr	Asp	Asn	Ser	Ser	Pro	Pro	Val	Val	Pro	Gln	Ser
			20					25					30		
Phe	Gln	Val	Ala	His	Leu	His	Ala	Pro	Thr	Gly	Ser	Gly	Lys	Ser	Thr
		35					40					45			
Lys	Val	Pro	Ala	Ala	Tyr	Ala	Ala	Gln	Gly	Tyr	Lys	Val	Leu	Val	Leu
	50					55					60				
Asn	Pro	Ser	Val	Ala	Ala	Thr	Leu	Gly	Phe	Gly	Ala	Tyr	Met	Ser	Lys
65					70					75					80
Ala	His	Gly	Ile	Asp	Pro	Asn	Ile	Arg	Thr	Gly	Val	Arg	Thr	Ile	Thr
				85					90					95	
Thr	Gly	Ser	Pro	Ile	Thr	Tyr	Ser	Thr	Tyr	Gly	Lys	Phe	Leu	Ala	Asp
			100					105					110		
Gly	Gly	Cys	Ser	Gly	Gly	Ala	Tyr	Asp	Ile	Ile	Ile	Cys	Asp	Glu	Cys
		115					120					125			
His	Ser	Thr	Asp	Ala	Thr	Ser	Ile	Leu	Gly	Ile	Gly	Thr	Val	Leu	Asp
		130					135					140			
Gln	Ala	Glu	Thr	Ala	Gly	Ala	Arg	Leu	Val	Val	Leu	Ala	Thr	Ala	Thr
145					150					155					160
Pro	Pro	Gly	Ser	Val	Thr	Val	Pro	His	Pro	Asn	Ile	Glu	Glu	Val	Ala
				165					170					175	
Leu	Ser	Thr	Thr	Gly	Glu	Ile	Pro	Phe	Tyr	Gly	Lys	Ala	Ile	Pro	Leu
			180					185					190		
Glu	Val	Ile	Lys	Gly	Gly	Arg	His	Leu	Ile	Phe	Cys	His	Ser	Lys	Lys
			195				200					205			
Lys	Cys	Asp	Glu	Leu	Ala	Ala	Lys	Leu	Val	Ala	Leu	Gly	Ile	Asn	Ala
	210					215					220				
Val	Ala	Tyr	Tyr	Arg	Gly	Leu	Asp	Val	Ser	Val	Ile	Pro	Thr	Ser	Gly
225					230					235					240
Asp	Val	Val	Val	Val	Ala	Thr	Asp	Ala	Leu	Met	Thr	Gly	Tyr	Thr	Gly
				245					250					255	
Asp	Phe	Asp	Ser	Val	Ile	Asp	Cys	Asn	Thr	Cys	Asn	Ser	Met	Gln	Lys
			260					265					270		
Lys	Asn	Lys	Arg	Asn	Thr	Asn	Arg	Arg	Pro	Gln	Asp	Val	Lys	Phe	Pro
		275					280					285			
Gly	Gly	Gly	Gln	Ile	Val	Lys	Thr	Ser	Glu	Arg	Ser	Gln	Pro	Arg	Gly
	290					295					300				
Arg	Arg	Gln	Pro	Ile	Pro	Lys	Ala	Arg	Arg	Pro	Glu	Gly	Arg	Thr	Trp
305					310					315					320
Ala	Gln	Pro	Gly	Tyr	Pro	Trp	Pro	Leu	Tyr	Gly	Asn	Glu	Gly	Cys	Gly
				325					330					335	

-continued

Trp Ala Gly Trp Leu Leu Ser Pro
340

<210> SEQ ID NO 21
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer

<400> SEQUENCE: 21

tatagaattc catggctggt gactttatcc 30

<210> SEQ ID NO 22
<211> LENGTH: 27
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer

<400> SEQUENCE: 22

ggaattgcaa gtgttgcaat cgataac 27

<210> SEQ ID NO 23
<211> LENGTH: 51
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer

<400> SEQUENCE: 23

gttatcgatt gcaacacttg caattccatg cagaaaaaaaa acaaacgtaa c 51

<210> SEQ ID NO 24
<211> LENGTH: 1066
<212> TYPE: DNA
<213> ORGANISM: Hepatitis C Virus
<220> FEATURE:
<223> OTHER INFORMATION: p9MB-20

<400> SEQUENCE: 24

gaattccatg gctgttgact ttatcccggg tgaaaatctc gagactacta tgcgtttctc 60
ggttttcact gacaactctt ctccgccggg tgttccgcag tctttccagg ttgctcaact 120
gcatgctccg actggttctg gtaaattctac taaagtcca gctgcttacg ctgctcaggg 180
ttacaaagt ctgggtctga acccgtctgt tgctgctact ctgggtttcg ggcctacat 240
gtctaaagct cacggatcg acccgaacat tcgtactggt gtacgtacta tcaactactg 300
ttctccgatc acttactcta cttacggtaa attcctggct gacggtggtt gctctgggtg 360
tgcttacgat atcatcatct gcgacgaatg ccaactctact gacgctactt ctatcctggg 420
tatcgggtacc gttctggacc aggctgaaac tgcaggtgct cgtctgggtg ttctggctac 480
tgctactccg ccgggttctg ttactgttcc gcacccgaac atcgaagaag ttgctctgtc 540
gactactggt gaaatcccgt tctacggtaa agctatcccg ctcgagggta tcaaagggtg 600
tcgctacctg attttctgcc actctaaaaa aaaaatgagc gaactggctg ctaagcttgt 660
tgctctgggt atcaacgctg ttgcttacta ccgtgggtctg gacgtttctg ttatcccgac 720
ttctgggtgac gttgttggtg tggccactga cgctctgatg actggttaca ctgggtgactt 780
cgactctggt atcgattgca acaattgcaa ttccgggtggg ggtggttcta tgcagaaaaa 840

-continued

```

aaacaaacgt aacaccaacc gtcgtccgca ggacgttaaa ttcccgggtg gtggtcagat    900
cgttaaaacc tctgaacgtt ctcagcccgcg tgggcgctcgt cagccgatcc cgaagctcg    960
tcgtccggaa ggtcgtacct gggctcagcc gggttaccgg tggccgctgt acggtaacga    1020
aggttgcggt tgggctgggt ggctgctgtc tccgtaataa ggatcc                    1066
    
```

```

<210> SEQ ID NO 25
<211> LENGTH: 349
<212> TYPE: PRT
<213> ORGANISM: Hepatitis C Virus
<220> FEATURE:
<223> OTHER INFORMATION: p9MB-20
    
```

<400> SEQUENCE: 25

```

Met Ala Val Asp Phe Ile Pro Val Glu Asn Leu Glu Thr Thr Met Arg
 1           5           10           15
Ser Pro Val Phe Thr Asp Asn Ser Ser Pro Pro Val Val Pro Gln Ser
          20           25           30
Phe Gln Val Ala His Leu His Ala Pro Thr Gly Ser Gly Lys Ser Thr
          35           40           45
Lys Val Pro Ala Ala Tyr Ala Gln Gly Tyr Lys Val Leu Val Leu
 50           55           60
Asn Pro Ser Val Ala Ala Thr Leu Gly Phe Gly Ala Tyr Met Ser Lys
65           70           75           80
Ala His Gly Ile Asp Pro Asn Ile Arg Thr Gly Val Arg Thr Ile Thr
          85           90           95
Thr Gly Ser Pro Ile Thr Tyr Ser Thr Tyr Gly Lys Phe Leu Ala Asp
          100          105          110
Gly Gly Cys Ser Gly Gly Ala Tyr Asp Ile Ile Ile Cys Asp Glu Cys
          115          120          125
His Ser Thr Asp Ala Thr Ser Ile Leu Gly Ile Gly Thr Val Leu Asp
130          135          140
Gln Ala Glu Thr Ala Gly Ala Arg Leu Val Val Leu Ala Thr Ala Thr
145          150          155          160
Pro Pro Gly Ser Val Thr Val Pro His Pro Asn Ile Glu Glu Val Ala
          165          170          175
Leu Ser Thr Thr Gly Glu Ile Pro Phe Tyr Gly Lys Ala Ile Pro Leu
180          185          190
Glu Val Ile Lys Gly Gly Arg His Leu Ile Phe Cys His Ser Lys Lys
195          200          205
Lys Cys Asp Glu Leu Ala Ala Lys Leu Val Ala Leu Gly Ile Asn Ala
210          215          220
Val Ala Tyr Tyr Arg Gly Leu Asp Val Ser Val Ile Pro Thr Ser Gly
225          230          235          240
Asp Val Val Val Val Ala Thr Asp Ala Leu Met Thr Gly Tyr Thr Gly
          245          250          255
Asp Phe Asp Ser Val Ile Asp Cys Asn Thr Cys Asn Ser Gly Gly Gly
260          265          270
Gly Ser Met Gln Lys Lys Asn Lys Arg Asn Thr Asn Arg Arg Pro Gln
275          280          285
Asp Val Lys Phe Pro Gly Gly Gly Gln Ile Val Lys Thr Ser Glu Arg
290          295          300
    
```

-continued

Ser	Gln	Pro	Arg	Gly	Arg	Arg	Gln	Pro	Ile	Pro	Lys	Ala	Arg	Arg	Pro
305					310				315						320
Glu	Gly	Arg	Thr	Trp	Ala	Gln	Pro	Gly	Tyr	Pro	Trp	Pro	Leu	Tyr	Gly
				325				330						335	
Asn	Glu	Gly	Cys	Gly	Trp	Ala	Gly	Trp	Leu	Leu	Ser	Pro			
			340					345							

<210> SEQ ID NO 26
 <211> LENGTH: 66
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

<400> SEQUENCE: 26

gttatcgatt gcaacacttg caattccggt ggtggtggtt ctatgcagaa aaaaaacaaa 60
 cgtaac 66

<210> SEQ ID NO 27
 <211> LENGTH: 1293
 <212> TYPE: DNA
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: p9MB-22

<400> SEQUENCE: 27

gaattccatg gctgttgact ttatcccggg tgaaaatctc gagactacta tgcgttctcc 60
 ggttttcact gacaactcct ctccgccggg tgttcccgag tctttccagg ttgctcacct 120
 gcatgctccg actggttctg gtaaactctac taaagttcca gctgcttacg ctgctcaggg 180
 ttacaaagtt ctggttctga acccgctctgt tgctgctact ctgggtttcg gcgctacat 240
 gtctaaagct cacggtatcg acccgaacat tcgtactggg gtacgtacta tcaactactgg 300
 ttctccgac acttactcta cttacggtaa attcctggct gacggtgggt gctctgggtg 360
 tgcttacgat atcatcatct gcgacgaatg ccaactctact gacgctactt ctatcctggg 420
 tatcgggtacc gttctggacc aggctgaaac tgcagggtgct cgtctgggtt ttctggctac 480
 tgctactccg ccgggttctg ttactgttcc gcacccgaac atcgaagaag ttgctctgtc 540
 gactactggg gaaatcccg tctacggtaa agctatcccg ctcgagggta tcaaagggtg 600
 tcgctacctg attttctgcc actctaaaaa aaaatgacgac gaactggctg ctaagcttgt 660
 tgctctgggt atcaacgctg ttgcttacta ccgtgggtctg gacggttctg ttatcccggc 720
 ttctgggtgac gttgtgtgtg tggccactga cgctctgatg actggttaca ctgggtgactt 780
 cgactctggt atcgattgca acaactgcaa ttccgggtgg ggtggttcta tgtctaccaa 840
 cccgaaaccg cagaaaaaaa acaaacgtaa caccaaccgt cgtccgcagg acgttaaatt 900
 cccgggtggg ggtcagatcg ttggtgggtt ttacctgctg ccgctcgtg gtcgcgctct 960
 ggggttctct gctacgcgta aaacctctga acgttctcag ccgctggggc gtcgctcagcc 1020
 gatcccgaaa gctcgtcgtc cggaaggctg tacctgggct cagccgggtt acccgtggcc 1080
 gctgtacggt aacgaagggt gcggttgggc tgggtggctg ctgtctccgc gttgatctcg 1140
 tccgtcttgg gttccgacgg acccgcgctg tcgttctcgt aaccttggta aagttatoga 1200
 taccctgacc tgcggtttct ctgacctgat gggttacata ccgctgggtt gagctccgct 1260
 ggggtgggtg gctcgtgctt aaccatgga tcc 1293

-continued

```

<210> SEQ ID NO 28
<211> LENGTH: 424
<212> TYPE: PRT
<213> ORGANISM: Hepatitis C Virus
<220> FEATURE:
<223> OTHER INFORMATION: p9MB-22

<400> SEQUENCE: 28

Met Ala Val Asp Phe Ile Pro Val Glu Asn Leu Glu Thr Thr Met Arg
 1           5           10           15
Ser Pro Val Phe Thr Asp Asn Ser Ser Pro Pro Val Val Pro Gln Ser
 20           25           30
Phe Gln Val Ala His Leu His Ala Pro Thr Gly Ser Gly Lys Ser Thr
 35           40           45
Lys Val Pro Ala Ala Tyr Ala Ala Gln Gly Tyr Lys Val Leu Val Leu
 50           55           60
Asn Pro Ser Val Ala Ala Thr Leu Gly Phe Gly Ala Tyr Met Ser Lys
 65           70           75           80
Ala His Gly Ile Asp Pro Asn Ile Arg Thr Gly Val Arg Thr Ile Thr
 85           90           95
Thr Gly Ser Pro Ile Thr Tyr Ser Thr Tyr Gly Lys Phe Leu Ala Asp
 100          105          110
Gly Gly Cys Ser Gly Gly Ala Tyr Asp Ile Ile Ile Cys Asp Glu Cys
 115          120          125
His Ser Thr Asp Ala Thr Ser Ile Leu Gly Ile Gly Thr Val Leu Asp
 130          135          140
Gln Ala Glu Thr Ala Gly Ala Arg Leu Val Val Leu Ala Thr Ala Thr
 145          150          155          160
Pro Pro Gly Ser Val Thr Val Pro His Pro Asn Ile Glu Glu Val Ala
 165          170          175
Leu Ser Thr Thr Gly Glu Ile Pro Phe Tyr Gly Lys Ala Ile Pro Leu
 180          185          190
Glu Val Ile Lys Gly Gly Arg His Leu Ile Phe Cys His Ser Lys Lys
 195          200          205
Lys Cys Asp Glu Leu Ala Ala Lys Leu Val Ala Leu Gly Ile Asn Ala
 210          215          220
Val Ala Tyr Tyr Arg Gly Leu Asp Val Ser Val Ile Pro Thr Ser Gly
 225          230          235          240
Asp Val Val Val Val Ala Thr Asp Ala Leu Met Thr Gly Tyr Thr Gly
 245          250          255
Asp Phe Asp Ser Val Ile Asp Cys Asn Thr Cys Asn Ser Gly Gly Gly
 260          265          270
Gly Ser Met Ser Thr Asn Pro Lys Pro Gln Lys Lys Asn Lys Arg Asn
 275          280          285
Thr Asn Arg Arg Pro Gln Asp Val Lys Phe Pro Gly Gly Gly Gln Ile
 290          295          300
Val Gly Gly Val Tyr Leu Leu Pro Arg Arg Gly Pro Arg Leu Gly Val
 305          310          315          320
Arg Ala Thr Arg Lys Thr Ser Glu Arg Ser Gln Pro Arg Gly Arg Arg
 325          330          335
Gln Pro Ile Pro Lys Ala Arg Arg Pro Glu Gly Arg Thr Trp Ala Gln
 340          345          350

```

-continued

Pro Gly Tyr Pro Trp Pro Leu Tyr Gly Asn Glu Gly Cys Gly Trp Ala
 355 360 365
 Gly Trp Leu Leu Ser Pro Arg Gly Ser Arg Pro Ser Trp Gly Pro Thr
 370 375 380
 Asp Pro Arg Arg Arg Ser Arg Asn Leu Gly Lys Val Ile Asp Thr Leu
 385 390 395 400
 Thr Cys Gly Phe Ala Asp Leu Met Gly Tyr Ile Pro Leu Val Gly Ala
 405 410 415
 Pro Leu Gly Gly Ala Ala Arg Ala
 420

<210> SEQ ID NO 29
 <211> LENGTH: 66
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

<400> SEQUENCE: 29

gttatcgatt gcaacacttg caattccggt ggtggtggtt ctatgtctac caaccgaaa 60
 ccgcag 66

<210> SEQ ID NO 30
 <211> LENGTH: 27
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

<400> SEQUENCE: 30

tataggatcc atgggttaag cacgagc 27

<210> SEQ ID NO 31
 <211> LENGTH: 1111
 <212> TYPE: DNA
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: p9MB-31

<400> SEQUENCE: 31

gaattccatg gctgttgact ttatcccggt tgaaaatctc gagactacta tgcgttctcc 60
 ggttttcact gacaactctt ctccgccggt tgttccgcag tctttccagg ttgctcaact 120
 gcatgctccg actggttctg gtaaactctac taaagtcca gctgcttacg ctgctcaggg 180
 ttacaaagt ctggttctga acccgtctgt tgctgctact ctgggtttcg ggcctacat 240
 gtctaaagct cacggtatcg acccgaacat tcgtactggt gtacgtacta tcaactactg 300
 ttctccgac acttactcta cttacggtaa attcctggct gacggtggtt gctctggtg 360
 tgcttacgat atcatcatct gcgacgaatg ccaactctact gacgctactt ctatcctggg 420
 tatcgggtacc gttctggacc aggctgaaac tgcaggtgct cgtctgggtt ttctggctac 480
 tgctactccg ccgggttctg ttactgttcc gcacccgaac atcgaagaag ttgctctgtc 540
 gactactggt gaaatcccg tctacggtaa agctatcccg ctcgagggta tcaaagggtg 600
 tcgtcacctg attttctgcc actctaaaaa aaaatgagc gaactggctg ctaagcttgt 660
 tgctctgggt atcaacgctg ttgcttacta ccgtggtctg gacgtttctg ttatcccgac 720

-continued

```

ttctggtgac gttggtggtg tggccactga cgctctgatg actggttaca ctggtgactt    780
cgactctggt atcgattgca acacttgcaa ttccatgtct accaaccoga aaccgcagaa    840
aaaaaacaaa cgtaacacca accgctcgcc gcaggacggt aaattcccgg gtggtggtca    900
gatcgtttac ctgctgccgc gtcgtggtcc gcgtctgggt gttacgcgta aaacctctga    960
acgttctcag ccgctggggc gtcgtcagcc gatcccgaaa gctcgtcgtc cgaaggtcg   1020
tacctgggct cagccggggt acccgtggcc gctgtacggt aacgaagggt gcggttgggc   1080
tggttggcta ctgtctccgt aataaggatc c                                     1111
    
```

```

<210> SEQ ID NO 32
<211> LENGTH: 364
<212> TYPE: PRT
<213> ORGANISM: Hepatitis C Virus
<220> FEATURE:
<223> OTHER INFORMATION: p9MB-31
    
```

<400> SEQUENCE: 32

```

Met Ala Val Asp Phe Ile Pro Val Glu Asn Leu Glu Thr Thr Met Arg
 1           5           10           15
Ser Pro Val Phe Thr Asp Asn Ser Ser Pro Pro Val Val Pro Gln Ser
          20           25           30
Phe Gln Val Ala His Leu His Ala Pro Thr Gly Ser Gly Lys Ser Thr
          35           40           45
Lys Val Pro Ala Ala Tyr Ala Ala Gln Gly Tyr Lys Val Leu Val Leu
 50           55           60
Asn Pro Ser Val Ala Ala Thr Leu Gly Phe Gly Ala Tyr Met Ser Lys
65           70           75           80
Ala His Gly Ile Asp Pro Asn Ile Arg Thr Gly Val Arg Thr Ile Thr
          85           90           95
Thr Gly Ser Pro Ile Thr Tyr Ser Thr Tyr Gly Lys Phe Leu Ala Asp
          100          105          110
Gly Gly Cys Ser Gly Gly Ala Tyr Asp Ile Ile Ile Cys Asp Glu Cys
          115          120          125
His Ser Thr Asp Ala Thr Ser Ile Leu Gly Ile Gly Thr Val Leu Asp
          130          135          140
Gln Ala Glu Thr Ala Gly Ala Arg Leu Val Val Leu Ala Thr Ala Thr
          145          150          155          160
Pro Pro Gly Ser Val Thr Val Pro His Pro Asn Ile Glu Glu Val Ala
          165          170          175
Leu Ser Thr Thr Gly Glu Ile Pro Phe Tyr Gly Lys Ala Ile Pro Leu
          180          185          190
Glu Val Ile Lys Gly Gly Arg His Leu Ile Phe Cys His Ser Lys Lys
          195          200          205
Lys Cys Asp Glu Leu Ala Ala Lys Leu Val Ala Leu Gly Ile Asn Ala
          210          215          220
Val Ala Tyr Tyr Arg Gly Leu Asp Val Ser Val Ile Pro Thr Ser Gly
          225          230          235          240
Asp Val Val Val Val Ala Thr Asp Ala Leu Met Thr Gly Tyr Thr Gly
          245          250          255
Asp Phe Asp Ser Val Ile Asp Cys Asn Thr Cys Asn Ser Met Ser Thr
          260          265          270
Asn Pro Lys Pro Gln Lys Lys Asn Lys Arg Asn Thr Asn Arg Arg Pro
    
```

-continued

275	280	285	
Gln Asp Val Lys Phe Pro Gly Gly Gly Gln Ile Val Tyr Leu Leu Pro			
290	295	300	
Arg Arg Gly Pro Arg Leu Gly Val Thr Arg Lys Thr Ser Glu Arg Ser			
305	310	315	320
Gln Pro Arg Gly Arg Arg Gln Pro Ile Pro Lys Ala Arg Arg Pro Glu			
	325	330	335
Gly Arg Thr Trp Ala Gln Pro Gly Tyr Pro Trp Pro Leu Tyr Gly Asn			
	340	345	350
Glu Gly Cys Gly Trp Ala Gly Trp Leu Leu Ser Pro			
	355	360	

<210> SEQ ID NO 33
 <211> LENGTH: 52
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

<400> SEQUENCE: 33

accgagacgc ggaccacgac gcggcagcag gtaaacgatc tgaccaccac cc 52

<210> SEQ ID NO 34
 <211> LENGTH: 54
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

<400> SEQUENCE: 34

ccgcgtcgtg gtccgcgtct ggggtgttacg cgtaaacct ctgaacgttc tcag 54

<210> SEQ ID NO 35
 <211> LENGTH: 1105
 <212> TYPE: DNA
 <213> ORGANISM: Hepatitis A Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: p9MB-24

<400> SEQUENCE: 35

gaattccatg gctgttgact ttatcccggg tgaaaatctc gagactacta tgcgttctcc 60
 ggttttcact gacaactctt ctccgccggg tgttccgcag tctttccagg ttgctcacct 120
 gcatgctccg actggttctg gtaaactctac taaagttcca gctgcttacg ctgctcaggg 180
 ttacaaagt ctggttctga acccgtctgt tgctgtact ctgggtttcg ggcctacat 240
 gtctaaagct cacggtatcg accogaacat tcgtactggt gtacgtacta tcaactactg 300
 ttctccgac acttactcta cttacggtaa attcctggct gacggtggtt gctctgggtg 360
 tgcttacgat atcatcatct gcgacgaatg ccaactctact gacgctactt ctatcctggg 420
 tatcgggtacc gttctggacc aggctgaaac tgcaggtgct cgtctgggtg ttctggctac 480
 tgctactccg ccgggttctg ttactgttcc gcaccogaac atcgaagaag ttgctctgtc 540
 gactactggt gaaatcccgt tctacggtaa agctatcccg ctcgagggta tcaaagggtg 600
 tcgtcacctg attttctgcc actctaaaaa aaaatgcgac gaactggctg ctaagcttgt 660
 tgctctgggt atcaacgctg ttgcttacta ccgtgggtctg gacgtttctg ttatcccgac 720
 ttctgggtgac gttgtgtgtg tggccactga cgctctgatg actggttaca ctggtgactt 780

-continued

```

cgactctgtt atcgattgca acacttgcaa ttccatgtct accaaccgga aaccgcagaa      840
aaaaaaciaa cgtaacacca accgctgtcc gcaggacgtt aaattcccgg gtgggtgtca      900
gatcgtttgt ctgctgccgc gtcgtgttcc gcgtctgggt cgtaaacct ctgaacgttc      960
tcagccgcgt gggcgtcgtc agccgatccc gaaagctcgt cgtccggaag gtcgtacctg     1020
ggctcagccg gttaccctt ggcgctgta cggtaacgaa gtttgcggtt gggctggtt      1080
gctgctgtct ccgtaataag gatcc                                           1105
    
```

```

<210> SEQ ID NO 36
<211> LENGTH: 358
<212> TYPE: PRT
<213> ORGANISM: Hepatitis C Virus
<220> FEATURE:
<223> OTHER INFORMATION: p9MB-24
    
```

<400> SEQUENCE: 36

```

Met Ala Val Asp Phe Ile Pro Val Glu Asn Leu Glu Thr Thr Met Arg
 1           5           10           15
Ser Pro Val Phe Thr Asp Asn Ser Ser Pro Pro Val Val Pro Gln Ser
           20           25           30
Phe Gln Val Ala His Leu His Ala Pro Thr Gly Ser Gly Lys Ser Thr
           35           40           45
Lys Val Pro Ala Ala Tyr Ala Ala Gln Gly Tyr Lys Val Leu Val Leu
 50           55           60
Asn Pro Ser Val Ala Ala Thr Leu Gly Phe Gly Ala Tyr Met Ser Lys
 65           70           75           80
Ala His Gly Ile Asp Pro Asn Ile Arg Thr Gly Val Arg Thr Ile Thr
           85           90           95
Thr Gly Ser Pro Ile Thr Tyr Ser Thr Tyr Gly Lys Phe Leu Ala Asp
          100          105          110
Gly Gly Cys Ser Gly Gly Ala Tyr Asp Ile Ile Ile Cys Asp Glu Cys
          115          120          125
His Ser Thr Asp Ala Thr Ser Ile Leu Gly Ile Gly Thr Val Leu Asp
          130          135          140
Gln Ala Glu Thr Ala Gly Ala Arg Leu Val Val Leu Ala Thr Ala Thr
          145          150          155          160
Pro Pro Gly Ser Val Thr Val Pro His Pro Asn Ile Glu Glu Val Ala
          165          170          175
Leu Ser Thr Thr Gly Glu Ile Pro Phe Tyr Gly Lys Ala Ile Pro Leu
          180          185          190
Glu Val Ile Lys Gly Gly Arg His Leu Ile Phe Cys His Ser Lys Lys
          195          200          205
Lys Cys Asp Glu Leu Ala Ala Lys Leu Val Ala Leu Gly Ile Asn Ala
          210          215          220
Val Ala Tyr Tyr Arg Gly Leu Asp Val Ser Val Ile Pro Thr Ser Gly
          225          230          235          240
Asp Val Val Val Val Ala Thr Asp Ala Leu Met Thr Gly Tyr Thr Gly
          245          250          255
Asp Phe Asp Ser Val Ile Asp Cys Asn Thr Cys Asn Ser Met Ser Thr
          260          265          270
Asn Pro Lys Pro Gln Lys Lys Asn Lys Arg Asn Thr Asn Arg Arg Pro
          275          280          285
    
```

-continued

Gln Asp Val Lys Phe Pro Gly Gly Gly Gln Ile Val Leu Leu Pro Arg
 290 295 300
 Arg Gly Pro Arg Leu Gly Arg Lys Thr Ser Glu Arg Ser Gln Pro Arg
 305 310 315 320
 Gly Arg Arg Gln Pro Ile Pro Lys Ala Arg Arg Pro Glu Gly Arg Thr
 325 330 335
 Trp Ala Gln Pro Gly Tyr Pro Trp Pro Leu Tyr Gly Asn Gly Trp Ala
 340 345 350
 Gly Trp Leu Leu Ser Pro
 355

<210> SEQ ID NO 37
 <211> LENGTH: 1120
 <212> TYPE: DNA
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: p9MB-25

<400> SEQUENCE: 37

gaattccatg gctgttgact ttatcccggg tgaaaaatctc gagactacta tgcgttctcc 60
 ggttttcact gacaactctt ctccgccggg tgttccgcag tctttccagg ttgctcaact 120
 gcatgctccg actggttctg gtaaactctac taaagttcca gctgcttacg ctgctcaggg 180
 ttacaaagtt ctggttctga acccgctctgt tgctgctact ctgggtttcg gcgctacat 240
 gtctaaagct cacggtatcg acccgaacat tcgtactggg gtacgtacta tcaactactgg 300
 ttctccgatc acttactcta cttacggtaa attcctggct gacggtggtt gctctgggtg 360
 tgcttacgat atcatcatct gcgacgaatg ccaactctact gacgctactt ctatcctggg 420
 tatcgggtacc gttctggacc aggctgaaac tgcagggtgct cgtctgggtg ttctggctac 480
 tgctactccg ccgggttctg ttaactgttc gcacccgaac atcgaagaag ttgctctgtc 540
 gactactggg gaaatcccgt tctacggtaa agctatcccg ctcgagggta tcaaagggtg 600
 tcgtcacctg attttctgcc actctaaaaa aaaatgagac gaactggctg ctaagcttgt 660
 tgctctgggt atcaacgctg ttgcttacta ccgtgggtctg gacgtttctg ttatcccgcg 720
 ttctgggtgac gttgttggtg tggccactga cgctctgatg actggttaca ctgggtgactt 780
 cgactctggt atcgattgca acaactgcaa ttccgggtgg ggtggttcta tgtctaccaa 840
 cccgaaaccg cagaaaaaaa acaaacgtaa caccaaccgt cgtccgcagg acgttaaatt 900
 cccgggtggg ggtcagatcg ttggtctgct gccgcgctgt ggtccgcgctc tgggtcgtaa 960
 aacctctgaa cgttctcagc cgcgtggggc tcgctcagccg atcccgaag ctgctcgtcc 1020
 ggaaggtcgt acctgggctc agccgggtta cccgtggccg ctgtacggta acgaaggttg 1080
 cgggtgggct ggttggtctg tgtctccgta ataaggatcc 1120

<210> SEQ ID NO 38
 <211> LENGTH: 367
 <212> TYPE: PRT
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: p9MB-25

<400> SEQUENCE: 38

Met Ala Val Asp Phe Ile Pro Val Glu Asn Leu Glu Thr Thr Met Arg
 1 5 10 15

-continued

Ser Pro Val Phe Thr Asp Asn Ser Ser Pro Pro Val Val Pro Gln Ser
 20 25 30

Phe Gln Val Ala His Leu His Ala Pro Thr Gly Ser Gly Lys Ser Thr
 35 40 45

Lys Val Pro Ala Ala Tyr Ala Ala Gln Gly Tyr Lys Val Leu Val Leu
 50 55 60

Asn Pro Ser Val Ala Ala Thr Leu Gly Phe Gly Ala Tyr Met Ser Lys
 65 70 75 80

Ala His Gly Ile Asp Pro Asn Ile Arg Thr Gly Val Arg Thr Ile Thr
 85 90 95

Thr Gly Ser Pro Ile Thr Tyr Ser Thr Tyr Gly Lys Phe Leu Ala Asp
 100 105 110

Gly Gly Cys Ser Gly Gly Ala Tyr Asp Ile Ile Ile Cys Asp Glu Cys
 115 120 125

His Ser Thr Asp Ala Thr Ser Ile Leu Gly Ile Gly Thr Val Leu Asp
 130 135 140

Gln Ala Glu Thr Ala Gly Ala Arg Leu Val Val Leu Ala Thr Ala Thr
 145 150 155 160

Pro Pro Gly Ser Val Thr Val Pro His Pro Asn Ile Glu Glu Val Ala
 165 170 175

Leu Ser Thr Thr Gly Glu Ile Pro Phe Tyr Gly Lys Ala Ile Pro Leu
 180 185 190

Glu Val Ile Lys Gly Gly Arg His Leu Ile Phe Cys His Ser Lys Lys
 195 200 205

Lys Cys Asp Glu Leu Ala Ala Lys Leu Val Ala Leu Gly Ile Asn Ala
 210 215 220

Val Ala Tyr Tyr Arg Gly Leu Asp Val Ser Val Ile Pro Thr Ser Gly
 225 230 235 240

Asp Val Val Val Val Ala Thr Asp Ala Leu Met Thr Gly Tyr Thr Gly
 245 250 255

Asp Phe Asp Ser Val Ile Asp Cys Asn Thr Cys Asn Ser Gly Gly Gly
 260 265 270

Gly Ser Met Ser Thr Asn Pro Lys Pro Gln Lys Lys Asn Lys Arg Asn
 275 280 285

Thr Asn Arg Arg Pro Gln Asp Val Lys Phe Pro Gly Gly Gly Gln Ile
 290 295 300

Val Gly Leu Leu Pro Arg Arg Gly Pro Arg Leu Gly Arg Lys Thr Ser
 305 310 315 320

Glu Arg Ser Gln Pro Arg Gly Arg Arg Gln Pro Ile Pro Lys Ala Arg
 325 330 335

Arg Pro Glu Gly Arg Thr Trp Ala Gln Pro Gly Tyr Pro Trp Pro Leu
 340 345 350

Tyr Gly Asn Glu Gly Cys Gly Trp Ala Gly Trp Leu Leu Ser Pro
 355 360 365

<210> SEQ ID NO 39
 <211> LENGTH: 1084
 <212> TYPE: DNA
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: p9MB-26
 <400> SEQUENCE: 39

-continued

```

gaattccatg gctgttgact ttatcccggg tgaaaatctc gagactacta tgcgttctcc    60
ggttttcact gacaactctt ctccgccggg tgttccgcag tctttccagg ttgctcacct    120
gcatgctccg actggttctg gtaaattctac taaagttcca gctgcttacg ctgctcaggg    180
ttacaaagt  ctgggttctga acccgtctgt tgctgtact  ctgggttctg ggcctacat    240
gtctaaagct cacggatcg  accggaacat tcgtactggg gtacgtacta tcaactactgg    300
ttctccgatc acttactcta cttacggtaa attcctggct gacggtggtt gctctggtg    360
tgcttacgat atcatcatct gcgacgaatg ccaacttact gacgctactt ctatcctggg    420
tatcgggtacc gttctggacc aggctgaaac tgcagggtgct cgtctggttg ttctggctac    480
tgctactccg ccgggttctg ttactgttcc gcacccgaac atcgaagaag ttgctctgtc    540
gactactggt gaaatcccgt tctacggtaa agctatcccg ctcgaggtta tcaaaggtgg    600
tcgtcacctg attttctgcc actctaaaaa aaaatgcgac gaactggctg ctaagcttgt    660
tgctctgggt atcaacgctg ttgcttacta ccgtggctg  gacgttctg  ttatcccgc    720
ttctggtgac gttgttgttg tggccactga cgctctgatg actggttaca ctggtgactt    780
cgactctgtt atcgattgca aacttgcaa  ttccgggtgg  ggtggttcta tgtctaccaa    840
cccgaaccg  cagaaaaaaaa acaaacgtaa caccaaccgt cgtccgcagg acgttaaatt    900
cccgggtggt ggtcagatcg ttaaaacctc tgaacgttct cagccgcgtg ggcgtcgtca    960
gccgatcccc aaagctcgtc gtccggaagg tcgtacctgg gtcagccgg gttaccctg    1020
gccgctgtac ggtaacgaag gttgcggttg ggctggttg  ctgctgtctc cgtaataagg    1080
atcc                                          1084
    
```

```

<210> SEQ ID NO 40
<211> LENGTH: 355
<212> TYPE: PRT
<213> ORGANISM: Hepatitis C Virus
<220> FEATURE:
<223> OTHER INFORMATION: p9MB-26
    
```

<400> SEQUENCE: 40

```

Met Ala Val Asp Phe Ile Pro Val Glu Asn Leu Glu Thr Thr Met Arg
 1           5           10           15
Ser Pro Val Phe Thr Asp Asn Ser Ser Pro Pro Val Val Pro Gln Ser
 20          25          30
Phe Gln Val Ala His Leu His Ala Pro Thr Gly Ser Gly Lys Ser Thr
 35          40          45
Lys Val Pro Ala Ala Tyr Ala Ala Gln Gly Tyr Lys Val Leu Val Leu
 50          55          60
Asn Pro Ser Val Ala Ala Thr Leu Gly Phe Gly Ala Tyr Met Ser Lys
 65          70          75          80
Ala His Gly Ile Asp Pro Asn Ile Arg Thr Gly Val Arg Thr Ile Thr
 85          90          95
Thr Gly Ser Pro Ile Thr Tyr Ser Thr Tyr Gly Lys Phe Leu Ala Asp
100         105         110
Gly Gly Cys Ser Gly Gly Ala Tyr Asp Ile Ile Ile Cys Asp Glu Cys
115         120         125
His Ser Thr Asp Ala Thr Ser Ile Leu Gly Ile Gly Thr Val Leu Asp
130         135         140
    
```

-continued

Gln Ala Glu Thr Ala Gly Ala Arg Leu Val Val Leu Ala Thr Ala Thr
 145 150 155 160

Pro Pro Gly Ser Val Thr Val Pro His Pro Asn Ile Glu Glu Val Ala
 165 170 175

Leu Ser Thr Thr Gly Glu Ile Pro Phe Tyr Gly Lys Ala Ile Pro Leu
 180 185 190

Glu Val Ile Lys Gly Gly Arg His Leu Ile Phe Cys His Ser Lys Lys
 195 200 205

Lys Cys Asp Glu Leu Ala Ala Lys Leu Val Ala Leu Gly Ile Asn Ala
 210 215 220

Val Ala Tyr Tyr Arg Gly Leu Asp Val Ser Val Ile Pro Thr Ser Gly
 225 230 235 240

Asp Val Val Val Val Ala Thr Asp Ala Leu Met Thr Gly Tyr Thr Gly
 245 250 255

Asp Phe Asp Ser Val Ile Asp Cys Asn Thr Cys Asn Ser Gly Gly Gly
 260 265 270

Gly Ser Met Ser Thr Asn Pro Lys Pro Gln Lys Lys Asn Lys Arg Asn
 275 280 285

Thr Asn Arg Arg Pro Gln Asp Val Lys Phe Pro Gly Gly Gly Gln Ile
 290 295 300

Val Lys Thr Ser Glu Arg Ser Gln Pro Arg Gly Arg Arg Gln Pro Ile
 305 310 315 320

Pro Lys Ala Arg Arg Pro Glu Gly Arg Thr Trp Ala Gln Pro Gly Tyr
 325 330 335

Pro Trp Pro Leu Tyr Gly Asn Glu Gly Cys Gly Trp Ala Gly Trp Leu
 340 345 350

Leu Ser Pro
 355

<210> SEQ ID NO 41
 <211> LENGTH: 39
 <212> TYPE: PRT
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: ALAM-16

<400> SEQUENCE: 41

Lys Thr Lys Arg Asn Thr Asn Arg Arg Pro Gln Asp Val Lys Phe Pro
 1 5 10 15

Gly Gly Gly Gln Ile Val Tyr Leu Leu Pro Arg Arg Gly Pro Arg Leu
 20 25 30

Gly Val Thr Arg Lys Thr Ser
 35

<210> SEQ ID NO 42
 <211> LENGTH: 44
 <212> TYPE: PRT
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: ALAM-17

<400> SEQUENCE: 42

Lys Thr Lys Arg Asn Thr Asn Arg Arg Pro Gln Asp Val Lys Phe Pro
 1 5 10 15

Gly Gly Gly Gln Ile Val Gly Gly Val Tyr Leu Leu Pro Arg Arg Gly
 20 25 30

-continued

Pro Arg Leu Gly Val Arg Ala Thr Arg Lys Thr Ser
 35 40

<210> SEQ ID NO 43
 <211> LENGTH: 25
 <212> TYPE: PRT
 <213> ORGANISM: Hepatitis C Virus
 <220> FEATURE:
 <223> OTHER INFORMATION: ALAM-18

<400> SEQUENCE: 43

Lys Thr Lys Arg Asn Thr Asn Arg Arg Pro Gln Asp Val Lys Phe Pro
 1 5 10 15

Gly Gly Gly Gln Ile Val Lys Thr Ser
 20 25

<210> SEQ ID NO 44
 <211> LENGTH: 18
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

<400> SEQUENCE: 44

gatcgctcga attcctcg 18

<210> SEQ ID NO 45
 <211> LENGTH: 20
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Primer

<400> SEQUENCE: 45

cgaggaattc gagcgatcctt 20

<210> SEQ ID NO 46
 <211> LENGTH: 18
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 46

Met Ser Thr Asn Pro Lys Pro Gln Lys Lys Asn Lys Arg Asn Thr Asn
 1 5 10 15

Arg Arg

<210> SEQ ID NO 47
 <211> LENGTH: 18
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 47

Asn Lys Arg Asn Thr Asn Arg Arg Pro Gln Asp Val Lys Phe Pro Gly
 1 5 10 15

Gly Gly

<210> SEQ ID NO 48

-continued

<211> LENGTH: 18
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 48

Asp Val Lys Phe Pro Gly Gly Gly Gln Ile Val Gly Gly Val Tyr Leu
 1 5 10 15

Leu Pro

<210> SEQ ID NO 49
 <211> LENGTH: 18
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 49

Val Gly Gly Val Tyr Leu Leu Pro Arg Arg Gly Pro Arg Leu Gly Val
 1 5 10 15

Arg Ala

<210> SEQ ID NO 50
 <211> LENGTH: 18
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 50

Gly Pro Arg Leu Gly Val Arg Ala Thr Arg Lys Thr Ser Glu Arg Ser
 1 5 10 15

Gln Pro

<210> SEQ ID NO 51
 <211> LENGTH: 18
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 51

Lys Thr Ser Glu Arg Ser Gln Pro Arg Gly Arg Arg Gln Pro Ile Pro
 1 5 10 15

Lys Ala

<210> SEQ ID NO 52
 <211> LENGTH: 18
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 52

Arg Arg Gln Pro Ile Pro Lys Ala Arg Arg Pro Glu Gly Arg Thr Trp
 1 5 10 15

Ala Gln

<210> SEQ ID NO 53
 <211> LENGTH: 18

-continued

<212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 53

Pro Glu Gly Arg Thr Trp Ala Gln Pro Gly Tyr Pro Trp Pro Leu Tyr
 1 5 10 15

Gly Asn

<210> SEQ ID NO 54
 <211> LENGTH: 19
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 54

Gln Tyr Pro Trp Pro Leu Tyr Gly Asn Glu Gly Cys Gly Trp Ala Gly
 1 5 10 15

Trp Leu Leu

<210> SEQ ID NO 55
 <211> LENGTH: 17
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 55

Cys Gly Trp Ala Gly Trp Leu Leu Ser Pro Arg Gly Ser Arg Pro Ser
 1 5 10 15

Trp

<210> SEQ ID NO 56
 <211> LENGTH: 25
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 56

Trp Leu Leu Ser Pro Arg Gly Ser Arg Pro Ser Trp Gly Pro Thr Asp
 1 5 10 15

Pro Arg Arg Arg Ser Arg Asn Leu Gly
 20 25

<210> SEQ ID NO 57
 <211> LENGTH: 25
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 57

Ser Trp Gly Pro Thr Asp Pro Arg Arg Arg Ser Arg Asn Leu Gly Lys
 1 5 10 15

Val Ile Asp Thr Leu Thr Cys Gly Phe
 20 25

<210> SEQ ID NO 58

-continued

<211> LENGTH: 25
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 58

Ser Arg Asn Leu Gly Lys Val Ile Asp Thr Leu Thr Cys Gly Phe Ala
 1 5 10 15

Asp Leu Met Gly Tyr Ile Pro Leu Val
 20 25

<210> SEQ ID NO 59
 <211> LENGTH: 25
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 59

Leu Thr Cys Gly Phe Ala Asp Leu Met Gly Tyr Ile Pro Leu Val Gly
 1 5 10 15

Ala Pro Leu Gly Gly Ala Ala Arg Ala
 20 25

<210> SEQ ID NO 60
 <211> LENGTH: 25
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 60

Tyr Ile Pro Leu Val Gly Ala Pro Leu Gly Gly Ala Ala Arg Ala Leu
 1 5 10 15

Ala His Gly Val Arg Val Leu Glu Asp
 20 25

<210> SEQ ID NO 61
 <211> LENGTH: 25
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 61

Gly Ala Ala Arg Ala Leu Ala His Gly Val Arg Val Leu Glu Asp Gly
 1 5 10 15

Val Asn Tyr Ala Thr Gly Asn Leu Pro
 20 25

<210> SEQ ID NO 62
 <211> LENGTH: 23
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 62

Leu Glu Asp Gly Val Asn Tyr Ala Thr Gly Asn Leu Pro Gly Cys Ser
 1 5 10 15

Phe Ser Ile Phe Leu Leu Ala

-continued

20

```

<210> SEQ ID NO 63
<211> LENGTH: 23
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: HCV-Core derived peptides

<400> SEQUENCE: 63

Leu Pro Gly Cys Ser Phe Ser Ile Phe Leu Leu Ala Leu Leu Ser Cys
 1             5             10             15

Leu Thr Val Pro Ala Ser Ala
 20

```

1. A method of simultaneously detecting at least one Hepatitis C Virus (HCV) antigen and at least one HCV antibody in a test sample comprising the steps of:

- a) contacting said test sample with: 1) at least one HCV viral antigen or portion thereof coated on a solid phase, for a time and under conditions sufficient for the formation of antibody/antigen complexes and 2) at least one antibody to HCV or portion thereof coated on said solid phase, for a time and under conditions sufficient for the formation of antigen/antibody complexes;
- b) detecting said antibody/antigen complexes, presence of said complexes indicating presence of HCV antigen in said test sample; and
- c) detecting said antigen/antibody complexes, presence of said complexes indicating presence of HCV antibody in said test sample.

2. The method of claim 1 wherein said at least one HCV antigen coated on the solid phase is selected from the group consisting of core antigen, NS3, NS4, NS5, and portions thereof.

3. The method of claim 2 wherein said at least one antibody coated on said solid phase is a monoclonal antibody selected from the group consisting of 13-959-270, 14-1269-281, 14-1287-252, 14-153-234, 14-153-462, 14-1705-225, 14-1708-269, 14-1708-403, 14-178-125, 14-188-104, 14-283-112, 14-635-225, 14-726-217, 14-886-216, 14-947-104, 14-945-218, 107-35-54, 110-81-17, 13-975-157, 14-1350-210, C11-3, C11-7, C11-10, C11-14 and C11-15.

4. The method of claim 3 wherein said at least one antibody coated on the solid phase is not reactive with said at least one antigen coated on the solid phase.

5. The method of claim 1 wherein said at least one antibody is a HCV anti-core monoclonal antibody and said at least one antigen is a recombinant HCV core protein.

6. The method of claim 5 wherein said recombinant core protein does not contain epitopes to which said anti-core monoclonal antibody binds.

7. The method of claim 1 wherein said solid phase is a microparticle.

8. A method for simultaneously detecting the presence of at least one HCV antigen and at least one HCV antibody in a test sample comprising the steps of:

- a) contacting said test sample with: 1) at least one HCV viral antigen or portion thereof coated on a solid phase, for a time and under conditions sufficient for the formation of antibody/antigen complexes and 2) at least one HCV antibody or portion thereof coated on said solid phase, for a time and under conditions sufficient for the formation of antigen/antibody complexes;
- b) adding a conjugate to the resulting antibody/antigen complexes for a time and under conditions sufficient to allow said conjugate to bind to the bound antibody in (a)(1), wherein said conjugate comprises a second antibody attached to a chemiluminescent compound capable of generating a detectable signal and simultaneously adding a second conjugate to the resulting antigen/antibody complexes for a time and under conditions sufficient to allow said conjugate to bind to the bound antigen in (a)(2), wherein said conjugate comprises a third antibody attached to said chemiluminescent compound capable of generating a detectable signal; and
- c) detecting said generated signal, presence of said signal indicating presence of at least one antigen in said test sample selected from the group consisting of HCV antigen and HCV antibody.

9. The method of claim 8 wherein said at least one HCV antigen coated on the solid phase is selected from the group consisting of core antigen, NS3, NS4, NS5, and portions thereof.

10. The method of claim 9 wherein said at least one antibody coated on said solid phase is a monoclonal antibody selected from the group consisting of 13-959-270, 14-1269-281, 14-1287-25, 14-153-234, 14-153-462, 14-1705-225, 14-1708-269, 14-1708-403, 14-178-125, 14-188-104, 14-283-112, 14-635-225, 14-726-217, 14-886-216, 14-947-104, 14-945-218, 13-975-157 and 14-1350-210, 107-35-54, 110-81-17, C11-3, C11-7, C11-10, C11-14 and C11-15.

11. The method of claim 10 wherein said at least one monoclonal antibody coated on the solid phase is not reactive with said at least one antigen coated on the solid phase.

12. A kit comprising:

- a) a container containing at least one HCV antigen coated on a solid phase; and

- b) a container containing at least one HCV antibody coated on a solid phase.
- 13.** A kit comprising:
- a container containing: 1) at least one HCV antigen coated on a solid phase and 2) at least one HCV antibody, coated on said solid phase.
- 14.** The kit of claim 12 or claim 13 further comprising at least one conjugate comprising a signal-generating compound attached.
- 15.** The kit of claim 14 wherein said signal-generating compound is acridinium.
- 16.** A method of detecting at least one HCV antigen in a test sample comprising the steps of:
- a) contacting said test sample with at least one HCV antibody coated on a solid phase, for a time and under conditions sufficient for the formation of antibody/antigen complexes; and
- b) detecting the presence of antibody/antigen complexes, presence of said complexes indicating presence of said at least one HCV antigen in said test sample.
- 17.** A method of detecting at least one HCV antigen in a test sample comprising the steps of:
- a) contacting said test sample with at least one HCV antibody coated on a solid phase, for a time and under conditions sufficient for the formation of antibody/antigen complexes;
- b) adding a conjugate to the resulting antibody/antigen complexes for a time and under conditions sufficient to allow said conjugate to bind to the bound at least one antibody, wherein said conjugate comprises a second antibody attached to a chemiluminescent compound capable of generating a detectable signal; and
- c) detecting said signal generated by said chemiluminescent compound, a signal generated by said chemiluminescent compound indicating the presence of at least one HCV antigen in said test sample.
- 18.** A recombinant protein comprising an amino acid sequence selected from the group consisting of SEQ ID NO: 6, SEQ ID NO: 8, SEQ ID NO: 12, SEQ ID NO: 16, and conservative amino acid substitutions thereof.
- 19.** A recombinant protein comprising an amino acid sequence encoded by a nucleotide sequence selected from the group consisting of, for example, SEQ ID NO: 5, SEQ ID NO: 7, SEQ ID NO: 11 and SEQ ID NO: 15.
- 20.** A vector or construct comprising a nucleotide sequence selected from the group consisting of SEQ ID NO: 5, SEQ ID NO: 7, SEQ ID NO: 11 and SEQ ID NO: 15.
- 21.** A host cell comprising said vector or construct of claim 20.
- 22.** An immunoassay which simultaneously detects at least one HCV antigen and at least one HCV antibody in a test sample.

* * * * *

专利名称(译)	同时检测HCV抗原和HCV抗体的方法		
公开(公告)号	US20030152948A1	公开(公告)日	2003-08-14
申请号	US10/173480	申请日	2002-06-17
[标]申请(专利权)人(译)	SHAH DINESHØ DAWSON GEORGE J MUERHOFF斯科特 蒋LILY 古铁雷斯罗宾 LEARY THOMAS P 德赛Suresh STEWART詹姆斯·L·		
申请(专利权)人(译)	SHAH DINESH O. DAWSON GEORGE J. MUERHOFF SCOTT A. 蒋LILY 古铁雷斯ROBIN A. LEARY THOMAS P. 德赛Suresh 斯图尔特, 杰姆斯L.		
当前申请(专利权)人(译)	亚培		
[标]发明人	SHAH DINESH O DAWSON GEORGE J MUERHOFF A SCOTT JIANG LILY GUTIERREZ ROBIN A LEARY THOMAS P DESAI SURESH STEWART JAMES L		
发明人	SHAH, DINESH O. DAWSON, GEORGE J. MUERHOFF, A. SCOTT JIANG, LILY GUTIERREZ, ROBIN A. LEARY, THOMAS P. DESAI, SURESH STEWART, JAMES L.		
IPC分类号	G01N33/543 C07K14/18 C12N1/15 C12N1/19 C12N1/21 C12N5/10 C12N15/09 C12P21/04 G01N21/76 G01N33/576 C12Q1/70 C12Q1/68 G01N33/53 C12N5/06 C12N5/16		
CPC分类号	C07K14/005 C12N2770/24222 G01N2469/20 G01N2469/10 G01N33/5767		
其他公开文献	US6727092		
外部链接	Espacenet USPTO		

摘要(译)

本发明涉及同时检测丙型肝炎病毒 (HCV) 抗原以及响应HCV抗原产生的抗体的方法。此外, 本发明允许人们甚至在抗体发展之前检测感染的早期急性期中的抗原, 从而允许早期检测受感染的血液和血液制品, 从而提高血液供应的安全性。

A PROTOTYPE ABBOTT PRISM[®] HCV Ag/Ab COMBO
ASSAY FORMAT

