



US 20160047807A1

(19) **United States**

(12) **Patent Application Publication**  
**Koulchin et al.**

(10) **Pub. No.: US 2016/0047807 A1**  
(43) **Pub. Date: Feb. 18, 2016**

(54) **METHOD FOR DETECTION OF  
LEGIONELLA BACTERIA EMPLOYING  
PURIFIED ANTIGEN-SPECIFIC ANTIBODIES**

(71) Applicant: **Alere Scarborough, Inc.**, Scarborough,  
ME (US)

(72) Inventors: **Vladimir A. Koulchin**, Portland, ME  
(US); **Norman J. Moore**, North  
Berwick, ME (US); **Elena V. Molokova**,  
Portland, ME (US)

(21) Appl. No.: **14/828,823**

(22) Filed: **Aug. 18, 2015**

**Related U.S. Application Data**

(63) Continuation of application No. 09/139,720, filed on  
Aug. 25, 1998, now Pat. No. 9,134,303.

**Publication Classification**

(51) **Int. Cl.**  
**G01N 33/569** (2006.01)  
**G01N 21/78** (2006.01)

(52) **U.S. Cl.**  
CPC ..... **G01N 33/56911** (2013.01); **G01N 21/78**  
(2013.01); **G01N 2400/10** (2013.01); **G01N**  
**2469/10** (2013.01); **G01N 2021/752** (2013.01)

(57) **ABSTRACT**

The present invention involves extracting from *Legionella* bacteria, particularly *L. pneumophila* bacteria, an essentially protein-free O-polysaccharide or carbohydrate antigen, coupling this antigen to an activated chromatographic column through a protein space molecule which is first conjugated to the antigen, utilizing the column thus prepared for the affinity purification of raw polyvalent antibodies to the same *Legionella* bacterium from which the O-polysaccharide or carbohydrate antigen was separated—thereby obtaining antigen-specific antibodies which are useful for the rapid detection of the corresponding *Legionella* bacterium or its antibody in human bodily fluids such as urine, sputum, blood and the like or in environmental samples suspected of harboring the *Legionella* bacterium.




FIG. 1




FIG. 1A





FIG. 1B

FIG. 1C





---

1 2 3 4 5 6 7 8 9

FIG. 2

**METHOD FOR DETECTION OF  
LEGIONELLA BACTERIA EMPLOYING  
PURIFIED ANTIGEN-SPECIFIC ANTIBODIES**

CROSS-REFERENCE TO RELATED  
APPLICATIONS

**[0001]** This application is a continuation application of U.S. application Ser. No. 09/139720, filed on Aug. 25, 1998; the contents of which is specifically incorporated by reference herein in its entirety.

**[0002]** This invention relates to essentially protein-free carbohydrate including polysaccharide, antigens separated from bacteria of the genus *Legionella*, and especially from serogroups and/or strains of *Legionella pneumophila*, including the O-polysaccharide antigen of *L. pneumophila* serogroup 1, and to the use of these antigens in the affinity purification of poly-valent antibodies to corresponding *Legionella* organisms. More particularly, the invention encompasses coupling the carbohydrate or polysaccharide antigen separated from a *Legionella* bacterium to an activated chromatographic column and using that column for affinity purification of the polyclonal antibodies to the same species, or the same serogroup of a species, of *Legionella*. The invention further encompasses the use of the affinity-purified polyvalent antibodies produced in immunochemical assays for the detection of *Legionella*-caused diseases such as Legionnaires disease and Pontiac fever in human patients and for the detection of environmental sources of *Legionella* infectious agents.

BACKGROUND

**[0003]** Legionnaires' disease, a pneumonia-like human infection caused by gram-negative bacteria of the genus *Legionella*, is virtually impossible to differentiate on a reliable clinical basis (involving assessment of patient symptoms without laboratory tests) from pneumonia and other similar lung infections. Because the disease may produce lung abscesses, infections in other bodily organs or bacteremia, and its mortality rate is significantly increased by delay in commencing appropriate therapy, there is a need for rapid and reliable diagnostic tests which has to date been only partially met. Stout, J. E. and Yu, V. L., "Legionellosis", 337, N. Eng. J. Med. 682-687 (1997). Efforts to develop such tests have been hampered by the fact that there are a number of *Legionella* species, at least some of which are known to have a number of distinct serogroups.

**[0004]** Legionnaires' disease was first recognized in the summer of 1976 and a number of techniques for detecting *Legionella* ("*L.*") *pneumophila*, which is now known to account for some 90% of cases (Stout and Yu, supra), have been developed in the interim. In general, these tests are time-consuming and incapable of identifying more than one serogroup of the 15 serogroups known so far to fall within the *L. pneumophila* species. It should be noted, however, that, as reported by Stout and Yu, supra, more than 80% of reported cases of Legionnaires' disease are attributable to *L. pneumophila* serogroup 1—a fact which makes the development of a rapid, reliable immunoassay for that entity of particular importance and has led researchers to focus on this as a priority.

**[0005]** The early efforts to establish the identity of the causative agent of Legionnaires' disease depended largely upon culturing the bacteria for 5 to 6 days and examining the culture microscopically. Efforts to speed up the identifica-

tion process led to numerous immunochemical tests of varying sensitivity and specificity including, inter alia, the presently commercially available EQUATE' radioimmunoassay ("RIA"), and the Binax enzyme immunoassay ("EIA"), both of which are sold in kit form by applicants' assignee, Binax, Inc. As indicated by Hackman, B. A. et al. in "Comparison of Binax *Legionella* Urinary Antigen EIA Kit with Binax RIA Urinary Antigen Kit for Detection of *Legionella pneumophila* Serogroup I Antigen", 34 *J. Clin. Microbiol.* 1579-1580 (1996), both of these assays were found to be specific for *L. pneumophila* Serogroup I antigen. The reported sensitivity of the EIA was 77%; that of the RIA was higher. The article indicates that each can be performed within "less than 3 h[ours] from beginning to end"; in Binax's own tests each requires at least 2 1/2 hours to perform. More information about this EIA assay appears in Kazandjian, D. et al., "Rapid Diagnosis of *Legionella pneumophila* Serogroup I Infection with Binax Enzyme Immunoassay Urinary Antigen Test", 35 *J. Clin. Immunobiol.* 954-956 (1997).

BRIEF DESCRIPTION OF THE INVENTION

**[0006]** Perhaps the most significant advantage of the immuno-chromatographic test ("ICT") described herein is its ability to give a test result within a 15-minute time span for the presence or absence of *L. pneumophila* Serogroup I (or its antigen), which result is of high specificity and sensitivity. The speed with which this test can be reliably conducted to yield a result of high specificity and sensitivity is believed to be due to the strongly reactive nature of the affinity purified antibodies prepared in accordance with this invention. The superior reactive properties of these antibodies, in turn, is believed to be attributable to the use for affinity purification of the novel, essentially protein-free O-polysaccharide antigen of *L. pneumophila* serogroup 1 which is also a part of this invention.

**[0007]** Another ICT test has been made possible by the separation, according to this invention, of an antigen of *L. pneumophila* serogroup 5 which is of non-proteinaceous, carbohydrate nature and is common to multiple serogroups of *L. pneumophila*. The presence of such a common antigen has been suggested in the scientific literature; see e.g., Nolte, F. S. et al., "Electrophoretic and Serological Characterization of the Lipopolysaccharides of *Legionella pneumophila*", 52 *Infection and Immunity* 676-681 (1986); Otten, S. et al., "Serospecific Antigens of *Legionella pneumophila*", 167 *J. Bacteriol.* 893-904 (1986); Barthe, C. et al., "Common Epitope on the Lipopolysaccharide of *Legionella pneumophila* Recognized by a Monoclonal Antibody", 26 *J. Clin. Microbiol.* 1016-1023 (1988); Knirel, Y. A. et al., "The Structure of the O-specific Chain of *Legionella pneumophila* Serogroup I", 227 *Eur. J. Biochem.* 239-245 (1994). Heretofore there has been no clear report of any separation of such an antigen for use in developing a useful broad spectrum immuno-chemical assay.

**[0008]** This invention comprises the extraction from any *Legionella* bacterium, and especially from a bacterium of any of the serogroups of *L. pneumophila* of an essentially protein-free carbohydrate antigen, the preparation of a conjugate of this antigen with a spacer protein, the coupling of the conjugate to an affinity column, the use of that column for the purification of polyvalent antibodies to the corresponding *Legionella* bacterium, such as a bacterium of a serogroup of *L. pneumophila* and the use of the antibodies thus purified in an ICT immunoassay for detecting *Legionella* bacteria or

their antigens, including antigens of *L. pneumophila* serogroups or their corresponding bacteria, as more specifically described hereinafter.

[0009] In particular, the invention includes the separation of an essentially protein-free O-polysaccharide antigen specific to *L. pneumophila* serogroup 1, its use in the affinity purification of the polyvalent antibody specific to the same microorganism and the use of that affinity purified antibody in an ICT immunoassay of high specificity and high sensitivity that is performable within 15 minutes.

[0010] In another embodiment, the invention includes the separation from *L. pneumophila* serogroup 5 of an essentially protein-free carbohydrate antigen, and its use in the affinity purification of the polyclonal antibody specific to the *L. pneumophila* serogroup 5 antigen (i.e., an antibody that was obtained from a rabbit immunized with the said antigen); this antibody showed an ability to cross-react with antigens of the *L. pneumophila* serotypes 1, 2 and 4 in addition to the antigen of serotype 5.

#### DESCRIPTION OF THE DRAWINGS

[0011] FIG. 1 and its related FIGS. 1A, 1B and 1C hereof show the structure of a typical ICT device which is suitably adapted to perform the *L. pneumophila* serogroup 1 specific assay, as described in Examples VII, VIII and IX hereof

[0012] FIG. 2 hereof shows the results of Western blot analyses of phosphate-buffered saline extracts of *L. pneumophila* serogroups 1, 2, 4 and 5 with *L. pneumophila* serogroup 5, essentially protein-free carbohydrate antigen-affinity-purified polyclonal antibody specific to serogroup 5 antigen, which analysis is described in Example X hereof.

#### DETAILED DESCRIPTION OF THE INVENTION

[0013] Previous experience in the art, including experience with the Binax RIA assay for *L. pneumophila* Serogroup 1 antigen sold under the trademark EQUATE and the Binax EIA assay for the same antigen, has shown that antigens of the *Legionella* species, including antigens of various *L. pneumophila* serogroups, are more conveniently detectable in specimens from the urine of patients infected with a *Legionella* microorganism than in specimens of blood, sputum or other fluids. This is in part because the *Legionella* antigens usually appear in urine within 1-3 days after infection of a human patient whereas their appearance at a detectable level in blood may occur later, and also in part because patients infected with Legionnaires' disease often do not produce much sputum. The ICT test which forms a part of this invention can be configured to run on blood, sputum or some other fluid such as cerebrospinal fluid, or on aqueous samples of environmental origin. It is noted that urine is generally the preferred sample fluid for diagnosis of human patients because it can be obtained non-invasively and easily, even in the doctor's office, and it is not as readily contaminated with other microorganisms, e.g., oral microflora present in sputum, which are innocuous but may affect the results obtained.

[0014] Broadly speaking, the ICT test for *L. pneumophila* serogroup 1 antigen that is specifically a part of the present invention may be designed to be run in any known disposable ICT device disclosed in the art. Preferably, the test is conducted using an ICT device of the type disclosed in copending U.S. application Ser. No. 07/706,639 of Howard Chandler, or one of its continuation-in-part applications, all assigned to Smith-Kline Diagnostics, Inc. but exclusively licensed to

Binax, Inc. in a wide area of applications that includes the present diagnostic field—i.e., diseases of the respiratory system. The device is suitably impregnated with the affinity purified polyvalent antibodies herein disclosed which are specific to the antigen of *L. pneumophila* serogroup 1 antigen. Positive results of the assay are shown by the appearance of a color upon reaction of suitably labeled antibodies with the antigen. Suitable labels may be any of those known in the art to produce visible color when antibodies conjugated thereto react with antigen, including finely divided metallics and various other labeling materials. Colloidal gold is the preferred label.

[0015] The invention contemplates that essentially protein free carbohydrate antigens may similarly be separated from each of the *Legionella* bacteria, including from bacteria of other *Legionella* species and from bacteria of other serogroups of *L. pneumophila*, that such antigens may similarly be utilized in the affinity purification of polyvalent antibodies to the corresponding *Legionella* species or serogroup bacteria and its antigen, and that these affinity-purified antibodies may be utilized in ICT and other immunoassays as specifically disclosed herein for affinity-purified antibodies to *L. pneumophila* serogroup 1.

[0016] Preliminary to the preparation of the device is obtaining the antibodies and effecting their affinity purification. The antibodies that may be used in this invention are conventional polyvalent (also called "polyclonal") antibodies obtained by the well-known process of immunizing a rabbit, goat or other animal to a *Legionella* antigen, e.g., *L. pneumophila* antigen of known serogroup and bleeding the animal after the passage of an appropriate time period to obtain serum containing the desired antibodies. See, e.g., Cherry, U. B. and McKinney, R. M., pp. 91-104 in Jones, G. L. and Herbert, G. A. (Eds.), "*Legionnaires' the disease, the bacterium and the methodology*", (Center for Disease Control, Atlanta, 1979). In this invention, relative to the ICT immunoassay for *L. pneumophila* serogroup 1 disclosed herein, the rabbit or other animal is immunized to *L. pneumophila* serogroup 1 antigen and the polyvalent antibodies recovered from its blood are antibodies to *L. pneumophila* serogroup 1.

[0017] The polyvalent antibodies to be used in this invention are further subjected to affinity purification on a specially prepared chromatographic column which employs the essentially protein-free carbohydrate, including polysaccharide, antigen of the same bacterium against which the antibodies are reactive as the purifying agent for them.

[0018] Preliminary to the affinity purification of the antibodies, it is therefore necessary according to this invention to prepare an essentially protein free purified carbohydrate antigen from a culture of known *Legionella* bacteria of the desired species or serogroup of a species.

[0019] The following examples I and II explain how the purified protein free carbohydrate antigen is obtained.

#### EXAMPLE I

[0020] Bacteria of *L. pneumophila* Serogroup 1 (strain Philadelphia-1) were obtained from Centers for Disease Control and Prevention (Atlanta, Ga.) and cultured on charcoal-yeast extract agar plates obtained from Northeast Laboratory (Waterville, Me.) for a period of 72 hours at 37° C. Cells were harvested with phosphate-buffered saline ("PBS") at pH 7.2 containing 0.2% of NaN<sub>3</sub> and collected by centrifugation at 8000 rpm for 25 minutes. The resulting cell pellet was stored at -20° C. until used.

[0021] Wet cells from this pellet were suspended in 20 ml. 0.1 M NaOH per gram of wet cells and stirred at room temperature for 45 minutes. The pH of the solution was then adjusted with concentrated acetic acid to 3.0 and the solution was subjected to centrifugation at 8000 rpm for 20 minutes. The supernatant from this step was then neutralized with aqueous NaOH and dialyzed against distilled water.

[0022] The resulting dialyzate was concentrated 10 times on a rotary vacuum evaporator and then sonicated for 5 minutes in an ultrasonic bath.

[0023] Proteinase K, in a concentration of 0.2 mg. per ml. of the concentrated product, was added to digest the remaining proteins and the mixture was incubated at 40° C. overnight. The next step was the addition of further Proteinase K, in a concentration of 0.1 mg. per ml., to the mixture, followed by further overnight incubation at 40° C. This second incubation was followed by concentration of the product on a rotary evaporator to a small volume, adjustment of its pH to 10-11 with 0.2% triethylamine and application of the thus-treated mixture to a column of Sephacryl S-200 from Pharmacia, equilibrated with 0.02% triethylamine. Material eluted in the first peak was pooled, adjusted with 0.1 N HCl to approximately neutral pH and dialyzed against distilled water for 18 hours, followed by lyophilization.

[0024] The yield of O-polysaccharide antigen from 16.5 grams of wet cells of *L. pneumophila* serogroup 1 strain Philadelphia-1 was 62 mg.

[0025] It should be noted that in place of Proteinase K, any broad spectrum protease enzyme preparation may be used. What is important in this procedure is the elimination of protein from the antigen to the maximum possible extent within reasonable limits of time feasibility.

#### EXAMPLE II

[0026] Example I was repeated using *L. pneumophila* serogroup 5 strain Dallas IE as the bacteria in the culture step. 11.5 wet grams yielded 21 mg. of a carbohydrate antigen. At a later time, the procedure was again repeated, this time using *L. pneumophila* serogroup 5 (strain U8W) as the bacteria in the culture step.

[0027] It is within the scope of this invention to repeat the culturing and extraction steps as herein described on any other Legionella bacterium and especially on bacteria of any other of the serogroups of *L. pneumophila*, i.e., on any of serogroups 2, 3, 4 and 6-15 inclusive, to obtain a carbohydrate antigen essentially free of proteinaceous material.

[0028] In order to couple the essentially protein free O-polysaccharide product of Example I or the carbohydrate of Example II to a chromatographic column for use in affinity purification of polyvalent antibodies to the corresponding bacterium from which the carbohydrate antigen was extracted, it is necessary to complex it with a protein spacer molecule suitably prepared to ensure that it will stably bond to both the carbohydrate antigen and the column and will itself remain inert during the affinity purification step. A modified bovine serum albumin ("BSA") was chosen as the preferred spacer molecule and was prepared as in Example III:

#### EXAMPLE III

##### Preparation of Protein Spacer Molecule

[0029] An 0.5 M aqueous solution of hydrazine dihydrochloride from Aldrich Chemical Co., Inc. (Milwaukee, Wis.)

was prepared and its pH was adjusted to 5.2 with dry NaOH. The solution was then mixed with dry BSA from Sigma Chemical Co. to a final concentration of 25 mg. per ml. of BSA. After the BSA was completely dissolved, N-(dimethylaminopropyl)-N<sup>1</sup>-ethylcarbodiimide hydrochloride obtained from Fluka Chemical Co. (St. Louis, Mo.) was added to a final concentration of 2.5 mg. per ml. The mixture was stirred overnight at ambient temperature and then dialyzed for a period of about five days against distilled water at 4° C., with daily changes of the water.

[0030] In lieu of this modified BSA, it is contemplated that other appropriate spacer molecules may be used.

[0031] The conjugation of the O-polysaccharide antigen from serogroup 1 of *L. pneumophila* to the spacer molecule was conducted as follows:

#### EXAMPLE IV

##### Conjugation of O-polysaccharide Antigen to Spacer Molecules

[0032] *L. pneumophila* serogroup 1 essentially protein-free O-polysaccharide antigen was dissolved in distilled water in a concentration of 4-5 mg/ml and the pH was adjusted to 5.0 with 1 M HCl. The modified BSA solution prepared as in Example III was adjusted to pH 5.0 with 1 M HCl and slowly added, in a weight ratio of 4-to-1 of the O-polysaccharide antigen solution, to the latter. After 5 minutes of stirring, about 100-200 mcl of distilled water containing the N-(dimethylamino-propyl)-N<sup>1</sup>-ethylcarbodiimide hydrochloride referred to in Example III was added in a weight ratio of 1:2 to the O-polysaccharide antigen/BSA solution. The resulting mixture was stirred at ambient temperature overnight.

[0033] To separate the conjugate of the O-polysaccharide antigen with modified BSA from any unreacted materials present, the reaction mixture was chromatographed on a Sepharose CL-4B column equilibrated with a buffer of 0.1 M NaH<sub>2</sub>PO<sub>4</sub>:0.5 M NaCl. All of the chromatography fractions were subjected to testing for antigen activity using the commercial Binax EIA test referred to above. All fractions that showed antigen activity in the tested were pooled in and used for affinity column preparation as shown in Example V.

#### EXAMPLE V

##### Coupling of Conjugate to Activated Chromatographic Column

[0034] An immunoadsorbent gel was prepared by conventionally activating Sepharose CL-4B with cyanogen bromide. The ligand of O-polysaccharide antigen with modified BSA, prepared as in Example IV, was coupled to it using procedures known in the art, e.g., as described in Hermanson, G. T., et al. in *Immobilized Affinity Ligand Techniques*, 53-56 (Academic Press, Inc. 1992). The gel was then packed in a column and washed successively with 5-10 volumes per volume of gel with PBS of pH 7.2, triple strength PBS of pH 7.2 and 0.2 M glycine-HCl of pH 2.5. The resulting activated column was used, as described below, for affinity purification of polyvalent antibodies to *L. pneumophila* serogroup 1 antigen. The column after elution should be stored in PBS or another neutral buffer until used.

[0035] Instead of cyanogen-bromide activated Sepharose, spherulose and various commercially available activated columns could be used in this step.

## EXAMPLE VI

## Affinity Purification of Antibodies

**[0036]** The activated column described in Example V was used for the affinity chromatography of rabbit-anti-*L. pneumophila* serogroup 1 polyvalent antibodies to *L. pneumophila* serotype 1 antigen according to the method described by Harlow E. and Lane, D. in *Antibodies: A Laboratory Manual* at 313-315 (Cold Spring Harbor Laboratory, 1988). Elution of the antibodies from the column was effected with 0.2M glycine-HCL buffer of pH 2.5. Alternative eluants such as 3 M NaSCN or 0.1 M Et<sub>3</sub>N could be substituted.

**[0037]** These affinity-purified antibodies were utilized in an ICT test specific to *L. pneumophila* serogroup 1 as described in the following example.

## EXAMPLE VII

## ICT Device and Its Preparation

**[0038]** A. Preparation of Test Device:

**[0039]** A test device comprising a hinged cardboard housing equipped with a window to allow the viewing of both the test results and control results was prepared as shown in FIG. 1. The device has a recess into which is placed a preformed plastic swab well for receiving the sample-wetted swab on the right-hand (labeled 1 in the drawing). An overlable shown in FIG. 1A is then placed over the entire right-hand side of the device. The overlable has been equipped with two holes—a lower one (marked B on FIG. 1A) into which the saturated swab is to be inserted and an upper one (marked A on FIG. 1A) toward which the swab will be pushed after insertion thereof into the hole B. The position of the overlable with its holes A and B, and the swab well cooperate to hold the swab in a proper position during the assay and to promote the expulsion of sorbed liquid from the swab.

**[0040]** A preassembled test strip (marked C on FIG. 1) described below, is inserted into the recess (labeled 2 on FIG. 1) and held in place by an adhesive applied to the bottom thereof. An overlable shown in FIG. 1B is placed atop the left-hand side. It has been equipped with a single hole (marked D in FIG. 1B) which mates to the right-hand side hole A when the device is closed for performance of the assay.

**[0041]** The assembled device is stored in a sealed pouch with desiccant until it is used. Prior to sealing the pouch and storing, a lightly adhesive tape is placed on the outer edge of the right-hand half of the device.

**[0042]** B. Construction and Preparation of the Preassembled Test Strip

**[0043]** FIG. 1C shows the construction of the preassembled strip. It is comprised of a conjugate pad of sorbent material in which a conjugate of gold particles and the affinity-purified rabbit anti-*Legionella pneumophila* serogroup 1 antibodies described above have been impregnated. In contact with this pad is a nitrocellulose pad onto which a capture line for the sample which reacts with the conjugate has been established by embedding a stripe of affinity-purified rabbit anti-*L. pneumophila* serogroup 1 antibodies, prepared as described above. The nitrocellulose pad also has a downstream control line established by striping the pad with goat anti-rabbit immunoglobulin (IgG). Following the nitrocellulose pad, the strip is ended by an absorbent pad which serves as a reservoir for liquid. All of these pads are backed by an adhesive strip when the device is ready to ship.

**[0044]** The conjugate pad is normally made from non-woven polyester or extruded cellulose acetate. To prepare this pad for use in the assay, gold particles of 50 nm. diameter are conjugated to affinity-purified rabbit anti-*Legionella pneumophila* serotype 1 antibodies prepared as described above. The conjugation is effected using a known method such as that described by DeMay' in Polak, J. M. and Van Norden, S. (Eds.), *Immunochemistry: Modern Methods and Application*, (Wright, Bristol, England, 1986). The gold conjugate particles are mixed with a drying agent consisting of aqueous 5 mM sodium tetraborate of pH 8.0 containing 1.0% BSA, 0.1% Triton X-100, 2.0% Tween 20, 6.0% sucrose and 0.02% sodium azide. The pad is heated sufficiently to remove all of the liquid present and stored in a low-humidity environment pending, assembly of the test strip. These pads and their treatment are especially chosen so that the pads will hold the dry conjugate and will release it only when later wetted by sample.

**[0045]** The nitrocellulose pad is first treated by embedding a stripe of affinity purified rabbit anti-*L. pneumophila* serotype 1 antibodies in a first portion thereof, using a carrier solution of phosphate buffered saline. These antibodies act as the capture line. In a second portion of the pad downstream of the first one in the assembled test device, the control line is established by striping goat anti-rabbit IgG in the same carrier solution on the surface of the pad. The nitro-cellulose pad is then subjected to desiccation at 18-25° C. to promote permanent absorption of the protein stripes thereto.

**[0046]** The absorbent pad used is of a commercially available cellulosic material sold under the name Ahlstrom 939. This pad requires no special treatment.

**[0047]** C. Kit Preparation

**[0048]** As sold in commerce, the test device containing the finished test strip is assembled. In practice, a number of devices are packaged with a commensurate number of swabs fashioned from fibrous Dacron and a bottle of "Reagent A" equipped with a top adapted to deliver Reagent A dropwise. "Reagent A" is a solution of 2.0% Tween 20, 0.05% sodium azide and 0.5% sodium dodecyl sulfate in a 0.05 M sodium citrate-sodium phosphate buffer of pH 6.5. Positive and negative controls are also included in each kit.

**[0049]** The use of the finished test devices to identify *L. pneumophila* serogroup 1 O-polysaccharide antigen is illustrated in the following example VIII:

## EXAMPLE VIII

Conducting the ICT Test for *L. pneumophila* Serogroup 1 O-Polysaccharide Antigen

**[0050]** In practice, the swab furnished with each device is dipped into the liquid sample, completely immersing the swab head. The use of the swab to act as a filter for undissolved solids, semisolids and colloids present in liquid biological samples such as urine, blood, lymph, etc. and also in liquid environmental samples is the subject of copending application Ser. No. 09/044,677 of Norman Moore and Vincent Sy filed Mar. 19, 1998. The swab is inserted into the hole at the bottom of the device (hole B of FIG. 1A) and gently pushed upward so that the swab tip is visible in the top hole (hole A of FIG. 1A). The Reagent A vial is held vertically above hole B and two drops of Reagent A are slowly added. The adhesive liner is then immediately peeled from the right edge of the device and the device is closed and securely sealed, thus pressing the swab in the swab well against the

gold conjugate pad. After 15 minutes, the result can be read in the window of the device. A negative sample—i.e., one containing no *L. pneumophila* serogroup 1 O-polysaccharide antigen—will exhibit only the control line in the top half of the window. A positive sample containing the target antigen will show two lines, the lower one of which is the patient (or sample) line; even a faint sample line indicates the presence of the target antigen in the sample. If no line appears in the window after 15 minutes, or only a sample line appears in the lower part of the window, the test is invalid and must be repeated.

**[0051]** Using the procedure described above, the devices prepared as described in Example VII were tested in the ICT procedure just described against 300 patient urine samples, 100 of which had been previously diagnosed as having *L. pneumophila* serogroup 1 infection.

**[0052]** These ICT tests according to this invention were conducted under circumstances such that the previous diagnoses were unknown to personnel performing the ICT tests. Overall, 98% of the ICT tests agreed with the previous positive diagnoses. Also overall, 98% of the urine samples previously diagnosed as negative for *L. pneumophila* serogroup 1 O-polysaccharide antigen gave results in agreement therewith when tested by the ICT procedure described herein, using the ICT device described in Example VII.

#### EXAMPLE IX

##### Use of the ICT to Test Environmental Samples

**[0053]** Applicability of this same test to environmental samples suspected of containing *L. pneumophila* serogroup 1 was also investigated as follows:

**[0054]** Water was seeded with *L. pneumophila* serogroup 1 bacteria obtained from a commercial source. The mixture was concentrated by filtering through a 0.22  $\mu\text{m}$  filter. A swab dipped in the sample was applied to the device, the device was closed and the assay was allowed to proceed. A positive result was observed within less than 15 minutes.

#### EXAMPLE X

##### Western Blot Immunoassay For Detection of Cross-Reactive Carbohydrate Antigens of *L. pneumophila* Serogroups 1, 2, 4 and 5

**[0055]** In order to perform the Western Blot immunoassay using a kit purchased from Bio-Rad Laboratories, *L. pneumophila* serogroup 5 cells were cultured as in Example II. A suspension of these cells was solubilized with 1% sodium dodecylsulfate in the presence of 10 mM mercaptoethanol at 100° C. for 5 minutes. The solubilized cells were treated with protease K and then subjected to electrophoretic separation of protein according to standard procedures provided by Bio-Rad.

**[0056]** The carbohydrate antigen from *L. pneumophila* serogroup 5 was conjugated to the spacer molecule described in Example III hereof in the manner described in Example IV and applied to an activated Sepharose column as described in Example V. This column was then used for the affinity purification of polyvalent rabbit antibodies specific to the carbohydrate antigen of *L. pneumophila* serogroup 5 (which were conventionally obtained from serum of a rabbit previously injected with the protein-containing of *L. pneumophila* serogroup 5) using the procedure of Example VI.

**[0057]** The Western immunoblot analysis was performed using a reagent kit from Bio-Rad and according to directions from this manufacturer. Briefly, the PBS extract of cells of *L. pneumophila* antigens 1, 2, 4 and 5 was subjected to the SDS-PAGE in 12% polyacrylamide gel blocked with 1% BSA with PBS transferred onto a nitrocellulose membrane. After this step, the membrane was incubated with affinity purified antibodies specific to carbohydrate of *L. pneumophila* serogroup 5. The membrane, washed as recommended by the manufacturer, and incubated with horseradish peroxidase conjugated to goat-anti-rabbit antibodies provided by Bio-Rad. After washing, the membrane was developed with a substrate system of 0.022 M 4-chloro-1 naphthol and 0.0012 M N; N-dimethyl-p-phenylene-diamine monohydrochloride in 0.1 M sodium citrate buffer of pH 6.9 containing 2.9 mM of hydrogen peroxide. FIG. 2 hereof shows the Western blot assay results compared with that of the prestained SDS-PAGE standard (in Lane 5) for the affinity purified antibodies of serogroup 5 of *L. pneumophila* against PBS extracts containing antigens of *L. pneumophila* as follows:

**[0058]** Lanes 1 and 7—*L. pneumophila* serogroup 2 (strain Togus-1)

**[0059]** Lanes 2 and 8—*L. pneumophila* serogroup 4 (strain Los Angeles-1)

**[0060]** Lanes 3 and 6—*L. pneumophila* serogroup 1 (strain Philadelphia-1)

**[0061]** Lanes 4 and 9—*L. pneumophila* serogroup 5 (strain U8W).

**[0062]** It is pointed out that the affinity purified antibodies for Lanes 1-4 were affinity purified on a column to which carbohydrate antigen from *L. pneumophila* serogroup 5 (strain U8W) was attached while those for Lanes 6-9 were affinity purified in the same manner on a column having attached carbohydrate antigen of *L. pneumophila* serogroup 5 (strain Dallas IE).

**[0063]** FIG. 2 clearly demonstrates that affinity purified antibodies as herein disclosed of *L. pneumophila* serogroup 5 react with antigens of *L. pneumophila* serogroups 1, 2, and 4 in addition to those of serogroup 5.

**[0064]** An ICT assay as described above in which affinity purified antibodies from *L. pneumophila* serogroup 5 are substituted for affinity purified antibodies from *L. pneumophila* serogroup 1 is contemplated.

**[0065]** Those skilled in the art of immunochemistry generally, and especially those skilled in immunoassays, will recognize that other materials and ingredients and at times, other procedural steps, can readily be substituted for those specifically recommended herein. A vast array of literature, both patent and non-patent, discusses the design and use of reliable, one-time-use, disposable immunoassay test devices that could be substituted for the preferred ICT device described and recommended herein. It is not intended that the present invention should be limited with respect to substitutable assay devices, materials, ingredients or process steps except insofar as the following claims may so limit it.

1-14. (canceled)

15. A method for determining the presence of *Legionella pneumophila* serogroup 1 in a liquid sample, the method comprising:

- (a) applying the liquid sample to a porous test strip, wherein the porous test strip defines a liquid flow path, the porous test strip comprising:
  - (i) a sample receiving zone for receiving the liquid sample;

- (ii) multiple binding agents disposed in a dried state at a binding zone along the flow path of the test strip downstream from the sample receiving zone, each binding agent comprising:  
antigen specific polyclonal antibodies conjugated to a detectable particle, wherein the antigen specific polyclonal antibodies specifically bind a *Legionella pneumophila* serogroup 1 O-polysaccharide antigen,
- (iii) multiple capture agents disposed in a dried state at a capture zone; and
- (iv) said capture zone downstream from the sample receiving zone; and
- (b) allowing said sample to flow laterally along the flow path of the test strip, mobilizing the binding agents;
- (c) binding, in the presence of *Legionella pneumophila* serogroup 1 O-polysaccharide antigen, a complex comprising one of the binding agents bound to the antigen in the capture zone, wherein the formation of a detectable line is indicative of the presence of *Legionella pneumophila* serogroup 1 in the liquid sample.
- 16.** The method of claim 15, wherein the capture agents are antibodies.
- 17.** The method of claim 16, wherein the capture agents are polyclonal antibodies.
- 18.** The method of claim 15, wherein the detectable particle is a color-producing particle.
- 19.** The method of claim 18, wherein the detectable particle is colloidal gold.
- 20.** The method of claim 15, wherein the detectable particle is a color-producing particle and the capture agents are antibodies.
- 21.** The method of claim 15, wherein the presence of *Legionella pneumophila* serogroup 1 in a liquid sample is determined within fifteen minutes of applying the liquid sample to the device.
- 22.** The method of claim 15, wherein the liquid sample is urine, serum, sputum or cerebrospinal fluid.
- 23.** The method of claim 15, wherein the antigen specific polyclonal antibodies are affinity-purified.
- 24.** The method of claim 23, wherein the antigen-specific polyclonal antibodies are affinity-purified with isolated and purified *Legionella pneumophila* serogroup 1 O-polysaccharide antigen.
- 25.** The method of claim 15, wherein the binding agents are essentially free of antibodies that bind proteins of *Legionella pneumophila* serogroup 1 O-polysaccharide antigen.
- 26.** A method for determining the presence of *Legionella pneumophila* serogroup 1 in a liquid sample, the method comprising:
- (a) contacting the liquid sample with a binding agent comprising an antigen specific antibody conjugated to a detectable particle, wherein the binding agent specifically binds to a *Legionella pneumophila* serogroup 1 O-polysaccharide antigen, in conditions under which the binding agent binds to the *Legionella pneumophila* serogroup 1 O-polysaccharide antigen if present in the sample, and
- (b) detecting the *Legionella pneumophila* serogroup 1 O-polysaccharide antigen, if present in the sample, bound by binding agent, wherein binding of the binding agent indicates that the *Legionella pneumophila* serogroup 1 O-polysaccharide antigen is present in the liquid sample.
- 27.** The method of claim 26, wherein binding of the binding agent to the *Legionella pneumophila* serogroup 1 O-polysaccharide antigen, if present in the sample, results in the formation of a detectable line that is indicative of the presence of *Legionella pneumophila* serogroup 1 in the liquid sample.
- 28.** The method of claim 26, wherein the liquid sample is applied to a porous test strip comprising (i) a sample receiving zone for receiving the liquid sample; and (ii) multiple capture agents disposed at a capture zone downstream from the sample receiving zone, wherein each of the capture agents bind, in the presence of *Legionella pneumophila* serogroup 1 O-polysaccharide antigen in the sample, a complex comprising one of the binding agents bound to the antigen.
- 29.** The method of claim 28, wherein the liquid sample is applied to a porous test strip that further comprises (iii) multiple binding agents disposed at a binding zone along the flow path of the test strip downstream from the sample receiving zone.
- 30.** The method of claim 26, wherein the antigen specific antibody is a polyclonal antibody.
- 31.** The method of claim 26, wherein the antigen specific antibody is affinity-purified.
- 32.** The method of claim 31, wherein the antigen specific antibody is affinity-purified with isolated and purified *Legionella pneumophila* serogroup 1 O-polysaccharide antigen.
- 33.** The method of claim 26, wherein the binding agents are essentially free of antibodies that bind proteins of *Legionella pneumophila* serogroup 1 O-polysaccharide antigen.
- 34.** The method of claim 26, wherein the liquid sample is urine, serum, sputum or cerebrospinal fluid.

\* \* \* \* \*

专利名称(译)	用纯化的抗原特异性抗体检测军团菌的方法		
公开(公告)号	<a href="#">US20160047807A1</a>	公开(公告)日	2016-02-18
申请号	US14/828823	申请日	2015-08-18
[标]申请(专利权)人(译)	了Alere SCARBOROUGH		
申请(专利权)人(译)	了Alere SCARBOROUGH INC.		
当前申请(专利权)人(译)	了Alere SCARBOROUGH INC.		
[标]发明人	KOULCHIN VLADIMIR A MOORE NORMAN J MOLOKOVA ELENA V		
发明人	KOULCHIN, VLADIMIR A. MOORE, NORMAN J. MOLOKOVA, ELENA V.		
IPC分类号	G01N33/531 G01N33/569 C07K1/22 C07K16/12 C12P19/04 C12R1/01 G01N21/78 G01N30/06 G01N30/88 G01N33/543		
CPC分类号	G01N33/56911 G01N21/78 G01N2400/10 G01N2021/7759 G01N2021/752 G01N2021/757 G01N2469 /10 C07K16/126 G01N33/53		
其他公开文献	US9989529		
外部链接	<a href="#">Espacenet</a> <a href="#">USPTO</a>		

摘要(译)

本发明涉及从军团菌中提取细菌，特别是嗜肺军团菌，一种基本上不含蛋白质的O-多糖或碳水化合物抗原，通过蛋白质空间分子将该抗原与活化的色谱柱偶联，该蛋白质空间分子首先与抗原偶联，利用由此制备的柱用于亲和纯化来自相同的军团菌细菌的原始多价抗体，其中O-多糖或碳水化合物抗原被分离 - 从而获得抗原特异性抗体，其可用于快速检测相应的军团菌人体液中的细菌或其抗体，例如尿液，痰液，血液等，或怀疑携带嗜肺军团菌细菌的环境样品中的细菌或其抗体。

